
- GB** Instructions for use/Technical description
CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm
- USA** Instructions for use/Technical description
CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm
- D** Gebrauchsanweisung/Technische Beschreibung
CRANIOFIX2-Titanklammern Ø 11 mm, Ø 16 mm, Ø 20 mm
- F** Mode d'emploi/Description technique
Clamps en titane CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- E** Instrucciones de manejo/Descripción técnica
Pinzas de titanio CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- I** Istruzioni per l'uso/Descrizione tecnica
Clamp in titanio CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- P** Instruções de utilização/Descrição técnica
Grampos de titânio CRANIOFIX2 de Ø 11 mm, Ø 16 mm, Ø 20 mm
- NL** Gebruiksaanwijzing/Technische beschrijving
CRANIOFIX2-titanklemmen Ø 11 mm, Ø 16 mm, Ø 20 mm
- S** Bruksanvisning/Teknisk beskrivning
CRANIOFIX2-titanklamrar Ø 11 mm, Ø 16 mm, Ø 20 mm
- RUS** Инструкция по применению/Техническое описание
Титановые клипсы CRANIOFIX2 Ø 11 мм, Ø 16 мм, Ø 20 мм
- CZ** Návod k použití/Technický popis
Titanové svorky CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- PL** Instrukcja użytkowania/Opis techniczny
Zaciski tytanowe CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- SK** Návod na používanie/Technický opis
Titánové svorky CRANIOFIX2 Ø 11 mm, Ø 16 mm, Ø 20 mm
- TR** Kullanım Kılavuzu/Teknik açıklama
CRANIOFIX2 titanyum klemensler Ø 11 mm, Ø 16 mm, Ø 20 mm
- KR** 사용 설명서 / 기술 설명
두개성형판 고정제 (CRANIOFIX, 형명별 개별기재)

B | BRAUN
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 011200 04/13 V6 Änd.-Nr. 46074

CE 0123 - DIR 93/42/EEC

Technical alterations reserved

Aesculap®

CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm

Legend

- Fig. 1
- 1 CRANIOFIX2 titanium clamps FF490T/FF490T-UNI, Ø 11 mm
 - CRANIOFIX2 titanium clamps FF491T/FF491T-UNI, Ø 16 mm
 - CRANIOFIX2 titanium clamps FF492T/FF492T-UNI, Ø 20 mm
- Fig. 2
- 1 CRANIOFIX2 titanium clamps FF491T/FF491T-UNI, Ø 16 mm
 - CRANIOFIX2 titanium clamps FF492T/FF492T-UNI, Ø 20 mm
 - 2 CRANIOFIX2 titanium clamps FF490T/FF490T-UNI, Ø 11 mm
- Fig. 3
- 1 Correct positioning of the CRANIOFIX2 titanium clamps FF491T/ FF491T-UNI/FF492T/FF492T-UNI: Trepanation hole Ø 12/15 mm
 - 2 Incorrect positioning of the CRANIOFIX2 titanium clamps FF491T/ FF491T-UNI/FF492T/FF492T-UNI: Craniotomy gap
- Fig. 4
- 1 CRANIOFIX2 titanium clamps
 - 2 New craniotomy line
 - 3 Old craniotomy line

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Product description

The CRANIOFIX2 titanium clamp system from Aesculap consists of sterile CRANIOFIX2 titanium clamps for single use and reusable neurosurgical instruments for the application of titanium clamps.

Intended use

The CRANIOFIX2 titanium clamp system from Aesculap is used for the fixation of craniotomized bone flaps and fractures to the neurocranium.

The CRANIOFIX2 titanium clamp FF490T/FF490T-UNI, Ø 11 mm, is inserted in the craniotomy gap (see Fig. 2) or in the fracture gap.

The CRANIOFIX2 titanium clamp FF491T/FF491T-UNI, Ø 16 mm, is inserted in a burr hole with Ø 12 mm (see Fig. 3) or is used to bridge defects in case of fractures.

The CRANIOFIX2 titanium clamp FF492T/FF492T-UNI, Ø 20 mm, is inserted in a burr hole with Ø 15 mm (see Fig. 3) or is used to bridge defects in case of fractures.

Materials

The materials used in the implant are listed on the packaging:

- Ti6Al4V acc. to ISO 5832-3

Indications

The CRANIOFIX2 titanium clamp system from Aesculap is used for the fixation of craniotomized bone flaps and fractures to the neurocranium.

Cautionary Notes

The following list shows some, but not all, conceivable contraindications.

- Systemic illnesses and metabolic disorders, hypersensitivity to metals or allergies to the implant materials, inflammations in the region of the implant site, insufficient compliance and co-operation of the patient
- Bone conditions that rule out the application of CRANIOFIX2 titanium clamps
- Artificial cranial bone flaps
- Bone tumors in the area supporting the implant
- Degenerative bone diseases
- Missing dura mater
- Application in the facial skull (viscerocranium) and in the orbita or skull-base region

Side effects and interactions

General perioperative risks include infections, damage to nervous tissue, hematomas and wound healing disturbances. The patient must be informed of these risks and other risks associated with neurosurgery, general surgery and the administration of anesthetics.

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or lack of asepsis.

- The user must follow the instructions that accompany each of the Aesculap implant components.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically removed components.
- Implants already used must not be reused.

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- Implants produce moderate MRI artifacts.

- Inappropriate fixation or any combination with other fixation methods (sutures, mini plates or clamp systems from other manufacturers) can result in the formation of steps, shifts of the calvarium and loosening and breakage of implant components. For maximum stability of a craniotomized bone flap, the use of three CRANIOFIX2 titanium clamps in a triangular formation are recommended. Larger bone flaps or a skull fracture with multiple fragments may require a larger number of CRANIOFIX2 titanium clamps to ensure sufficient stability.
- During application of CRANIOFIX2 titanium clamps, friction between the upper titanium disc and the pin indentation can lead to spark formation. Appropriate safety precautions must be taken to prevent this. Due to the risk of fire, surgical interventions with this product must not be carried out in rooms where flammable anesthetics or other flammable gases, liquids, objects or oxidizing agents are stored.
- After the operation, exceptionally heavy mechanical strains on the fixated bone flap must be avoided.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants.

Application

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information about product training, please contact your national B. Braun/ Aesculap agency.

- The operating surgeon shall devise an operation plan that specifies and accurately documents the following:
- Selection of the implant components and their dimensions
 - Positioning of the implant components in the bone
 - Location of intraoperative landmarks
- The following conditions must be fulfilled prior to application:
- All requisite implant components are ready to hand
 - Operating conditions are highly aseptic
 - The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
 - The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
 - The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
 - The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

Preparing and handling the implants

Damage to/destruction of the titanium clamps due to use an incompatible applicator!

- Use CRANIOFIX2 applicator (FF494R) only for CRANIOFIX2 titanium clamps (FF490T/FF490T-UNI–FF492T/FF492T-UNI).
- Use CRANIOFIX applicator (FF102R, FF106R and FF107R) only for CRANIOFIX titanium clamps (FF099T–FF101T).

Prior to the surgical intervention, the surgeon must prepare an operation plan to ensure that there is a sufficient supply of sterile CRANIOFIX2 titanium clamps available, that the operation will be performed in aseptic conditions, that the required CRANIOFIX2 applicator and all other CRANIOFIX instruments are complete and ready for use, and that the surgeon and operating team are familiar with the instruments and the surgical technique. Attention must be paid to an osteoplastic incision for the craniotomy.

The correct handling of the CRANIOFIX2 titanium clamps before and during the operation is an important prerequisite for the success of the bone fixation.

- The CRANIOFIX2 titanium clamps are individually packed in protective packaging and sterilized by gamma radiation.
- The CRANIOFIX2 titanium clamps must be kept in their original packaging until immediately before use. Prior to use, check the expiry date and make sure that the sterile packaging is intact.
 - Check to make certain that the surfaces of the CRANIOFIX2 titanium clamps are not bent or damaged.
 - CRANIOFIX2 titanium clamps are intended for single use only and must not be reused under any circumstances.
 - Damaged or surgically removed CRANIOFIX2 titanium clamps must not be reused.
 - To avoid spark formation, immerse the CRANIOFIX2 titanium clamps in a sterile liquid prior to application.

Application of the CRANIOFIX2 titanium clamps

- Caution
- The CRANIOFIX2 titanium clamps Ø 16 mm from Aesculap may only be applied in the burr hole Ø 12 mm (see Fig. 3) or used to bridge defects in case of fractures.
- The CRANIOFIX2 titanium clamps Ø 20 mm from Aesculap may only be applied in the burr hole Ø 15 mm (see Fig. 3) or used to bridge defects in case of fractures.
- The operating surgeon is responsible for the proper application of the CRANIOFIX2 titanium clamps. For maximum stability of a bone flap approx. 20 cm² in size, it is recommended to apply three CRANIOFIX2 titanium clamps at equal distance from each other in the craniotomy gap or burr hole. Figure 2 shows a typical application of a craniotomized bone flap.
- Brief description of operative procedure:
- Position the three CRANIOFIX2 titanium clamps at an equal distance from each other on the bone flap.
 - Slide the lower disc between the dura mater and skull cap.
 - Insert the bone flap in its original position.
 - Grasp the first CRANIOFIX2 titanium clamp with the holding forceps. To ensure a secure fixation and to prevent damage to the pin indentation, the pin of the CRANIOFIX2 titanium clamp must be held in the recess in the jaw of the holding forceps intended for this purpose.
 - Position the CRANIOFIX2 applicator over the pin of the CRANIOFIX2 titanium clamp.

- Compress the handle of the **CRANIOFIX2** applicator to move the upper disc along the pin indentation, closer to the lower disc. The spring leaves will audibly engage in the pin indentation of the **CRANIOFIX2** titanium clamp.
- Only remove the holding forceps immediately before closing the **CRANIOFIX2** titanium clamp.
- The **CRANIOFIX2** applicator has an automatic strain relief mechanism, which ensures that the **CRANIOFIX2** applicator disengages when the maximum allowable application force is reached. After the **CRANIOFIX2** applicator has disengaged, the handle should be pressed down completely several times until no further audible engaging of the spring reeds in the pin indentation of the **CRANIOFIX2** titanium clamp can be heard.

The upper disc can loosen after application if the pin is not cut off properly!

- After tightening the upper disc, cut the excess pin off at a right angle to the pin's axis using the **CRANIOFIX** cutting forceps.
- Be sure that the indentation above the upper disc remains completely intact.

The same procedure is used for the other **CRANIOFIX2** titanium clamps.

- Remove the remaining excess part of the pin with the cutting forceps. When doing so, the excess pin of the **CRANIOFIX2** titanium clamp must be held onto, to prevent it springing into the operating site in an uncontrolled manner.

Removing the CRANIOFIX2 titanium clamps in a repeat intervention

In rare cases in which a repeat intervention is necessary, the surgeon can leave the **CRANIOFIX2** titanium clamps in position on the bone flap. Aesculap recommends that the surgeon makes a second incision around the original clamp positions. If this is not possible (e.g. for anatomical reasons), Aesculap recommends a second incision within the original clamp positions. The incision must in no way be made through a **CRANIOFIX2** titanium clamp. Decisions regarding the approach for the repeat intervention fall under the responsibility of the surgeon (see Fig. 4). The operative procedure for the application of **CRANIOFIX2** titanium clamps is described in the **CRANIOFIX2** brochure or in the list of literature (www.aesculap.de). Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact us directly:

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-0
Fax: +49 7461 95-2600
Website: www.aesculap.de

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:
Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Patient ID

The Patient ID enables a written record of important information concerning the implanted **CRANIOFIX2** titanium clamps and the surgical intervention to be kept:

- Patient data
- Hospital data
- Date of provision
- Attending physician and hospital
- Art. no. and batch no. of the **CRANIOFIX2** titanium clamps implanted

To facilitate post-operative radiological examinations, the patient ID should be given to every patient. It can be requested from Aesculap separately. Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact Aesculap directly.

Accessories

Art. no.	Designation
FF490T	CRANIOFIX2 titanium clamps Ø 11 mm
FF490T-UNI	CRANIOFIX2 titanium clamps Ø 11 mm
FF491T	CRANIOFIX2 titanium clamps Ø 16 mm
FF491T-UNI	CRANIOFIX2 titanium clamps Ø 16 mm
FF492T	CRANIOFIX2 titanium clamps Ø 20 mm
FF492T-UNI	CRANIOFIX2 titanium clamps Ø 20 mm
FF494R*	CRANIOFIX2 applicator, cannot be disassembled
FF105R*	CRANIOFIX holding forceps
FF103R*	CRANIOFIX cutting forceps
FF104R*	CRANIOFIX removal forceps
FF094P	CRANIOFIX storage tray

* Regarding the operation, function and reprocessing of this product, follow instructions for use TA011161.

Aesculap®

CRANIOFIX2 titanium clamps Ø 11 mm, Ø 16 mm, Ø 20 mm

Legend

- Fig. 1
- 1 **CRANIOFIX2** titanium clamps FF490T/FF490T-UNI, Ø 11 mm
 - CRANIOFIX2** titanium clamps FF491T/FF491T-UNI, Ø 16 mm
 - CRANIOFIX2** titanium clamps FF492T/FF492T-UNI, Ø 20 mm
- Fig. 2
- 1 **CRANIOFIX2** titanium clamps FF491T/FF491T-UNI, Ø 16 mm
 - CRANIOFIX2** titanium clamps FF492T/FF492T-UNI, Ø 20 mm
- Fig. 3
- 1 Correct position of **CRANIOFIX2** titanium clamps FF491T/ FF491T-UNI/FF492T/FF492T-UNI: Trepanation hole Ø 12/15 mm
 - 2 Incorrect position of **CRANIOFIX2** titanium clamps FF491T/ FF491T-UNI/FF492T/FF492T-UNI: Craniotomy gap
- Fig. 4
- 1 **CRANIOFIX2** titanium clamps
 - 2 New craniotomy line
 - 3 Old craniotomy line

Symbols on product and packages

	Sterilization using irradiation
	This product is not for reuse in terms of its proper applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Product description

The Aesculap **CRANIOFIX2** titanium clamp system comprises sterile **CRANIOFIX2** titanium clamps for single use and diverse reusable neurosurgical instruments for applying titanium clamps.

Intended use

The Aesculap **CRANIOFIX2** titanium clamp system is used for fixating craniotomized bone flaps and fractures at the brain skull (neurocranium).

The **CRANIOFIX2** titanium clamp FF490T/FF490T-UNI, Ø 11 mm, is applied in the craniotomy gap (see Fig. 2) or in the fracture gap.

The **CRANIOFIX2** titanium clamp FF491T/FF491T-UNI, Ø 16 mm, is applied in a trepanation hole of Ø 12 mm (see Fig. 3) or for bridging the defect in case of a fracture.

The **CRANIOFIX2** titanium clamp FF492T/FF492T-UNI, Ø 20 mm, is applied in a trepanation hole of Ø 15 mm (see Fig. 3) or for bridging the defect in case of a fracture.

Materials

The **CRANIOFIX2** titanium clamps are made of titanium alloy (Ti6Al4V) acc. to ISO 5832-3. This titanium alloy is a biocompatible material.

The **CRANIOFIX2** titanium clamps are non-ferromagnetic and MR-Conditional up to 3.0 Tesla. Consequently, the magnetic fields applied during NMR tomography do not present any risk.

Stronger magnetic fields or an enlargement of the imaging field can lead to a dramatically increased incidence of positioning errors and artifacts.

The size of the artifacts can vary considerably, depending on the NMR-pulse frequency, and compromise the significance of the NMR image if the area concerned is in the immediate vicinity (within a distance of some millimeters) of the clamp position.

Contraindications

- The following list shows some, but not all conceivable contraindications.
- Patients with systemic illnesses and metabolic disorders, oversensitivity to metals or allergies against the implant materials, inflammations in the region of the implant site, and patients that show insufficient compliance and co-operation
 - Bone conditions that rule out the application of **CRANIOFIX2** titanium clamps
 - Artificial cranial bone flaps
 - Bone tumors in the area on which the implant is supported
 - Degenerative bone diseases
 - Missing dura mater
 - Application in the facial skull (viscerocranium) and in the orbita or skull-base regions

Notes

- The following list shows some, but not all, conceivable contraindications.
- Improper fixation or combination with other fixation methods (sutures, mini plates or clamp systems from other manufacturers) can lead to step formation, position shifts of the skull calotte and loosening or breakage of implant components. For maximum stability of a craniotomized bone flap, we recommend using three **CRANIOFIX2** titanium clamps in a triangular arrangement. For larger bone flaps or multiple fragments of a skull fracture, additional **CRANIOFIX2** titanium clamps may be required to ensure adequate stability.
 - During application of **CRANIOFIX2** titanium clamps, friction between the upper titanium disc and the pin indentation can lead to spark formation. Appropriate safety measures must be taken in order to prevent this. Because of the fire hazard, surgical interventions of this kind must not be carried out in rooms where flammable anesthetics or other flammable gases, liquids, objects or oxidizing agents are stored.
 - After the operation, exceptionally high mechanical loads on the fixated bone flap must be avoided.

Safety measures

The operating surgeon should be thoroughly familiar with the surgical techniques required for using this product. The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. The operating surgeon is also responsible for complications resulting from misdiagnosis and flawed operation techniques.

Possible side effects

Infections, damage to nervous tissue, hematomas and wound dehiscence are among the general perioperative risks. The patient should be made aware of these risks as well as other risks of neurosurgery, general surgery and the administration of anesthetics.

MRI Information

MR-Conditional.

Non-clinical testing demonstrated that the **CRANIOFIX2** System is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 3,000 Gauss/cm or less
- First Level Controlled Operating Mode, with a maximum whole-body averaged SAR of 4-W/kg

MRI-Related Heating

In non-clinical testing, the **CRANIOFIX2** System produced the following temperature rises during MRI performed for 15 min of scanning (i.e., per pulse sequence) in 1.5 Tesla/64 MHz (Magnetom, Siemens Medical Solutions, Malvern, PA. Software Numaris/4, Version Syngo MR 2002B DHHS Active-shielded, horizontal field scanner) and 3 Tesla (3 Tesla/128 MHz, Excite, HDx, Software 14X, M5, General Electric Healthcare, Milwaukee, WI) MR Systems:

	1.5 Tesla	3 Tesla
First Level Controlled Operating Mode, scaled to whole body averaged SAR of 4 W/kg, highest temperature change	+2.6 °C	+2.8 °C

These temperature changes will not pose a hazard to a human subject under the conditions indicated above.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the **CRANIOFIX2** System. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary. The maximum artifact size (i.e., as seen on the gradient echo pulse sequence) extends approximately 5 mm relative to the size and shape of this implant.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	503 mm ²	196 mm ²	701 mm ²	563 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Delivery format

Aesculap **CRANIOFIX2** titanium clamps are delivered in sterile condition and are designed for single use only. The **CRANIOFIX2** titanium clamps must not be resterilized.

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product.

Provision and handling of the implants

CAUTION

Federal law restricts this device to sale by or on order of a physician

Damage to/destruction of the titanium clamps caused by using a non-compatible applier!

- Use the **CRANIOFIX2** applier (FF494R) only for **CRANIOFIX2**-titanium clamps (FF490T/FF490T-UNI-FF492T/FF492T-UNI).
- Use the **CRANIOFIX** applier (FF102R, FF106R and FF107R) only for the **CRANIOFIX** titanium clamps (FF099T-FF101T).

Prior to the surgical intervention, the operating surgeon must prepare an operation plan in order to ensure that a sufficient number of sterile **CRANIOFIX2** titanium clamps are available, aseptic operation conditions are in place, the **CRANIOFIX2** applier required and other **CRANIOFIX** instruments are complete and ready for operation, and the operating surgeon and the O.R. team are familiar with the instruments and operation technique. The craniotomy incisions must be performed in an osteoplastic fashion.

Correct handling of the **CRANIOFIX2** titanium clamps before and during the operation is essential for successful bone flap fixation.

- The **CRANIOFIX2** titanium clamps are individually wrapped in protective packaging and gamma-ray sterilized.
- The **CRANIOFIX2** titanium clamps have to be kept in their original packaging until immediately before use. Prior to use, check the expiry date and make sure that the sterile packaging is intact.
 - Make certain that the surfaces of the **CRANIOFIX2** titanium clamps are not bent or damaged.
 - The **CRANIOFIX2** titanium clamps have been designed as a disposable product and must not be reused.
 - Damaged or surgically excised **CRANIOFIX2** titanium clamps must not be used again.
 - To avoid spark formation, immerse the **CRANIOFIX2** titanium clamps in a sterile liquid prior to application.

Application of CRANIOFIX2 titanium clamps

Caution

The Aesculap **CRANIOFIX2** titanium clamp of Ø 16 mm may only be applied in a trepanation hole of Ø 12 mm (see Fig. 3) or for bridging the defect in case of a fracture.

The Aesculap **CRANIOFIX2** titanium clamp of Ø 20 mm may only be applied in a trepanation hole of Ø 15 mm (see Fig. 3) or for bridging the defect in case of a fracture.

The operating surgeon is responsible for the proper application of the **CRANIOFIX2** titanium clamps. For maximum stability of one bone flap approx. 20 cm² in size, we recommend applying three **CRANIOFIX2** titanium clamps in the craniotomy gap or the trepanation hole, respectively, at equal distances to each other. Figure 2 shows a typical application of a craniotomized bone flap.

The following is a brief description of the surgical procedure:

- 1 Position the three **CRANIOFIX2** titanium clamps in equal distances to each other at the bone flap.
- 2 The lower discs are pushed in between the dura and the calvarium.
- 3 Replace the bone flap in its original position.
- 4 Apply the holding forceps on the first **CRANIOFIX2** titanium clamp. To ensure a safe fixation and avoid damage to the pin indentation, the pin of the **CRANIOFIX2** titanium clamp must be held in the recess intended for this purpose.

- 5 Position the **CRANIOFIX2** applier over the pin of the **CRANIOFIX2** titanium clamps.
- 6 Compress the handle of the **CRANIOFIX2** applier so that the upper disk approaches the lower disc along the pin indentation.
The spring reeds engage in the pin indentation of the **CRANIOFIX2** titanium clamp, with an audible "click".
- 7 Do not remove the holding forceps until immediately before closing the **CRANIOFIX2** titanium clamp.
- 8 The **CRANIOFIX2** applier is equipped with an automatic strain relief mechanism, which ensures that the **CRANIOFIX2** applier is disengaged as soon as the maximum allowable force is applied. Once the **CRANIOFIX2** applier has disengaged, the handle should be fully compressed several times until there is no more audible clicking of the spring reeds as they engage in the indentation of the pin of the **CRANIOFIX2** titanium clamp.

CAUTION

The top disk can loosen after application if the pin is not cut off properly!

- After tightening the top disk, cut the excess pin off at a right angle to the pin's axis using the **CRANIOFIX** cutting forceps.
- Be sure that the indentation above the top disk remains completely intact.

- 9 Use the cutting forceps to remove the remaining, outward projecting section of the pin. When doing so, the protruding pin of the **CRANIOFIX2** titanium clamp must be held securely to prevent the pin from snapping off into the operation field in an uncontrolled manner.
- 10 The same procedure is carried out with the other **CRANIOFIX2** titanium clamps.
- The surgical procedure for applying **CRANIOFIX2** titanium clamps is described in the **CRANIOFIX2** brochure and in the literature list (www.aesculap.de). Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact us directly:

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-0
Fax: +49 7461 95-2600
Website: www.aesculap.de

Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Patient ID

The Patient ID serves to keep a written record of important information on the implanted **CRANIOFIX2** titanium clamps and the surgical intervention that was carried out:

- Patient data
- Hospital data
- Date of provision
- Physician in charge and hospital
- Art. no. and lot no. of the **CRANIOFIX2** titanium clamps implanted

To facilitate the postoperative radiological examination, each patient should be issued with a Patient ID, which can be ordered separately from Aesculap. Please address any inquiries in this respect to your Aesculap Medical Products Consultant, or contact Aesculap directly.

Removing the **CRANIOFIX2** titanium clamps in a revision operation

In the rare case that the operation has to be repeated, the operating surgeon can leave the **CRANIOFIX2** titanium clamps in place on the bone flap. Aesculap recommends applying a second incision around the original clamp position. If this is not feasible (e.g. for anatomic reasons), Aesculap recommends a second incision within the original clamp site. Under no circumstances must the incision cut through a **CRANIOFIX2** titanium clamp. The operating surgeon is responsible for any decision about how to approach the revision intervention (see Fig. 4).

Accessories

Art. no.	Designation
FF490T	CRANIOFIX2 titanium clamps Ø 11 mm
FF490T-UNI	CRANIOFIX2 titanium clamps Ø 11 mm
FF491T	CRANIOFIX2 titanium clamps Ø 16 mm
FF491T-UNI	CRANIOFIX2 titanium clamps Ø 16 mm
FF492T	CRANIOFIX2 titanium clamps Ø 20 mm
FF492T-UNI	CRANIOFIX2 titanium clamps Ø 20 mm
FF494R*	CRANIOFIX2 applicator, cannot be disassembled
FF105R*	CRANIOFIX holding forceps
FF103R*	CRANIOFIX cutting forceps
FF104R*	CRANIOFIX removal forceps
FF094P	CRANIOFIX storage tray

* Concerning the handling, functioning and reprocessing of this product, the instructions for use TA011161 must be observed.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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