Aesculap® activ C

Aesculap Spine

- GB

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
Technical alterations reserved

Aesculap® activ C

Distraction set with lockable distraction screws

Legend

- 1 activ C distractor for lockable distraction screws, right
- 2 activ C distractor for lockable distraction screws, left
- 3 activ C screwdriver for lockable distraction screws
- 4 activ C lockable offset distraction screw
- 5 Hexagon
- 6 Shaft
- 7 activ C lockable distraction screw
- 8 Thumb screw
- 9 Ratchet lever
- 10 Toothed rack
- 11 Guide
- 12 Distraction arm (at guide)
- 13 Buttons (for locking and release)
- 14 Distraction arm (at toothed rack)

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

The activ C distractor for lockable distraction screws is used in combination with the corresponding activ C lockable distraction screws for the parallel distraction of adjacent vertebral bodies in the cranial section for the purpose of extending the height of the intervertebral disk space.

Available sizes

Art. no.	Designation
FW848R	activ C distractor for lockable distraction screws, right
FW849R	activ C distractor for lockable distraction screws, left
FW850R	activ C distraction screwdriver for lockable distraction screws
FW851SU	activ C lockable offset distraction screw, 12 mm
FW852SU	activ C lockable offset distraction screw, 14 mm
FW853SU	activ C lockable offset distraction screw, 16 mm
FW854SU	activ C lockable offset distraction screw, 18 mm
FW855SU	activ C lockable distraction screw, 16 mm
FW856SU	activ C lockable distraction screw, 18 mm
FW861SU	activ C lockable distraction screw, 12 mm
FW862SU	activ C lockable distraction screw, 14 mm

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Observe the instructions for use of the activ C intervertebral disk prosthesis, TA011995.
- Follow the operating manual 031302.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).
- Do not use products that are in open or damaged sterile packaging.

DANGER

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!
► Do not reprocess the product.

The product is gamma-sterilized and supplied in sterile packaging.
The product must not be reused.

- Do not use the product after its use-by date.

Safe operation

In the text below, short designations are used for the following components:

- activ C distractor for lockable distraction screws is referred to as "distractor"
- activ C distraction screwdriver for lockable distraction screws is referred to as "distraction screwdriver"
- activ C lockable distraction screws are referred to as "distraction screws"
- activ C lockable offset distraction screws are referred to as "offset distraction screws"

Note

Any application of the distraction screw includes the application of the offset distraction screw, even if the latter is not mentioned explicitly!

Applying the distraction screw, see Fig. A

WARNING

Compression of or damage to the spinal canal due to selection of a distraction screw of excessive length!
► Use the radiographic image to determine the correct length of the distraction screw according to the vertebral body.
► Always maintain x-ray control when applying the distraction screw.

Note

For very small vertebral bodies, offset distraction screws must be used to ensure that there is sufficient space for the implantation instruments!

- Choose distraction screw 7 or offset distraction screw 4.
- Insert distraction screw 7 or offset distraction screw 4 in distraction screwdriver 3.
- Place distraction screwdriver 3 on hexagon 5 of distraction screw 7 and make sure it engages.
- Check the self-holding function of distraction screw 7 in distraction screwdriver 3.
 - or -
- Place distraction screwdriver 3 on hexagon 5 of offset distraction screw 4 and make sure it engages.

Inserting the distraction screw into the vertebral body under X-ray control, see Fig. B

Note

For very small vertebral bodies, offset distraction screws must be used to ensure that there is sufficient space for the implantation instruments!

- Position and screw in distraction screw 7/offset distraction screw 4 at the center of the vertebral body.
- Remove distraction screwdriver 3 from distraction screw 7/offset distraction screw 4.
- Apply second distraction screw 7/offset distraction screw 4, ensuring that the distraction screw 7/offset distraction screw 4 is positioned in parallel.

Mounting the distractor and carrying out the distraction, see Fig. C/D

WARNING

Risk of loosening/tear-out of the distraction screws when distracting with the distractor and the distraction screws!
► Use the distraction forceps for the distraction process.
► Do not use the distractor and the distraction screws for the actual distraction process.
► Lock the distractor in the extended condition.

WARNING

Irritation of the facet joints/strain on spinal cord and nerve roots due to excessive distraction!
► Extend the distraction forceps slowly and carefully under X-ray control.
► Observe the tactile response.
► Use the heights of adjacent (healthy) disks as a guide for the required degree of distraction.

- Select distractor right 1 or distractor left 2.
- Turn thumb screw 8 to transfer the distance of parallel distraction screws 4/7 to distraction arms 12/14. The distance between distraction arms 12/14 corresponds to the distance between distraction screws 4/7.
- Mount distractor right 1 or distractor left 2 on parallel distraction screws 4/7 so that they click into position.
- Distract with the distraction forceps.
- While distracting with the distraction forceps, extend distractor right 1 or distractor left 2 by turning thumb screw 8.
- Once the distraction process carried out with the distraction forceps is complete, lock distractor right 1 or distractor left 2 and hold the vertebrae in the distracted position.

Removing the distractor, see Fig. E

- To release the distraction after the implantation in the intervertebral disk space, press ratchet lever 9 at the distractor.
Distraction arms 12/14 move toward each other.
The distractor is fully released from any load.
- Press buttons 13 at the side of distraction arm 12/14 and slide off the distractor from distraction screws 4/7.
- Remove distraction screws 4/7 from the vertebrae, using distraction screwdriver 3.

Disassembling

- Pull distraction arm 12 (labeled "RIGHT/LEFT") along toothed rack 10 until toothed rack 10 is not seated in guide 11 anymore.

Assembling

Note

Combine distraction arm 12 labeled "RIGHT" only with distraction arm 14 labeled "FW848R"!
Combine distraction arm 12 labeled "LEFT" only with distraction arm 14 labeled "FW849R"!

- Push the ratchet lever 9.
- Mount distraction arm 12 (labeled "RIGHT/LEFT") with guide 11 on the toothed rack, ensuring the correct orientation of the distraction arms 12/14.

Validated reprocessing procedure

General safety instructions

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this. The recommended chemistry was used for validation.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Single-use products

WARNING

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

► **Do not reprocess the product!**

The following products must not be reprocessed:

Art. no.	Designation
FW851SU	activ C lockable offset distraction screw, 12 mm
FW852SU	activ C lockable offset distraction screw, 14 mm
FW853SU	activ C lockable offset distraction screw, 16 mm
FW854SU	activ C lockable offset distraction screw, 18 mm
FW855SU	activ C lockable distraction screw, 16 mm
FW856SU	activ C lockable distraction screw, 18 mm
FW861SU	activ C lockable distraction screw, 12 mm
FW862SU	activ C lockable distraction screw, 14 mm

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary. Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.

► Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

► Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for the material in question (e.g., aluminum, plastics, high-grade steel),
 - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum permitted cleaning temperature of 60 °C.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfection.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfection.
 - as an integrated mechanical support measure for mechanical cleaning/disinfection.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfection.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FW848R-FW850R	■ Brush length: 30 mm/Ø: 5.5 mm, e.g., TA006874 ■ 20 ml disposable syringe ■ When cleaning instruments with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with immersion disinfection
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW848R-FW850R	■ Brush length: 30 mm/Ø: 5.5 mm, e.g., TA006874 ■ 20 ml disposable syringe ■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Place instruments in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfection with manual pre-cleaning

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

- Phase I
- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
 - Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
 - If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
 - Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
 - Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

- Phase II
- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
 - Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the instrument is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in sterile Aesculap containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

The product can be sterilized either in disassembled or in assembled condition.

To avoid breakage due to stress crack corrosion, sterilize the instruments with the lock fully open or locked no further than on the first ratchet tooth.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Disassemble the instrument
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Or in the US:
Aesculap Implant Systems LLC
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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