Aesculap Spine activ C

(GB)

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


A

B

Distraction set with lockable distraction screws

Legend

- 1 activ C distractor for lockable distraction screws, right
- 2 activ C distractor for lockable distraction screws, left
- 3 activ C screwdriver for lockable distraction screws
- 4 activ C lockable offset distraction screw
- 5 Hexagon
- 6 Sheath
- 7 activ C lockable distraction screw
- 8 Flap screw
- 9 Ratchet lever
- 10 Toothed rack
- 11 Guide
- 12 Distraction arm (at guide)
- 13 Buttons (for locking and release)
- 14 Distraction arm (at toothed rack)

Figure captions

- A Illustration on fold-out page
- B Illustration on fold-out page
- C Illustration on fold-out page

Symbols on product

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Attention, see instructions for use

Intended use

The activ C distractor for lockable distraction screws is used in combination with the corresponding activ C lockable distraction screws for parallel distraction of adjacent vertebral bodies in the cranial section for the purpose of extending the height of the intervertebral space.

Available sizes

Art.-no.	Designation
FW848R	activ C distractor for lockable distraction screws, right
FW849R	activ C distractor for lockable distraction screws, left
FW850R	activ C distraction screwdriver for lockable distraction screws
FW851SU	activ C lockable offset distraction screw, 12 mm
FW852SU	activ C lockable offset distraction screw, 14 mm
FW853SU	activ C lockable offset distraction screw, 16 mm
FW854SU	activ C lockable offset distraction screw, 18 mm
FW855SU	activ C lockable distraction screw, 16 mm
FW856SU	activ C lockable distraction screw, 18 mm
FW861SU	activ C lockable distraction screw, 12 mm
FW862SU	activ C lockable distraction screw, 14 mm

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Read, follow and keep the instructions for use.
- Follow the instructions for use of activ C intervertebral disk prosthesis TA011995.
- Follow the operating manual (O31302).
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to the initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working tip: Exercise caution when passing products through the working channel (e.g. trocars).

The product is gamma-sterilized and supplied in sterile packaging.

The product must not be reused.

- Do not use products from open or damaged sterile packaging.
- Do not re-sterilize the product.
- Do not use the product after its use-by date.

Distraction set with lockable distraction screws

Safe operation

In the text below, the following short designations are used for the product components:

- activ C distractor for lockable distraction screws is referred to as "distractor"
- activ C distraction screwdriver for lockable distraction screws is referred to as "distraction screwdriver"
- activ C lockable distraction screws are referred to as "distraction screws"
- activ C lockable offset distraction screws are referred to as "offset distraction screws"

Note

Any application of the distraction screw includes the application of the offset distraction screw, if the latter is not mentioned explicitly!

Applying the distraction screw

Fig. 1 Distraction screwdriver with distraction screw and offset distraction screw

WARNING

Compression of, or damage to the spinal canal due to selection of a distraction screw of excessive length!

- Use the radiographic image to determine the correct length of the distraction screw according to the vertebral body.
- Always maintain x-ray control when applying the distraction screw.

Note

For very small vertebral bodies, offset distraction screws must be used to ensure that there is sufficient space for the implantation instruments!

- Choose distraction screw **7** or offset distraction screw **4**.
- Insert distraction screw **7** or offset distraction screw **4** in distraction screwdriver **3**.
- Place distraction screwdriver **3** on hexagon **5** of distraction screw **7** and make sure it engages.
- Check the self-holding function of distraction screw **7** in distraction screwdriver **3**.
 - or -
- Install distraction screwdriver **3** on hexagon **5** of offset distraction screw **4**, making sure that it engages.

Inserting the distraction screw into the vertebral body under X-ray control


Fig. 2

Note

For very small vertebral bodies, offset distraction screws must be used to ensure that there is sufficient space for the implantation instruments!

- Position and screw in distraction screw **7**/offset distraction screw **4** at the center of the vertebral body.
- Remove distraction screwdriver **3** from distraction screw **7**/offset distraction screw **4**.

Distraction set with lockable distraction screws

- Apply second distraction screw **7**/offset distraction screw **4**, ensuring the parallel positioning of distraction screw **7**/offset distraction screw **4**.

Mounting the distractor and applying the distraction, see Fig. A/B

WARNING

Risk of loosening/tear-out of the distraction screws when distracting with the distractor and the distraction screws!

- Use the distraction forceps for the distraction process.
- Do not use the distractor and the distraction screws for the actual distraction process.
- Lock the distractor in extended condition.

WARNING

Irritation of the facet joints/strain on spinal cord and nerve roots due to excessive distraction!

- Extend the distraction forceps slowly and carefully under X-ray control.
- Observe the tactile response.
- Use the heights of adjacent (healthy) disks as a guide for the required degree of distraction.

- Select distractor right **1** or distractor left **2**.
- Turn flap screw **8** to transfer the distance of parallel distraction screws **4/7** to distraction arms **12/14**.

The distance of distraction arms **12/14** corresponds to the distance of distraction screws **4/7**.

- Mount distractor right **1** or distractor left **2** on parallel distraction screws **4/7** so that they click into position.
- Distract with the distraction forceps.

- While distracting with the distraction forceps, extend distractor right **1** or distractor left **2** by turning flap screw **8**.
- Once the distraction process carried out with the distraction forceps is complete, lock distractor right **1** or distractor left **2** and hold the vertebrae in the distracted position.

Removing the distractor, see Fig. C

- To release the distraction after the implantation in the intervertebral space, press ratchet lever **9** at the distractor.
Distraction arms **12/14** move towards each other. The distractor is fully released from any load.
- Press buttons **13** at the side of distraction arm **12/14** and slide off the distractor from distraction screws **4/7**.
- Remove distraction screws **4/7** from the vertebrae, using distraction screwdriver **3**.

Disassembling

Disassembling the distractor, see Fig. C

- Pull distraction arm **12** (labeled "RIGHT/LEFT") along toothed rack **10** until toothed rack **10** is not seated in guide **11** anymore.

Assembling

Assembling the distractor, see Fig. C

Note

*Combine distraction arm **12** labeled "RIGHT" only with distraction arm **14** labeled "FW848R"!*

*Combine distraction arm **12** labeled "LEFT" only with distraction arm **14** labeled "FW849R"!*

- Push the ratchet lever **9**.
- Mount distraction arm **12** (labeled "RIGHT/LEFT") with guide **11** on the toothed rack, ensuring the correct orientation of distraction arms **12/14**.

Validated processing procedure

Note

Observe all relevant national regulations and standards with regard to processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, manufacturer's specifications can only serve as guide values for the processing procedures applied by the individual user.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

Single-use products

WARNING

Breakage of distraction screws/offset distraction screws due to reuse!

- The distraction screws/offset distraction screws are intended for single use only.

CAUTION

Possible damage to or destruction of products, due to processing!

Intended use: These products have been designed for single use only.

- Do not reprocess these products.

General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of "stainless" steels. Consequently the time interval between application and processing should not exceed 6 h, and neither potentially fixating pre-cleaning temperatures > 45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) should be used.

Excessive doses of neutralizers or basic detergents can cause chemical deterioration and/or fading of the laser engraving on stainless steel.

Residues containing chlorine or chlorides, as found e.g. in surgical residues, tinctures, medicines, saline solutions, service water used for cleaning, and in cleaning/disinfecting agents, cause corrosion damage (pitting, stress corrosion) and thereby destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals recommended by the chemicals' manufacturers as effective cleaning/disinfecting agents compatible with the materials used in the respective product. Strictly observe all application instructions, e.g. regarding temperatures, concentrations, exposure times, etc. Failure to do so can result in the following problems:

- Optical changes to the material, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH > 8 to cause visible surface changes or
- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use oxidizing process chemicals, which could cause bleaching/layer loss of the product.

Distraction set with lockable distraction screws

- Do not use process chemicals that cause stress cracking of plastics such as PPSU, or attack and cause brittleness e.g. of silicone softeners.
- Clean the product immediately after use.
Encrusted debris on HF instruments is softened most effectively and conservatively by a 5 minutes immersion treatment in a 3 % H₂O₂ solution. Encrustations can be removed by hand, with a soft brush and/or in an ultrasound bath. After that treatment, proceed through the normal steps of processing.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, Publications, Red Brochure – Proper maintenance of instruments.
- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove visible residues as completely as possible, using a lint-free, single-use cleaning tissue.

Preparation prior to cleaning

- Have the product dry in a disposal container and ready for immediate cleaning and disinfecting within 30 min after use.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use a cleaning and disinfecting agent according to the manufacturer's instructions. The cleaning and disinfecting agent – must be approved for the product materials (e.g. aluminum, plastics, high-grade steel) – must not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 60 °C.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Manual cleaning/disinfecting

- Clean hinged or jointed products in open and closed positions.
- Check visible surface for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfecting and cleaning with brush

FW848R–FW850R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Final rinse	RT (cold)	0.5	–	FD–W	–
IV	Drying	RT	–	–	–	–

D–W: Drinking water

FD–W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product under running tap water, using a suitable cleaning brush if necessary, until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives, which would damage the surfaces, for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage III

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage IV

- Dry the product completely with compressed air.

Distraction set with lockable distraction screws

Mechanical cleaning/disinfecting with manual pre-cleaning

FW848R–FW850R

Manual pre-cleaning with brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	–	D–W	–

D–W: Drinking water

RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives that would damage the surfaces for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and channels directly to the injector carriage attachment.
- Store products with hinges or joints on the tray in such a way that the joints are kept open.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	< 25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	- Concentrate, alkaline: pH = 10.9 < 5 % anionic tensides - 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD-W	- Concentrate, acid: pH = 2.6 Basis: Citric acid - 0.2 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD-W	-
V	Thermal disinfecting	94/201	10	FD-W	-
VI	Drying	90/194	40	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

Distraction set with lockable distraction screws

Control, care and inspection

- Allow the product to cool down to room temperature.
- Lightly lubricate moving parts (e.g. joints and links) with a maintenance oil suitable for the respective sterilization process (e.g. Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- Assemble the separable product, see Assembling.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Check for compatibility with associated products.
- Set aside the product if it is damaged.

Packaging

- Appropriately protect products with fine working tips.
- Secure product with retractor in the first notch.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters

Note

The product can be sterilized either in disassembled or in assembled condition.

Note

To prevent breakage due to stress crack corrosion, sterilize the instruments with the locking mechanism open, or locked on the first ratchet tooth.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilization acc. to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST 46-1993, validated acc. to EN ISO 17665 or EN 554/ISO 13683
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270— 275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

Storage

- Store processed products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

Technical Service

For service, maintenance and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:
Aesculap Implant Systems Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420
Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG