


## Surecan® Safety II

**B|BRAUN**

ENGLISH P 1-2

FRANÇAIS P 3-4

DEUTSCH P 5-6

ESPAÑOL P 7-8

ITALIANO P 9-10

SVENSKA P 11-12

PORTUGUÊS P 13-14

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NEDERLANDS P 19-20

NORSK P 21-22

DANSK P 23-24

SUOMI P 25-26

POLSKI P 27-28

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TÜRKÇE P 31-32

CESKY JAZYK P 33-34

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TIẾNG VIỆT P 39-40

ROMÂNĂ P 41-42

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










(B)

## Surecan® Safety II

Power-injectable safety non-coring needle for Access Ports

### Intended Use:

The Surecan® Safety II power-injectable safety non-coring needle is a device intended for insertion into the septum of a subcutaneously implanted port for the infusion of fluids into the port. The Surecan® Safety II safety feature is manually activated during needle removal, and is designed to aid in the prevention of accidental needle-sticks. When used with ports that are indicated for power injection of contrast media into the central venous system, the Surecan® Safety II needle is also indicated for power injection of contrast media.

For power injection of contrast media, the maximum recommended pressure is 325psi (22.4 bars).

### MR Information

#### General Information

According to IEC Standard 60601-2-33(2008), the scanner must be operated in Normal Operating Mode (defined as the mode of operation of the MR system in which none of the outputs have a value that causes physiological stress to the patient):

- The whole body averaged specific absorption rate (SAR) must be  $\leq 2.0 \text{ W/kg}$
- The head SAR must be  $< 3.2 \text{ W/kg}$

#### Device Information

The Surecan® Safety II was determined to be **MR-Conditional** according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-08. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

Non-clinical testing demonstrated that the **Surecan® Safety II** is MR Conditional. A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla and 1.5-Tesla
- Maximum spatial gradient magnetic field of 710 Gauss/cm or less
- Maximum whole body averaged specific absorption rate (SAR) of 2.9 W/kg for 15 minutes of scanning

#### MRI-Related Heating

In non-clinical testing, the **Surecan® Safety II** produced the maximum temperature rise during MRI performed for 15-min (i.e., per pulse sequence) in 3-Tesla (Excite, Software 14X.M5, General Electric Healthcare, Milwaukee, WI) MR systems, as:

MRI Condition	MR System Reported, Maximum Whole Body Average SAR (W/kg)	Calorimetry Value (W/kg)	Highest Temperature Change	Time for MRI (per pulse sequence)
3-T / 128-MHz	2.9	2.7	2.0 °C	15-min

#### Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the **Surecan® Safety II**. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

#### Contraindications:

- DO NOT USE, if the presence of an access port-related infection, bacteremia, or septicemia is known or suspected.
- DO NOT USE, if local tissue factors will prevent proper Surecan® Safety II device stabilisation and/or access.

#### Warnings:

- Verify needle length is correct based on port reservoir depth, tissue thickness, and the thickness of any dressing beneath the base of the needle; if too long, the needle and/or port may be damaged at insertion; if too short, the needle may not completely pierce the port septum, and medication may be delivered into surrounding tissue and/or the needle may be blocked.
- Fully tighten all connections or needle-less connectors before use (Fig.1). Failure to attach an end cap after removing a male Luer locking end cap can result in an embolism or bleeding.
- Failure to use the device correctly when removing the needle from the port site could result in the needle tip re-emerging from the base, resulting in an accidental needle-stick with a contaminated needle. A needle-stick with a contaminated needle may cause infectious disease.

- Handle and discard in accordance with accepted medical practice and all applicable regulations. After use, this product is a potential biohazard.
- The device and its accessories are not re-usable nor designed to be re-used. Any re-use would definitely compromise the performance and safety of the device.
- Do not alter this device.

**Cautions:**

- Store at room temperature, protected from light.
- Single use only, do not re-sterilise the product, and destroy after use.
- The product is sterilized using ethylene oxide and is sterile and non-pyrogenic, undamaged, individual packages.
- Check the expiry date before use and never use a device past the expiry date.
- Do not use this product if it appears damaged or if the package has been previously opened or damaged.
- Use only under aseptic conditions.
- Care must be taken to avoid accidental needle-sticks.
- Only qualified healthcare practitioners should insert, manipulate and remove these devices.
- Carefully read and follow all instructions prior to use.
- Follow all instructions, contraindications, warnings and cautions for all infusates, ports, IV pumps, IV sets and needle-less systems, as specified by the respective manufacturer.
- Follow local standard infection control precautions.
- Confirm correct needle placement in the port reservoir by aspiration of blood before infusion of any substance. If there is doubt regarding proper needle placement, perform a radiographic dye procedure to confirm placement per local protocol.
- Do not remove and reinsert the needle into the port.
- Avoid excessive manipulation once the needle is in the port.
- Do not access the needle-less Y-injection site valve with a needle. Valve puncture may result in an embolism.

**Warning for high flow rate/high pressure injection:**

- Always verify that the access port implanted is power injectable.
- Always verify that the port and catheter are functional by aspirating 2mL of blood into a syringe and injecting 5mL of sodium chloride (NaCl) 0.9% into the port/catheter before attempting to start an infusion of medication.
- Do not exceed the recommended pressure (325psi – 22.4 bars) and access port maximum flow rate recommended by the manufacturer as device failure may occur.
- Do not exceed the maximum pressure indicated on the needle (325psi). (Fig. 6)
- Contrast media should be warmed to 37°C (98.6°F) before use. Failure to do so will result in up to 50% lower flow rates and/or device failure.
- Depending on the technical characteristics of the injector system, the target flow rate might not be attained.
- Ensure that the needle is correctly placed in the port, securely taped to the skin and covered with an adhesive dressing before starting high pressure injection.
- Vigorously flush the Surecan® Safety II needle using a 10 ml or larger syringe and sterile sodium chloride (NaCl) 0.9% prior to and immediately following the completion of power injection studies.
- The Surecan® Safety II indication for power injection of contrast media implies the ability of the system to withstand the procedure, but does not imply appropriateness of the procedure for a particular patient. A suitably trained clinician is responsible for evaluating the health status of a patient as it pertains to a power injection procedure.

**Instructions for Use:**

Use aseptic technique. Refer to the numbered illustrations when using this device.

**Needle insertion:**

- Prepare port site for sterile needle insertion. (fig. 2)
- Attach a 10 ml syringe containing sodium chloride (NaCl) 0.9% to the proximal Luer lock connection of the safety non-coring needle. Grasp the wings and remove the needle cover. Prime and flush the infusion set. (fig. 3)
- Insert the needle perpendicular into the port septum. (fig. 4) You may position the base on the skin as needed.
- Verify patency. Dress and secure the site per institutional protocol. (fig. 5) Begin injection, infusion, or blood aspiration per institutional procedure.

**Needle removal**

- After completion of treatment, flush the port per institutional protocol. Stabilise the port by securely holding the base down. (fig. A)
- Firmly pull the wings up until you feel a firm stop (fig. B) and the needle is locked in the safety position. A green dot on the base of the needle confirms safety device deployment.
- Dispose of set in a sharps container

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