Aesculap Neurosurgery

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KAIRison pneumatic bone punches

DISASSEMBLY


I Pull back the locking mechanism to release the shaft.

II Turn the shaft in the direction of the arrow.

III Unhook the shaft and lift it off the handle.

IV Pull back the upper shaft (1) and lift it off (2).

Safe handling and preparation

- Read and follow the instructions for use, keeping them in a safe place. Use the product only in accordance with its intended purpose (Instructions for Use TA011795).
- Clean the new instrument thoroughly, either manually or mechanically, prior to initial sterilization.
- Visually inspect the product before each use for loose, bent, broken, cracked, worn or missing parts.
- Do not use the product if it is damaged or defective. Set aside damaged products immediately.
- Replace damaged components immediately with original spare parts.

Cleaning, Disinfection and Sterilization

- For general information on reprocessing, please also refer to the Aesculap Extranet: www.aesculap-extranet.net
- For further detailed instructions on safe, hygienic reprocessing to preserve the substance and value of your instruments, see www.a-k-i.org (red brochure: Proper Maintenance of Instruments).
- The shafts of the pneumatic bone punches can be sterilized either assembled or disassembled.


- Close the opening with silicone caps before washing.
- Do not clean the handle and reducer in an ultra-sonic tub.

After each use

- Reprocess contaminated instruments as quickly as possible. Clean manually or mechanically – for cleaning parameters, please refer to the Aesculap Extranet at www.aesculap-extranet.net or to the instructions for use.
- Disassemble the shaft into its individual components (see Disassembly). Dry disposal is recommended. For wet disposal, use an active cleaning disinfectant (follow the manufacturer's instructions for use). Make sure that all exterior and interior surfaces, lumina and openings are covered. Prior to mechanical cleaning and disinfection, rinse the instrument thoroughly in clear running water.
- After chemical disinfection, rinse again thoroughly and intensively using clear running water. Follow the disinfectant manufacturer's instructions.
- Remove encrusted dirt with a soft synthetic brush. Do not use abrasive cleaning agents or metal brushes. Clean lumina and channels with a soft round synthetic brush of the appropriate diameter.

- Use demineralised water for the final rinse.
- Dry the instrument completely with compressed air.
- Oil moving parts (joint, inner surfaces of sliding component, guiding mechanisms) well with steam permeable Aesculap Sterilit maintenance oil JG598 (see illustration).
- Sterilization must be carried out in accordance with a validated steam sterilization procedure (e.g. in a sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI/ST46-1993 and validated in conformance with EN 554/ISO 13683). If the fractionated vacuum procedure is used, sterilization must be performed using the 134°C/2 bar program with a minimum holding time of 5 minutes.

Width	Punch shaft 235 mm	Punch shaft 285 mm
2 mm	FK882R	FK892R
3 mm	FK883R	FK893R
4 mm	FK884R	FK894R
5 mm	FK885R	FK895R
6 mm	FK886R	FK896R


The instrument areas indicated should be adequately lubricated after each cleaning cycle with Aesculap Sterilit oil, e.g. JG598.

- When sterilizing several products simultaneously in a steam sterilizer, ensure that the maximum permitted loading of the steam sterilizer as specified in the manufacturer's instructions is not exceeded.


Storage

- Keep the instrument in a suitable storage tray under low contamination conditions in a dry, dark and cool place.


Individual components

- FK899R Pneumatic punch handle
- FK898R Pressure reducer
- FK887R Punch bracket pressure reducer
- FK888SU Disposable tubes
- GA466 Pneumatic hose with Aesculap-Dräger connection
- GA468 Pneumatic hose with Aesculap-Schrader connection

ASSEMBLY

I Insert the ejector (1), push forward until the ball clicks into place (2).

II Check shaft position, hook shaft into handle.

III Turn shaft towards handle.

IV Check correct shaft position before use.

Aesculap AG & Co. KG
Am Aesculap-Platz, 78532 Tuttlingen, Germany
Phone +49 7461 95-0, Fax +49 7461 95-2600
www.aesculap.de