GB	Instructions for use
USA	Metha rasp handle
D	Gebrauchsanweisung
	Metha-Raspelhandgriff
F	Mode d'emploi
	Porte-râpe Metha
E	Instrucciones de manejo
	Mango de raspador Metha
I	Istruzioni per l'uso
	Impugnatura per raspa Metha
P	Instruções de utilização
	Cabo de grossa Metha
DK	Brugsanvisning
	Metha-raspehåndgreb
FIN	Käyttöohjeet
	Metha-kaavinkahva
GR	Οδηγίες χρήσης
	Metha-χειρολαβή ράσπας
CZ	Návod k použití
	Rukojet' Metha-rašple

B | BRAUN

SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 011592 04/11 V5 Änd.-Nr. 43048

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Símbolo CE, em conformidade com a Directiva 93/42/CEE
CE-markering iht. retningslinje 93/42/EEC
93/42/EEC-standardin mukainen CE-hyväksyntä
Ενδειξη CE σύμφωνα με την Οδηγία 93/42/EEC
CE-označení podle směrnice 93/42/ES

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Retten til tekniske ændringer forbeholdes
Oikeusteknisistä syistä johtuviin muutoksiin pidätetään
Επιφυλασσόμεθα για τεχνικές αλλαγές
Technické změny vyhrazeny


Fig. 1


Fig. 2


Fig. 3

Legend

- 1Rasp handle
- 2Hammer plate
- 3Handle part
- 4Pushbutton
- 5Spring
- 6Locking screw
- 7Rasp

Symbols on product and packages

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The rasp handle is used in hip-endoprosthetic interventions. It holds prosthesis-specific form rasps, which are driven in and out during the preparation of the femoral marrow cavity

Available sizes

Different variants of the rasp handle are available for different operating techniques, e.g.:

Designation	Art. no.
Metha rasp handle	
lateral access	
- straight	NF180R
- left angled	NF141R
- right angled	NF142R
posterior access	
- straight	NF144R
- left angled	NF142R
- right angled	NF141R
anterior access	
- straight	NF140R
- left angled	NF138R
- right angled	NF139R

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Note
Use of the Metha rasp handle requires precise knowledge in hip endoprosthetics and the biomechanical situation at the hip!

- Make certain that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only for its intended purpose, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

Risk of injury and/or malfunction!
➤ Always carry out a function check prior to using the product.

WARNING

Inserting a rasp into the rasp handle

- Push and hold down pushbutton 4.
- Slide rasp 7 into handle part 3, as shown in Fig. 1, making certain that the notch in the rasp connector points towards the pushbutton.
- Release pushbutton 4.
- To check that the rasp is firmly seated, make certain that
 - there is no gap between the endplates,
 - the rasp cannot be rotated against the handle.

Driving the rasp in and out

Risk of injury and/or damage to the product due to incorrect application of the hammer!
➤ When driving the rasp in and out, using a hammer, only hit the hammer plate of the Metha rasp handle.

WARNING

Disassembling

- Push in pushbutton 4 at least halfway down.
- Rotate pushbutton 4 by approx. 90° in any direction and release the pushbutton.
The locking mechanism becomes detached.
- Remove pushbutton 4 and spring 5.

Note
Never unscrew lateral locking screw 6 on the rasp handle. This screw has been sealed with a special adhesive during the manufacturing process; it must not be removed.

Assembling

- Put spring 5 onto pushbutton 4.
- Slide pushbutton 4 into handle part 3, making sure that the notch on the pushbutton is aligned as shown in Fig. 3.
- Push in pushbutton 4 at least halfway down.
- Rotate pushbutton 4 clockwise by approx. 90°, or turn the notch on pushbutton 4 to the front.

Validated reprocessing procedure

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note
Mechanical processing should be favored over manual cleaning as it gives better and more reliable results.

Note
It should be noted that successful processing of this medical device can only be guaranteed following prior validation of the processing method. The operator/sterile processing technician is responsible for this.

Note
Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General information

To prevent increased contamination of loaded instrument trays during use, please ensure that contaminated instruments are collected separately and not returned to the instrument tray.

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion of stainless steel. You should therefore leave no more than 6 hours between use and cleaning, not use pre-clean fixing temperatures >45 °C, and ensure that non-fixing disinfectants (aldehyde-free/alcohol-free) are used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chloride - e.g. in surgical residues, drugs, saline solutions and water for cleaning, disinfection and sterilization - may cause corrosion damage to stainless steel (pinholing, stress corrosion), thus rendering the products unusable. These must be removed by rinsing thoroughly with demineralized water and then drying.

Only use process chemicals which have been tested and approved (e.g. VAH/DGHM or FDA-certified or CE marking) and recommended by the chemical manufacturer as being compatible with the materials. All the chemical manufacturer's application specifications regarding temperature, concentration and contact time should be strictly observed. Failure to do so can result in the following problems:

• Visible changes to materials e.g. fading or changes in the color of titanium or aluminum. As regards aluminum, visible changes to the surface may occur at a pH level of >8 in the application/working solution.

• Material damage such as corrosion, cracks, fractures, premature deterioration or swelling.

Preparations at the place of use

- Disassemble the product immediately after use, as described in these instructions for use, see Disassembling.
- Rinse non-visible surfaces such as those in instruments with concealed crevices, lumens or complex geometries, preferably with distilled water, using a disposable syringe for instance.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Cleaning/disinfection

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!
➤ Use cleaning and disinfecting agents approved for aluminum, plastics, and high-grade steel, for example, according to the manufacturer's instructions.
➤ Observe specifications regarding concentration, temperature and exposure time.
➤ Do not exceed the maximum permitted cleaning temperature of 120 °C .

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The disinfector must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note
For thermal disinfection, always use fully desalinated water (demineralized, low microbial contamination: drinking water quality at least) and ensure that Ao is >3 000 for the process.

Note
The disinfector used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	BBraun Stabimed; aldehyde-free, phenol-free and QUAT-free
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room temperature

Stage I

- Clean the product in the ultrasonic cleaning bath (frequency 35 kHz) for at least 15 minutes. Make certain that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden gaps, lumens or complex geometries, for at least 1 minute or until no more residues can be removed, using a suitable cleaning brush (e.g. TA006874).
- Mobilize non-rigid components such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).

- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	BBRAUN HELIMATIC CLEANER alkaline with tensides; application solution 0.5%
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfection	90/194	5	FD–W	–
V	Drying	–	–	–	According to disinfectant program

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

Inspection, maintenance and checks

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks , lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting, and drying cycle, check that the instrument is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Make certain that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable instruments, see Assembling.
- Check for compatibility with associated products.

Packaging

- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization

- Make certain that the sterilant can reach all external and internal surfaces (e.g. by opening valves and stopcocks).
- Validated sterilization process
 - Disassemble the instrument.
 - Steam sterilization using the fractionated vacuum process
 - Steam sterilizer specified in DIN EN 285 and validated in accordance with DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer: ensure that the maximum permitted load specified by the manufacturer for the steam sterilizer is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre–vacuum	270 °F–275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the ineffectiveness of reprocessing and sterilization to eliminate the risk of cross-contamination!

Technical Service

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technical Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Aesculap Implant Systems, LLC
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems, LLC
3773 Corporate Parkway
Center Valley, PA 18034
USA

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