

- (GB)** Instructions for use/Technical description
- (USA)** MINOP micro-instruments
- (D)** Gebrauchsanweisung/Technische Beschreibung
MINOP-Mikro-Instrumente
- (F)** Mode d'emploi/Description technique
Micro-instruments MINOP
- (E)** Instrucciones de manejo/Descripción técnica
Microinstrumental MINOP
- (I)** Istruzioni per l'uso/Descrizione tecnica
Microstrumenti MINOP
- (P)** Instruções de utilização/Descrição técnica
Microinstrumentos MINOP
- (NL)** Gebruiksaanwijzing/Technische beschrijving
MINOP-micro-instrumenten
- (S)** Bruksanvisning/Teknisk beskrivning
MINOP-mikroinstrument
- (RUS)** Инструкция по применению/Техническое описание
Микроинструменты MINOP
- (CZ)** Návod k použití/Technický popis
MINOP mikro nástroje
- (PL)** Instrukcja użytkowania/Opis techniczny
Mikroinstrumenty MINOP
- (SK)** Návod na používanie/Technický opis
Mikronástroje MINOP
- (TR)** Kullanım Kılavuzu/Teknik açıklama
MINOP-mikro-enstrüman

B | BRAUN

SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 010327 07/15 V6 Änd.-Nr. 51280

CE 0123 – DIR 93/42/EEC

Aesculap®
MINOP micro-instruments

Legend

- 1 Locking knob
- 2 Handle with shaft
- 3 Jaw part connected (on handle)
- 4 Jaw parts with notched rod

Symbols on product and packages

Caution, general warning symbol
Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

Intended use

The micro-instruments are intended for use in neurosurgery, and are specially designed for use in narrow spaces. They can be used for cutting, grasping, dissecting, etc. of blood vessels and tissue.

Indications

Indications, see Intended use.

Note
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Contraindications

None known.

Risks and side effects

As part of the legal duty to inform, the following typical risks and side effects associated with the use of surgical instruments are referred to. These are predominantly procedure-specific, non product-specific and not limited to unwanted damage to surrounding tissue resulting in e.g. bleeding, infection, material incompatibilities or instrument parts remaining unnoticed in the patient, etc.

Available sizes

The micro-instruments are available as complete instruments and as separate, combinable components. Various jaw parts (e.g. straight scissors, curved scissors or forceps) can be used with the appropriate handles.

Note
The current product range is described in Aesculap brochure C 74511.

Safe handling and preparation

- CAUTION
- Federal law restricts this device to sale by or on order of a physician!
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
 - Read, follow, and keep the instructions for use.
 - Use the product only in accordance with its intended use, see Intended use.
 - Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
 - Store any new or unused products in a dry, clean, and safe place.
 - Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Replace any damaged components immediately with original spare parts.

Safe operation

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

Risk of injury when using the product beyond the field of view!

- Apply the product only under visual control.

Disassembling

- Turn the locking knob 1 by 180° so that it disengages.
- Hold the handle stems 2 with one hand and the jaw parts with notched rod 4 with the other hand.
- Pull to separate the jaw part mounting 3 and the jaw parts with notched rod 4.
- Turn the locking knob 1 back to its original position.

Assembling

- Ensure that the bridge of the jaw part 4 and the groove in the jaw part mounting 3 are in the same axis when the handle stems 2 are in closed position.
- Push the jaw part mounting 3 and jaw parts with notched rod 4 together until the notched rod engages with an audible click.
- Function check: Press the handle stems 2 together, all the way to the stop. Check that the jaw parts operate properly.

Validated reprocessing procedure

General safety notes

- Note*
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.
- Note*
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.
- Note*
It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this. The specified chemistry was used for validation.
- Note*
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at [www.extranet.bbraun.com](https://extranet.bbraun.com).
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

- Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.
- Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result damage to metallic products. These must be removed by rinsing thoroughly with demineralized water and then drying.
- Perform additional drying, if necessary.
- Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:
- Optical changes in materials, e.g., fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
 - Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.
 - Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
 - For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-l.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open up products with hinges.

Preparations at the place of use

- If applicable, rinse surfaces that are not accessible to visible inspection (preferably with demineralized water), using a disposable syringe, for example.
- Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

Danger to the patient!

- Reprocess the product only with manual pre-cleaning followed by mechanical cleaning.

Risk to patient due to cross contamination!

- Do not clean contaminated products together with uncontaminated products in a tray.

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for (e.g. aluminum, plastics, high-grade steel),
 - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum cleaning temperature of 55 °C.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfecting machine must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Ultrasonic cleaning	RT (cold)	>15	1	D–W	Dr. Weigert neodisher® SeptoClean®
II	Rinsing	RT (cold)	1	–	D–W	–

D–W: Drinking water
RT: Room temperature
*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

- Phase I
- ▶ Clean the product in an ultrasound cleaning bath (frequency of 35 kHz) for at least 15 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
 - ▶ Clean the product with a suitable cleaning brush (e.g. FM32800) in the solution until all visible residues have been removed from the surface.
 - ▶ If applicable, brush non-visible surfaces for at least 1 minute with a suitable cleaning brush (e.g. FM32800).
 - ▶ Mobilize non-rigid components, such as set screws and hinges, during cleaning.
 - ▶ Then flush these areas thoroughly at least five times with the cleaning disinfectant solution using a disposable syringe (20 ml).

- Phase II
- ▶ Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
 - ▶ Mobilize non-rigid components, such as set screws and hinges, during rinsing.

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	Dr. Weigert neodisher® SeptoClean 1 %* working solution
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfection	90/194	5	FD–W	–
V	Drying	–	–	–	In accordance with the program for the cleaning and disinfecting machine

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

- ▶ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots)
- ▶ Place instruments in the tray with their hinges open.
- ▶ Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ Dry the product if it is wet or damp.
- ▶ Repeat cleaning and disinfection of products that still show impurities or contamination.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Assemble dismountable products, see Assembling.
- ▶ Check for compatibility with associated products.

Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note
The product can be sterilized either in disassembled or in assembled condition.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractional vacuum process at 134 °C, holding time 18 minutes for prion inactivation
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{–6}, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- ▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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