GB
USA

Instructions for use

MINOP®/Paediscscope Disposable Introducer for Neuroendoscopy

D

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
TA-Nr. 012165 08/10 Änd.-Nr. 41794

0123

CE marking according to directive 93/42/EEC

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MINOP®/Paediscopes Disposable Introducer for Neuroendoscopy

Legend

- 1 Obturator
- 2 Sheath
- 3 Obturator and sheath locked assembly
- 4 Locking flange
- 5 Locking knob
- 6 Sheath handles
- 7 Locking threads
- 8 Handles

Symbols on product and packages

	Sterilization using ethylene oxide
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The purpose of the Disposable Introducer is to obtain and maintain a temporary pathway to the ventricular system of the patient. The Disposable Introducer is designed to be split lengthwise and peeled down to the skull level of the patient, accommodating different depth requirements.

The Disposable Introducer consists of two main parts: Obturator and sheath, which are locked to an assembly.

Available sizes

Disposable Introducer	Size (Ø inside of sheath)
FH603SU	10 French (3.2 mm)
FH604SU	19 French (6.1 mm)

Disposable Introducer	Typical usage with
FH603SU	Paediscopes PF010A
FH604SU	MINOP® 6 mm trocar FF399R

Note

These products DO NOT contain latex.

Safe handling and preparation

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product.

The product is EO sterilized and sterile packed.

The product must not be reused.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only for its intended purpose., see Intended use.
- Do not use products from open or damaged sterile packaging.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Do not re-sterilize the product.
- Do not use the product after its use-by date.

Safe operation

- Ensure that the obturator 1 is locked into the sheath 2 by twisting the locking flange 4 of the obturator, thereby locking it into the locking threads 7 of the sheath. Compare with picture 3.
The tip of the obturator extends from the sheath.
There is no gap between the tip of the sheath and the outer circumference of the obturator.
- Determine the depth the introducer has been inserted by means of the graduation on the side of the obturator.
- Insert the Disposable Introducer to the desired depth, unlock and remove the obturator, being sure to maintain its sterility.
- Insert the scope carefully, making sure that the sheath does not advance.
The scope allows for visual confirmation that the sheath is in the desired location.
- If the sheath is not in the desired location, carefully remove the scope and sheath.

- Reinsert and lock the obturator into the sheath.
- Repeat ventricular insertion. The tip of the sheath should be maintained inside the ventricular system after removal of the obturator.
This reduces the likelihood of the tip to cause hemorrhaging due to the movement or advancement of the scope.

The circular tip of the sheath could cause haemorrhaging!

- Never insert or advance the sheath into the brain without having the obturator locked into the place.
- Extend the obturator completely and lock it into the sheath.

- Two sterile persons are required to affix the sheath.
- One person has to hold the depth of the sheath constant.
- The other firmly grips the sheath handles 6.
- Pulling apart the handles 8 will result in the splitting of the sheath lengthwise.
- Split the sheath down to a level just above the skull.
- Firmly affix the split sheath in place by suturing or stapling down the two split halves of the sheath.
This has now established the pathway to the ventricle.
- When you now insert anything through the sheath, always be cognizant that the sheath must not move, so tight fitting objects should be avoided.

Single-use products

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- Do not reprocess the product.

Technical Service

For service, maintenance or repairs, please contact your national B. Braun/ Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-mail: ats@aesculap.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA
Phone: +1 800 282-9000
www.aesculapusa.com

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