
|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|-----|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| en  | <b>Instructions for use/Technical description</b><br>MINOP®/PaediScope single-use introducer for neuro-endoscopy                                                                                                                                                                                                                                                                                                                                                                                                                     |
| USA | Note for U.S. users<br>This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at <a href="http://www.aesculapusa.com">www.aesculapusa.com</a> . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost. |
| de  | <b>Gebrauchsanweisung/Technische Beschreibung</b><br>MINOP®/PaediScope Einmal-Introducer für Neuro-Endoskopie                                                                                                                                                                                                                                                                                                                                                                                                                        |
| fr  | <b>Mode d'emploi/Description technique</b><br>Intubateur à usage unique MINOP®/PaediScope pour neuro-endoscopie                                                                                                                                                                                                                                                                                                                                                                                                                      |
| es  | <b>Instrucciones de manejo/Descripción técnica</b><br>Introductor desechable MINOP®/PaediScope para neuroendoscopia                                                                                                                                                                                                                                                                                                                                                                                                                  |
| it  | <b>Istruzioni per l'uso/Descrizione tecnica</b><br>Introduttore monouso per neuro-endoscopia MINOP®/PaediScope                                                                                                                                                                                                                                                                                                                                                                                                                       |
| pt  | <b>Instruções de utilização/Descrição técnica</b><br>MINOP®/introdutor descartável PaediScope para endoscopia neurológica                                                                                                                                                                                                                                                                                                                                                                                                            |
| nl  | <b>Gebruiksaanwijzing/Technische beschrijving</b><br>MINOP®/PaediScope Wegwerp-introducer voor neuro-endoscopie                                                                                                                                                                                                                                                                                                                                                                                                                      |
| da  | <b>Brugsanvisning/Teknisk beskrivelse</b><br>MINOP® /PaediScope engangsintroducer til neuroendoskopi                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| sv  | <b>Bruksanvisning/Teknisk beskrivning</b><br>MINOP®/PaediScope engångsintroducer för neuro-endoskopi                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| fi  | <b>Käyttöohje/Tekninen kuvaus</b><br>Neurokirurginen MINOP®/PaediScope -enteroskooppi                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| lv  | <b>Lietošanas instrukcijas/tehniskais apraksts</b><br>MINOP®/PaediScope vienreiz lietojams ievadītājs neiroendoskopijai                                                                                                                                                                                                                                                                                                                                                                                                              |
| lt  | <b>Naudojimo instrukcija/techninis aprašas</b><br>MINOP®/PaediScope vienkartinis introduseris neuroendoskopijai                                                                                                                                                                                                                                                                                                                                                                                                                      |
| ru  | <b>Инструкция по применению/Техническое описание</b><br>Одноразовый интродьюсер для нейроэндоскопии MINOP®/PaediScope                                                                                                                                                                                                                                                                                                                                                                                                                |
| cs  | <b>Návod k použití/Technický popis</b><br>MINOP®/PaediScope – zavaděč k jednomu použití pro neuroendoskoppii                                                                                                                                                                                                                                                                                                                                                                                                                         |
| pl  | <b>Instrukcja użytkowania/Opis techniczny</b><br>Jednorazowy instrument do wprowadzania MINOP®/PaediScope do neuroendoskopii                                                                                                                                                                                                                                                                                                                                                                                                         |
| sk  | <b>Návod na použitie/Technický opis</b><br>MINOP®/PaediScope jednorazový vkladáč pre neuroendoskopi                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| hu  | <b>Használati útmutató/Műszaki leírás</b><br>MINOP®/PaediScope eldobható introducer neuroendoszkópiához                                                                                                                                                                                                                                                                                                                                                                                                                              |
| sl  | <b>Navodila za uporabo/Tehnični opis</b><br>MINOP®/ PaediScope za enkratno uporabo za nevroendoskopijo                                                                                                                                                                                                                                                                                                                                                                                                                               |
| hr  | <b>Upute za uporabu/Tehnički opis</b><br>MINOP®/PaediScope jednokratna uvodnica za neuroendoskopiju                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| ro  | <b>Manual de utilizare/Descriere tehnică</b><br>Introducător de unică folosință MINOP®/PaediScope pentru neuroendoscopie                                                                                                                                                                                                                                                                                                                                                                                                             |
| bg  | <b>Упътване за употреба/Техническо описание</b><br>MINOP®/PaediScope инструмент за въвеждане за невроендоскопия за еднократна употреба                                                                                                                                                                                                                                                                                                                                                                                               |
| tr  | <b>Kullanım Kılavuzu/Teknik açıklama</b><br>MINOP®/PaediScope endoskopi için tek kullanımlık uygulayıcı                                                                                                                                                                                                                                                                                                                                                                                                                              |
| el  | <b>Οδηγίες χρήσης/Τεχνική περιγραφή</b><br>Διάταξη εισαγωγής μίας χρήσης MINOP®/PaediScope για νευροενδοσκόπηση                                                                                                                                                                                                                                                                                                                                                                                                                      |

Aesculap®  
MINOP®/PaediScope single-use introducer for neuro-endoscopy

Legend

- 1 Obturator
- 2 Sheath
- 3 Obturator and sheath assembled and interlocked
- 4 Interlock flange
- 5 Release locking knob
- 6 Sheath handles
- 7 Set screws
- 8 Handles

Symbols on product and packages

|                                                                                  |                                                                                         |
|----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
|  | Sterilization using ethylene oxide                                                      |
|  | Not for reuse in intended applications as defined by the manufacturer                   |
|   | Use by                                                                                  |
|  | Caution, general warning symbol<br>Caution, see documentation supplied with the product |

Intended use

The single-use introducer serves the creation and preservation of a temporary access canal to the patient's ventricular system. For the adjustment to different entry depths, the single-use introduce is constructed in such a way that the sheath covering can be cleaved axially and the sheath halves can be hinged away right up to the cranial surface. The single-use introducer consists of two main components: Obturator and sheath, which are interlocked.

Indications

Indications, see Intended use.

*Note*  
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Contraindications

Absolute contraindications

None known.

Relative contraindications

Based on our current knowledge, there are no product-specific relative contraindications. However, there are medical or surgical conditions that may interfere with endoscopic technique, such as strong bleeding or high protein levels within the ventricular system limiting the view in the surgical field. In case of relative contraindications, the user shall decide on the use of the product on an individual basis.

Risks and side effects

As part of the legal obligation to provide information, the possible risks and side effects in connection with the use of surgical instruments known to the manufacturer are highlighted below. These are predominantly process-specific, not product-specific, and include unwanted damage to surrounding tissue, such as bleeding, infection, incompatibility of materials, or parts of instruments left in the patient unperceived.

Available sizes

| Art. no. Single-use introducer | Color of sheath handle | Size (interior part of sheath Ø) | For use in conjunction with (typical combination) |
|--------------------------------|------------------------|----------------------------------|---------------------------------------------------|
| FH603SU                        | Purple                 | 10 French (3.2 mm)               | PaediScope PF010A                                 |
| FH604SU                        | Pink                   | 19 French (6.1 mm)               | MINOP® 6 mm trocar FF399R                         |

*Note*  
These products are LATEX-FREE.

Safe handling and preparation

**Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!**

► Do not reprocess the product!

The product is EO sterilized and sterile packed.  
The product must not be reused.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Do not use products from open or damaged sterile packaging.
- Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Do not re-sterilize the product.
- Do not use the product after its use-by date.

Safe operation

- Turn the interlock flange 4 of the obturator, thereby interlocking the obturator into the set screws 7 of the sheath, in order to ensure that the obturator 1 is interlocked into the sheath 2. Compare image 3. The tip of the obturator protrudes from the sheath. At the end of the sheath, there is no visible gap between the obturator and the shaft.
- With the aid of the scale on the side of the obturator, determine the depth to which the introducer should be inserted.
- Insert the single-use introducer to the desired depth and unlock and remove the obturator. When doing this, ensure that the obturator remains sterile.
- Carefully insert the endoscope with pushing forward the sheath. It can be visually confirmed with the endoscope that the sheath is positioned correctly.
- If the sheath is not in the desired position, carefully remove the endoscope and the sheath.
- Insert the obturator into the sheath again and interlock.
- Repeat the ventricular insertion. After removing the obturator, the sheath end should remain in the ventricular system. This reduced the risk of the tip causing hemorrhage when the endoscope is pushed forwards.

**The round sheath end may cause hemorrhages!**

► Under no circumstances push the sheath forwards into the brain without the obturator being interlocked in the correct position.

► Lower the obturator completely and interlock it in the sheath.

*Note*  
Two sterile individuals are required for fixing the sheath.

- Person 1: Hold the sheath depth steady.
- Person 2: Hold sheath handles 6 securely.
- Separate handles 8. The sheath splits lengthways.
- Split the sheath to just above the cranial bone.
- Fix the split sheath securely to the cranium by sewing or tacking both sheath halves. The access canal to the ventricle is created.
- When introducing objects through the sheath, attention must be paid to ensuring that the sheath does not move. Do not insert any objects that barely fit through the sheath.

Validated reprocessing procedure

Single-use products

**Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!**

► Do not reprocess the product!

Storage

- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Technical Service

**Risk of injury and/or malfunction!**

► Do not modify the product.

► For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

**Service addresses**  
Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 (7461) 95-1602  
Fax: +49 (7461) 16-5621  
E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging! Incineration of the product, its components and its packaging produces harmless residues.
- Comply with national regulations for incineration under all circumstances!

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