- (GB)** Instructions for use/Technical description
MINOP TREND (TRansnasal ENDoscopy system)
- (USA)** MINOP TREND (TRansnasal ENDoscopy system)
- (D)** Gebrauchsanweisung/Technische Beschreibung
MINOP TREND (TRansnasales ENDoskopiesystem)
- (F)** Mode d'emploi/Description technique
MINOP TREND (TRansnasal ENDoscopy system)
- (E)** Instrucciones de manejo/Descripción técnica
MINOP TREND (sistema TRansnasal de ENDoscopia)
- (I)** Istruzioni per l'uso/Descrizione tecnica
MINOP TREND (sistema per ENDoscopia TRansnasale)
- (P)** Instruções de utilização/Descrição técnica
MINOP TREND (TRansnasal ENDoscopy system)
- (NL)** Gebruiksaanwijzing/Technische beschrijving
MINOP TREND (TRansnasaal ENDoscopiesysteem)
- (S)** Bruksanvisning/Teknisk beskrivning
MINOP TREND (TRansnasalt ENDoskopisystem)
- (RUS)** Инструкция по применению/Техническое описание
MINOP TREND (Трансназальная ЭНДоскопическая система)
- (CZ)** Návod k použití/Technický popis
MINOP TREND (TRansnasální ENDoskopický systém)
- (PL)** Instrukcja użytkowania/Opis techniczny
MINOP TREND (TRansnasal ENDoscopy System)
- (SK)** Návod na používanie/Technický opis
MINOP TREND (TRansnazálny ENDoskopický systém)
- (TR)** Kullanım Kılavuzu/Teknik açıklama
MINOP TREND (TRansnazal ENDoskopi sistemi)

B | BRAUN

SHARING EXPERTISE



Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company


TA-Nr. 012407 11/12 V6 Änd.-Nr. 45488


Technical alterations reserved


1

2


3

4

5

6

Aesculap®

MINOP TREND (TRansnasal ENDoscopy system)

Legend

- 1 MINOP TREND (TRansnasal ENDoscopy system)
- 2 MINOP TREND suction/irrigation sheath, rotatable
- 3 MINOP TREND handle part
- 4 Irrigation button
- 5 MINOP TREND endoscope
- 6 Adapter for holding arm
- 7 Conical sheath connection

Symbols on product and packages

Caution, general warning symbol

Caution, see documentation supplied with the product

Scope

- MINOP TREND
- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com

Intended use

The MINOP TREND (TRansnasal ENDoscopy System) is used for operating sites that can be reached through a transnasal approach. The MINOP TREND (TRansnasal ENDoscopy System) keeps the field of view of the adapted MINOP TREND endoscope clear from any blood and secretion, and optimum viewing conditions are guaranteed during transnasal access in the operational area.

The MINOP TREND (TRansnasal ENDoscopy System) can be operated freehand or via the adapter for holding arm RT099R on the corresponding holding arm, and used with a rinsing liquid container, suction pump as well as the instrument set for transnasal access.

Components

Art. no.	Designation
FH615	MINOP TREND handle (with irrigation button)
FH610R	MINOP TREND suction/irrigation sheath 0°
FH611R	MINOP TREND suction/irrigation sheath 30°
PE487A	MINOP TREND endoscope 0°, Ø 4 mm
PE507A	MINOP TREND endoscope 30°, Ø 4 mm
RT099R	Adapter for holding arm

Safe handling and preparation

- CAUTION**
- Federal law restricts this device to sale by, or on order of a physician!**
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
 - Read, follow, and keep the instructions for use.
 - Use the product only in accordance with its intended use, see Intended use.
 - Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
 - Store any new or unused products in a dry, clean, and safe place.
 - Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Replace any damaged components immediately with original spare parts.

Safe operation

The MINOP TREND handle can be combined with the following suction/irrigation sheaths:

- MINOP TREND suction/irrigation sheath 0° and MINOP TREND endoscope 0°
- MINOP TREND suction/irrigation sheath 30° and MINOP TREND endoscope 30°

When irrigation button 4 is pressed, the field of view of the working MINOP TREND endoscope 5 is rinsed and thus cleansed of any blood and secretion.

- Rinse MINOP TREND endoscope 5: Push through irrigation button 4 to the positive stop.
Suction is interrupted while irrigation is active.
- Stop rinsing MINOP TREND endoscope 5: Release irrigation button 4.
The valve automatically switches to continuous suction and interrupts the irrigation.

Note

Suction is permanently active, except for the short periods when irrigation is activated.

Disassembling

Loosening MINOP TREND suction/irrigation sheath from MINOP TREND handle

- Turn the MINOP TREND suction/irrigation sheath 2 with the MINOP TREND endoscope used 5 and, at the same time, pull backwards out of the MINOP TREND handle 3.

Disconnecting the tubes and the irrigation button

- Remove the tubes from the tube connectors.
- Unscrew irrigation button 4, see Fig. 1.

Disconnecting adapter for holding arm RT099R

- Undo the cap nut.
- Remove adapter for holding arm RT099R.

Removing MINOP TREND endoscope out of the MINOP TREND suction/irrigation sheath

- Undo the light cable from the MINOP TREND endoscope 5.
- Take the camera from the eyepiece of the MINOP TREND endoscope 5.
- Unlatch the locking lever of the MINOP TREND suction/irrigation sheath 2 and remove the MINOP TREND endoscope 5.

Assembling

Inserting MINOP TREND endoscope into the MINOP TREND suction/irrigation sheath

Damage to components caused by excessive force!

- Handle the components with appropriate care.

- Insert MINOP TREND endoscope 5 axially into the MINOP TREND suction/irrigation sheath 2, see Fig. 2. Whilst completing this step, ensure that the lever of the MINOP TREND suction/irrigation sheath 2 is open and the peg of fits into the groove on the bayonet of the MINOP TREND endoscope 5.
- Turn the locking lever of MINOP TREND suction/irrigation sheath 2 as far as it will go.
The MINOP TREND endoscope 5 has been securely locked.

Connecting the tubes and the irrigation button

Restricted functionality of the MINOP TREND (TRansnasal ENDoscopy System) if the irrigation button is not screwed in completely!

- Make certain that irrigation button is screwed into MINOP TREND handle down to the stop.

Restricted functionality of the MINOP TREND (TRansnasal ENDoscopy System) due to incorrect connecting direction of the tubes!

- Note the arrow on the MINOP TREND handle.
- Carry out function checks before using the product.

- Verify that irrigation button 4 does not show any damage (valve seals, surfaces, etc.)
- Screw the irrigation button 4 all the way into the MINOP TREND handle 3, see Fig. 3.
- Connect the suction tube to the suction connector, making certain that the arrow points away from MINOP TREND handle 3.
- Connect the irrigation tube to the irrigation connector, making certain that the arrow points towards the MINOP TREND handle 3.

Note

In order to facilitate connection, a MINOP TREND suction/irrigation tube is available, the irrigation tube of which is designated along its entire length with a blue line running through it. The suction tube is not marked.

Connecting the MINOP TREND suction/irrigation sheath with the MINOP TREND handle

- Insert the MINOP TREND suction/irrigation sheath 2 with the MINOP TREND endoscope 5 from the rear into the MINOP TREND handle, until the MINOP TREND suction/irrigation shaft 2 engages, see Fig. 4.
MINOP TREND suction/irrigation sheath 2 is rotatable.
- Check the irrigation and suction tubes and the conical seating for any leaks.
- Press irrigation button 4 to check that it functions properly.

Mounting adapter for holding arm RT099R

- Undo the cap nut for adapter for holding arm 6.
- Insert adapter for holding arm 6 from either side (left/right) into the borehole until the flange engages in the appropriate groove, see Fig. 5.
- Screw on the cap nut and tighten it until adapter for holding arm 6 is fixed, see Fig. 6.
The MINOP TREND (TRansnasal ENDoscopy System) is ready for connection to the the holding arm.

Validated reprocessing procedure

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.extranet.bbraun.com

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- ▶ Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- ▶ If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- ▶ Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- ▶ Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- ▶ Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

DANGER

Risk to patients!

- ▶ The product must only be cleaned mechanically!

WARNING

Inadequate cleaning result if the irrigation button remains screwed in!

- ▶ Remove the irrigation button prior to cleaning.

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - be approved for plastic material and high-grade steel),
 - do not attack softeners (e.g. in silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum allowable temperature of 55 °C.

Note
The endoscopes have a separate set of instructions for use, where the applicable validated cleaning procedure can be found.

- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- ▶ Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- ▶ Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Mechanical alkaline cleaning and thermal disinfection ■ MINOP TREND handle part ■ MINOP TREND suction/irrigation shafts	■ Cleaning brush PM995200 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfecting and sub-chapter: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ Adapter for holding arm ■ Irrigation button	■ Cleaning brush PM995200 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Mechanical cleaning/disinfecting

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	■ Concentrate, alkaline: – pH = 13 – <5 % anionic surfactant ■ 0.5 % working solution – pH = 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disin- fecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	–	D–W	–

D–W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	–	D–W	–

D–W: Drinking water
RT: Room temperature
*Recommended: BBraun Stabimed

- ▶ Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- ▶ If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- ▶ Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- ▶ Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- ▶ Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Concentrate, alkaline: - pH = 13 <5 % anionic surfactant 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disin- fecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water
FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

► Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note
The product can be sterilized either in disassembled or in assembled condition.

► Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).

► Validated sterilization process

- Steam sterilization using fractional vacuum process
- Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
- Sterilization using fractional vacuum process at 134 °C/holding time 5 min

► When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

► Do not modify the product.

► For service and repairs, please contact your national B. Braun/Aesculap agency.
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.
Modifications carried out on medical technical equipment may result in loss of warranty rights and applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Or in the US:
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Art. no.	Designation
FF357R	MINOP TREND tray for trocars and optics
FH605SU	MINOP TREND irrigation/suction tube
FH615801	Irrigation button
PM995200	Cleaning brush

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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