GB Instructions for use/Technical description

MINOP Intraventricular trocars

USA Instructions for use/Technical description

MINOP Intraventricular trocars

Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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


Aesculap®
MINOP Intraventricular trocars

Legend

- 1 MINOP trocar (operative), FF399R
- 2 Single-obturator (for working channel)
- 3 MINOP trocar (diagnostic), FF398R
- 4 MINOP trocar (diagnostic), FF397R
- 5 Locking cap (obturator)
- 6 Obturator (for irrigation/drainage channels)
- 7 Obturator (for irrigation/drainage channels)
- 8 Guideway
- 9 Sealing caps
- 10 Obturator (optics and working channel)
- 11 Release button
- 12 Irrigation or drainage channel
- 13 Obturator (for optics channel)
- 14 Endoscope
- 15 Electrode
- 16 Instrument

Symbols on product and packages

Caution, general warning symbol
Caution, see documentation supplied with the product

Intended use

The MINOP trocars are used for endoscopic examinations and therapies in the central nervous system.

Indications

Indications, see Intended use.

Note
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Contraindications

None known.

Risks and side effects

As part of the legal duty to inform, the following typical risks and side effects associated with the use of surgical instruments are referred to. These are predominantly procedure-specific, non product-specific and not limited to unwanted damage to surrounding tissue resulting in e.g. bleeding, infection, material incompatibilities or instrument parts remaining unnoticed in the patient, etc.

Available sizes

Designation	Art. no.
MINOP trocar (diagnostic)	FF397R
MINOP trocar (diagnostic)	FF398R
MINOP trocar (operative)	FF399R

Safe handling and preparation

CAUTION
Federal law restricts this device to sale by, or on order of a physician!

Note
For the MINOP trocar system, there are specially developed Aesculap endoscopes available. For reasons of compatibility, the MINOP trocars may only be used with these special Aesculap endoscopes!

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.
- ▶ For monopolar contact coagulation, make certain that the output power at the HF device does not exceed 50 W.
- ▶ To avoid burns when using the ventriculoscope in combination with HF electrodes, ensure that HF current is activated only under visual control.
- ▶ When using the HF electrode in combination with a flexible instrument, ensure that the working tips are kept as far apart as possible.
- ▶ Set the HF device to the lowest possible output power required for the intended effect.
- ▶ Prior to each use, visually inspect the product for: damage or surface changes to the insulation.
- ▶ Never place the product on or next to the patient.
- ▶ Follow the instructions for use of the HF device.

Safe operation

WARNING

Risk of injury and/or malfunction!
▶ Always carry out a function check prior to using the product.

WARNING

Risk of injury when using the product beyond the field of view!
▶ Apply the product only under visual control.

WARNING

Risk of injury caused by protruding 30°-endoscope in locked position!
▶ Insert the endoscope only when the trocar is in its final position.

WARNING

Risk of burns due to high temperature of the instrument-side end of the optical cable of a light source!
▶ Apply proper care when operating the light source.

CAUTION

Damage to the endoscope due to incorrect handling or operation!
▶ Insert the endoscope in the MINOP trocar only if the trocar is not deformed, bent or kinked.

Note
The endoscopes offered for application with the MINOP system must be used with a light source equipped with a spare bulb!

Mounting the camera

Note
Mount the camera only with the endoscope locked in the working trocar!

- ▶ Lock the endoscope in the MINOP trocar.
- ▶ Mount the camera.

Demounting the camera

Note
Extract the endoscope from the working trocar only after the camera has been demounted from the endoscope!

- ▶ Demount the camera from the endoscope.
- ▶ Extract the endoscope from the MINOP trocar.

Disassembling

Removing the obturator (for optics and working channel)

- ▶ Press the release button and carefully remove the obturator from the MINOP trocar.

Removing the obturator (for irrigation/drainage channels)

- ▶ Turn the locking cap by 90° and carefully remove the obturator from the MINOP trocar.

Assembling

Inserting the endoscope/obturator

- ▶ Carefully insert the endoscope/obturator in the MINOP trocar, making sure that the locking pin of the endoscope/obturator is on the same side as the release button of the MINOP trocar.
- ▶ Insert the endoscope/obturator down to the positive stop.

Inserting the obturators (for irrigation/drainage channels)

Note
Be careful not to interchange the obturators. Mind the markings, "L" (= left-hand obturator) and "R" (= right-hand obturator)!

- ▶ Hold the obturators at the distal end and carefully slide them into the respective channels, down to the positive stop. Ensure that the guideway on the obturator is aligned with the pin in the main body of the MINOP trocar.
- ▶ Turn the locking cap by 90°.

Validated reprocessing procedure

General safety notes

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.
The specified chemistry was used for validation.

Note
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com.
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine. Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result damage to metallic products. These must be removed by rinsing thoroughly with demineralized water and then drying. Perform additional drying, if necessary. Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes in materials, e.g., fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see [www.a-k-i.org](#), link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove the sealing cap from the Luer lock connector.

Preparations at the place of use

- If applicable, rinse surfaces that are not accessible to visible inspection (preferably with demineralized water), using a disposable syringe, for example.
- Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

DANGER

Danger to the patient!

- Reprocess the product only with manual pre-cleaning followed by mechanical cleaning.

DANGER

Risk to patient due to cross contamination!

- Do not clean contaminated products together with uncontaminated products in a tray.

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
 - are approved for (e.g. aluminum, plastics, high-grade steel),
 - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable temperature of 55 °C.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

Manual cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfecting machine must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Ultrasonic cleaning	RT (cold)	>15	1	D-W	Dr. Weigert neodisher® SeptoClean®
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

Phase I

- Clean the product in an ultrasound cleaning bath (frequency of 35 kHz) for at least 15 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush (e.g. TE654202) in the solution until all visible residues have been removed from the surface.
- If applicable, brush non-visible surfaces for at least 1 minute with a suitable cleaning brush (e.g. TE654202).
- Mobilize non-rigid components, such as set screws and hinges, during cleaning.
- Then flush these areas thoroughly at least five times with the cleaning disinfectant solution using a disposable syringe (20 ml).

Phase II

- Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws and hinges, during rinsing.

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	Dr. Weigert neodisher® SeptoClean 1 %* working solution
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine

D-W:

Drinking water

FD-W:

Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

- Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots)
- Place instruments in the tray with their hinges open.
- Connect components with lumens and channels directly to the rinsing port of the injector carriage.
- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note
The product can be sterilized either in disassembled or in assembled condition.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Disassemble the product
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractional vacuum process at 134 °C, holding time 18 minutes for prion inactivation
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

Risk of injury and/or malfunction!
► Do not modify the product.

► For service and repairs, please contact your national B. Braun/Aesculap agency.
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042

Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Art. no.	Designation
EJ751200	MINOP sealing caps for side channels (black)
PE184A	0°-optics Ø 2.7 mm/0°/180 mm
PE204A	30°-optics Ø 2.7 mm/30°/180 mm
PF893800	MINOP cleaning brush

Disposal

► Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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