Neuro-Patch®

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Neuro-Patch®

Description

The new product is a finely fibrillar, microporous, nonwoven material (Neuro-Patch®) manufactured from a special high purity polyesterurethane.

The structure of the Neuro-Patch® is characterized by "open" micropores on the surface of the nonwoven material, that make the rapid immigration of cells appreciably easier. The Neuro-Patch® is extraordinarily well tolerated by tissue and is completely biostable.

Composition

The Neuro-Patch® consists of high purity aliphatic polyester urethane, ETO-sterilized.

Indications

Neuro-Patch® is employed in neurosurgery as dura mater replacement:

- for covering cerebral and cerebellar dura defects
- for cerebral decompression surgery when there is elevated intracranial pressure
- for covering spinal dura defects
- for spinal decompression surgery

Mode of action

The process of healing is primarily influenced by the porosity of the implant. When the structure is microporous the large surface area frequently leads to foreign body reactions with giant cells, which then prevent the inward growth of connective tissue.

The microporous structure present in Neuro-Patch® allows the immediate inward migration of fibroblasts and the patch is anchored into the tissue by the collagen that is secreted; there is no infiltration of foreign body giant cells. At the implant there is no collection of cells, such as lymphocytes, eosinophile and neutrophile granulocytes or mononuclear macrophages, which could indicate a chronic inflammation reaction or a hyperallergic reaction. Neomembranes are produced which surround the PUR patch (Neuro-Patch®) used as a permanent replacement. These neomembranes remain very thin and do not form any material-dependent adhesions with the brain.

Long-term stability

Long-term studies have revealed that there is no loss of substance and no changes in physical or chemical properties after 6 years of implantation. This result is confirmed by histological examination.

Contra-indications

- implantation in infected regions
- when there is open cerebrocranial traumata
- when there is open spina bifida

Precautionary measures in application

See under mode of application.

Interactions with other drugs

None known.

Warning notes

Neuro-Patch® should not be used together with bone cement as the patch material could be damaged depending on the application.

Dosage instructions and mode of application

If not otherwise prescribed a Neuro-Patch® is chosen to fit the application area and cut to fit the size of the defect.

In order to embed with as little stress as possible the Neuro-Patch® should be cut down to fit the defect as closely as possible. The Neuro-Patch® should be fixed with nonabsorbable suture material (polyester, polypropylene); then further sutures can be added with absorbable or nonabsorbable suture material. Then a seal can be made with fibrin adhesive.

Atraumatic round-bodied needles make it possible to suture without causing great damage to the implant.

Before implantation of the Neuro-Patch® the person using it must be acquainted with the surgical technique required, the specific applications and the properties of Neuro-Patch® in vivo.

Side effects

Possible complications include infection, cerebrospinal fluid leakage, adhesion and foreignbody reaction.

Notes and details of stability

Neuro-Patch® should not be employed after the expiry date given.

Neuro-Patch® should be stored at 25 ± 5 °C.

Other information

The Neuro-Patch® may only be employed if the packing is undamaged. Opened packs of Neuro-Patch® and pieces of implant, that are no longer required, must not be used later. Resterilization must not be carried out, because the structure of the Neuro-Patch® and, hence, its behaviour in vivo can be detrimentally affected.

Neuro-Patch® may not be resterilized.

Product range

Packs of 1 piece		Packs of 2 pieces	
Size	REF	Size	REF
12 x 14 cm	106 4002	4 x 5 cm	106 4045
6 x 14 cm	106 4010	2 x 10 cm	106 4053
4 x 10 cm	106 4037	1,5 x 3 cm	106 4061
6 x 8 cm	106 4029		
8 x 9 cm	106 4020		
5 x 6 cm	106 4040		

Symbols used on labelling

Do not reuse!

Use until Year + Month

Temperature limitation

Sterile unless package is opened or damaged.
Method of Sterilization - Ethylene oxide

CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.

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Neuro-Patch®

Product description

Neuro-Patch® is a microporous fabric, manufactured from specially polyester urethane and used as a dura substitute. The fine-fiber microstructure is characterized by intercommunicating pores and many openings on the surface. This allows immigration of the body's own connective tissue cells. The device is provided sterile, for single use only.

Composition

Neuro-Patch® consists of specially purified aliphatic polyester urethane. It is EtO sterilized.

Indications for use

Neuro-Patch® is indicated as a dura mater substitute in neurological procedures for soft tissue reconstruction of damaged, impaired, or missing tissue, such as for:

- Closing cerebra l and cerebellar dural defects, such as after excision of tumors, operations on the posterior cranial fossa and after the removal of basal meningiomas;
- Decompressive surgery in the presence of increased intracranial pressure;
- For closure of spinal dural defects after removal of spinal tumors and after spinal trauma; and
- Spinal decompressive surgery.

Contraindications

Do not implant in infected areas, open cranial trauma or open spina bifida.

How supplied

The product is supplied sterile, non-pyrogenic and packaged as single or double pieces. The products are available in various dimensions ranging from 1,5 x 3 cm sections to as large as 6 x 14 cm pieces. Specific sizes are noted by catalog number below.

Precautions

Before using Neuro-Patch®, the surgeon should be familiar with the surgical technique, special applications and the properties of Neuro-Patch® in vivo.

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Neuro-Patch® should only be used if the package is intact. Visually inspect the packaging; do not use product if the package is found opened, punctured, torn or tampered with as sterility may be compromised.

Do not resterilized, as the sterilization process may alter the structure of Neuro-Patch® and thus its properties in vivo.

U.S.A. Federal Law restricts this device to sale by or on the order of a physician.

Possible adverse reactions

Hypersensitivity or an immune response to Neuro-Patch® may occur. The physician should closely monitor the patient for such a reaction and treat accordingly.

Postoperative adjacent tissue adhesions in situ may occur.

Directions for use

The correct size of Neuro-Patch® should be chosen and cut as closely as possible to fit the defect.

The Neuro-Patch® should be fixed with non-absorbable suture material (polyester, polypropylene). To improve suture security, non-absorbable suture material (polyester, polypropylene) is recommended. The use of atraumatic round-bodied needles will allow suturing without significant damage to the implant.

Storage

Store the Neuro-Patch® box in a clean, dry and protected area with protection from extreme changes in temperature and humidity. Storage conditions should not exceed temperature of 25 ± 5 °C. (68 °F – 86 °F)

Product range

Packs of 1 piece		Packs of 2 pieces	
Dimensions	REF	Dimensions	REF
12 x 14 cm	106 4002	4 x 5 cm	106 4045
6 x 14 cm	106 4010	2 x 10 cm	106 4053
4 x 10 cm	106 4037	1,5 x 3 cm	106 4061
6 x 8 cm	106 4029		
8 x 9 cm	106 4020		
5 x 6 cm	106 4040		

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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0123 CE-mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.

Batch Number

See Instructions for Use

Cat. No.

DIM Size