
Aesculap Neurosurgery

-
- (GB)** **Instructions for use/Technical description**
PaediScope: Endoscope PF010A, Eyepiece PF011A
 - (USA)** **Instructions for use/Technical description**
PaediScope: Endoscope PF010A, Eyepiece PF011A
 - (D)** **Gebrauchsanweisung/Technische Beschreibung**
PaediScope: Endoskop PF010A, Okular PF011A
 - (F)** **Mode d'emploi/Description technique**
PaediScope: Endoscope PF010A, oculaire PF011A
 - (E)** **Instrucciones de manejo/Descripción técnica**
PaediScope: Endoscopio PF010A, ocular PF011A
 - (I)** **Istruzioni per l'uso/Descrizione tecnica**
PaediScope: Endoscopio PF010A, oculare PF011A

PaediScope: Endoscope PF010A, Eyepiece PF011A

Legend

- 1 Endoscope
- 2 Aesculap-Storz fiberoptic light connection
- 3 Eyepiece
- 4 Eyepiece window
- 5 Eyepiece cup
- 6 Focusing ring
- 7 Eyepiece cover class
- 8 Fiberoptic light guide end of endoscope
- 9 Image channel
- 10 Light bundle
- 11 Irrigation and overflow channel
- 12 Working channel
- 13 Endoscope tip
- 14 Luer lock connectors for irrigation and overflow channel
- 15 Silicone sealing cap

Symbols on product and packages

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product
	Classification Type BF
	Date of manufacture

Intended use

The PaediScope is used for intraventricular procedures in neurosurgery. Its small size makes it particularly suitable for pediatric indications.

With the endoscope **1** and the eyepiece **3** assembled, the operating environment at the endoscope tip is displayed in the eyepiece window **4**. Instruments can be inserted into the endoscope through the working channel.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information regarding such training, please contact your national B. Braun/Aesculap agency.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
 - Read, follow and keep the instructions for use.
 - Use the product only in accordance with its intended use, see Intended use.
 - Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
 - Store any new or unused products in a dry, clean and safe place.
 - Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Replace any damaged components immediately with original spare parts.
 - To avoid damage to the working end: Exercise caution when passing products through the working channel (e.g. trocar).
 - Replace the product if its image quality is degrading.
 - Do not, under any circumstances, put down the product on the patient or on the surgical drape covering the patient.
 - To avoid burns caused by heating up of the distal end of the product: Prior to putting down the product, turn down the power of the light source.
- To avoid HF burns:
- The product's working tip must be in the user's field of vision whenever HF power is activated.
 - Prior to activating the HF device, make certain that the working end of the product is not touching any electrically conductive accessories.
 - Prior to each use, visually check the product for: damage to, or surface changes on, the insulation.
 - Deactivate the automatic switch-on mode of the HF device when using accessories for endoscopic or laparoscopic applications.
 - Follow the instructions for use of the HF device.

Safe operation

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

Function checks

- Check that the image projected on the eyepiece window **4** is clear and distinct.
- Hold the endoscope with the Aesculap-Storz fiberoptic light connection **2** against a light source and check that the four points of the light bundle **10** at the endoscope tip **13** are evenly illuminated.
- Ensure that the eyepiece window **4**, the window of the image channel **9** at the endoscope tip **13** and the input surface at the Aesculap-Storz fiberoptic light connection **2** are not dull or dirty.

Use of the endoscope

WARNING

Danger to patients due to inappropriate use of the irrigation channel!

- Feed the irrigation channel of the endoscope only by a passive gravitational (drip) system.
- Set the height difference of the gravitational system to not more than 50 cm.
- Do not kink the flexible parts of the endoscope **1**.
- Do not put heavy objects on top of the endoscope **1**.
- Insert endoscope **1** with its distal tip, avoiding any strong impact.
- Turn the focusing ring **6** until the eyepiece window **4** shows a sharp image.
- For improved vision: Feed the irrigation channel **11** of the endoscope by a passive gravitational (drip) system without pressure sleeve.
 - Do not employ active pumps.
 - Allow liquid to drain through the overflow channel, against atmospheric pressure.
 - Do not employ suction devices or vacuum pumps.
 - Ensure that the overflow channel **11** is not closed during irrigation.
- Push additional instruments through the silicone sealing cap **15** in order to insert them into the working channel **12** of the endoscope **1**.
- To remove instruments, pull the latter out of the working channel **12** at the silicone sealing cap **15**.
- To avoid burns caused by heating up of the distal end of the product: Adjust the light source to the minimum required power for optimal illumination of the endoscopic image.

Using the endoscope with additional instruments

The working channel **12** allows instruments being inserted and applied with the endoscope **1**.

WARNING

Risk of burns from the hot tips of fiberoptic bundles on the instrument side!

- Make certain the instrument-side end of the fiberoptic bundle does not touch human tissue or highly combustible materials while the light source is active.
- Follow the user instructions of any additional instruments used.
- Only use instruments that can be easily inserted into the working channel **12** (maximum outer diameter of the instrument used: 1 mm or 3 french).
- When using the endoscope with an electromedical device, adhere to rules for the use of high-frequency devices.
- When using HF electrodes: only activate HF electrodes if their insulation is visible in the eyepiece window **4**.

Disassembling

- Remove the fiberoptic bundle from the Aesculap-Storz fiberoptic connection **2**.
- Remove the tube from the irrigation and overflow channel **11**.
- Rotate the eyepiece **3** counterclockwise and remove it from the endoscope **1**.
- Remove the silicone sealing cap **15** from the endoscope **1**.
- Remove the covers from the Luer lock connectors **14**.

Assembling

- Make certain that the tip of the fiberoptic light guide **8** on the endoscope and the cover glass **7** on the eyepiece **3** are not soiled.
- Rotate the eyepiece **3** clockwise to screw it onto the fiberoptic light guide end **8** of the endoscope until there is palpable, slight resistance. Then rotate the eyepiece **3** by an additional one-quarter clockwise turn.

The threaded connection between the endoscope and eyepiece is now watertight.
- Attach the irrigation and overflow tubes at the Luer lock connectors **14** of the irrigation and overflow channel.
- Attach the fiberoptic bundle to the Aesculap-Storz fiberoptic connection **2**.

PaediScope: Endoscope PF010A, Eyepiece PF011A

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Therefore the time interval between application and cleaning should not exceed 6 h, no potentially fixating pre-cleaning temperatures >45 °C be applied, and no fixating disinfectants (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfecting and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals that have been tested and approved (e.g. VAH/ DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process solution is sufficient to cause visible surface changes.
- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use process chemicals that would cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.

Any gross debris can be removed manually with a soft brush. After this pre-treatment, continue through the normal steps of the processing procedure.

For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

- Use appropriate cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.
- Make certain that any bending or coiling radius of the endoscope is not less than 11 cm.
- Avoid excessive pressure on the product.

Preparations at the place of use

- Disassemble the product immediately after use, as described in these instructions for use, see Disassembling.
- Remove the covers from the Luer lock connectors and the sealing cap from the working channel.
- Rinse surfaces that cannot be visually inspected, e.g. on products with hidden gaps or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.
- Aesculap recommends putting away the product in dry condition.

Preparation prior to cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents according to the manufacturer's instructions. The cleaning/disinfecting agents must
 - be approved for flexible endoscopes,
 - not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature, exposure time and use-by date.
- Do not exceed the maximum allowable cleaning temperature of 45 °C.

CAUTION

Damage to the optical system caused by loosening of connections during ultrasound cleaning!

- Do not clean the endoscope by ultrasonic treatment under any circumstances.

CAUTION

Damage to the endoscope due to inappropriate cleaning/disinfecting agents!

- Only apply manual cleaning of the endoscope.
- Follow the manufacturer's instructions.

WARNING

If the PaediScope was used on patients with suspected CJD or the new variant CJD, the endoscope cannot be reprocessed and sterilized anymore to the required safety standards, due to its materials and design, and has to be disposed of appropriately.

Manual cleaning/disinfecting

- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.
- Keep the product immersed in the cleaning/disinfecting solution for not longer than 1 hour.

Manual cleaning with immersion disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	15	1.0	D-W	BBraun Cleaner N
II	Intermediate rinse	RT (cold)	3 x 1	-	D-W	-
III	Disinfecting	RT (cold)	15	1.5	FD-W	BBraun Helipur H Plus N
IV	Final rinse	RT (cold)	3 x 2	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: max. 10 microorganisms/ml, low endotoxin: max. 0.25 endotoxin units/ml)

RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Check that all accessible surfaces are moistened and all lumens, channels and complex geometries are bubble-free and filled with cleaning solution.
- Clean the outer surfaces of the immersed product with a soft, lint-free, disposable wipe.
- Brush through all areas not accessible to visual inspection, e.g. in products with hidden crevices, lumens (e.g. open/working channel) or complex geometries, for at least 1 min or until no more residues can be removed. Mobilize non-rigid components, e.g. set screws, joints, etc. to the positive stop, 3 times in each direction.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Clean open/working channels with the cleaning brush supplied with the product (PF893800) for at least 1 min.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).

- To clean with the brush, insert the brush with its metal tip into the channel and slowly extract it from the other end as soon as it emerges there.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Completely rinse (through) the product 3 times for at least 1 min, using clean, fresh water for each rinse. Mobilize non-rigid components, e.g. set screws, joints, etc. to the positive stop, 3 times in each direction.
- For each rinsing cycle, rinse through lumens and channels with a single-use syringe (20ml) at least 5 times.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution. Check that all accessible surfaces are moistened and all lumens, channels and complex geometries are bubble-free and filled with cleaning solution. Mobilize non-rigid components, e.g. set screws, joints, etc. to the positive stop, 3 times in each direction.

PaediScope: Endoscope PF010A, Eyepiece PF011A

Stage IV

- Completely rinse (through) the product 3 times for at least 2 min, using clean, fresh water for each rinse. Mobilize non-rigid components, e.g. set screws, joints, etc. to the positive stop, 3 times in each direction.
- For each rinsing cycle, rinse through lumens and channels with a single-use syringe (50ml) at least 5 times.
- Allow water to drip off for a sufficient length of time.

Stage V

- Thoroughly dry the external surfaces of the product with a soft, lint-free disposable wipe or with medical compressed air.
- Thoroughly blow-dry lumens and channels with air from a single-use syringe (20 ml). Do not use compressed air for this purpose.

Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- To remove any residues of cleaning/disinfecting agents, wipe the eyepiece window **4**, the window of the image channel **9** at the endoscope tip **13** and at the input surface of the Aesculap-Storz fiberoptic light connection **2** with a swab moistened with alcohol.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Assemble the separable product, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product into the appropriate storage holders or on an appropriate tray.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Make certain that any bending or coiling radius of the endoscope is not less than 11 cm.
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

Note

The product may only be sterilized in disassembled condition.

CAUTION

Risk of damage to the optical system through loosening of connections as the result of hot air or flash sterilization!

- **Do not sterilize the endoscope with dry hot air or by flash sterilization.**

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Do not sterilize damaged products.
- Make certain the product is cleaned manually.
- Sterilize the product in the appropriate storage rack (FF373R).
- Avoid knocking or shaking the product.
- Avoid temperature fluctuations after sterilization (e.g. due to cooling with water).
- Aesculap recommends sterilizing the endoscope with steam.

Note

Aesculap guarantees steam sterilization of the product up to 15 autoclave cycles.

- Validated sterilization process
 - Disassemble the product, see Disassembling.
 - Steam sterilization through fractionated vacuum process.
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665.
 - Sterilization through fractionated vacuum process at 134 °C/2 bar, holding time 5 min.
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum drying time	
		Wrapped	In a sterile container system
Pre-vacuum	270 °F–275 °F	6 min	6 min (solid or perforated bottom)

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.
- Make certain that any bending or coiling radius of the endoscope is not less than 11 cm.
- Do not store the product in its transport packaging.

Maintenance

The endoscope requires no maintenance for the first 15 autoclave cycles. After that, the image quality of the endoscope can be diminished by normal wear due to steam sterilization.

Technical Service

Risk of injury and/or malfunction!

- Do not modify the product.

WARNING

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Repairs may be carried out only by persons authorized by AESCULAP FLEXIMED. Any repairs carried out by other persons will void all warranty and guarantee entitlements.

- Prior to dispatching the product for repairs:
 - Have the product cleaned and disinfected or sterilized, and mark it as "disinfected" or "sterilized", respectively.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

AESCULAP FLEXIMED GmbH
Robert-Bosch-Strasse 1
79211 Denzlingen / Germany
Phone: +49 (0) 7666 / 9321-0
Fax: +49 (0) 7666 / 9321-580
E-mail: techserv@aesculap-fleximed.de
www.aesculap-fleximed.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Art. no.	Designation
EJ670202	Silicone sealing cap for working channel
EJ751200	Covers for Luer lock connectors
FH603SU	PaediScope insertion aid 10 Fr.
PF893800	Cleaning brush

Technical specifications

Ø Distal end	3.0 mm
Ø Working channel	1.2 mm
Ø Irrigation/Overflow channel	0.8 mm
Length of shaft	150 mm
Long flexible component	1000 mm
Image angle	80°
Viewing direction	0°
Depth of field	3 mm to ∞

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!
- Put away the product in dry condition.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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Aesculap AG