Aesculap Spine

GB	Instructions for use PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
USA	Instructions for use PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK <i>Note for U.S. users</i> <i>U.S. cleared device instructions for use can be found at</i> <i>www.aesculapusa.com "Products Et Services".</i> <i>If you do not have internet access please call your local Aesculap</i> <i>representative or customer service at 1-800-258-1946.</i>
D	Gebruiksaanweisung PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
F	Mode d'emploi PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
E	Instrucciones de manejo PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
I	Istruzioni per l'uso PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
P	Instruções de utilização PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
RUS	Инструкция по применению PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
CZ	Návod k použití PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
PL	Instrukcja użytkowania PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
SK	Návod na používanie PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
CN	使用说明书 腰椎后路椎间融合器： A-Space PEEK， CeSpace PEEK， ProSpace PEEK， T-Space PEEK
KR	사용방법및사용시주의사항 추간체고정보형체 CESPAC PEEK CERVICAL IMPLANT, 형명개별기재)

B | BRAUN
SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 012095 11/10 Änd.-Nr. 32607/40431

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Símbolo CE, em conformidade com a Directiva 93/42/CEE
Маркировка CE согласно директиве 93/42/ЕЭС
CE-označení podle směrnice 93/42/EHS
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
Označenie CE podľa smernice 93/42/EHS
根据 93/42/EWG 指令的 CE 标志

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Мы оставляем за собой право вносить технические изменения
Technické zmeny vyhradeny
Zmiany techniczne zastrzeżone
Technické zmeny vyhradené
保留技术变动的权力

PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK

Intended use

PEEK cages are used for the following purposes:

- A-Space PEEK:
 - Stabilization of the lumbar spine through anterior, lateroanterior or lateral approach, mono- and multisegmental
- CeSpace PEEK:
 - Stabilization of the cervical spine C2 -T1 through anterior approach mono- and multisegmental
- ProSpace PEEK:
 - Stabilization of the lumbar and thoracic spine through posterior approach, mono- and multisegmental
- T-Space PEEK:
 - Stabilization of the lumbar and thoracic spine through transforaminal approach, mono- and multisegmental

Materials

The materials used in the implant are listed on the package. The PEEK cages are made of PEEK-OPTIMA®. They contain marker pins made of tantalum to ensure radiological visibility for inspecting the implant position. PEEK-OPTIMA® is a registered trademark of Invivo Ltd, Lancashire FY5 4QD, UK.

Indications

Surgically installed implants serve to support normal healing processes. They should neither replace normal structures of the body nor permanently bear the loads occurring in the case of incomplete healing.

Use for:

- Degenerative instability
- Spondylolisthesis
- Post-discectomy syndrome
- Post-traumatic instabilities

Contraindications

The operation should not be carried out against the following contraindications:

- Acute or chronic infections or severe defects of the osseous structures of the vertebral bodies, which need to be sound for the stable implantation of the PEEK cages
- Bone tumors in the region of the implant anchoring
- Unwillingness or inability of the patient to follow the instructions for postoperative treatment
- Any medical or surgical condition that could preclude the potential success of the implantation
- Pregnancy
- Osteoporosis or similar bone density loss
- Systemic or metabolic illnesses
- Drug abuse or alcoholism
- Generally poor condition of the patient
- Adiposity
- Psychosocial issues; lack of co-operation by the patient
- All cases that are not listed under indications

Risks

The surgical intervention involves the following potential risks:

- Neurological complications caused by overdistraction or trauma of the nerve roots or dura
- Loss of intervertebral disk height due to removal of healthy bone material

Complications that can generally occur in connection with intervertebral surgery:

- Pseudarthrosis
- Incorrect implant position
- Spondylolisthesis
- Loss of fixation; dislocation or migration

Side-effects and adverse interactions

None known.

Safety notes

- The implant components were tested and approved for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Only combine Aesculap modular implant components with each other.
- Do not use damaged or surgically excised components under any circumstances.
- Implants that have already been used must not be reused.
- Delayed healing can result in implant rupture due to material fatigue.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation (min. dose 25 kGy).
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiry date or whose packaging is damaged.

Implants supplied in sterile condition must not be resterilized under any circumstances. If the present instructions for use do not offer explicit information in this respect, the product user bears all responsibility for any resterilization.

- Danger to the patient and possible loss of implant functionality due to resterilization!**
- Do not resterilize and reuse PEEK implants under any circumstances.
 - Never resterilize used, explanted or damaged implants.

Application

The operating surgeon draws up an operating plan that specifies and appropriately documents the following steps:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating team are aware of information concerning the operating technique and the range of implants and associated instruments; this information is complete and ready to hand.
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific publications by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The intervention has been explained to the patient, whose consent concerning the following information has been documented:

- In the case of delayed or incomplete fusion, the implants can break and loosen due to high loads.
- The life-span of the implant depends on the patient's body weight.
- Corrective surgery may become necessary if the implant loosens.
- The patient must undergo regular check-ups of the implant components, performed by a physician.

Implanting the PEEK cages

- Increased risk of migration due to over-preparation of the vertebral body endplates!**
- When preparing the implant bed, make certain that the base and cover plates of the adjacent vertebral bodies are not weakened.

- Select the appropriate PEEK implant size according to the indication, the preoperative planning and the bone situation found intraoperatively.
- Correctly apply the preparation instruments (rasps, curettes and chisels) for preparing the implant bed, as well as the implantation instrument.
- Apply appropriate care when inserting the implant.

A-Space PEEK

- Check implant height and angle using the trial implants.

- Damage to the implant thread!**
- Keep to the thread axis when screwing the implant onto the insertion instrument.
 - Screw in the implant as far as it will go so that the dihedron of the instrument rests in the groove on the implant.

CeSpace PEEK

- Determine the appropriate implant size, using the trial implants.
- Apply the insertion instrument correctly. Pay attention to any inscriptions/labels.

Note

If not enough stability is produced, additional fixation may be necessary.

ProSpace PEEK

- Always implant two implants per layer (PLIF technique).

T-Space PEEK

- Determine the implant height, using the trial impacts.

Note

ProSpace PEEK and T-Space PEEK should always be used with internal fixation devices!

Further information on B. Braun/Aesculap implant systems is always available from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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