- GB

Instructions for use
- USA

PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
- D

Gebrauchsanweisung
PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
- F

Mode d'emploi
PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
- E

Instrucciones de manejo
PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
- I

Istruzioni per l'uso
PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK
- P

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PEEK cages: A-Space PEEK, CeSpace PEEK, ProSpace PEEK, T-Space PEEK

Intended use

PEEK cages are used for the following purposes:

- A-Space PEEK: stabilization of the lumbar spine through anterior, lateroanterior, or lateral approach, monosegmental and multisegmental
- CeSpace PEEK: stabilization of the cervical spine C2-T1 through anterior approach, monosegmental and multisegmental
- ProSpace PEEK: stabilization of the lumbar and thoracic spine through posterior approach, monosegmental and multisegmental
- T-Space PEEK: stabilization of the lumbar and thoracic spine through transforaminal approach, monosegmental and multisegmental

Materials

The materials used in the implant are listed on the packaging:

- PEEK-OPTIMA®
 - Tantal for the marker pins used to ensure radiological visibility for inspecting the implant position
- PEEK-OPTIMA® is a registered trademark of Invivio, Ltd Lancashire FY5 4QD / UK.

Indications

Surgically installed implants serve to support normal healing processes. They are not intended for use either as replacements for natural body parts or to bear loads over the long term if healing does not occur.

Use for:

- Degenerative instability
- Spondylolisthesis
- Post-discectomy syndrome
- Post-traumatic instabilities

Contraindications

Do not use in the presence of:

- Fever
- Acute or chronic infections or severe defects of the osseous structures of the vertebral bodies, which need to be sound for the stable implantation of the PEEK cages
- Bone tumors in the region of implant fixation
- Unwillingness or inability of the patient to follow the instructions for postoperative treatment
- Any medical or surgical condition that could impair the potential success of the implantation
- Pregnancy
- Osteoporosis or similar bone density loss
- Systemic or metabolic diseases
- Drug abuse or alcoholism
- Generally poor condition of the patient
- Adiposity
- Psychosocial issues; lack of cooperation by the patient
- Any conditions not listed under Indications

Side effects and interactions

The application or improper use of this system entails the following risks:

- Neurological complications caused by overdistraction or trauma of the nerve roots or dura
- Loss of intervertebral disk height due to removal of healthy bone material

Complications that can generally occur in connection with intervertebral surgery:

- Pseudarthrosis
- Incorrect implant position
- Spondylolisthesis
- Loss of fixation; dislocation, or migration

Safety notes

CAUTION

Federal law restricts this device to purchase by, or on instruction by a physician!

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Delayed healing can cause implant breakage due to metal fatigue.
- The attending physician shall make any decision with regard to the removal of implant components that have been used.
- Damage to the load-bearing structures of the implant can lead to loosening of components, dislocation, migration, and other severe complications.

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.

WARNING

Damage to implants caused by processing and resterilization!

- Do not reprocess or resterilize the implants.

Application

WARNING

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information about product training, please contact your national B. Braun/Aesculap agency.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- In the case of delayed or incomplete fusion, the implants can break and loosen due to high loads.
- The life-span of the implant depends on the patient's body weight.
- The implant components must not be overloaded by extreme strains, hard physical labor or sports.
- Corrective surgery may be necessitated by implant loosening, fracture or loss of correction.
- Smokers present an increased risk of bone fusion failure.
- The patient must undergo medical check-ups of the implant components at regular intervals.

Implantation of the PEEK cages

WARNING

Increased risk of migration due to over-preparation of the vertebral body endplates!

- When preparing the implant bed, make certain that the endplates of the adjacent vertebral bodies are not weakened.

- Select the appropriate PEEK implant size according to the indication, the preoperative planning, and the bone situation found intraoperatively.
- Correctly apply the preparation instruments (rasps, curettes, and chisels) for preparing the implant bed, as well as the implantation instrument.
- Take appropriate care when inserting the implant.

A-Space PEEK

- Check implant height and angle using the trial implants.

CAUTION

Damage to the implant thread!

- Keep to the thread axis when screwing the implant onto the insertion instrument.
- Screw in the implant as far as it will go so that the dihedral of the instrument rests in the groove on the implant.

CeSpace PEEK

- Determine the appropriate implant size, using the trial implants.
- Apply the insertion instrument correctly. Pay attention to any inscriptions/labels.
- Use a cervical plate for additional stabilization.

ProSpace PEEK

CAUTION

Damage to the connector on the implant!

- When inserting the implant into the intervertebral space, avoid canting and take care to maintain an alignment parallel to the endplate.

- Always implant two implants per layer (PLIF technique).
- Always use ProSpace PEEK in conjunction with an internal fixator.

T-Space PEEK

- Determine the implant height, using the trial implants.
- Always use T-Space PEEK in conjunction with an internal fixator.

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the relevant B. Braun/Aesculap office.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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