en	Instructions for use/Technical description Plasmafit® Dual Mobility hip endoprosthesis components
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung Plasmafit® Dual Mobility Hüftendoprothesenkomponenten
fr	Mode d'emploi/Description technique Composants d'endoprothèse de hanche humaine Plasmafit® Dual Mobility
es	Instrucciones de manejo/Descripción técnica Componentes endoprotésicos de cadera Plasmafit® Dual Mobility
it	Istruzioni per l'uso/Descrizione tecnica Componenti per protesi d'anca Plasmafit® Dual Mobility
pt	Instruções de utilização/Descrição técnica Componentes de endopróteses de quadril Plasmafit® Dual Mobility
nl	Gebruiksaanwijzing/Technische beschrijving Plasmafit® Dual Mobility endoprothetische heupcomponenten
sv	Bruksanvisning/Teknisk beskrivning Plasmafit® Dual Mobility höftendoproteskomponenter
ru	Инструкция по применению/Техническое описание Plasmafit® Dual Mobility компоненты эндопротеза тазобедренного сустава
cs	Návod k použití/Technický popis Komponenty kyčelních endoprotéz Plasmafit® Dual Mobility
pl	Instrukcja użytkowania/Opis techniczny Komponenty endoprotezy biodrowej Plasmafit® Dual Mobility
sk	Návod na použitie/Technický opis Komponenty endoprotéz bedrového kĺbu Plasmafit® Dual Mobility
tr	Kullanım Kılavuzu/Teknik açıklama Plasmafit® Dual Mobility Kalça Endoprotez Bileşenleri

Aesculap® Plasmafit® Dual Mobility hip endoprosthesis components

1. Intended use

- The implant is used
- As a subcomponent of a human hip endoprosthesis: Dual Mobility insert and Dual Mobility head
 - to combine with Aesculap hip endoprosthesis components
 - Use with Aesculap Plasmafit® Plus or Plasmafit® Revision cup systems
- The implant components are packed individually; each has a separate article number.

2. Combination possibilities

12/14 cone

Head type	Ø	Neck length	CoreHip®	Prevision®	Excia® T	Trilliance®
Standard prosthesis head	22.2 mm	M	✓	✗	✗	✗
	28 mm	S	✓	✓	✓	✓
		M	✓	✓	✓	✓
		L	✓	✗	✗	✗
Dual Mobility prosthesis head	28 mm	L	✓	✓	✓	✓

✓

Combination possible

✗

Combination not possible

For all stems with 8/10 or 12/14 cone that are not listed, a Dual Mobility prosthesis head must be used.

► Follow the appropriate instructions for use.

3. Materials

- The materials used in the implant are listed on the packaging:
- Dual Mobility insert:
 - ISODUR®; cobalt-chromium forged alloy CoCrMo according to ISO 5832-12
 - Ceramic multi-layer coating made from chromium nitride-chromium-carbonitride-zirconium nitride
 - Dual Mobility head:
 - Vitelene® UHMWPE-XE (Vitamin-E-stabilized, highly crosslinked polyethylene)
- ISODUR® and Vitelene® are trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Note
The ceramic multilayer coating can undergo changes to the surface color. Such changes in color have no effect on the function or properties of the coating.

4. Indications

- The indications fundamentally depend on the Aesculap cup system used.
- Follow the appropriate instructions for use.
- Can furthermore be used for:
- chronic joint instabilities

5. Containdications

5.1 Absolute contraindications

- The absolute contraindications fundamentally depend on the Aesculap cup system used.
- Follow the appropriate instructions for use.
- Furthermore not usable in the case of:
- breakage of ceramic implant components

⚠ WARNING
Risk of wear and tear due to ceramic fragments inside the joint following revision surgery that involved the breakage of ceramic components!
► Do not implant any Dual Mobility components (danger of third-body wear).

5.2 Relative contraindications

- The relative contraindications fundamentally depend on the Aesculap cup system used.
- Follow the appropriate instructions for use.

6. Side effects and interactions

- The side effects and interactions fundamentally depend on the Aesculap cup system used.
- Follow the appropriate instructions for use.

7. Safety notes

⚠ WARNING
The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Combine modular implant components only with suitable Aesculap hip endoprostheses.

- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

8. Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The Dual Mobility inserts are radiation-sterilized.
- The Dual Mobility heads are Eo-sterilized (ethylene oxide).
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.

⚠ WARNING
Damage to implants caused by processing and resterilization!
► Do not process or resterilize the implants.

9. Application

- The operating surgeon shall devise an operation plan that specifies and accurately documents the following:
- Selection of the implant components and their dimensions
 - Location of intraoperative landmarks
- The following conditions must be fulfilled prior to application:
- All requisite implant components are ready to hand.
 - All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
 - The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
 - The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
 - The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

- The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:
- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
 - The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
 - The prosthetic joint can come loose due to excessive load, wear and tear, or infection.
 - The lifespan of the prosthetic joint depends on the patient's body weight and the strain to which the joint is subjected.
 - The prosthetic joint must not be overstrained by extreme loads, hard physical labor or sporting activities.
 - Corrective surgery may become necessary if the implant loosens.
 - In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
 - The revision of a hip endoprosthesis stem is a complex joint replacement operation.
 - Generally, the revision joint replacement will be inferior in its performance to the primary joint replacement.
 - The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.

9.1 Implanting the socket and stem components

- Cup and stem components are implanted according to the respective Aesculap instructions for use and the appropriate operation room manual.
- Before deciding on the implantation of the Dual Mobility components: install a plastic trial insert into the cup and implant the stem component.
 - Before inserting the implants, carry out a test reduction and examine joint mobility, joint stability and leg length.
 - Ensure that the connection dimensions of the Dual Mobility components, identified by letters, correspond with the cup.
 - Ensure that the diameter of the prosthesis head corresponds with the inside diameter of the Vitelene® Dual Mobility head.

⚠ WARNING
Damage to implant components due to loose anchoring screws or central sealing cap!
► Ensure that, if you are using anchoring screws or a central sealing cap, they are securely tightened.

- Before inserting the Dual Mobility insert, carefully clean and dry the internal cone of the cup.
- Carefully remove all foreign matter such as pieces of tissue, or bone or cement chips, especially any such material found around the cone clamp.
- Insert Dual Mobility insert into the cup and joint with plastic element with a hammer strike on the insertion instrument.
- Insert Dual Mobility insert into the cup only once.
- Examine the seating of the Dual Mobility insert in the cup. The insert must connect flush to the cup edge.

Note
Apply proper care and attention when inserting the Dual Mobility insert into the cup. After inserting, examine the exact seating of the insert using a finger.

- Perform joining process:
 - Connect Vitelene® Dual Mobility head with a compatible prosthesis head using the joining instrument.
 - After the joining process, examine whether the prosthesis head can move freely in the Vitelene® Dual Mobility head.
- Before inserting the joined implant components, carefully clean and dry the outer cone of the stem.
- Place joined implant components on the cone of the stem and join to plastic element using a hammer strike on the insertion instrument.
- Reduce Vitelene® Dual Mobility head with insertion instrument with plastic element in the Dual Mobility insert.
- Check the mobility and stability of the joint and the leg length.

9.2 Intraoperative replacement of a Dual Mobility insert

A surgical replacement of a Dual Mobility insert can become necessary during a revision intervention. A Dual Mobility insert punched into the metal cup can be replaced with a new insert, provided any damage to the inner cone of the cup during the removal of the insert can be excluded. Specialized Aesculap instruments must be used to remove the Dual Mobility insert which release the Dual Mobility insert by multiple, pulse-like strikes on the metal cup edge.

10. Disposal

⚠ WARNING
Risk of infection due to contaminated products!
► Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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