proGAV 2.0®

ⓓ Gebrauchsanweisung | ⓖ Instructions for use
① Istruzioni per l'uso | ⓔ Instrucciones de manejo

Ⓢ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000.

A paper copy will be provided to you upon request at no additional cost.

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CONTENT

INDICATION	18
TECHNICAL DESCRIPTION	18
FUNCTION OF THE VALVE	19
PHYSICS BACKGROUND	20
SELECTING THE APPROPRIATE SHUNT	21
USING THE INSTRUMENTS	21
ADJUSTING THE ADJUSTABLE DP-UNIT	24
READING THE PRESSURE SETTING FROM AN X-RAY IMAGE	24
POSSIBLE SHUNT COMPONENTS	25
TUBE SYSTEMS	26
SURGICAL PROCEDURE	26
PREOPERATIVE VALVE TEST	27
PRESSURE-FLOW CHARACTERISTICS	28
CONTRAINDICATIONS	29
INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS	29
RE-IMPLANTATION	29
SAFETY MEASURES	29
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	29
POSTOPERATIVE VALVE TEST	29
FUNCTIONAL SAFETY	29
ADVERSE REACTION	29
STERILIZATION	28
RESTERILIZATION	29
MEDICAL PRODUCTS CONSULTANTS	30
REQUIREMENTS OF THE MDD 93/42/EEC	30
NOTE ON THE INSTRUCTIONS FOR USE	30
VARIATIONS	31

INDICATION

The *proGAV 2.0* is used for cerebrospinal fluid (CSF) drainage in the treatment of hydrocephalus.

TECHNICAL DESCRIPTION

The *proGAV 2.0* is a posture dependent hydrocephalus valve. It comprises an adjustable differential pressure unit and a fixed gravitational unit (fig. 1 and 2).

Fig. 1: *proGAV 2.0* - side view

The adjustable DP-unit is composed of a solid titanium body with a well-ried ball-cone valve (1) integrated in its proximal part. A bow spring (2) defines the opening pressure of the ball-cone valve. The pretensioning of the spring, and thus the valve opening pressure, can be adjusted by turning a rotor (3), with the valve implanted under the patient's skin. The gravitational unit contains a tantalum ball (4), which defines the opening pressure of this valve, and a sapphire ball (5), which ensures the precise closure of the valve. Underneath the silicone catheters (6), a connector is provided.


Fig. 2: Schematic cross section of the *proGAV 2.0*

FUNCTION OF THE VALVE

The opening pressure of the *proGAV 2.0* is defined by the opening pressure of the adjustable DP-unit and the opening pressure of the gravitational unit.

Horizontal position

When the patient is lying down, the gravitational unit is always open and therefore does not present any resistance to the fluid flow. (fig. 3).


Fig. 3: Gravitational unit in horizontal position

Hence, the opening pressure of the *proGAV 2.0* is defined by the adjustable DP-unit. The operational principle of the adjustable DP-unit is illustrated in figures 3 and 4.

In fig. 4a, the ball-cone valve is closed. The drainage is blocked. In fig. 4b, the adjustable DP-unit is shown in the open condition. The patient's IVP is increased and the spring force, which otherwise keeps closed the ball-cone valve, is overcome. The closing ball moves out of the cone and a gap opens to allow drainage.


Fig. 4: Adjustable DP-unit
a) closed b) open

Vertical position

When the patient moves into an upright position, in that moment the gravitational unit closes (fig. 5a). Now, additionally to the opening pressure of the adjustable DP-unit, the weight of the tantalum ball has to be exceeded (opening pressure of the gravitational unit), thus the opening pressure of the *proGAV 2.0* is significantly increased. Only when the sum of the IVP and the hydrostatic pressure exceeds the opening pressure of the *proGAV 2.0*, drainage will be possible again (fig. 5b).


Fig. 5: Gravitational unit in vertical position
a) closed, b) open

Note: During physical activity which is associated with shocks (e.g.: running), laboratory results have shown that the *proGAV 2.0*'s opening pressure can decrease temporarily around 25% to 35%.

This concerns the single valve as well as the combination with the gravitational unit. Principally, the reliability of the *proGAV 2.0* is preserved. Once physical activity is stopped the original opening pressure will be re-established.

PHYSICS BACKGROUND

The intraventricular pressure is positive in a healthy human in a horizontal position. To adjust this pressure through shunt drainage, one has to choose the appropriate pressure range, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (fig. 6).

The ventricular pressure in a healthy human in a vertical position becomes slightly negative. To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure.


Fig. 6: Calculating the intraventricular pressure for horizontal and vertical body position

SELECTING THE APPROPRIATE SHUNT

The *proGAV 2.0* is a position-dependent shunt, meaning the opening pressure changes depending on the position of the patient.

To choose the suitable *proGAV 2.0* for an individual patient, one opening pressure is set for the horizontal position (patient lying down), and one for the vertical position (patient standing upright).

Horizontal position

The opening pressure for the horizontal position is defined by the adjustable DP-unit. The pressure level should be chosen according to the clinical situation and indication. The unit can be adjusted to a pressure setting between 0 cmH₂O and 20 cmH₂O. The valve is preset to 5 cmH₂O.

Vertical position

The opening pressure of the *proGAV 2.0* for the vertical position is calculated by the sum of the opening pressure of both the adjustable unit and of the gravitational unit. The selection of the gravitational unit depends on the activity and the abdominal pressure (adiposity).

USING THE INSTRUMENTS

With the *proGAV 2.0 Tool Set* the selected opening pressure of the *proGAV 2.0* can be determined, varied and controlled.

The *proGAV 2.0 Compass* is used to locate and verify the DP adjustable unit.


Fig. 7: *proGAV 2.0 Compass*
a) open b) closed

The *proGAV 2.0 Adjustment Tool* is used for adjusting the valve opening pressure of the *proGAV 2.0* from 0 to 20 cmH₂O.

Fig. 8: *proGAV 2.0 Adjustment Tool*

Each *proGAV 2.0* is calibrated under strict quality control procedures. The presetting of the adjustable DP-unit is 5 cmH₂O, but it must be checked before implantation. The setting is changed in the following steps:

1. Locating the valve

If the instrument is opened a template is visible (fig. 9). Then the valve can be located on the patient's head with the forefinger.

Fig. 9 Locating the valve with the *proGAV 2.0 Compass*

The *proGAV 2.0 Compass* must be positioned centrally on the valve. The markings on the instrument "proximal" and "distal" show the flow direction.

2. Verifying the opening pressure

When the compass is closed, the pressure setting is indicated automatically. (fig. 10)

Fig. 10: Verifying the pressure setting with the *proGAV 2.0 Compass*

Caution: Placing the *proGAV 2.0 Compass* in a non-central position on the valve can lead to erroneous readings!


The *proGAV 2.0 Compass* is sensitive to external magnetic fields. To exclude undesirable interactions the *proGAV 2.0 Adjustment Tool* should not be in the immediate vicinity of the *proGAV 2.0 Compass* while determining the opening pressure. We recommend a distance of about 30 cm.

3. Adjusting the opening pressure

The *proGAV 2.0 Adjustment Tool* must be positioned centrally on the valve. For a correct placement the valve should be palpated with the forefinger through the opening in the middle of the instrument. The desired pressure setting must point on the scale in direction of the inlet connector and the ventricular catheter. By applying light pressure the rotorbrake will be released and the pressure of the *proGAV 2.0* can be changed.

a)

b)

Fig. 11: a) and b)

Adjustment with the *proGAV 2.0 Adjustment Tool*

The *proGAV 2.0* is equipped with a feedback-mechanism. When using the *proGAV 2.0 Adjustment Tool*, pressure on the housing of the valve is created and a resulting acoustic signal (a clicking sound) is produced due to the unique construction of the valve housing. This clicking sound indicates that the rotorbrake is released. Now the rotor can rotate freely. Once the pressure on the valve is released, a clicking sound is heard and the rotorbrake is again locked safely so that the valve is safe against spontaneous re-adjustments.

The clicking sound is well recognizable before implantation. However after implantation, once the valve is filled up, depending on place and

texture of the surrounding area of the implant, the acoustic signal could be considerably muted. The clicking sound should generally be audible by the patient itself or via a stethoscope.

Caution: The new opening pressure setting of the valve must not differ from the measured opening pressure by more than 8 cmH₂O in any one setting (see chapter 4 “verifying the adjustment”).

Example: Opening pressure is to be changed from 3 to 18 cmH₂O. With only one adjustment procedure the rotor would turn in the wrong direction (short way) and would stop at the position 0 cmH₂O. The correct adjustment is in 2 steps: Adjustment from 3 to 11, and from 11 to 18 cmH₂O. The rotor turns correctly.


Fig. 12: Rotor rotation during adjustment
a) false b) correct

Caution: From proGAV 2.0 Adjustment Tool a magnetic field emanates. Metallic objects and magnetic media storages should have a sufficient safety margin.

4. Verifying the adjustment

After adjusting the valve by using the *proGAV 2.0 Adjustment Tool*, it must be verified using the *proGAV 2.0 Compass* as described in step 2. If the measured pressure now differs from the intended pressure level, the adjustment procedure has to be repeated from step 3.

Due to postoperative swelling of the skin the adjustment of the valve setting may be difficult within the first few days.

If the pressure configuration of the valve cannot be determined with complete certainty by the *proGAV 2.0 Compass*, the use of imaging techniques is recommended.

proGAV Check-mate

The *proGAV Check-mate* is delivered sterile and is intended to be re-sterilised. It is possible to change and to verify an applied pressure setting on the valve directly. To verify the actual pressure setting the *proGAV Check-mate* has to be put centrally over the valve. The *proGAV Check-mate* will immediately start to move. If it remains stable, the pressure setting can be read in alignment to the inlet connector.

To adjust a new pressure setting, the *proGAV Check-mate* has to be placed centrally over the valve. The new pressure setting has to point towards the proximal catheter (leading to the ventricle). By pressing down slightly the *proGAV Check-mate*, the brake of the valve is decoupled, the rotor turns and the opening pressure of the *proGAV 2.0* is changed.

Please be aware that the steps for changing the pressure setting should not be more than 8 cmH₂O per step.

Fig. 13: *proGAV Check-mate*

Caution: Due to magnets inside the proGAV 2.0 Tools, do not use the proGAV 2.0 Tools nearby pacemakers. Further more do not use the proGAV 2.0 Tools near-by MRI scanner, since there is a danger of damaging the MRI-scanner.

ADJUSTING THE ADJUSTABLE DP-UNIT

Please verify specifically before using any tool for verifying or adjusting the opening pressure:

For the adjustable DP-unit use either the *proGAV 2.0 Tools*


Fig. 14:

- a) *proGAV 2.0 Compass*
- b) *proGAV 2.0 Adjustment Tool*

and with a combination of *proGAV 2.0* with an **adjustable gravitational-unit, use only *proSA Tools***


Fig. 15: *proSA Adjustment Tool*

READING THE PRESSURE SETTING FROM AN X-RAY IMAGE

The pressure setting of the *proGAV 2.0* should be measured either with the *proGAV 2.0 Compass* or via X-ray.

The position of the rotor determines and indicates the pressure setting in the X-ray.

The rotor includes four magnets, arranged in pairs and face to face on top of the rotor, which are recognizable in the X-ray as four white dots. Two additional drill holes on one side of the rotor, arranged on the left respectively right side of the magnets, serve as an additional support for orientation: the two drill holes are recognizable as black dots in the X-ray.

The rotor side with the two drill holes can be described as the back side of the rotor. The gap of the two magnets opposite to those on the back side of the rotor is regarded as the triangle tip. A triangle is formed by visually connecting the magnets and the drill holes. The sharp corner of the triangle indicates the position of the triangle tip. The direction of the triangle tip indicates the pressure setting of the valve (see fig. 16).


Fig. 16: Schematic X-ray image

The triangle tip can occupy any position outside the region 'range not adjustable' (fig. 16). Thus, the opening pressure of the *proGAV 2.0* can be adjusted in increments of 1 cmH₂O between 0 and 20 cmH₂O.

In order to avoid misidentification of the adjusted opening pressure in the X-ray image, the valve is marked with a orientation indicator on one side (recognizable as a black cut-out in the X-ray see fig. 17). On the schematic top view as in fig. 16 the valve indication is visible on the right hand side of the valve's housing.


Fig. 17: Radiograph pressure setting 14 cmH₂O

The following opening pressure ranges for the gravitational unit are possible, the pressure range selected can be checked postoperatively on X-ray image:

Opening pressure for vertical posture	Coding of gravitational unit
10 cmH ₂ O	no ring, small (fig. 18b)
15 cmH ₂ O	no ring, large
20 cmH ₂ O	1 ring, large (fig. 18a)
25 cmH ₂ O	2 rings, large
30 cmH ₂ O	3 rings, large
35 cmH ₂ O	4 rings, large


Fig. 18: X-ray image of the gravitational unit
a) large, 1 ring = 20 cmH₂O,
b) small, no ring = 10 cmH₂O


Fig. 19: Schematic representation of the proGAV 2.0 with its external dimensions
a) with large gravitational unit
b) with small gravitational unit

POSSIBLE SHUNT COMPONENTS

The *proGAV 2.0* is available in different configurations. These shunt variations (for pediatric hydrocephalus and adult hydrocephalus) comprise of a variety of components, which are briefly described below.

The *Burrhole reservoir*, the *SPRUNG RESERVOIR* or the *McLANAHAN RESERVOIR* are positioned in the cranial burrhole. They allow flushing and extraction of CSF, measurement of the intraventricular pressure and the injection drugs. Their solid titanium base is highly puncture-resistant.

With the *SPRUNG RESERVOIR* it is possible to flush CSF one-way in the direction of the valve, this is achieved due to a non-return valve in the bottom of the reservoir. With the aid of this mechanism a flow in direction of the ventricular catheter is avoided during the pumping procedure. A control of the distal part of the shunt system (the reservoir is hard to pump) and whether the ventricular catheter is occluded (the reservoir does not refill after pumping) can be carried out.

The *CONTROL RESERVOIR* or the (pediatric) *prechamber* are positioned on the cranium. They allow the measurement of intraventricular pressure, injection of drugs, extraction of CSF and palpatory ventricle inspection. Similarly to the *SPRUNG RESERVOIR* the *CONTROL RESERVOIR* contains a non-return valve.

Its solid titanium base is highly puncture-resistant. A puncture of the *prechamber* or the *CONTROL RESERVOIR* should be performed with a cannula with a maximum diameter of 0,9 mm as perpendicular to the reservoir surface as possible. 30 punctures are possible without any restrictions. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR* or the *CONTROL RESERVOIR*.

Warning note: Frequent pumping can lead to overdrainage and thus to pressure conditions outside the normal physiological range. The patient should discuss the risks (involved) with their surgeon.

Tight tolerancing of the deflector ensures a good fit with the ventricular catheter. By adjusting the deflector (prior to implantation) the length of catheter penetrating into the skull can be optimised. The ventricular catheter is "deflected" at a right angle in the burrhole (see chapter "Variations").

TUBE SYSTEMS

The *proGAV 2.0* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows using catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. Regardless, the catheters must be carefully fixed, with a ligature, to the valve connectors. It is essential that kinks in the catheter are avoided. The included catheters have virtually no effect on the Pressure-flow characteristics.

SURGICAL PROCEDURE

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible

after applying the burrhole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The *proGAV 2.0* is available in different shunt variations:

When using a *proGAV 2.0 SHUNTSYSTEM with Burrhole reservoir, SPRUNG RESERVOIR* or *McLANAHAN RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the *burrhole reservoir* is connected, with the connection secured with a ligature. The skin incision should not be located directly above the reservoir.

The *proGAV 2.0 SHUNTSYSTEM with (pediatric) prechamber* or *CONTROL RESERVOIR* comes with a *Ventricular Catheter with deflector*. This deflector is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the prechamber is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the valve

The adjustable DP-unit of the *proGAV 2.0* is supplied with a factory setting of 5 cmH₂O. This opening pressure can be set to a different value prior to implantation (see chapter "Adjusting the *proGAV 2.0*"). The gravitational unit of the *proGAV 2.0* is a posture-dependent valve. Therefore, care must be taken that the unit is implanted parallel to the body axis. A suitable implantation site is behind the ear.

After the skin incision and tunneling under the skin, the catheter is pushed forward, from the burrhole to the intended shunt implantation site. The catheter is shortened, if necessary, and secured at the *proGAV 2.0* with a ligature. The shunt should not be located directly under the skin incision. The valve is marked with arrows pointing in the direction of flow (arrows pointing to distal or downward).

Warning note: The adjustable DP-unit must be placed over a hard bony surface and should not be implanted within an area, which makes finding and feeling the valve more difficult (e. g. under a scar).

Frequent pumping can lead to overdrainage and

thus to unphysiological pressure conditions. The patient should be informed about the risk.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied e. g. para-umbilically in a horizontal direction or transrectally at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to the intended position of the catheter. The peritoneal catheter, which is usually securely attached to the *proGAV 2.0*, has an open distal end, but no wall slits. Following the exposure of, and the entry into, the peritoneum by means of a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

PREOPERATIVE VALVE TEST


Fig. 20: Patency test

The *proGAV 2.0* can be filled by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see fig. 20).

Caution: Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end. Contaminations in the solution used for the test can impair the product's performance.

PRESSURE-FLOW CHARACTERISTICS

Horizontal valve position

The following diagrams show the pressure-flow characteristics of the adjustable DP-unit of the *proGAV 2.0* for the pressure settings 0, 10 and 20 cmH₂O in the horizontal position of the valve.


Fig. 21: Pressure-flow characteristics for some pressure settings of the adjustable DP-unit

Vertical valve position

The opening pressure of the *proGAV 2.0* in the vertical position is the sum of the opening pressure of the adjustable DP-unit and the gravitational unit.

The following diagrams (see fig. 21) show the pressure-flow characteristics for some pressure settings in the vertical body position. The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher.







Fig. 22: Pressure-flow characteristics for some pressure settings of the adjustable DP-unit (0, 10, 20 cmH₂O setting) and the gravitational unit of a) 10 cmH₂O, b) 15 cmH₂O, c) 20 cmH₂O, d) 25 cmH₂O, e) 30 cmH₂O, f) 35 cmH₂O

CONTRAINDICATIONS

The adjustable DP-unit should not be implanted within an area which makes locating and sensing the valve more difficult (e. g. under a scar). The valve should lie on the perist or the bone to make an adjustment after implantation possible.

INTERACTIONS WITH PRODUCTS FROM OTHER MANUFACTURERS

The *proGAV 2.0* with gravitational unit should not be used under any circumstances in conjunction with any additional hydrostatic valves, as this can bring about abnormally high ventricular pressure outside of the normal physiological range. Hydrostatic valves allow for changes in hydrostatic pressure in the drainage system caused by changes in position. If in doubt, please contact the medical product consultants at Christoph Miethke GmbH & Co. KG.

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

SAFETY MEASURES

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible, or of the entire shunt system.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried without endangering or impairing the functionality of the shunt. The *proGAV 2.0* is MR Conditional (ASTM-F2503-13). The *proGAV 2.0* will not change when subjected to an MRI of 1.5T or 3.0T. All components are visible via X-ray. The provided catheters are MRI Safe. Reservoirs, deflectors and connectors are MR Conditional.

Warning note: When using a magnetic field and simultaneous pressing on the valve an adjustment of the valve cannot be excluded. The *proGAV 2.0* will produce artifacts or signal-intensity voids in MR images larger than the physical size of the device.

POSTOPERATIVE VALVE TEST

The *proGAV 2.0* has been designed as a safe and reliable unit even without the provision of a pumping device. However, there are ways of testing the unit if a shunt system with a prechamber or a *burrhole reservoir* is used. Valve tests can be carried out by flushing or pressure measurements.

FUNCTIONAL SAFETY

The valves have been designed for long-term reliable and precise operation. Still, it cannot be excluded that the shunt system needs to be replaced for technical or medical reasons. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation. The valve should be stored dry and clean.

Warning note for carriers of pacemakers: Due to the implantation of a *proGAV 2.0* the function of a pacemaker can be affected.

ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

STERILIZATION

The products are sterilized with steam under closely monitored conditions. The double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

RESTERILIZATION

The functional safety and reliability of resterilized products cannot be guaranteed, therefore resterilisation is not recommended.

MEDICAL PRODUCTS CONSULTANTS

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH & Co. KG names medical products consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz
Michaela Funk-Neubarth

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
D-14469 Potsdam · Germany
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:
AESCULAP AG
Am Aesculap Platz
D-78532 Tuttlingen · Germany
Phone: +49 (0) 7461 95-0
Fax: +49 (0) 7461 95-26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042

AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

VARIATIONS

The *proGAV 2.0* is available as a single valve or as a shunt system comprising various components.


proGAV 2.0

proGAV 2.0 SHUNTSYSTEM

proGAV 2.0 SHUNTSYSTEM with SPRUNG RESERVOIR or McLANAHAN RESERVOIR or (pediatric) Burrhole reservoir

proGAV 2.0 SHUNTSYSTEM with CONTROL RESERVOIR or with (pediatric) prechamber

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany