Aesculap Spine

(GB) Instructions for use/Technical description

S⁴ Cervical – Navigated Instruments

(USA) Instructions for use/Technical description

S⁴ Cervical – Navigated Instruments

Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.

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







S⁴ Cervical – Navigated Instruments

Legend

- A** Navigated tapper Ø 3.5 mm (FW655R)
- B** Reference star unit (pre-calibrated) (Brainlab Art. no. 55830-20A)
- C** Handle (FW165R with ratchet or FW067R without ratchet)
- D** Navigated screwdriver for polyaxial screw (FW656R)
- E** Reference star unit ML (calibration with ICM4) (Brainlab Art. no. 55830-25A)
- F** Navigated drill guide Ø 3.5 mm, (FW664R)
- G** Navigated drill guide Ø 4.0 mm, (FW665R)
- H** Aesculap star unit navigation attachment (FW652R)
- I** Navigated guide sleeve Ø 4.0 mm for smooth-shaft screws (FW658R)
- J** Reduction sleeve Ø 13 mm for Brainlab instrument calibration matrix version 4 (FW657R)
- K** Navigated guide sleeve Ø 3.5 mm/4.0 mm (FW660R)
- L** C1/C2 inner drill sleeve Ø 4 mm (FJ985R)
- M** C1/C2 inner drill sleeve Ø 3.5 mm (FW661R)
- N** C1/C2 drill Ø 2.4 mm (FW662SU)
- O** C1/C2 tapper Ø 3.5 mm (FW663R)
- P** Centering tool for navigated drill guide Ø 3.5 mm FW664R (FW688R)
- Q** Centering tool for navigated drill guide Ø 4.0 mm FW665R (FW689R)
- R** Drill Ø 2.4 mm for FW664R (FW666SU)
- S** Drill Ø 2.9 mm for FW665R (FW667SU)

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product
	Date of manufacture

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

Intended use

The S⁴ Cervical system (S⁴C system) is used for posterior occipital cervical and thoracic stabilization and fusion. Indications and contraindications are specified in the instructions for use for the implants (TA011796). The instruments listed in the key form part of this system. They are used for the individual adaptation, positioning and insertion of S⁴C implants in patients. For safe usage, follow instructions for use for S⁴C System Instruments (TA011984) and Operating Techniques (O34202/O83002).

The instruments listed in the key may only be used with the Brainlab navigation system. For safe handling prior to the operation, read Spine & Trauma user manual for Brainlab instruments and the corresponding software manual for the Brainlab spine application used.

Combination specifications

Note

Aesculap and Brainlab accept absolutely no responsibility if instruments, centering tools or drills other than those named below are used with the corresponding drill guides and guide sleeves.

- Only combine S⁴C navigated drill guide **F** with:
 - S⁴C centering tool **P** for S⁴C drill guide **F**
 - S⁴C drill **R** Ø 2.4 mm for Ø 3.5 mm screws
- Only combine S⁴C navigated drill guide **G** with:
 - S⁴C centering tool **Q** for drill guide **G**
 - S⁴C drill **S** Ø 2.9 mm for Ø 4.0 mm screws
- Only combine S⁴C navigated guide sleeve Ø 4.0 mm **I** for smoothshank screws with:
 - Smooth-shaft screw bone awl (FW085R)
 - Drill for Ø 2.9 mm smooth-shaft screws (FW086SU)
 - Smooth-shaft (Shank) screw tapper (Tap) (FW087R)
 - Polyaxial screwdriver (FW070R/FW131R)
 - Navigated screwdriver for polyaxial screw (FW656R)
 - Apfelbaum ball head screwdriver, short (FJ968R)
- Only combine S⁴C navigated guide sleeve Ø 3.5 mm/4.0 mm **K** with:
 - Apfelbaum C1/C2 obturator (FJ983R)
 - Apfelbaum trocar (FJ984R)
 - Favored Angle screw drill Ø 2.9 mm for Ø 4 mm screws (FW088SU)
 - Favored Angle tapper (Tap) Ø 4 mm (FW089R)
 - C1/C2 inner drill sleeve Ø 4 mm **L** (FJ985R)
 - S⁴C Favored Angle screwdriver (FW069R/FW132R)
 - C1/C2 drill Ø 2.4 mm for Ø 3.5 mm screws (FW662SU) **N**
 - C1/C2 tapper (Tap) Ø 3.5 mm (FW663R) **O**
 - C1/C2 inner drill sleeve Ø 3.5 mm (FW661R) **M**
 - Apfelbaum ball head screwdriver (FJ988R)

Note

K-wires generally may not be used with the S⁴C system.

Safe handling and preparation

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.
- ▶ Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Replace any damaged components immediately with original spare parts.
- ▶ To avoid damage to the working end: Carefully insert the product through the working channel (e.g. trocar).

Single-use products

Art. no.	Designation
FW086SU	Drill Ø 2.9 mm for smooth-shaft screws
FW088SU	Favored Angle screw drill Ø 2.9 mm for Ø 4 mm screws
FW662SU	S ⁴ C1/C2 drill Ø 2.4 mm
FW666SU	Drill Ø 2.4 mm
FW667SU	Drill Ø 2.9 mm

DANGER

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

- ▶ Do not reprocess the product.

The product is gamma-sterilized and supplied in sterile packaging.

The product must not be reused.

- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Read, follow, and keep the instructions for use.
- ▶ Use the product only in accordance with its intended use, see Intended use.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Prior to each use, inspect the product for loose, bent, broken, cracked or fractured components.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- ▶ Do not use the product after its use-by date.

Safe operation

WARNING

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

WARNING

Risk of injury to patients due to damaged instruments!

S⁴C instruments are extremely precise and very sensitive.

- ▶ Handle S⁴C instruments with the greatest of care.
- ▶ Check that dropped or damaged S⁴C instruments are accurately calibrated, or send them to the Aesculap Technical Service.

WARNING

Risk of injury to patients due to interruption of the navigation!

- ▶ Prior to the operation, plan the configuration of the operating room, the assembly of the instruments and the alignment of the reference star.
- ▶ Ensure that the navigation camera has an unrestricted view of the reflective marker spheres of the instruments.

WARNING

Risk of injury to patients due to imprecise instruments!

- ▶ Ensure that the instruments used are not bent or damaged.
- ▶ Before use, check the precision of the instruments, particularly that of fine instruments. For this, hold the instrument tip in the pivot point of the Brainlab instrument calibration matrix version 4.

Aesculap®

S⁴ Cervical – Navigated Instruments

WARNING

Risk of injury to patients due to incorrect navigation!

- Use navigated S⁴C instruments only with Brainlab disposable reflective marker spheres.
- For further information on the correct handling of the marker spheres, see corresponding Brainlab user manual.

Preparing holes for S⁴C screws with navigation

WARNING

Risk of injury to patients due to incorrect navigation!

- Observe the combination specifications

To center the screw holes in the cortical bone layer for the self-tapping S⁴C screws Ø 3.5 mm with navigation, use the S⁴C navigated drill guide Ø 3.5 mm **F** only together with S⁴C centering tool for S⁴C navigated drill guide **P**.

To center the screw holes for the self-tapping S⁴C screws Ø 4.0 mm with navigation, use the S⁴C navigated drill guide Ø 4.0 mm **G** only together with S⁴C centering tool for S⁴C navigated drill guide **Q**.

S⁴C navigated drill guide **G** and centering tool **Q** have a purple marking on the shaft for the self-tapping S⁴C screws Ø 4.0 mm, see Fig. 1.


Fig. 1 Purple markings on S⁴C navigated drill guide and centering tool

To center the screw holes in the cortical layer without navigation, see instructions for use of S⁴ Cervical instruments (TA011984).

WARNING

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue when inserting the centering tool without drill guide!

- Use the centering tool only with drill guide FW664R/FW665R.

CAUTION

Restriction of the function of the centering tool or inaccurate reading of the depth gage!

- Do not repaint the centering tool.
- Replace blunt centering tools.

- For mounting the reflective Brainlab marker spheres **4** onto Aesculap star unit **H**, see Brainlab user manual.
- Retract and hold locking sleeve **1** of the Aesculap star unit **H** against the spring pressure in the direction of the arrow.
- Push the Aesculap star unit **H** onto adapter **2** of the S⁴C drill guide **F/G**. When doing so, ensure that the pin of the adapter **3** engages in the recess on the star unit.
- Release locking sleeve **1**.


Fig. 2 Mounting the Aesculap star unit on S⁴C drill guide FW664R/FW665R

- Prior to each use, check and/or validate the instrument with the Brainlab instrument calibration matrix version 4, see Brainlab user manual. When doing so, ensure that the centering tool for navigated drill guide **P/Q** is removed.
- To ensure that the camera has an unrestricted view of the reflective marker spheres, loosen knob **6** on the S⁴C handle of the S⁴C drill guide **F/G** and rotate Aesculap star unit **H** to the desired position, see Fig. 2.
- Once the desired position is reached, tighten knob **6** again.
- Adjust the center punch depth by turning the depth stop **5** on the S⁴C drill guide **F/G**. The maximum depth is 6 mm.
- Insert S⁴C centering tool **P/Q** into S⁴C drill guide **F/G**.
- Check the pre-set center punch depth with a caliper (e.g. AA845R).
- To open the cortical bone, press S⁴C centering tool to the preset depth under control with the Brainlab navigation system.

Drilling holes for the S⁴C screws

DANGER

Risk of injury from an incorrectly placed hole or a hole that is too deep!

- Do not sharpen the drill, as this would cause imprecise or incorrect readings on the depth gauge.
- Replace blunt drills with new ones.

The drill is applied with a S⁴C drill guide and drilled in either manually with the drill handle (FJ839R) or with a motor system with the Aesculap Intra-handpiece (e.g. GD450R/GD456R).

Assembling the drill and drill handle (for manual drilling only)

DANGER

Risk of damage to the spinal cord, nerve roots, adjacent intervertebral space or soft tissue through incorrect drilling!

- Use only the correct S⁴C drill guides to drill holes. Insert drill only with the correct drill guide.
- Before drilling, the pre-set drill length must be checked with a caliper (e.g. AA845R, Caspar instrument for anterior cervical fusion).

DANGER

Injury to spinal cord and nerve roots caused by application of a drill that is too long!

- Use the X-ray image to select an appropriate drill length prior to the operation.
- Drill may only be aligned and inserted under radiographic control and/or with the aid of a navigation system.
- Select a drill of a length equivalent to the intended drill hole depth.

- Insert drill in the drill handle (FJ839R), see Fig. 3.
- Retract and hold locking sleeve against the spring pressure in the direction of the arrow.

- Push the drill into the receptacle of the drill handle as far as it will go.
- Lightly rotate the drill and release locking sleeve. The drill engages audibly.


Fig. 3 Assembling the drill

DANGER

Risk of injury and/or damage to the drill if the drill rotation speed is too high!

- Use the lowest drilling speed possible, so that you can control the drilling depth.
- Do not bend the drill during the drilling process.

Drilling holes for S⁴C screws Ø 3.5 mm/4.0 mm

WARNING

Risk of injury to patients due to incorrect navigation!

- Observe the combination specifications

For controlled drilling of the holes for S⁴C screws Ø 3.5 mm with drill Ø 2.4 mm R (FW666SU), the S⁴C navigated drill guide Ø 3.5 mm F (FW664R) must always be used.

For controlled drilling of the holes for S⁴C screws Ø 4.0 mm with drill Ø 2.9 mm S (FW667SU), the S⁴C navigated drill guide Ø 4.0 mm G (FW665R) must always be used.

Aesculap®

S⁴ Cervical – Navigated Instruments

S⁴C navigated drill guide **G** and drill Ø 2.9 mm **S** have a purple marking on the shaft for the self-tapping S⁴C screws Ø 4.0 mm screws, see Fig. 4.


Fig. 4 Purple marking on S⁴C navigated drill guide and drill

Note

The drill Ø 2.4 mm (FW051SU) and Ø 2.9 mm (FW052SU) must not be used with S⁴C navigated drill guide F/G.

- For mounting the reflective Brainlab marker spheres **4** onto Aesculap star unit **H**, see Brainlab user manual.
- Retract and hold locking sleeve **1** of the Aesculap star unit **H** against the spring pressure in the direction of the arrow.
- Aesculap star unit **H** is then pushed onto adapter **2** of the S⁴C drill guide F/G. When doing so, ensure that the pin **3** of the adapter engages in the recess on the star unit.
- Release locking sleeve **1**.


Fig. 5 Mounting the Aesculap star unit on S⁴C drill guide FW664R/FW665R

- Prior to each use, check and/or validate the instrument with the Brainlab instrument calibration matrix version 4, see Brainlab user manual. When doing so, ensure that the drill **R/S** is removed.
- To ensure that the camera has an unrestricted view of the reflective marker spheres, loosen knob **6** on the handle of the S⁴C drill guide **F/G** and rotate Aesculap star unit **H** to the desired position, see Fig. 2.
- Once the desired position is reached, tighten knob **6** again.
- Adjust the drill depth by turning the depth stop **5** on the S⁴C drill guide **F/G**.
The scale is in mm. Each half turn on the depth stop moves a distance of 0.5 mm.
- Insert the drill with mounted handle or Intra-handpiece into S⁴C drill guide **F/G**.
- Before drilling, the pre-set drill length must be checked with a caliper (e.g. AA845R, Caspar instrument for anterior cervical fusion).
- Drill to the pre-set depth under control with the Brainlab navigation system.

Tapping (optional)

S⁴C screws are self-tapping screws. However, if the bone quality is found to be hard during the operation, the surgeon can also pre-tap the thread with the S⁴C taper.

- For tapping the drill holes for screws Ø 3.5 mm with navigation, use the S⁴C navigated taper **A**.
- For tapping the drill holes for screws Ø 4 mm, use the standard taper without navigation (FW047R), see TA011984.

DANGER

Risk of tissue trauma when using the S⁴C taper (A) and damage to bone thread!

- Prior to using the S⁴C taper, ensure that the moveable sleeve of the taper retracts correctly.

- For mounting the reflective Brainlab marker spheres **4** onto the reference star unit **B** (pre-calibrated), see Brainlab user manual.
- Push the reference star unit onto the shaft **7** of the S⁴C taper (Tap) **A**. Ensure that the star unit is securely fitted onto the shaft of the S⁴C taper (Tap).

Note

The star unit can be rotated on the shaft of the S⁴C taper (Tap).

- ▶ Retract and hold locking sleeve **8** against the spring pressure.
- ▶ Push the S⁴C handle FW067R/FW165R onto the shaft of the S⁴C taper **A**.
- ▶ Release locking sleeve **8**, see Fig. 6.
Check that the S⁴C handle is engaged.


Fig. 6 Mounting the reference star unit ("pre-calibrated") and S⁴C handle onto the S⁴C taper

WARNING

Risk of injury to the patient!

- ▶ Before use, ensure that the selected instrument has been correctly assembled.
- ▶ Ensure that the arrow on the underside of the reference star unit (pre-calibrated) is pointing to the tip of the tool.

- ▶ Prior to each use, check and/or validate the instrument with the Brainlab instrument calibration matrix version 4, see Brainlab user manual.
- ▶ To tap the thread, hold the reference star unit **B** at the planned indentations with one hand, and with the other hand screw in S⁴C handle **C** slowly and steadily under control with the Brainlab navigation system, until the required depth is reached.
- ▶ Use the scale behind the retractable sleeve of the S⁴C taper to read the depth during the tapping process, see Fig. 7.

Fig. 7 Taper with readable thread depth

Position S⁴C screw under navigation and temporarily fix it in place

DANGER

Risk of injury to the patient!

- ▶ Use S⁴C screwdriver only with reference star unit ML for manual calibration.
- ▶ If you change screws, perform the calibration again.

- ▶ To position S⁴C screws $\varnothing$ 3.5 mm and $\varnothing$ 4 mm with navigation, use S⁴C navigated screwdriver **D**.
- ▶ For mounting the reflective Brainlab marker spheres **4** onto the reference star unit ML **E**, see Brainlab user manual.
- ▶ Push the reference star unit ML **E** onto the shaft **10** of the S⁴C screwdriver **D**. Ensure that the star unit is securely fitted onto the shaft of the S⁴C screwdriver.

Note

The star unit can be rotated on the shaft of the S⁴C screwdriver.

- ▶ Retract and hold locking sleeve **8** against the spring pressure.
- ▶ Push the S⁴C handle FW067R/FW165R onto the shaft of the S⁴C screwdriver **D**.
- ▶ Release locking sleeve **8**, see Fig. 8.
Check that the S⁴C handle is engaged.


Fig. 8 Mounting the reference star unit ("calibration with ICM 4") and S⁴C handle (FW067R or FW165R) onto the S⁴C screwdriver

Note

The S⁴C screwdriver is fitted with a self-retaining function to prevent the S⁴C screw from falling out when it is passed to the surgeon.

- ▶ Retract and hold the holding sleeve **12** on the S⁴C screwdriver **D**, see Fig. 9.
- ▶ Insert the tip of the S⁴C screwdriver **D** fully into the hexagon of the screw **11**.

Aesculap®

S⁴ Cervical – Navigated Instruments

- Release the holding sleeve 12.

Ensure that the screw 11 is securely in place on the S⁴C screwdriver D and that the polyaxiality of the screw 11 is blocked.


Fig. 9 Pick up the S⁴C screw with S⁴C screwdriver

DANGER

Risk of injury to the patient!

- Before use, ensure that the selected instrument has been correctly assembled. Ensure that the arrow on the underside of the reference star unit ML is pointing to the tip of the tool.

- Prior to using the instrument with the correctly held screw, perform manual calibration with the Brainlab instrument calibration matrix version 4, see Brainlab user manual.
- Screw in the screw under control with the Brainlab navigation system. When doing so, hold reference star unit ML at the planned indentations with one hand, and turn the S⁴C handle C to screw in the screw 11 with the other hand.

Smooth-shaft screw instruments

WARNING

Risk of injury to patients due to incorrect navigation!

- Observe the combination specifications

Instruments for smooth-shaft screws are marked with a light-blue ring. They are used to center punch, drill and tap holes for smooth-shaft screws Ø 4 mm.

DANGER

Serious complications for the patient can be caused by incorrect positioning of instruments or implants!

- Carry out operative steps with radiographic visualization.
- When removing the smooth-shaft screw awl (FW085R) and during the further operating steps, ensure that the S⁴C navigated smooth-shaft screw guide sleeve remains securely fixed in place.
- Ensure that the window on the S⁴C navigated smooth-shaft screw guide sleeve is closed during the preparation of the screw hole and while the screw is being inserted, see laser marking on the inner sleeve.
- Take care that no tissue gets caught when opening and closing the window on the S⁴C navigated smooth-shaft screw guide sleeve, after the screw has been put in place.

- For mounting the reflective Brainlab marker spheres 4 onto Aesculap star unit H, see Brainlab user manual.
- Retract and hold locking sleeve 1 against the spring pressure in the direction of the arrow.
- Aesculap star unit H is then pushed onto adapter 2. When doing so, ensure that the recess on the star unit is seated over the pin 3 of the adapter.
- Release locking sleeve 1.


Fig. 10 Mounting the Aesculap star unit on the S⁴C guide sleeve for smooth-shaft screws.

Risk of injury to patients due to incorrect navigation!

- Slide reduction sleeve into the Brainlab instrument calibration matrix version 4, until you hear an audible click.

- Prior to each use, check and/or validate the instrument with special reduction sleeve J with the Brainlab instrument calibration matrix version 4, see Brainlab user manual. Ensure that all other instruments (centering tool, drill, tapper etc.) are removed for the instrument validation.


Fig. 11 Reduction sleeve for validation/verification with the Brainlab instrument calibration matrix version 4

- To ensure that the camera has an unrestricted view of the reflective marker spheres, loosen knob 6 on the S⁴C handle of the S⁴C guide sleeve I and rotate Aesculap star unit H to the desired position, see Fig. 10.
- Once the desired position is reached, tighten knob 6 again.
- Place S⁴C guide sleeve I for smoothshank screws in the operating field. When doing so, ensure that the window of the guide sleeve is closed during the preparation and insertion of the screw with the guide sleeve, see laser marking on the inner sleeve 13.
- Center punch the cortical layer of the vertebral body with the smooth-shank screw aw (FW085R), see TA011984 and Operating Technique (O34202).
- If necessary, insert the awl into the inner sleeve 13 and center punch the bone as far as the stop-position. The stop-position is indicated with a marking on the awl.
- Remove the awl from the operating field.
- Insert the smooth-shaft screw drill (FW086SU) 14 with mounted handle (FJ839R) or Intra-handpiece into the inner sleeve 13.
- Before drilling, the pre-set drill length must be checked with a caliper (e.g. AA845R, Caspar instrument for anterior cervical fusion).
- Under control with the Brainlab navigation system, drill a hole in the bone until the adjustable stop-position is reached. For further information on drilling, see instructions for use of the S⁴ Cervical instruments (TA011984) and Operating Technique (O34202).

- Remove the drill from the operating field.
- To prepare the drill holes for the screws, tap the thread with the smooth-shaft screw tapper (FW087R).
- Insert the smooth-shaft screw tapper into the inner sleeve and slowly and steadily tap in until the desired depth is reached. When doing so, read the thread depth on the tapper scale.
- Remove smooth-shaft screw tapper from the operating field.

Note

If the S⁴C screw is to be inserted under navigation with the S⁴C screwdriver FW656R using the guide sleeve I, remove the Aesculap star unit on the S⁴C guide sleeve for smoothshank screws.

Risk of injury to patients due to incorrect navigation!

- The navigated S⁴C screwdriver FW656R or other S⁴C screwdrivers are only intended for navigation in the non-navigated S⁴C guide sleeve FW658R.
- For navigating the S⁴C screwdriver FW656R, see section on positioning the S⁴C screw under navigation and temporarily fixing it in place.

Risk of injury to the patient through freely rotating screws!

- Do not screw in screw so far that the screw head comes into contact with the S⁴C guide sleeve.

Note

Use the navigated S⁴C screwdriver FW656R only with non-navigated S⁴C guide sleeve FW658R.

- Navigate S⁴C screwdriver, see Position S⁴C screw under navigation and temporarily fix it in place.
- Insert the screw through the S⁴C guide sleeve but do not screw it in completely (smooth shaft must remain free). Remove S⁴C screwdriver D from the operating field.
- Remove instrument from the screw:
 - Turn the blue inner sleeve 13 and open the window on the S⁴C guide sleeve.
 - Carefully push away S⁴C guide sleeve I laterally from the screw.
 - Remove S⁴C guidesleeve I from the operating field.

Instruments for Favored Angle screws

WARNING

Risk of injury to patients due to incorrect navigation!

- Observe the combination specifications

Favored Angle instruments are marked with a gold-colored ring.

DANGER

Serious complications for the patient can be caused by incorrect positioning of instruments or implants!

- Carry out operative steps with radiographic visualization.
- When removing the obturator (FJ983R) and during the further operating steps, ensure that the S⁴C guide sleeve remains securely fixed in place.

Drilling holes for Favored Angle screws

- For mounting the reflective Brainlab marker spheres **4** onto Aesculap star unit **H**, see Brainlab user manual.
- Aesculapstar unit **H** is mounted onto S⁴C navigated guide sleeve $\varnothing$ 3.5/4.0 mm **K**. When doing so, retract and hold the locking sleeve **1** against the spring pressure in the direction of the arrow.
- Push the Aesculap star unit **H** onto adapter **2** of the S⁴C guide sleeve **K**. When doing so, ensure that the recess is seated over the pin **3** of the adapter.
- Release locking sleeve **1**.


Fig. 12 Mounting the Aesculap star unit onto the S⁴C guide sleeve; Inserting the S⁴C inner drill sleeve

- Prior to each use, check and/or validate the instrument with the Brainlab instrument calibration matrix version 4, see Brainlab user manual. Ensure that the inner drill sleeve **L** is mounted prior to this check.
- To ensure that the camera has an unrestricted view of the reflective marker spheres, loosen knob **6** on the S⁴C handle of the S⁴C guide sleeve **K** and rotate Aesculap star unit **H** to the desired position, see Fig. 12.
- Once the desired position is reached, tighten knob **6** again.
- Remove inner drill sleeve **L** from the S⁴C guide sleeve **K**.
- Insert obturator (FJ983R) into the inner sleeve **15** of the S⁴C guide sleeve **K**.
The obturator engages in the inner sleeve and can still be rotated.
- Bring S⁴C guide sleeve **K** with mounted obturator into the operating field through the stab incision and position it in place.
- Press the button **16** on the obturator (FJ983R) and withdraw the obturator from the inner sleeve **15**.
- If necessary, slide Apfelbaum trocar (FJ984R) into the inner sleeve **15** and insert into the bone to center punch the screw entry point.
- Remove the trocar from the operating field.
- Slide the inner drill sleeve **L** onto the inner sleeve **15**.
The inner drill sleeve engages on the inner sleeve and can still be rotated, see Fig. 12.
- Insert drill for Favored Angle screws (FW088SU) with mounted handle (FJ839R) or Intra-handpiece into the inner drill sleeve **L**.
- Drill to the pre-set depth under control with the Brainlab navigation system. The drill depth can be read on the scale on the inner drill sleeve **L**.

Note

*So as not to lose the entry opening, it can be helpful to keep the drill in the drill hole, press the button **16** on the inner drill sleeve **L** and push down the guide sleeve until it reaches the stop on the bone surface. Then remove the drill and inner drill sleeve **L** from the S⁴C guide sleeve **K**.*

- Press the button **16** and remove the inner drill sleeve **L** from the inner sleeve **15**.

S⁴C screws are self-tapping screws. However, if the bone quality is found to be hard during the operation, the surgeon can also pre-tap the thread with the S⁴C taper.

- To prepare screw holes, insert the Favored Angle taper $\varnothing$ 4 mm (FW089R) into the inner sleeve **15** and tap. The drill depth can be read on the scale on the taper (Tap).
- Turn the taper (Tap) counterclockwise until it exits the bone.

Note

So as not to lose the entry opening, it can be helpful to turn the taper (Tap) counterclockwise until it almost exits the bone. Then turn the inner sleeve 15 counterclockwise and at the same time push down the S⁴C guide sleeve K until it reaches the stop on the bone surface. Then unscrew the Favored Angle taper Ø 4 mm (FW089R) completely from the bone and remove it together with the inner sleeve 15 from the S⁴C guide sleeve K.

Inserting the screw

- Ensure that the inner sleeve 15 has been removed from the S⁴C guide sleeve K by turning it counterclockwise.
- Pick up a suitable Favored Angle screw Ø 4.0 mm with the self-retaining S⁴C screwdriver (FW069R/FW132R). When doing so, retract and hold the black holding sleeve against the spring pressure.

Note

The self-retaining function of the instrument prevents the screw from falling out of the S⁴C screwdriver when it is being passed to the operating surgeon.

- Press the working end of the S⁴C screwdriver fully into the hexagon of the screw 11.
- Release the black holding sleeve.
- Screw in the screw under control with the Brainlab navigation system.
- Tighten the screw. When doing so, work through the S⁴C guide sleeve.
- Activate the black holding sleeve and release the S⁴C screwdriver from the screw.
- Remove S⁴C guide sleeve and S⁴C screwdriver from the operating field.

Standard screw (Ø 3.5 mm) instrumentation with Favored Angle instruments

WARNING

Risk of injury to patients due to incorrect navigation!

- Observe the combination specifications

To insert standard screws Ø 3.5 mm with the Favored Angle instruments, the instruments M, N and O must also be used. These instruments are marked with a black ring.

DANGER

Serious complications for the patient can be caused by incorrect positioning of instruments or implants!

- Carry out operative steps with radiographic visualization.
- When removing the obturator (FJ983R) and during the further operating steps, ensure that the guide sleeve remains securely fixed in place.

- For mounting the reflective Brainlab marker spheres 4 onto Aesculap star unit H, see Brainlab user manual.
- Aesculap star unit H is mounted onto S⁴C navigated guide sleeve Ø 3.5/4.0 mm K. When doing so, retract and hold the securing sleeve 1 against the spring pressure in the direction of the arrow.
- Aesculap star unit H is then pushed onto adapter 2 of the S⁴C guide sleeve K. When doing so, ensure that the recess is seated over the pin 3 of the adapter.
- Release locking sleeve 1.


Fig. 13 Mounting the Aesculap star unit onto the S⁴C guide sleeve; Inserting the S⁴C drill guide

- Prior to each use, check and/or validate the instrument with the Brainlab instrument calibration matrix version 4, see Brainlab user manual. Ensure that the inner drill sleeve M is mounted prior to this check.
- To ensure that the camera has an unrestricted view of the reflective marker spheres, loosen knob 6 on the S⁴C handle of the S⁴C guide sleeve K and rotate Aesculap star unit H to the desired position, see Fig. 13.
- Once the desired position is reached, tighten knob 6 again.
- Remove inner drill sleeve M from the guide sleeve K.
- Insert obturator (FJ983R) into the inner sleeve 15 of the S⁴C guide sleeve K. The obturator engages in the inner sleeve and can still be rotated.
- Bring S⁴C guide sleeve K with mounted obturator into the operating field through the stab incision and position it in place.
- Press the button 16 on the obturator (FJ983R) and withdraw the obturator from the inner sleeve 15.

Aesculap®

S⁴ Cervical – Navigated Instruments

- ▶ If necessary, slide Apfelbaum trocar (FJ984R) into the inner sleeve **15** and insert into the bone to center punch the screw entry point.
- ▶ Remove the trocar from the operating field.
- ▶ Slide the inner drill sleeve **M** onto the inner sleeve **15**. The inner guide sleeve engages on the inner sleeve and can still be rotated, see Fig. 13.
- ▶ Insert drill Ø 2.4 mm **N** with mounted handle (FJ839R) or Intra-hand-piece into the inner drill sleeve **M**.
- ▶ Drill to the pre-set depth under control with the Brainlab navigation system. The drill depth can be read on the scale on the inner drill sleeve **M**.

Note

So as not to lose the entry opening, it can be helpful to keep the drill in the drill hole, press the button **16** on the inner drill sleeve **M** and push down the guide sleeve until it reaches the stop on the bone surface. Then remove the drill and inner drill sleeve **M** from the S⁴C guide sleeve **K**.

- ▶ Press the button **16** and remove the inner drill sleeve **M** from the inner sleeve **15**.

S⁴C screws are self-tapping screws. However, if the bone quality is found to be hard during the operation, the surgeon can also pre-tap the thread with the S⁴C taper.

- ▶ To prepare screw holes, insert the C1/C2 taper (Tap) Ø 3.5 mm **O** into the inner sleeve **15** and tap. The drill depth can be read on the scale on the taper (Tap).
- ▶ Turn the taper (Tap) **O** counterclockwise until it exits the bone.

Note

So as not to lose the entry opening, it can be helpful to turn the taper (Tap) counterclockwise until it almost exits the bone. Then turn the inner sleeve **15** counterclockwise and at the same time push down the S⁴C guide sleeve **K** until it reaches the stop on the bone surface. Then unscrew the taper **O** completely from the bone and remove it together with the inner sleeve **15** from the S⁴C guide sleeve **K**.

Inserting the screw

- ▶ Ensure that the inner sleeve **15** has been removed from the S⁴C guide sleeve **K** by turning it counterclockwise.
- ▶ Pick up a suitable standard screw Ø 3.5 mm with the self-retaining S⁴C screwdriver (FW069R/FW132R). When doing so, retract and hold the black holding sleeve against the spring pressure.

Note

The self-retaining function of the instrument prevents the screw from falling out of the S⁴C screwdriver when it is being passed to the operating surgeon.

- ▶ Press the working end of the S⁴C screwdriver fully into the hexagon of the screw **11**.
- ▶ Release the black holding sleeve.

- ▶ Screw in the screw under control with the Brainlab navigation system.
- ▶ Tighten the screw. When doing so, work through the S⁴C guide sleeve.
- ▶ Activate the black holding sleeve and release the S⁴C screwdriver from the screw.
- ▶ Remove S⁴C guidesleeve S⁴C from the operating field.

Disassembling

S⁴C navigated taper Ø 3.5 mm, FW655R

- ▶ Loosen nut **19** and unscrew from the sleeve **17**.
- ▶ Remove sleeve **17** in the direction of the arrow.

Note

Spring **18** is fixed on shaft.


Fig. 14 Disassembling the S⁴C taper (Tap)

S⁴C navigated screwdriver, FW656R

- ▶ Retract and hold green holding sleeve **22** against the spring pressure.
- ▶ Unscrew screw sleeve **20** at the working end of the S⁴C screwdriver and remove from the shaft.
- ▶ Release the green holding sleeve **22**.
- ▶ Push the green holding sleeve **22** with the retaining tongues in the direction of the instrument's working end and remove from the screwdriver shaft.


Fig. 15 Disassembling the S⁴C screwdriver

S⁴C navigated drill sleeve Ø 3.5 mm (FW664R)/Ø 4.0 mm (FW665R)

- ▶ Turn guide sleeve **24** clockwise and remove. Be aware that it is a left-handed thread.
- ▶ Loosen knob **25** counterclockwise, unscrew and remove from the handle by pulling in the direction of the arrow.


Fig. 16 Disassembling the S⁴C drill guide

S⁴C navigated guide sleeve for smooth-shaft screws, FW658R

- ▶ Turn inner sleeve **26** to the "remove" position and remove from the outer sleeve in the direction of the arrow.
- ▶ Loosen knob **25** counterclockwise, unscrew and remove from the handle by pulling in the direction of the arrow.


Fig. 17 Disassembling the guide sleeve for smoothshank screws

S⁴C navigated guide sleeve Ø 3.5 mm/4.0 mm FW660R

- ▶ Loosen inner sleeve **27** counterclockwise, unscrew and remove from the handle by pulling in the direction of the arrow.
- ▶ Loosen knob **25** counterclockwise, unscrew and remove from the handle by pulling in the direction of the arrow.


Fig. 18 Disassembling the guide sleeve Ø 3.5 mm/4.0 mm

Assembling

S⁴C navigated tapper Ø 3.5 mm, FW655R

- ▶ Push guide sleeve **17** onto the S⁴C tapper **A** in the direction of the arrow and turn nut **19** clockwise to tighten.

Note

Spring **18** is fixed on shaft.


Fig. 19 Mounting the S⁴C tapper (Tap)

S⁴C navigated screwdriver FW656R

CAUTION

The S⁴C screwdriver will not function properly if the retaining tongues are bent or kinked!

- ▶ Do not bend or crimp the retaining tongues.

- ▶ Push green holding sleeve **22** with retaining tongues **21** onto the screwdriver shaft so that the retaining tongues engage in the grooves **23** of the screwdriver shaft.
- ▶ Retract and hold green holding sleeve **22** against the spring pressure.
- ▶ Screw the screw sleeve **20** onto the working end of the S⁴C screwdriver **D** and tighten.
- ▶ Release the green holding sleeve **22**.


Fig. 20 Mounting the S⁴C screwdriver

Aesculap®

S⁴ Cervical – Navigated Instruments

S⁴C navigated drill sleeve Ø 3.5 mm (FW664R)/Ø 4.0 mm (FW665R)

- Screw on guide sleeve **24** counterclockwise. Be aware that it is a left-handed thread. You will hear and feel the guide sleeve clicking into position every half turn.
- Push knob **25** into the handle of the instrument, screw in clockwise and tighten.


Fig. 21 Mounting the S⁴C drill guide

S⁴C navigated guide sleeve for smooth-shaft screws, FW658R

- Insert inner sleeve **26** into the outer sleeve in the direction of the arrow in the "remove" position.
- Then turn to "closed" position.
- Push knob **25** into the handle of the instrument, screw in clockwise and tighten.


Fig. 22 Mounting the guide sleeve for smoothshank screws

S⁴C navigated guide sleeve Ø 3.5 mm/4.0 mm, FW660R

- Push inner sleeve **27** into the outer sleeve in the direction of the arrow, screw in clockwise and tighten.
- Push knob **25** into the handle of the instrument, screw in clockwise and tighten.


Fig. 23 Mounting the guide sleeve Ø 3.5 mm/4.0 mm

Validated reprocessing procedure

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

Single-use products

Art. no.	Designation
FW086SU	Drill Ø 2.9 mm for smooth-shaft screws
FW088SU	Favored Angle screw drill Ø 2.9 mm for Ø 4 mm screws
FW662SU	S ⁴ C1/C2 drill Ø 2.4 mm
FW666SU	Drill, Ø 2.4 mm
FW667SU	Drill, Ø 2.9 mm

WARNING

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!

► Do not reprocess the product!

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Open up products with hinges.
- Open all valves/faucets.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Dismantle the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning agents and disinfectants according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable temperature of 55 °C.

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Note

For cleaning and disinfecting the reference star unit (pre-calibrated) and the reference star unit ML (calibration with ICM4), see Brainlab user manual.

Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ FW652R–FW658R ■ FW660R ■ FW663R–FW665R ■ FW688R–FW689R	■ Cleaning brush: 30 mm/Ø: 4.5 mm, e.g. TA011944 and cleaning brush: 20 mm/Ø: 2.5 mm, e.g. TE654202 and cleaning brush: 50 mm/Ø: 10 mm, e.g. TA007747 ■ 20 ml disposable syringe ■ Keep working ends open for cleaning. ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with immersion disinfection
Manual cleaning with ultrasound and immersion disinfection ■ FJ985R ■ FW067R ■ FW165R ■ FW661R	■ Cleaning brush: 30 mm/Ø: 4.5 mm, e.g. TA011944 and cleaning brush: 20 mm/Ø: 2.5 mm, e.g. TE654202 and cleaning brush: 50 mm/Ø: 10 mm, e.g. TA011327 ■ 20 ml disposable syringe ■ Keep working ends open for cleaning. ■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning. ■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and sub-chapter: ■ Chapter Manual cleaning with ultrasound and immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ FW657R ■ FW663R ■ FW688R–FW689R	■ Connect components with lumens and channels directly to the rinsing port of the injector carriage. ■ Keep working ends open for cleaning. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfecting and sub-chapter: ■ Chapter Mechanical alkaline cleaning and thermal disinfecting
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ FW067R ■ FW165R ■ FW652R ■ FW655R–FW656R ■ FW658R ■ FW660R ■ FW664R–FW665R	■ Cleaning brush: 30 mm/Ø: 4.5 mm, e.g. TA011944 and cleaning brush: 20 mm/Ø: 2.5 mm, e.g. TE654202 and cleaning brush: 50 mm/Ø: 10 mm, e.g. TA007747 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Keep working ends open for cleaning. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Validated procedure	Specific requirements	Reference
Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection ■ FJ985R ■ FW661R	■ Cleaning brush: 30 mm/Ø: 4.5 mm, e.g. TA011944 and cleaning brush: 20 mm/Ø: 2.5 mm, e.g. TE654202 and cleaning brush: 50 mm/Ø: 10 mm, e.g. TA007747 ■ 20 ml disposable syringe ■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots). ■ Keep working ends open for cleaning. ■ Place products in the tray with their hinges open.	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: ■ Chapter Manual pre-cleaning with ultrasound and brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces).
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Manual cleaning with ultrasound and immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfection	RT (cold)	15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Drain any remaining water fully.

Phase III

- Fully immerse the product in the disinfectant solution.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
- Rinse lumens at least five times at the beginning of the exposure time with an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
- Rinse lumens with an appropriate disposable syringe at least five times.
- Drain any remaining water fully.

Phase V

- Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfecting

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Aesculap®

S⁴ Cervical – Navigated Instruments

Mechanical cleaning/disinfection with manual pre-cleaning

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfectant cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	>15	2	D-W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I

- Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 min. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
- If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.

- Mobilize non-rigid components, such as set screws, links, etc. during cleaning.

- Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II

- Rinse/flush the product thoroughly (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Aesculap®

S⁴ Cervical – Navigated Instruments

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

► Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

Note

Separable products must only be sterilized after they have been dismantled.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
 - Disassemble the product
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.
- Store sterile packed single-use products dust-protected in a dry, dark and temperature-controlled room.

Technical Service

Risk of injury and/or malfunction!

► **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95-1602

Fax: +49 (7461) 16-5621

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

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