- (GB)** Instructions for use/Technical description  
S<sup>4</sup>® Cervical Spinal System – Occipital Torque Wrench FW103R
- (USA)** S<sup>4</sup>® Cervical Spinal System – Okziput Drehmomentschlüssel FW103R
- (D)** Gebrauchsanweisung/Technische Beschreibung  
S<sup>4</sup>® Cervical Spinal System – Okziput Drehmomentschlüssel FW103R
- (F)** Mode d'emploi/Description technique  
Système spinal cervical S<sup>4</sup>® – clé dynamométrique occipitale FW103R
- (E)** Instrucciones de manejo/Descripción técnica  
Sistema espinal cervical S<sup>4</sup>® – Llave de torque occipital FW103R
- (I)** Istruzioni per l'uso/Descrizione tecnica  
Sistema spinale cervicale S<sup>4</sup>® – Chiave dinamometrica occipitale FW103R
- (P)** Instruções de utilização/Descrição técnica  
S<sup>4</sup>® Cervical Spinal System – Chave dinamométrica occipital FW103R
- (NL)** Gebruiksaanwijzing/Technische beschrijving  
S<sup>4</sup>®-cervicaal wervelkollomsysteem – Occipitale momentsleutel FW103R
- (S)** Bruksanvisning/Teknisk beskrivning  
S<sup>4</sup>® Cervikalt spinalsystem – Occipital momentnyckel FW103R
- (RUS)** Инструкция по применению/Техническое описание  
S<sup>4</sup>® Cervical Spinal System — окципитальный динамометрический ключ FW103R
- (CZ)** Návod k použití/Technický popis  
S<sup>4</sup>® Krční páteřní systém – okcipitální klíč s měřením krouticího momentu FW103R
- (PL)** Instrukcja użytkowania/Opis techniczny  
System do kręgosłupa szyjnego S<sup>4</sup>® – Potyliczny klucz manometryczny FW103R
- (SK)** Návod na používanie/Technický opis  
Cervikálny spinálny systém S<sup>4</sup>® – okcipitálny momentový kľúč FW103R
- (TR)** Kullanım Kılavuzu/Teknik açıklama  
S<sup>4</sup>® Servikal Spinal Sistem – Okcipital tork anahtarı FW103R

**B | BRAUN**  
SHARING EXPERTISE


Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

TA-Nr. 012367 04/13 V6 Änd.-Nr. 45682

**CE** 0123 – DIR 93/42/EEC

Technical alterations reserved

Aesculap®

S<sup>4</sup>® Cervical Spinal System – Occipital Torque Wrench FW103R

Legend

- 1 Thumb nut
- 2 Hexagonal bit

Symbols on product and packages

Caution, general warning symbol

Caution, see documentation supplied with the product

Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at [www.extranet.bbraun.com](http://www.extranet.bbraun.com)

Intended use

The S<sup>4</sup> Cervical Spinal System torque wrenches with click function are used to tighten the S<sup>4</sup> Cervical System set screws with a defined torque:

- 2.8 Nm for tightening the set screw in the screw head/hook head/occipital plate/transverse connector and other connectors.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

WARNING

**Risk of injury caused by incorrect operation of the product!**

- Attend appropriate product training before using the product.
- For information about product training, please contact your national B. Braun/ Aesculap agency.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Store any new or unused products in a dry, clean, and safe place.
- Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

- Make sure that the instrument is correctly mounted and that the hexagonal bit is firmly fixed.
- Tighten the set screws, turning the torque wrench evenly in a clockwise direction. In doing so, make sure that the torque wrench is turned without jerking.
- Instrument will click when the defined torque is obtained.

WARNING

**Risk of injury/correction loss due to improper mounting of the set screw!**

- Set the set screw in place correctly.
- Make certain that the rods are correctly positioned in the body of the polyaxial screw head.
- When tightening the set screw always apply counter torque using the countering instrument intended for this purpose.

WARNING

**Risk of injury due to improper mounting of the transverse connector!**

- When tightening the set screw make sure that the transverse connector is completely seated on the rods.
- Fully tighten the set screw on the transverse connector using the torque wrench.

CAUTION

**Risk of damage to the implant by applying excess torque to the set screw!**

- Slowly tighten the set screw using the torque wrench, until the torque wrench clicks.

CAUTION

**Risk of damage to the torque wrench due to improper usage!**

- To undo a tightened set screw always use the set screw removal screwdriver FW064R.

CAUTION

**Damage to the countering instrument due to improper usage!**

- Always slide the countering instrument over the rod receptacle head as far as it will go. In doing so, make sure that the implant entirely engages into the slot at the working tip.
- Make sure that the rod protrudes on both sides of the countering instrument.

Disassembling

- Remove the hexagonal bit 2 by unscrewing the thumb nut 1 from the rear of the wrench head. Take care not to lose the hexagonal bit and the thumb screw.

Assembling

- After cleaning re-assemble the hexagonal bit 2 in the torque wrench by inserting it into the bit cavity on the torque wrench head.
- Engage the thumb nut 1 on the rear of the wrench head and securely tighten.

Validated reprocessing procedure

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this. The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at [www.extranet.bbraun.com](http://www.extranet.bbraun.com)

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Disassemble the product immediately after use, as described in the respective instructions for use.

Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

Preparation before cleaning

- Disassemble the product prior to cleaning, see Disassembling.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

CAUTION

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
  - are approved for (e.g. aluminum, plastics, high-grade steel),
  - do not attack softeners (e.g. in silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable washing temperature of 55 °C.

- Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- Carry out ultrasound cleaning:
  - as an effective mechanical supplement to manual cleaning/disinfecting.
  - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
  - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
  - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Validated cleaning and disinfection procedure

| Validated procedure                                                                                 | Specific requirements                                                                                                                                                                                                                                                | Reference                                                                                                                                                                                                                                             |
|-----------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Manual cleaning with immersion disinfection                                                         | <ul style="list-style-type: none"><li>■ Cleaning brush: e.g. TA011944, TE654202, GK469200</li><li>■ When cleaning products with movable hinges, ensure that these are in an open position and, if applicable, move the joint while cleaning.</li></ul>               | Chapter Manual cleaning/disinfection and sub-chapter: <ul style="list-style-type: none"><li>■ Chapter Manual cleaning with immersion disinfection</li></ul>                                                                                           |
| Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection | <ul style="list-style-type: none"><li>■ Cleaning brush: e.g. TA011944, TE654202, GK469200</li><li>■ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots).</li><li>■ Place products in the tray with their hinges open.</li></ul> | Chapter Mechanical cleaning/disinfection with manual pre-cleaning and sub-chapter: <ul style="list-style-type: none"><li>■ Chapter Manual pre-cleaning with a brush</li><li>■ Chapter Mechanical alkaline cleaning and thermal disinfecting</li></ul> |

Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning /disinfection process if necessary.

Manual cleaning with immersion disinfection

| Phase | Step                     | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                             |
|-------|--------------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------------|
| I     | Disinfecting<br>cleaning | RT (cold)    | >15        | 2            | D–W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| II    | Intermediate<br>rinse    | RT (cold)    | 1          | –            | D–W              | –                                                                    |
| III   | Disinfection             | RT (cold)    | 15         | 2            | D–W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| IV    | Final rinse              | RT (cold)    | 1          | –            | FD–W             | –                                                                    |
| V     | Drying                   | RT           | –          | –            | –                | –                                                                    |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
RT: Room temperature  
\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

**Phase I**  
► Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.  
► Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.  
► If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.  
► Mobilize non-rigid components, such as set screws, links, etc. during cleaning.  
► Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

**Phase II**  
► Rinse/flush the product thoroughly (all accessible surfaces) under running water.  
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.  
► Drain any remaining water fully.

**Phase III**  
► Fully immerse the product in the disinfectant solution.  
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.  
► Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

**Phase IV**  
► Rinse/flush the product thoroughly (all accessible surfaces).  
► Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.  
► Rinse lumens with an appropriate disposable syringe at least five times.  
► Drain any remaining water fully.

**Phase V**  
► Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

Mechanical cleaning/disinfection with manual pre-cleaning

*Note*  
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

*Note*  
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

| Phase | Step                     | T<br>[°C/°F] | t<br>[min] | Conc.<br>[%] | Water<br>quality | Chemical                                                             |
|-------|--------------------------|--------------|------------|--------------|------------------|----------------------------------------------------------------------|
| I     | Disinfectant<br>cleaning | RT (cold)    | >15        | 2            | D–W              | Aldehyde-free, phenol-free,<br>and QUAT-free concentrate,<br>pH ~ 9* |
| II    | Rinsing                  | RT (cold)    | 1          | –            | D–W              | –                                                                    |

D–W: Drinking water  
RT: Room temperature  
\*Recommended: BBraun Stabimed

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

**Phase I**  
► Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.  
► Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.  
► If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.  
► Mobilize non-rigid components, such as set screws, links, etc. during cleaning.  
► Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

**Phase II**  
► Rinse/flush the product thoroughly (all accessible surfaces) under running water.  
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

| Phase | Step                  | T<br>[°C/°F] | t<br>[min] | Water<br>quality | Chemical                                                                                                     |
|-------|-----------------------|--------------|------------|------------------|--------------------------------------------------------------------------------------------------------------|
| I     | Prerinse              | <25/77       | 3          | D–W              | –                                                                                                            |
| II    | Cleaning              | 55/131       | 10         | FD–W             | ■ Concentrate, alkaline:<br>– pH = 13<br>– <5 % anionic surfactant<br>■ 0.5 % working solution<br>– pH = 11* |
| III   | Intermediate<br>rinse | >10/50       | 1          | FD–W             | –                                                                                                            |
| IV    | Thermal disinfecting  | 90/194       | 5          | FD–W             | –                                                                                                            |
| V     | Drying                | –            | –          | –                | According to the program for cleaning and disinfection device                                                |

D–W: Drinking water  
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)  
\*Recommended: BBraun Helimatic Cleaner alkaline

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

**Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!**

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or damp.
- Repeat cleaning and disinfection of products that still show impurities or contamination.
- Check that the product functions correctly.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- Assemble dismountable products, see Assembling.
- Check for compatibility with associated products.

Packaging

- Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

*Note*  
The product can be sterilized either in disassembled or in assembled condition.

- Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- Validated sterilization process
  - Steam sterilization through fractionated vacuum process
  - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
  - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10<sup>-6</sup>, Aesculap recommends the following parameters:

| Aesculap Orga Tray/Sterile container (perforated bottom)<br>Minimum cycle parameters* |               |       |                     |
|---------------------------------------------------------------------------------------|---------------|-------|---------------------|
| Sterilization method                                                                  | Temp.         | Time  | Minimum drying time |
| Prevacuum                                                                             | 270 °F/275 °F | 4 min | 20 min              |

\*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

**WARNING for the US market**  
If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Maintenance

The torque wrenches are calibrated by the manufacturer to 2.8 Nm.  
► After a period of no more than 6 months send the torque wrench to Aesculap for re-calibration.

## Technical Service

**Risk of injury and/or malfunction!**  
► Do not modify the product.

► For service and repairs, please contact your national B. Braun/Aesculap agency.  
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

### Service addresses

Aesculap Technischer Service  
Am Aesculap-Platz  
78532 Tuttlingen / Germany  
Phone: +49 (7461) 95-1602  
Fax: +49 (7461) 16-5621  
E-Mail: ats@aesculap.de

Or in the US:  
Aesculap Implant Systems LLC  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042  
Aesculap Repair Hotline  
Phone: +1 (800) 214-3392  
Fax: +1 (314) 895-4420  
Other service addresses can be obtained from the address indicated above.

## Disposal

► Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

## Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC  
3773 Corporate Parkway  
Center Valley, PA, 18034,  
USA

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