


Aesculap Spine


GB	Instructions for use/Technical description
USA	S⁴ instruments
D	Gebrauchsanweisung/Technische Beschreibung S⁴-Instrumente
F	Mode d'emploi/Description technique Instruments S⁴
E	Instrucciones de manejo/Descripción técnica Instrumental S⁴
I	Istruzioni per l'uso/Descrizione tecnica Strumenti S⁴
P	Instruções de utilização/Descrição técnica Instrumentos S⁴
NL	Gebruiksaanwijzing/Technische beschrijving S⁴-instrumenten
CZ	Návod k použití/Technický popis S⁴-instrumenty

11


12

13


14

15


16

17


18

19


20

21

22

23

24

25

26

S⁴ instruments

Legend

- 1 Set screw starter FW177R
- 2 Counter torque FW178R (with L-handle)
- 3 Rod persuader FW208R
- 4 Reduction lever (consisting of lever threadpipe FW233R, screwdriver FW234R and outer sleeve FW235R)
- 5 Reduction lever percutaneous (consisting of lever threadpipe FW229R, screwdriver FW228R, outer sleeve with long slot FW353R and spacer FW141P)
- 6 Counter torque FW236R
- 7 Tommy bar FW237R
- 8 Distractor (consisting of distractor FW238R and two distraction arms FW239R)
- 9 Rod insertion instrument FW240R
- 10 Distraction spindle FW241R (consisting of distractor and extension)
- 11 Rod length measuring instrument FW242R
- 12 Screw length measuring instrument FW351R
- 13 Slotted hammer FW243R
- 14 Tissue protection sleeve FW244P
- 15 Tissue protection sleeve FW355P
- 16 K-wires FW257S and FW247S
- 17 Trocar FW271R
- 18 K-wire aiming device FW258R
- 19 Handle (for removal of K-wire aiming device) FW274R
- 20 K-wire protection sleeve FW352R
- 21 Cannulated monoaxial screwdriver FW262R
- 22 Pedicle dissector FW263R
- 23 Cannulated screw tabs FW264R, FW265R, FW266R, FW267R and FW268R
- 24 Cannulated polyaxial screwdriver FW270R
- 25 Dilation tube FW272R
- 26 Dilation tube FW354R

Symbols on product and packages

	Not for reuse in intended applications as defined by the manufacturer
	Caution: General warning symbol Caution: See documentation supplied with the product

Intended use

The S⁴ instruments are used for the implantation of the S⁴ implants with the following systems:

- S⁴ Spinal System
- S⁴ CS Cannulated System
- S⁴ FRI/S⁴ FRI percutaneous

The S⁴ fracture reduction instruments are used for the reduction and stabilization of fractures with the cannulated S⁴ implants of the S⁴ spinal system.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

DANGER

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product after use!

➤ Do not reuse the product.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.
- To avoid damage to the working end: Exercise caution when passing products through the working channel (e.g. trocar).
- Read, observe and keep the operating manuals S⁴ Spinal System Instrumentation, S⁴ FRI and S⁴ SRI.

Safe operation

DANGER

Danger of death of the patient caused by pushing the guide wire into the aorta!

- When inserting the S⁴ implants or S⁴ instruments with the guide wire, always observe the markings on the guide wire.
- Hold the guide wire with a forceps.
- Remove the guide wire at the right moment.

WARNING

Excessive trocar insertion depth due to missing end stop!

- Use trocar FW271R only in combination with K-wire aiming device FW258R.

WARNING

Risk of injury and/or malfunction!

- Always carry out a function check before using the product.

S⁴ instruments

Disassembling


Fig. 1

- Dismantle the set screw starter with an SW6 open-ended spanner, see Fig. 1.

Note

The set screw starter does not need to be dismantled for cleaning.


Fig. 2

Legend

A Button

- To dismount reduction lever **4**, see Fig. 2:
Firmly press button **A** and remove screwdriver FW234R.


Fig. 3

- Dismount reduction lever percutaneous **5**, see Fig. 3.
Disassemble it completely and remove the lock nut.


Fig. 4

Legend

B Slide

- To dismount distractor **8** with the distraction arms, see Fig. 4:
Actuate slide **B** and remove distraction arm FW239R.


Fig. 5

- To dismount rod insertion instrument 9, see Fig. 5: Unscrew the inner part.


Fig. 6

- Dismount distraction spindle 10, see Fig. 6.


Fig. 7

- Dismount rod length measuring instrument 11, see Fig. 7.


Fig. 8

- Dismount the handle, see Fig. 8.

S⁴ instruments

Assembling

Assembling the set screw starter


Fig. 9

- Screw together the set screw and the starter so that they are hand-tight, see Fig. 9.

Note

The set screw starter is not intended for loosening set screws. If the working end becomes loose, it can be tightened until hand-tight, without tools. Use the revision screwdriver FW193R to loosen the set screws.

Mounting the reduction lever


Fig. 10

- Mount reduction lever 4, see Fig. 10.

Mounting the reduction lever percutaneous


Fig. 11

Legend

C Lock nut

- Screw lock nut C onto the lever threadpipe, see Fig. 11.


Fig. 12

- Insert the screwdriver in the lever threadpipe, secure it with the spacer and install locking screw D, see Fig. 12.

Fig. 13

- Insert the subassembly into the FRI outer sleeve, see Fig. 13.

Mounting other product components


Fig. 14

- Mount distractor 8 with the distraction arms, see Fig. 14.


Fig. 15

- Mount rod insertion instrument 9, see Fig. 15.


Fig. 16

- Mount distraction spindle 10, see Fig. 16.


Fig. 17

- Mount rod length measuring instrument 11, see Fig. 17.


Fig. 18

- Mount handle, see Fig. 18.

S⁴ instruments

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, the manufacturer's specifications can only serve as an approximate guide for assessing the processing procedures applied by the individual operator/processors.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

Single-use products

K-wires FW247S and FW257S

DANGER

Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product after use!

➤ **Do not reuse the product.**

- Do not reprocess single-use products following their direct or indirect contamination with blood.
- Do not reprocess or reuse any used, damaged or contaminated single-use products.

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. To avoid this, the time interval between application and processing should not exceed 6 h, and neither fixating pre-cleaning temperatures >45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process can already cause visible surface changes.
- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use process chemicals that cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, Publications Red Brochure – Proper maintenance of instruments.
- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Rinse surfaces inaccessible to visual inspection, e.g. on products with hidden crevices or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.

Preparation prior to cleaning

- Have the dry product in a disposal container ready for cleaning and disinfecting within 30 min after use.
- Disassemble dismountable products, see Disassembling.

Cleaning/Disinfecting

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use a cleaning and disinfecting agent according to the manufacturer's instructions. The cleaning and disinfecting agent
 - must be approved for the product materials (e.g. aluminum, plastics, high-grade steel)
 - must not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.

S⁴ instruments

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.
- Clean and disinfect microsurgical products mechanically, provided they can be securely fixed in machines or storage devices in such a way that they will be thoroughly cleaned.

Manual cleaning/disinfecting

- Keep working ends open during cleaning.
- Clean products with movable links with the links opened or mobilized.
- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfection

FW141P, FW228R, FW229R, FW233R–FW235R, FW238R, FW239R, FW241R, FW244P, FW258R, FW262R–FW268R, FW270R–FW272R, FW274R, FW351R–355P

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Disinfecting	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	0.5	–	FD–W	–
V	Drying	RT	–	–	–	–

D–W: Drinking water

FD–W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make sure that all accessible surfaces are moistened.
- Clean the product under running tap water, if necessary with a suitable cleaning brush (TA011944/TE654202, GK469200), until all visible residues have been removed from the product surfaces.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometry, for at least 1 min or until no more residues can be removed. Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning/disinfecting solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution. Make sure that all accessible surfaces are moistened.

Stage IV

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

S⁴ instruments

Manual ultrasonic cleaning and immersion disinfection

FW177R, FW208R, FW240R, FW242R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	5	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	1st Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	2nd Intermediate rinse	RT	1	-	D-W	-
V	Final rinse	RT (cold)	0.5	-	FD-W	-
VI	Drying	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Remove the jaw protection and clean the product with an appropriate cleaning brush (TA no. 011944/TE654202, GK469200) until all visible residues have been removed from the product surfaces.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometry, for at least 1 min or until no more residues can be removed. Mobilize non-rigid components, such as set screws, links, etc. during cleaning.

- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution. Make sure that all accessible surfaces are moistened.

Stage IV

- Rinse the product completely (all accessible surfaces) under running water.

Stage V

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage VI

- Thoroughly dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHH or FDA approval/license or CE mark according to DIN EN ISO 15883).

Note

For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that A_0 is $>3\ 000$ for the process.

Note

The disinfectant used for processing must be serviced and checked at regular intervals.

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and canals directly to the special rinsing port of the injector cart.
- Keep working ends open during cleaning.
- Ensure that product links and joints are opened when placing the product on the tray.

S⁴ instruments

Mechanical alkaline cleaning and thermal disinfecting

FW141P, FW234R, FW238R, FW239R, FW241R, FW271R

Machine type: Single-chamber washer/disinfector without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER alkaline with tensides, applica- tion solution 0.5 %
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfectant pro- gram

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

Mechanical cleaning/disinfecting with manual pre-cleaning

Note

The disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval/license or CE mark according to DIN EN ISO 15883).

Note

The disinfectant used for processing must be serviced and inspected at regular intervals.

Note

For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that A_0 is $>3\,000$ for the process.

Manual pre-cleaning with brush

FW178R, FW208R, FW228R, FW229R, FW233R, FW235R, FW244P, FW258R, FW262R–FW268R, FW270R, FW272R, FW274R, FW351R–FW355P

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Irrigation	RT (cold)	1	–	D–W	–

D–W: Drinking water

RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make sure that all accessible surfaces are moistened.
- Clean the product with an appropriate cleaning brush (TA no. 011944/TE654202, GK469200) until all visible residues have been removed from the product surfaces.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometry, for at least 1 min or until no more residues can be removed. Mobilize non-rigid components, such as set screws, links, etc. during cleaning.

- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

S⁴ instruments

Manual pre-cleaning with ultrasound and brush

FW177R, FW240R, FW242R

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Irrigation	RT (cold)	1	-	D-W	-

D-W: Drinking water

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with an appropriate cleaning brush (TA no. 011944/TE654202, GK469200) until all visible residues have been removed from the product surfaces.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometry, for at least 1 min or until no more residues can be removed. Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly rinse through these components (at least five times) with the cleaning solution, using a disposable syringe (20 ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfector without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Connect components with lumens and canals directly to the special rinsing port of the injector cart.
- Keep working ends open during cleaning.
- Place the product on the tray with all product links and joints open.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER alkaline with tensides, applica- tion solution 0.5 %
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfector pro- gram

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

S⁴ instruments

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Assemble the separable product, see Assembling.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.

- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

Note

The product can be sterilized either in disassembled or in assembled condition.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterilcontainer (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F— 275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap Sterilcontainer cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.
- Store single-use products in sterile packaging in dust-protected conditions in a dry, dark and temperature-controlled room.

Technical Service

Risk of injury and/or malfunction!

- **Do not modify the product.**

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-Mail: ats@aesculap.de

Or in the US:

Aesculap Implant Systems Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Accessories/Spare parts

Designation	Art.-no.
Replacement nut (for lever threadpipe FW233R, distractor FW238R and distraction spindle FW241R)	FW233201

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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- S⁴ Spinal System
- S⁴ CS – Cannulated System
- S⁴ FRI/S⁴ FRI perkutan

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Aesculap AG