Aesculap Spine

- GB

USA

D

F

E

I

P

NL

CZ
- Instructions for use

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S⁴ Spinal System – augmentation screw

Note
Generally, the S⁴ Spinal System – augmentation screw is covered by the instructions for use of the S⁴ Spinal System – Lumbar/Deformity TA011187! The present information on the augmentation screw is intended to complement the respective information contained in the instructions for use of the S⁴ Spinal System – Lumbar/Deformity.

Symbols on product

Symbol	Explanation
	Sterilization using ethylene oxide
	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Attention, see instructions for use

Intended use

The S⁴ Spinal System implants are used for dorsal monosegmental and multisegmental stabilization of the lumbar and thoracic spine. The injection cannula is inserted in the S⁴ Spinal System – augmentation screw to apply bone cement.

System components

- The S⁴ Spinal System – augmentation screw comprises:
- Monoaxial/ polyaxial screws (augmentation screw), supplied in sterile condition
 - Cement injections cannula (sterile)
 - Extension tube for cement injections cannula (sterile)

Note
There are special S⁴ instruments provided for the implantation of these system components and for the augmentation, distraction, compression and reduction of the lumbar and thoracic spine.

Materials

The materials used in the implant are listed on the packaging:

Materials
ISOTAN®F titanium forged alloy Ti6Al4V according to ISO 5832-3
ISOTAN®P pure titanium according to ISO 5832-2

The titanium implants are coated with a colored oxide. Possible, slight color changes do not affect the quality of the implant. ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Safe handling and preparation

CAUTION
Federal law restricts this device to sale by or on order of a physician!

Risk of lasting damage caused by non-compliance with the instructions for use of the S⁴ Spinal System – Lumbar/Deformity TA011187!
➤ **Always adhere to the instructions for use of the S⁴ Spinal System – Lumbar/Deformity TA011187.**

- The product has been radiation/EO-sterilized and is supplied in sterile packaging.
- The product must not be reused.
- Read, follow and keep the instructions for use.
 - Read, follow and keep the instructions for use of all components applied.
 - Use the product only in accordance with its intended use, see Intended use.
 - Do not use products from open or damaged sterile packaging.
 - Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Do not reprocess the product.
 - Do not use the product after its use-by date.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Single-use products

Possible damage to/or destruction of products, due to processing!
Intended use: These products have been designed for single use only.
➤ **Do not reprocess these products.**

Indications

Surgically installed implants serve to support normal healing processes. They are not supposed to replace normal body structures or to support permanently loads occurring in cases where healing did not take place. Use for:

- Fractures
- Pseudarthrosis or delayed healing
- Dislocation
- Degenerative instability
- Post-discectomy syndrome
- Spondylolisthesis
- Kyphosis
- Degenerative scoliosis
- Stenosis
- Post-trauma instability
- Tumors

Risk of screw breakage caused by pedicle screws used for spondylolisthesis!
➤ **Support the stabilization measure with an intercorporeal fusion (ALIF, PLIF or TLIF).**

Contraindications

- Do not use in the presence of:
- Fever
 - Acute or chronic vertebral infections of a local or systemic nature
 - Severe osteoporosis or osteopenia:
 - The operating surgeon has to weigh possible risks affecting the success of treatment with the S⁴ Spinal System – augmentation screw.
 - In cases of foreseeable reduced stability of the standard S⁴ Spinal System – augmentation screw anchoring, additional stability can be achieved by fixation with bone cement.
 - Pregnancy
 - Infection
 - systemic
 - in the spine
 - local
 - Medical or surgical conditions that could negatively affect the success of the implantation
 - Mental illness
 - Dependency on pharmaceutical drugs, drug abuse, or alcoholism
 - Adiposity
 - Severely damaged bone structures that could prevent stable implantation of implant components
 - Neuromuscular disorders or illnesses
 - Bone tumors in the region of implant fixation
 - Wound healing disorders
 - Inadequate patient compliance
 - Foreign body sensitivity to the implant materials
 - Cases not listed under indications

Side-effects and adverse interactions

- Implant failure due to overloading:
- Bending
 - Loosening
 - Breakage
- Inadequate fixation
 - Failed or delayed fusion
 - Cement leakages (cement emerging from the vertebral body)
 - Infection
 - Fractured vertebral body or discs
- Injuries to:
- nerve roots
 - spinal cord
 - blood vessels
 - organs
- Changes of the normal spine lordosis
 - Impairment of the gastrointestinal, urological and/or reproductive systems
 - Pain or indisposition
 - Bursitis
 - Decreased bone density due to load avoidance
 - Bone atrophy/fracture above or below the spine section provided for
 - Limited performance
 - Persistence of symptoms that were to be treated by the implantation
 - Death

Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present instructions for use.
- The operating surgeon must have a thorough command of all hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with the bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.

- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Only combine Aesculap modular implant components with each other.
- Do not use damaged or surgically excised components under any circumstances.
- Implants that have already been used must not be reused.
- The operating surgeon is responsible for combining and implanting the implant components, with or without bone cement.
- Delayed healing can cause implant breakage due to metal fatigue.
- The physician in charge decides whether an implant component implanted should be removed.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Damage to the weight-bearing bone cement and/or bone structures can give rise to loosening of the components, fracture of the bone or implant, as well as to other grave complications.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation (min. dose 25 kGy).
- The implant components have been sterilized with ethylene oxide.
 - Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
 - Prior to use, check the product use-by date and verify the integrity of the sterile packaging. Do not use implant components that are past their use-by date or whose packaging is damaged.

Implants supplied in sterile condition must not be resterilized.
If the present instructions for use do not offer explicit information in this respect, the product user carries all responsibility for any resterilization!

Danger to the patient and possible loss of implant functionality due to resterilization!
➤ **Do not resterilize and reuse, under any circumstances, surface-coated implants (e.g. PLASMAPORE®), implants made of ceramics or polyethylene, or implants containing plastic components.**
➤ **Never resterilize used, explanted or damaged implants.**

Application

- The operating surgeon draws up an operating plan that specifies and appropriately documents the following steps:
- Selection of the implant components and their dimensions
 - Positioning of the implant components in the bone
 - Location of intraoperative landmarks
- The following conditions must be fulfilled prior to application:
- All requisite implant components are ready to hand
 - Operating conditions are highly aseptic
 - The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
 - The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
 - The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
 - The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The intervention has been explained to the patient, whose consent concerning the following information has been documented:

- Delayed healing or incomplete fusion can cause the implant to fracture or loosen as a result of the extreme load to which it is subjected.
- Smokers present an increased risk of fusion failure.
- The life-span of the implant depends on the patient's body weight.
- Implant components must not be subject to extreme load arising from sports activities or strenuous physical work.
- If the implant becomes loose, fractures or deviates from its correction position, a revision operation must be performed.
- The patient must undergo regular check-ups of the implant components, performed by a physician.

- The implantation of the S⁴ Spinal System implant (augmentation screw) involves the following steps:
- Select the appropriate S⁴ Spinal System variant and anchoring implants according to indication, preoperative planning and intraoperative bone situation.
 - When choosing the S⁴ Spinal System implant, be careful to select the correct diameter, length and canal orientation.

Incorrect application can result in trauma to the spinal cord and nerve roots!
➤ **Align and screw in the instruments, the pedicle screw and the polyaxial screw only under X-ray control or using a navigation system.**

- To avoid internal stress on, and weakening of the implant: avoid scoring or scratching of the implant components.
- Only allow axial load on the screws during screw insertion. Avoid induction of lateral forces during screw insertion.

Damage to the hexagon socket in the screw head due to incorrect application of the Allen key!

- Make certain the Allen key is fully inserted in the hexagon socket of the screw head when turning the screw.

Note

Cementing is carried out only after all screws have been aligned and the connecting rods have been adapted, and before locking the S⁴ Spinal System - augmentation screw!

Risk of breaking the injection cannula by inducing excessive shear forces!

- Do not apply excessive transverse load to the injection cannula.
- Do not bend the injection cannula.

Note

Once the S⁴ Spinal System - augmentation screw has been cemented, it will be impossible to change the implant position, insert the reduction lever or apply corrections!

Note

With the S⁴ Spinal System - augmentation screw finally cemented, explanting the screw can become extremely difficult, especially in cases of osteopenia.

- Prior to applying the bone cement, fully screw in the injection cannula at the thread in the bone screw head, without applying perceptible force.
- Check the bone cement for proper texture and viscosity before positioning the applicator on the injection cannula. Observe the manufacturer's advisories for the bone cement!
- Cement the S⁴ Spinal System - augmentation screw only with cement that is suitable for screw augmentation.

Risk of injury from cement leakages!

- Always apply imaging techniques for control during the preparation and positioning of the implants.
- Make certain the slots of the S⁴ Spinal System - augmentation screw are positioned fully inside the bone.
- Make sure that the pedicle cortex is not perforated.
- Do not penetrate the opposite cortex and ensure that the opposite cortex is undamaged.
- Do not insert or cement the S⁴ Spinal System - augmentation screw in a fractured or partly resected vertebral body. Check by preoperative thin-layer CT.
- Screw the injection cannula into the screw head without applying any perceptible force. Make certain the axes of the pedicle screw and the injection cannula are aligned with each other.
- Apply cement only under high-quality radiographic control (max. volume of cement per polyaxial screw: 2 ml).

Risk of injury due to material incompatibility!

- Only use application systems that have been approved for use with bone cement.

Contamination of the locking thread/insufficient fixation when removing the application system, due to incorrect viscosity of the cement!

- Make certain to be using cement of the correct viscosity.

- Do not allow the cement to cure completely before removing the injection cannula. Remove the cement application system incl. cannula as soon as the cement shows a doughy texture (allow for a certain amount of cement for reference purposes when preparing the cement).

Wound contamination due to insufficient curing of cement! Risk of caulking of the locking screw thread, resulting in system instability!

- Observe the correct curing period after removing the injection cannula.
- For monitoring the cement texture, prepare a reference volume of cement together with the cement to be applied.
- Make certain the system is not exposed to any load before the cement has reached its final strength.

Risk of connecting rod breakage!

- Do not bend back or bend excessively the connecting rods and cross connectors.

- If the screw head loosens during revision of the larger polyaxial screws (Ø 7.5 mm), remove the remaining screw component with the SW 3.5 screwdriver.

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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