en	Instructions for use/Technical description S ^{4®} Spinal System – Sterile-packaged
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapImplantsystems.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung S ^{4®} Spinal System – Steril verpackt
fr	Mode d'emploi/Description technique S ^{4®} Spinal System – Sous emballage stérile
es	Instrucciones de manejo/Descripción técnica S ^{4®} Spinal System en embalaje estéril
it	Istruzioni per l'uso/Descrizione tecnica S ^{4®} Spinal System – In confezione sterile
pt	Instruções de utilização/Descrição técnica S ^{4®} Spinal System – Em embalagem esterilizada
nl	Gebruiksaanwijzing/Technische beschrijving S ^{4®} Spinal System – Steriel verpakt
sv	Bruksanvisning/Teknisk beskrivning S ^{4®} Spinal System – sterilt förpackat
ru	Инструкция по применению/Техническое описание S ^{4®} Spinal System в стерильной упаковке
cs	Návod k použití/Technický popis S ^{4®} Spinal System – sterilně balený
pl	Instrukcja użytkowania/Opis techniczny S ^{4®} Spinal System – pakowany sterylnie
sk	Návod na použitie/Technický opis S ^{4®} Spinal System – sterilne balené
tr	Kullanım Kılavuzu/Teknik açıklama S ^{4®} Spinal System – Steril Paketlenmiş
zh	使用说明 / 技术说明 S ^{4®} 胸腰椎后路固定系统组件 – 无菌产品

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap® – a B. Braun brand

TA013579 2019-08 V6 Change No. 60104

en	
Aesculap®	
S ⁴ Spinal System – Sterile-packaged	
Symbols on product and packages	
STERILE R	Sterilization using irradiation
2	Not for reuse as per the manufacturer-defined intended use
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

Note
The S⁴ Spinal System – in sterile condition is addressed in general in the operating instructions for the S⁴ Spinal System – Lumbar/Deformity TA011187. This information on the sterile-packaged S⁴ implants supplements the respective information in the instructions for use of the S⁴ Spinal System – Lumbar/Deformity.

The S⁴ Spinal System implants are used for dorsal monosegmental and multisegmental stabilization of the lumbar and thoracic spine.
The parallel (closed and open) and axial rod connectors are connected to S⁴ Spinal System rods in order to connect a rod parallel or in a line with another rod.
The lateral offset connectors are connected to the S⁴ Spinal System rods in order to place a screw offset. The rod connectors thus extend the rod to the adjacent spinal column segments.
The S⁴ Spinal System – sterile-packaged comprises:
■ Rod connector – parallel (closed and open), axial and lateral offset connectors

Note
Special S⁴ instruments must be used for implanting these components, as well as for the distraction, compression and reduction of the lumbar and thoracic spine.

Materials

The materials used in the implant are listed on the packaging:
■ ISOTAN®F titanium forged alloy Ti6Al4V acc. to ISO 5832-3
The titanium implants are coated with a colored oxide. Slight changes in coloration may occur, but do not affect the implant quality.
ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen/Germany.

Indications

Surgically installed implants are designed to support the normal healing process. They are not intended for use either as replacements for natural body parts or to bear loads over the long term if healing does not occur.
Use for acute disorders of the spine that cannot be treated by other therapies::
■ Fractures
■ Pseudarthrosis or delayed healing
■ Dislocation
■ Degenerative instability
■ Post-discectomy syndrome
■ Spondylolisthesis
■ Kyphosis
■ Scoliosis
■ Stenosis
■ Post-trauma instability
■ Tumors
■ Osteoporosis (with the simultaneous use of bone cement and fenestrated screws)

WARNING

Risk of screw breakage caused by pedicle screws used for spondylolisthesis!
► Support the stabilization measure with interbody fusion (ALIF, PLIF or TLIF).

Absolute contraindications

Do not use in the presence of:
■ Severe damage to the bone structures of the spine that could prevent the stable implantation of the implant components; for example, osteopenia, severe osteoporosis, Paget's disease, bone tumors etc.
■ Metabolic or degenerative metabolic bone diseases that could compromise the stable anchoring of the implant system
■ Suspected allergy or sensitivity to the implant materials
■ Acute or chronic vertebral infections of a local or systemic nature
■ Poor patient compliance or limited ability to follow medical instructions, particularly in the post-op phase, including with regard to the restrictions on range of movement in terms of physical exercise and occupational activity
■ Cases not listed under indications

Relative contraindications

In the following circumstances, use of the implant system could represent an increased clinical risk and therefore requires precise, individual assessment by the surgeon:
■ Medical or surgical conditions that could negatively impact the success of the implantation, including wound healing disorders
■ Conditions that could subject the spine and implants to excessive pressure; for example, pregnancy, obesity, neuromuscular diseases or disorders
■ Generally poor condition of the patient; for example, drug or alcohol addiction

Side effects and interactions

The general risks of surgery are assumed known and are therefore not described in these instructions for use. Potential risks and clinical consequences associated with use of the implant system include, but are not limited to:
■ Implant failure resulting from excessive load
– Bending
– Loosening
– Breakage
■ Poor stabilization, requiring subsequent revision surgery or surgery to remove the implant system
■ Failed or delayed fusion
■ Infection
■ Vertebral fractures
■ Injuries to
– Nerve roots
– Spinal cord
– Blood vessels
– Organs
■ Changes of the normal spine lordosis
■ Neurological disorders, including impairment of the gastrointestinal, urological and/or reproductive system
■ Pain or indisposition
■ Bursitis
■ Decreased bone density due to load avoidance
■ Bone atrophy/fracture above or below the spine section provided with the implant
■ Limited physical ability
■ Persistence of symptoms that were to be treated by the implantation
■ Paresis
■ Subsequent surgery to remove the implant system
■ Metal/foreign body allergy or hypersensitivity
■ Dural tears
■ Meningitis

Safety notes

CAUTION

According to US federal law, this device may only be sold by or on the orders of a doctor!

WARNING

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with the bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user must follow the instructions that accompany each of the Aesculap implant components.
- The implant components have been tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically removed components.
- Implants already used must not be used again.
- Only combine Aesculap modular implant components with each other.
- Delayed healing can cause implant breakage due to metal fatigue.
- The surgeon in charge of the therapy shall make any decision with regard to the removal of implant components that have been used.
- Damage to the weight-bearing implant and/or bone structures can cause loosening of the components, fracture of the bone or implant, dislocation and migration and other grave complications.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged

WARNING

Damage to implants caused by processing and resterilization!
► Do not reprocess or resterilize the implants.

Applications

The operating surgeon shall devise an operation plan that specifies and precisely documents the following:
■ Selection of the implant components by type and size
■ Positioning of the implant components in the bone
■ Definition of the intraoperative orientation points
The following conditions must be fulfilled prior to application:
■ All requisite implant components are ready to hand
■ Highly aseptic operating conditions (not contaminated)
■ All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems
■ The surgeon and the operating team are thoroughly conversant with the surgical technique and the available range of implants and instruments; appropriate information materials are complete and ready to hand.
■ The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
■ The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- In cases of delayed or incomplete fusion, the implant may break or loosen as a result of excessive strains.
- The life span of the implant depends on the patient's body weight.
- Implant components must not be subject to extreme load arising from sports activities or strenuous physical work.
- Corrective surgery may be necessitated by implant loosening, fracture, or loss of correction.
- Smokers present an increased risk of bone fusion failure.
- The patient must undergo regular medical follow-up examinations of the implant components.
- The patient must be advised that being overweight places such excessive pressure on the implant that the prospect of failure becomes likely.
- The patient must be advised that strenuous physical activity in the form of heavy lifting, excessive muscle strain, exaggerated twisting, repeated bending and bending over or running should be avoided until the bones have healed fully.
- The patient must be instructed to keep pressure off the implant post surgery in order to prevent the onset of clinical problems associated with implant failure.

The implantation of the S⁴ Spinal System rod connectors requires the following steps:

- ▶ Select the appropriate S⁴ Spinal System variant and anchoring implants according to the individual indication, preoperative planning and bone situation found intraoperatively.
- ▶ To prevent internal stresses and weakening of the implants, avoid scoring or scratching any components.
- ▶ Position the rod connector on the rod at the site planned by the surgeon using the provided insertion/holding instrument.
- ▶ Do not unscrew the locking screws from the rod connector construct.
- ▶ Do not assemble the locking screws with the other rod connector components.
- ▶ Do not implant the rod connector if any locking screws are missing.

CAUTION

If the locking screws were removed or if they are missing, the rod connector cannot be used!

- ▶ In this case, choose a different rod connector for implantation.

- ▶ Before tightening up the locking screws, ensure that the rod connector is positioned correctly.
- ▶ For tightening and loosening the locking screws, always use the screwdriver and countering instrument intended for this purpose.
- ▶ Always tighten the locking screws with the 4 Nm torque wrench FW207R intended for this purpose.

WARNING

Risk of injury if the locking screws are not fixed sufficiently!

- ▶ Make certain that there is no patient tissue caught in the rod connector construct.
- ▶ Position the rod connector correctly.
- ▶ Ensure that the rods have been completely inserted into the rod connector.
- ▶ Always tighten the locking screws with the 4 Nm torque wrench FW207R intended for this purpose.

WARNING

Risk of danger if the locking screws are tightened too firmly or not firmly enough!

- ▶ Always tighten up the locking screws with the defined torque of exactly 4 Nm.

WARNING

Inadequate fixation due to incorrect rod position!

- ▶ Always position the rod connector in such a way that the rod is completely inserted into the rod connector.
- ▶ Always position rods in such a way that the hexagon or the tip is located outside the clamping area of the rod connector.

WARNING

If patient tissue is caught in the rod connector, there is insufficient fixing!

- ▶ Make certain that there is no patient tissue caught in the rod connector construct.
- ▶ Remove patient tissue from the rod connector construct or select another rod connector for implantation.

- ▶ When tightening or loosening, only subject the locking screws to axial pressure. Avoid exerting lateral forces when tightening or loosening the locking screws.

WARNING

If the screwdriver or torque wrench are used incorrectly, there is damage to the hexagon socket in the locking screws.

- ▶ Ensure that the tip of the hexagon of the screwdriver or torque wrench is positioned entirely within the hexagon socket of the locking screw when tightening or loosening.

More information on Aesculap implant systems is available from B. Braun/Aesculap or from the corresponding B. Braun/Aesculap office.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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