


- GB** Instructions for use/Technical description
USA Slim-line aneurysm clip applicator forceps for Aesculap YASARGIL aneurysm clips
D Gebrauchsanweisung/Technische Beschreibung
 Schlanke Aneurysmen-Clip-Anlegezange für Aesculap YASARGIL Aneurysmenclips
F Mode d'emploi/Description technique
 Pince de pose pour clips d'anévrisme fine Aesculap YASARGIL
E Instrucciones de manejo/Descripción técnica
 Pinza fina aplicadora de clips de aneurisma para clips de aneurismas YASARGIL de Aesculap
I Istruzioni per l'uso/Descrizione tecnica
 Pinza applicatrice sottile per clip da aneurisma Aesculap YASARGIL
P Instruções de utilização/Descrição técnica
 Pinça estreita de aplicação de cliques de aneurisma, para cliques de aneurisma da Aesculap YASARGIL
NL Gebruiksaanwijzing/Technische beschrijving
 Slanke aneurysmaclip-applicatietang voor Aesculap YASARGIL aneurysmaclips
S Bruksanvisning/Teknisk beskrivning
 Smal Aneurysm clip appliceringspincett för Aesculap YASARGIL aneurysmclip
RUS Инструкция по применению/Техническое описание
 Клиппаппликаторы Слим Лайн (Slim Line) для наложения клипсов сосудистых нейрохирургических Aesculap, далее именуемые клиппаппликаторы YASARGIL
CZ Návod k použití/Technický popis
 Štíhlé kleště k nakládání aneurysmatických svorek pro aneurysmatické svorky YASARGIL
PL Instrukcja użytkowania/Opis techniczny
 Wąskie kleszczyki do aplikacji zacisków tętniakowych Aesculap YASARGIL
SK Návod na použitie/Technický opis
 Úzke kliešte na zakladanie svoriek pre aneurizmy Aesculap YASARGIL
TR Kullanım Kılavuzu/Teknik açıklama
 Aesculap YASARGIL anevrizma klipsleri için ince anevrizma klipsi yerleştirme forsepsi
KR 사용 설명서 / 기술 설명
 Aesculap YASARGIL 동맥류 클립을 위한 슬림형 동맥류 클립 적용 포셉

B | BRAUN
 SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
 Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

TA-Nr. 015002 08/16 V6

CE 0123 - DIR 93/42/EEC

Aesculap®

Slim-line aneurysm clip applier forceps for Aesculap YASARGIL aneurysm clips

Legend

- 1 Jaw part
- 2 Branches left/right
- 3 Color-coding button
- 4 Lock with fastening screws
- 5 Clip applier forceps
- 6 Rotary plate
- 7 Vario clip applier forceps
- 8 Adjusting plates

Symbols on product and packages

Caution, general warning symbol
Caution, see documentation supplied with the product

Intended use

The slim-line aneurysm clip applier forceps for Aesculap YASARGIL aneurysm clips are used for opening, closing and applying Aesculap YASARGIL titanium aneurysm clips.

Note
In the following text, Aesculap YASARGIL titanium aneurysm clips are referred to as aneurysm clips.

Indications

Indications, see Intended use.

Note
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Contraindications

None known.

Risks and side effects

As part of the legal duty to inform, the following typical risks and side effects associated with the use of surgical instruments are referred to. These are predominantly procedure-specific, non product-specific and not limited to unwanted damage to surrounding tissue resulting in e.g. bleeding, infection, material incompatibilities or instrument parts remaining unnoticed in the patient, etc.

Available sizes

All slim-line aneurysm clip applier forceps are marked according to size (Mini or Standard) and clip material (titanium) so that their correct application with aneurysm clips of the appropriate size and material is ensured. For detailed information on the available sizes, see the YASARGIL Aneurysm Clip System brochure.

Safe handling and preparation

- CAUTION**
Federal law restricts this device to sale by, or on order of a physician!
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
 - Read, follow, and keep the instructions for use.
 - Use the product only in accordance with its intended use, see Intended use.
 - Clean and disinfect the new product mechanically after removing its transport packaging and prior to its initial sterilization.
 - Store any new or unused products in a dry, clean, and safe place.
 - Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
 - Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
 - Replace any damaged components immediately with original spare parts.

Safe operation

Damage to, imprecise functioning and incorrect closing force of the aneurysm clips caused by using the wrong applier forceps!
► Use aneurysm clips only with the matching applier forceps; note the respective color of the color-coding button.

- The instruments can be distinguished by means of the color-coding button 3 and labeling of the applier forceps.

Size	Color-coding button
Standard	Blue
Mini	Red

- Mini and Standard aneurysm clips must only be applied with the clip applier forceps of the appropriate size (Mini or Standard).
- Slim-line aneurysm clip applier forceps must only be used with YASARGIL titanium aneurysm clips from Aesculap.

Locking function

Malfunction or damage to the product by changing the lock!
► Do not re-mount a lock that has already been dismantled once.
► Do not manipulate the lock by bending.

Note
The clip applier forceps can be used with or without locking function.

Use with locking function

- Insert the aneurysm clip into jaw piece 1.
- Compress branches left/right 2 until the lock 4 engages.
The clip is in transfer position.
- Move clip applier forceps 5 with clip into the operating field.
- To release the lock, compress branches left/right 2.
The clip opens.
The lock automatically returns to the unlocked position.
- Fully release branches left/right 2.
The aneurysm clip closes.
- Detach the clip applier forceps 5 from the aneurysm clip.

Use without locking function

- Dismount the lock 4 by undoing the fastening screws.
- Insert the aneurysm clip into jaws 1.
- Compress branches left/right 2.
The inserted aneurysm clip opens.
- Fully release branches left/right 2.
The aneurysm clip closes.
- Detach the clip applier forceps 5 from the aneurysm clip.

Vario clip applier forceps

The Aesculap Vario clip applier forceps 7 offer an additional rotary function for positioning of the clip in the clip applier forceps.
The adjusting plates 8 of the Aesculap Vario clip applier forceps can be used for the adjustment in parallel of the rotary plates 6 in the jaw part. The adjusting plates (FT402800) can be used to turn the rotary plates 6.

Validated reprocessing procedure

General safety notes

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.
The specified chemistry was used for validation.

Note
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at www.extranet.bbraun.com.
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and reprocessing should not exceed 1 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.
Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result damage to metallic products. These must be removed by rinsing thoroughly with demineralized water and then drying.
Perform additional drying, if necessary.

- All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:
- Optical changes in materials, e.g., fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
 - Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.
 - Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
 - For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-l.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Disassembling the product before carrying out the reprocessing procedure

- Open up products with hinges.

Preparations at the place of use

- If applicable, rinse surfaces that are not accessible to visible inspection (preferably with demineralized water), using a disposable syringe, for example.
- Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- Put the wet product into a closed disposal container and have it transferred to cleaning and disinfecting within 1 h.

Cleaning/disinfection

Product-specific safety notes on the reprocessing procedure

DANGER

Danger to the patient!

- ▶ Reprocess the product only with manual pre-cleaning followed by mechanical cleaning.

DANGER

Risk to patient due to cross contamination!

- ▶ Do not clean contaminated products together with uncontaminated products in a tray.

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning agents and disinfectants according to the manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 55 °C.

- ▶ Use suitable cleaning/disinfecting agents if the product is put away in a wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.
- ▶ Clean and disinfect microsurgical products mechanically if they can be placed securely in the machine or on the positioning aids.

Mechanical cleaning/disinfection with manual pre-cleaning

Note
The cleaning and disinfecting machine must be of tested and approved effectiveness (e.g. meets the requirements of EN ISO 15883).

Note
The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Disinfectant cleaning	RT (cold)	>15	1	D–W	Dr. Weigert neodisher® SeptoClean*
II	Rinsing	RT (cold)	1	-	D–W	-

D–W: Drinking water
RT: Room temperature
*Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

- Phase I
- ▶ Fully immerse the product in the cleaning/disinfecting solution for at least 15 minutes. Ensure that all accessible surfaces are moistened.
 - ▶ Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed.
 - ▶ If applicable, brush non-visible surfaces for at least 1 minute with a suitable cleaning brush.
 - ▶ Mobilize non-rigid components, such as set screws and hinges, during cleaning.
 - ▶ Then flush these areas thoroughly at least five times with the cleaning disinfectant solution using a disposable syringe (20 ml).
- Phase II
- ▶ Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
 - ▶ Mobilize non-rigid components, such as set screws and hinges, during rinsing.

Mechanical alkaline cleaning and thermal disinfection

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound ¹⁾

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D–W	-
II	Cleaning	55/131	10	FD–W	Dr. Weigert neodisher® SeptoClean 1 % ²⁾ working solution
III	Intermediate rinse	>10/50	1	FD–W	-
IV	Thermal disinfection	90/194	5	FD–W	-
V	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine

D–W: Drinking water
FD–W: Fully desalinated water
¹⁾ To demonstrate the cleanability, the following cleaning and disinfecting machine was used: Miele 7836 CD
²⁾ Prion-inactivating detergent (see Technical information Dr. Weigert neodisher® SeptoClean)

- ▶ Place the product in a tray that is suitable for cleaning (avoiding rinsing blind spots)
- ▶ Place instruments in the tray with their hinges open.
- ▶ Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- ▶ Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- ▶ Allow the product to cool down to room temperature.
- ▶ After each complete cleaning, disinfecting and drying cycle, check that the product is dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- ▶ If necessary, dry the product with a lint-free cloth.
- ▶ Check that the product functions correctly.
- ▶ Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical Service.
- ▶ Check for compatibility with associated products.

Packaging

- ▶ Appropriately protect products with fine working tips.
- ▶ Place the product in its holder or on a suitable tray. Ensure that all cutting edges are protected.
- ▶ Pack the trays suitable for the sterilization process in sterile packaging that complies with EN ISO 11607-1 (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging provides sufficient protection against recontamination of the product during storage.

Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to EN 285 and validated according to EN ISO 17665
 - Sterilization using fractional vacuum process at 134 °C, holding time 18 minutes for prion inactivation
- ▶ When sterilizing several products at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Technical Service

WARNING

Risk of injury and/or malfunction!

- ▶ Do not modify the product.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency. Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 (7461) 95-1602
Fax: +49 (7461) 16-5621
E-Mail: ats@aesculap.de
Or in the US:
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420
Other service addresses can be obtained from the address indicated above.

Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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