Aesculap Sterile Technology

Sterile container system
Edition 01/2008

-
- (GB)** **Instructions for use/Technical description**
Sterile container system
 - (D)** **Gebrauchsanweisung/Technische Beschreibung**
Sterilcontainer-System
 - (F)** **Mode d'emploi/Description technique**
Système de containers de stérilisation
 - (E)** **Instrucciones de manejo/Descripción técnica**
Sistema de contenedores estériles
 - (I)** **Istruzioni per l'uso/Descrizione tecnica**
Sterile container system

9

9


10

11

12

13

Sterile container system

Legend

- 1 Outer lid
- 2 Outer lid lock
- 3 Inner lid
- 4 Inner lid lock
- 5 Bottom
- 6 Indicator seal retainer
- 7 Indicator seal
- 8 Handle
- 9 Slot (for the indicator seal on the inner lid lock)
- 10 Universal filter retainer
- 11 Cap
- 12 Push buttons
- 13 Lid latch (for optional outer lid of BASIS container)
- 14 Plastic cover
- 15 Perforation field cover
- 16 Germ retention system
- 17 Adapter frame
- 18 Finger tray
- 19 Handling pin
- 20 Ribbed cover grid
- 21 Mounting handle
- 22 Plastic cover lock

Symbols on the product

	Attention, see instructions for use
---	-------------------------------------

Contents

1. Safe handling	3
2. Product description	3
2.1 Intended use	3
2.2 Operating principle	3
3. Preparation and setup	3
3.1 First use	3
4. Working with the sterile container	4
4.1 System set-up	4
4.2 Function checks	5
4.3 Safe operation	5
5. Maintenance	7
6. Care	7
7. Troubleshooting list	8
8. Technical Service	9
9. Accessories/Replacement parts	9
10. Technical specifications	9
11. Extracts from relevant standards	9
11.1 Quoted standards	9
11.2 Storage period for sterile materials — extract from DIN 58953-9, sterile supplies	9

1. Safe handling

WARNING

Risk of contamination of sterile supplies due to inadequately sealed sterile container!

The sealing of the sterile container and its functioning as a germ barrier will be compromised if the sterile container is combined with components from other manufacturers.

➤ **Only combine Aesculap sterile container products with each other.**

- Prior to use, check for proper condition and functioning of the product.
- To avoid damaging the product through improper setup or operation and mitigating the warranty and liability on the part of the product's manufacturer:
 - Use the product only according to these instructions for use.
 - Follow the safety instructions and maintenance advisories.
 - Never use damaged or faulty sterile containers. Replace any damaged components immediately with original spare parts.
 - If sterile container components are repaired in a way that could affect the container's germproof qualities: Inspect the sterile container thoroughly before use.
- Ensure that only persons with the requisite training, expertise or experience will handle the product and its accessories.
- Keep the instructions accessible for personnel.
- Adhere to applicable standards.
- Follow general guidelines and principles of hygiene when handling contaminated items that have undergone or are to undergo sterilization.

2. Product description

2.1 Intended use

The Aesculap sterile container system is intended for use as sterile packaging for instruments and textiles due for steam sterilization through a validated steam sterilization process (e.g. in a sterilizer complying with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated according to EN ISO 17665 or EN 554/ISO 13683).

After sterilization, the sterile materials are stored in the sterile container until they are used.

Note

Please contact your Aesculap representative if your Aesculap sterile containers are to be used in any other steam sterilization process.

2.2 Operating principle

The Aesculap sterile container system conforms to DIN 58953-9 and EN ISO 11607.

Sterile containers with a perforated lid and a closed bottom have been validated for steam sterilization in a sterilizer acc. to EN 285 through:

- fractionated vacuum process

Sterile containers with a perforated lid and a perforated bottom are also suitable for steam sterilization in a sterilizer acc. to EN 285 through:

- gravitational process

Note

The suitability of any specific process must be evidenced by a validation at the site of application.

3. Preparation and setup

3.1 First use

- Thoroughly clean the new container prior to first use.
- After cleaning, insert a suitable filter, see Filter change.

Sterile container system PrimeLine:

The germ retention system **16** is integrated with the system.

Sterile container system

4. Working with the sterile container

4.1 System set-up

Removing the outer lid/inner lid/plastic cover

If an outer lid is used, it can be removed to clean the sterile container and, if it is soiled, it can be separated from the inner lid.

VARIO container (by default with outer and inner lid):

- Push down the outer lid lock **2** and remove the outer lid **1** and the inner lid **3**.

BASIS container (with retrofitted outer lid):

- Dismount the combined outer lid **1** and inner lid **3** from the bottom **5**.
- Loosen lid latch **13** and remove the outer lid **1**.

Sterile container system PrimeLine:

- Remove the plastic cover **14** from the bottom **5**.

Cleaning and disinfecting the sterile container

Risk of corrosion and stress cracks due to inappropriate cleaning and disinfection!

- Follow the instructions of the cleaning and disinfecting agents' manufacturers when cleaning and disinfecting the product.

Additional instruction for PrimeLine:

- Do not use solvents such as acetone or benzene.
- Carry out final rinse with fully desalinated water only.
- Verify that there are no residues from the cleaning process left on the lid.

- Use fully desalinated water of a quality according to EN 285 Annex B for reprocessing the product.

Note

Literature on cleaning and disinfecting:

"Test Series and Statements; Working Group Instrument Preparation; C 75911", available through Aesculap.

Further advice for the careful reprocessing of reusable medical products can be found in the recommendations of the Working Group Instrument Preparation:

"Proper Maintenance of Instruments"; see "www.a-k-i.org".

Material-specific requirements for the cleaning and disinfection of sterile containers

	Disposal container, stainless steel	Standard container, anodized aluminum	PrimeLine container, plastic
Ultrasound treatment in manual and mechanical processing	√	√	√
Cleaning/disinfecting agent*	acid, pH-neutral, alkaline	pH-neutral, mildly alkaline	pH-neutral, mildly alkaline
Thermal disinfection with fully desalinated water, 5 min, at 93 °C	√	√	√
Mechanical drying with hot air at up to 120 °C	√	√	√
* Mildly alkaline cleaning/disinfecting may cause surface changes discoloration of the aluminum lids (fading, staining). Such changes do not affect the functionality of the product. For colored aluminum lids use, if possible, a process involving a neutral cleaning and disinfecting agent and fully desalinated water.			

Manual processing

- Use appropriate cleaning and disinfecting agents. Always adhere to the manufacturers' instructions.
- Wipe surfaces with a clean, lint-free cloth and a pH-neutral cleaning and disinfecting agent.
- Remove encrusted debris with a soft nylon brush. Do not use harsh cleaning agents or metal brushes.
- Carry out the final rinse with fully desalinated water.

Mechanical processing

- Use appropriate cleaning and disinfecting agents. Always adhere to the manufacturers' instructions.
- Use suitable storage racks for mechanical cleaning according the manufacturers' instructions for the respective devices.
- Carry out the final rinse with fully desalinated water.

Filter change

Change the filter at the following intervals, depending upon filter type:

- Disposable filter: prior to each sterilization cycle
- PTFE permanent filter: after a maximum of 1 000 sterilizations
- PrimeLine germ retention system: after a maximum of 5 000 sterilizations

VARIO container and BASIS container:

- Simultaneously press both push buttons **12** on universal filter retainer **10**.
- Remove universal filter retainer **10**.
- Insert a new filter and remount universal filter retainer **10**.
- Push down cap **11** on universal filter retainer **10** until you hear it click into place.

Sterile container system PrimeLine:

- Turn ribbed cover grid **20** anticlockwise until it is unlocked.
- Remove ribbed cover grid **20** of germ retention system **16**.
- Turn germ retention system **16** anticlockwise at finger tray **18** or mounting handle **21** until it disengages from adapter frame **17**.
- Lift germ retention system **16**, gripping it at handling pin **19**, and remove it.
- For installing the germ retention system **16**, retrace the above steps.

4.2 Function checks

- Prior to each use, carry out a visual inspection for damage and correct functioning of all components of the sterile container:
 - Metal parts not deformed
 - Aluminum lid not warped
 - Edges of universal filter retainer **10** sitting in full-surface contact
 - Seals intact (no cracks, ...)
 - Plastic parts not cracked
 - Permanent filter intact
 - Lock functions properly (engages)
- Use sterile containers only if they are in mint condition. Replace any damaged components immediately with original spare parts or have the affected components repaired.

4.3 Safe operation

Risk of contamination of sterile supplies due to inadequately sealed sterile container!

The sealing of the sterile container, and its effectiveness as a germ barrier will be compromised if the sterile container is combined with components from other manufacturers.

- Only combine Aesculap sterile container products with each other.

Loading a sterile container

Instruments

In accordance with DIN EN 868-8 and DIN 58953-9, we recommend the following maximum load of the container:

- Standard container: 10 kg
 - Half container: 5 kg
 - 3/4-container: 7 kg
- Put instruments into appropriate storage units in a tray, with hollow components, bowls, dishes, etc. with their openings at a downward angle.

VARIO container and BASIS container:

- Load sterile containers in such a way that their universal filter retainers **10** remain unobstructed. Max. loading height: to approx. 2 cm below the rim of the container base.
- Use inner lid lock **4** to lock inner lid **3** to container bottom **5**.
Inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile container system PrimeLine:

- Load the sterile container in such a way that ribbed cover grid **20** in plastic cover **14** remains unobstructed. Max. loading height: to approx. 2 cm below the rim of the container base.
- Lock plastic cover **14** with plastic cover lock **22** on bottom **5**.
Plastic cover lock **22** must clearly engage. If this is not the case: Have the sterile container repaired, see Technical Service.

Textiles

- Pack the textiles together in such a way that they fit upright in the sterile container.
- Make certain that, when the sterile container is fully loaded, it is still possible to insert a flat hand between the individual items without difficulty.
- Load sterile containers in such a way that their universal filter retainers **10** remain unobstructed.
- Use inner lid lock **4** to lock inner lid **3** to container bottom **5**.
Inner lid lock **4** has to engage perceptibly. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile container system PrimeLine:

- Load the sterile container in such a way that ribbed cover grid **20** in plastic cover **14** remains unobstructed.
- Lock plastic cover **14** with plastic cover lock **22** on bottom **5**.
plastic cover lock **22** must clearly engage. If this is not the case: Have the sterile container repaired, see Technical Service.

Sterile container system

Inserting the indicator seal

- Once the sterile container is loaded, enter the following data on indicator seal **7**: Date of sterilization, sterilization number, expiry date, name and signature.
- Slide the indicator seal **7** from outside into indicator seal retainer **6** in such a way that the indicator segment engages in the slot on the lid lock and seals the lock.

Loading the sterilizer

Observe the following rules when prepare the sterile container and the sterilizer for loading the sterilizer:

Risk of vacuum damage to the sterile container due to unbalanced pressurization!

- Do not use outer packaging for the sterile containers.
- Do not block the air flow through the perforation fields in the bottom and the inner lid under any circumstances (applicable to all container types).
- Do not place foil packaging directly on the sterile container.

- Follow the sterilizer manufacturer's instructions.
- Whenever possible, load the sterilizer with a mixture of textiles and instrument containers.
- Always place heavy sterile containers at the bottom of the sterilizer.
- VARIO containers must be sterilized with their outer lids attached.
- Always carry the sterile containers by their handles.
- Stack sterile containers in the sterilizer if necessary.
- Transport stacks of sterile containers in such a way that the stacks will not topple.

Sterilization method and parameters

Risk of sterilization failure!

- Sterilize the sterile containers though approved and validated sterilization processes only.
- Sterilize with steam and comply with the following rule:
Sterilization must be carried out through a validated steam sterilization process (e.g. in a sterilizer according to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, and validated according to EN ISO 17665 or EN 554/ISO 13683).

Unloading the sterilizer

Risk of burns due to hot sterile containers after sterilization!

- Always work with protective gloves.

Transporting the sterile container

Risk of unsterility of sterile supplies!

- Always carry the sterile containers by their handles.
- Never lift or carry a sterile container at the lid.
- Transport sterile containers in such a way that mechanical damage is prevented.

Storing sterile containers

Note

Sterile containers may be stored in stacks.

- Store sterile containers in a dry, clean and protected place.
- Observe storage period and storage conditions acc. to DIN 58953-9, see Storage period for sterile materials – extract from DIN 58953-9, sterile supplies.
Recommended guide figure for the storage period acc. to DIN 58953-9: 6 months.

Checking and setting up of sterile supplies

The contents of a sterile container are only deemed sterile if the container has been properly sterilized, stored and transported.

If this is not the case, the sterile supplies must be reprocessed.

Risk of contamination due to improperly sterilized sterile supplies!

- **Before setting up the sterile supplies, check to ensure that the sterilization was successful.**
- Before opening the sterile container, make certain that its indicator seal is not damaged.
- Make certain that the color of the indicator point has changed.
- Make certain that all container components, especially the lid locks, are intact.
- Make certain that there is no condensate in the container.
If any condensate is found, eliminate the cause, see Troubleshooting list.

5. Maintenance

Note

The sterile containers must only be tested and repaired by persons with the appropriate training, expertise or experience.

The seals in the container lids have been tested in conformance with DIN EN 868-8, Appendix G, for a minimum service life of 5,000 cycles.

- In case of visible damage, change the seal immediately.
- Have the lid sent to Aesculap, see Technical Service.

6. Care

- Lightly lubricate moving metal parts, e.g. latch hinges, with sterilizable and steam-permeable maintenance oil (e.g. Aesculap STERILIT® spray JG600 or maintenance oil JG598).

Sterile container system

7. Troubleshooting list

Malfunction	Cause	Remedy
Excessive condensate formation inside the sterile container	Temperature of sterile materials too low prior to sterilization	Pre-heat sterile supplies to room temperature (approx. 20 °C) before sterilization
	Textiles too damp	Sterilize dry textiles only
	Sterile container too heavy	Standard container with instruments: max. load 10 kg with textiles: max. load 8 kg
		Half container with instruments: max. load 5 kg
		3/4-container with instruments: max. load 7 kg
	Sterile supplies improperly positioned	Position hollow components, bowls, plates etc. face down at an angle.
		Arrange textiles in loose vertical piles; do not press them in
	Sterile container incorrectly positioned in sterilizer	Always place heavy sterile containers at the bottom
	Sterile containers set up for use immediately after sterilization	Allow sterile containers to cool down to room temperature before setting them up for use
	Sterile containers improperly positioned during cooling phase	Do not store sterile containers on a floor or in a drafty place. Store sterile containers in air-conditioned areas according to DIN 58953-9.
	Sterilizer properties do not meet DIN EN 285 standards	Have sterilizer serviced regularly Check the drying vacuum Check the drying time
		Check the steam quality and upgrade if necessary
	Empty-cycle and vacuum test not carried out daily prior to sterilization	Run empty-cycle and vacuum test daily prior to sterilization
Condensate in groove of (inner) lid	Unsuitable sterilizer program selected	Select a program according to the load
	Sterilizer door left open for too long; sterilizer cooled down	Load and unload sterilizer quickly
Stacking of several sterile containers or heavy sterile supplies		Subdivide sterilizer with shelving Use outer lid
Indicator seal color unchanged	Sterilization not completed correctly; sterilizer defective	Have sterilizer repaired by manufacturer
BASIS/VARIO sterile container or PrimeLine sterile container system deformed	Perforation field covered during sterilization	Never cover the perforation field
Lid lock defective	Outer cover used to carry sterile container	Always carry sterile containers by their handles

8. Technical Service

For service, maintenance or repairs, please contact your national B. Braun/Aesculap agency.

All guarantees and warranty rights as well as licenses will be void if modifications on the device are performed.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

9. Accessories/Replacement parts

Accessories and expendables are listed in brochure no. C40401.

10. Technical specifications

The range of models and sizes of sterile containers is listed in brochure no. C40401.

11. Extracts from relevant standards

11.1 Quoted standards

The following standards are quoted in connection with the sterile containers:

- DIN 58953-9
- DIN 58946-1, EN 285, EN 866-3
- EN 868-1, EN 868-8
- EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993, further applicable EN ISO 17665 or EN 554/ISO 13683
- EN ISO 11607

11.2 Storage period for sterile materials – extract from DIN 58953-9, sterile supplies

Loss of sterility is less likely caused by excessive storage periods as by external influences and impacts during storage, transport and handling. Therefore, there can be no general statement concerning appropriate storage periods. The appropriate length of storage periods is determined by the hygiene control committee. Responsibility for storage conditions and storage periods lies with the hospital's operator or its medical director.

Where appropriate, different storage periods for individual department should be laid down in writing, also taking account of the type of packaging and the storage conditions.

Recommended storage period for sterile materials in sterile containers

according to EN 868-1 and EN 868-8: 6 months

The storage period just quoted represents a recommendation that, if followed, would limit the risk of contamination during transport and opening of the sterile container. Storage for longer periods can lead to contamination of outer surfaces, depending upon storage conditions. Although such longer periods do not in and of themselves lead to recontamination of the contents of the package during storage, they do increase the risk of contamination when sterile containers are transported or opened. The recommended storage period is based on storage under dry and dust-protected conditions. Longer storage periods are possible under appropriate conditions at the storage site.

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Aesculap AG & Co. KG