- GB

Instructions for use
- USA

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Targon® femur and tibia nail system

Intended use

The Targon® femur and tibia nail system is used as a splint for fractures of long tubular bones, the femur and the tibia.

System configurations

- Targon® universal femur nail: various lengths and diameters
- Targon® universal tibia nail: various lengths and diameters
- Targon® solid titanium tibia nail: various lengths and diameters
- Targon® solid titanium femur nail: various lengths and diameters
- For further information regarding diameters and lengths, see operating manual O113 and brochure O111.
- Targon® retrograde femur nail: various lengths and diameters
- For further information regarding diameters and lengths, see brochure O190.

Associated implants:

- Fixation screws in various diameters and lengths
- Threaded sleeves for fixation screws
- Fixation screws for the end of the nail on the connecting side

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®: Titanium forged alloy Ti6Al4V acc. to ISO5832-3
- ISONIC® FeCrNiMoN forged alloy acc. to ISO 5832-9

The titanium implants are coated with a colored oxide. Slight changes in coloration may occur, but do not affect the implant quality.

ISOTAN® and ISONIC® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use Targon® universal femur nails, Targon® universal tibia nails, Targon® solid titanium tibia nails and Targon® solid titanium femur nails for the following indications:

- Open or closed fractures of the femur shaft or tibia shaft
- Repositioning osteotomies in the shaft region
- Reconstruction following tumor resection in the shaft region
- Osteotomies in the shaft region to lengthen or shorten the limb
- Pathological fractures
- Pseudarthroses
- Knee joint arthrodesis (Targon® F special length)

Use Targon® retrograde femur nails for the following indications:

- Femur shaft fractures
- Extra- and intra-articular fractures of the distal femur (these fractures also with ipsilateral femur shaft fractures)
- Patients with ipsilateral knee prosthesis, hip prosthesis or other proximal implants, e.g. dynamic hip screw
- Shaft fractures necessitating additional proximal implants, e.g. for ipsilateral femoral neck fractures
- Shaft fracture and additional nailing of the ipsilateral tibia (identical, extended approach incision)
- Pathological fractures
- Pseudarthroses
- Ankle arthrodesis

The operating surgeon is responsible for establishing the above indications and indications beyond those, taking into account the specific clinical, biological, and biomechanical situation.

Contraindications

Do not use Targon® universal femur nails, Targon® universal tibia nails, Targon® solid titanium tibia nails and Targon® solid titanium femur nails for the following:

- Open or closed fractures in joint regions
- Trochanteric fractures of the femur

Do not use Targon® retrograde femur nails for the following:

- Septic knee arthritis

Generally, do not use in cases of:

- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Anticipated overloading of the implant
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Side effects and interactions

- Dislocation, loosening and fracture of implant components
- Delayed healing, or non-union and development of pseudarthrosis
- Primary and secondary infections
- Venous thrombosis, lung embolism, cardiac arrest
- Tissue reaction to implant materials
- Trauma to nerves, tendons, or vessels
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Compartment syndrome
- Pain in the region of the nail insertion point and in the region of locking components

Safety notes

CAUTION

Federal law restricts this device to purchase by, or on instruction by a physician!

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.

- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Mobilization, dynamic locking or use of compression involve the risk of migration of the nail and/or the locking component into the joint region or out of the bone. Therefore the migration paths of the implants and those of the fragments must be considered in the implantation and taken into account when selecting the implant length.
- Considerably higher strain due to increased leverage is to be expected if there are locking holes at the height or in the region of the fracture line. This can lead to implant failure. Accordingly, postoperative strain has to be reduced, or full load may be envisaged only after complete callus formation.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudarthrosis develops, or if the implant is overloaded for an excessively long period. One way to achieve this is through mobilization.
- Always insert the Targon® implants in such a way that forces to be transmitted are slight and can be transmitted by the bone at an early stage.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Do not apply strong force when implanting or explanting the implants. Should a problem arise, check the position and orientation of the implants, instruments, and fragments and analyze the origin of the problem. If necessary, repeat previous operating steps and check the instruments (e.g. for blocked drill helix).

Interactions between MRI and implant components!

- MRI examinations using magnetic fields of 1.5 and 3.0 tesla do not present an additional risk to implant wearers.
- MRI causes non-critical, localized heat development.
- Implants produce moderate MRI artifacts.

- Targon® implants may only be implanted using the Aesculap implantation instruments specially designed for this purpose.
- Use implant components only in the nail holes intended for this purpose.
- Only use cutting instruments, e.g. drill bits, guide pins, etc., that are sharp.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.
- The attending physician shall make any decision with regard to the removal of implant components that have been used,
- Increase the load on the fracture and implant in accordance with the progress of healing.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- At the time of follow-up exams, the scope of postoperative load and aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. Also, the dimensions of the implant should be taken into account.
- Recommended follow-up intervals
 - prior to the patient's discharge from hospital
 - 10-12 weeks post operationem
 - 6 months post operationem
 - 12 months post operationem

For implanting the nail via a guide pin:

- Use the appropriate guide pin
- Avoid jamming
- Monitor the tip of the guide pin under X-ray control
- Remove the guide pin prior to locking

Aesculap and the Aesculap Academy regularly hold general and product-specific training events, offering users helpful application advice, especially for complicated or critical situations.

Sterility

- The implant components are supplied in an unsterile condition.
- The implant components are packaged individually.
- Store the implant components in their original packaging and only remove them from their original and protective packaging immediately prior to application.
- Use the implant system storage devices for processing, sterilization and sterile setup.
- Use a suitable tray for cleaning/disinfection.
- Use the system storage device only for sterilization and sterile provision.
- Ensure that the implant components in their implant system storage devices do not come into contact with each other or with instruments.
- Ensure that the implant components are not damaged in any way.

Prior to initial sterilization and subsequent resterilization, the implant components must be cleaned using the following validated reprocessing procedure:

Validated reprocessing procedure

General safety instructions

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For the latest information on reprocessing and material compatibility see also the Aesculap extranet at www.aesculap-extra.net

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying. Additional drying, if necessary. Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.ak-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

Validated cleaning and disinfecting procedure

For mechanical cleaning:

Note
Use appropriate (neutral, enzymatic and mild-alkaline) detergents on this product with its aluminum components.

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection	■ 20 ml disposable syringe ■ Drying phase: Use a lint-free cloth or medical compressed air.	Chapter Manual cleaning/disinfection and subsection: ■ Chapter Manual cleaning with immersion disinfection
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfecting	■ 20 ml disposable syringe ■ Place the instrument in a tray that is suitable for cleaning (avoiding rinsing blind spots).	Chapter Mechanical cleaning/disinfection with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning with a brush ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

Manual cleaning/disinfection

Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	D–W	-
III	Disinfection	RT (cold)	15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	RT	-	-	-	-
D–W:	Drinking water					
FD–W:	Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)					
RT:	Room temperature					
*Recommended:BBraun Stabimed						

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfecting procedure.

Phase I
► Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
► Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
► If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
► Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
► Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II
► Rinse/flush the product thoroughly (all accessible surfaces) under running water.
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
► Drain any remaining water fully.

Phase III
► Fully immerse the product in the disinfectant solution.
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
► Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV
► Rinse/flush the product thoroughly (all accessible surfaces).
► Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
► Rinse lumens with an appropriate disposable syringe at least five times.
► Drain any remaining water fully.

Phase V
► Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfecting procedure.

Mechanical cleaning/disinfection with manual pre-cleaning

Targon® implants not falling under the above restrictions can be resterilized after proper processing. For implant components that are to be resterilized:

WARNING

Intraoperative contamination with blood, secretions, and other fluids may render the affected component unsuitable for resterilization!

- Handle the implants with new gloves only.
- Keep the implant system storage devices covered or closed.
- Process implant system storage devices separately from instrument trays.
- Clean implants must not be processed together with contaminated implants.
- Process the implant components individually and separately if no implant system storage devices are available, ensuring that the implant components are not damaged in the process.
- Mechanically clean and disinfect the implant components.
- Do not reuse surgically contaminated implants!

WARNING

Direct or indirect contamination may render implants unsuitable for resterilization!

- Do not reprocess implants that have been directly or indirectly contaminated with blood.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Disinfectant cleaning	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Rinsing	RT (cold)	1	-	D–W	-
D–W: Drinking water RT: Room temperature *Recommended:BBraun Stabimed						

- Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfecting procedure.

Phase I
► Fully immerse the product in the cleaning/disinfecting solution for at least 15 minutes. Ensure that all accessible surfaces are moistened.
► Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
► If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 minute.
► Mobilize non-rigid components, such as set screws and hinges, during cleaning.
► Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a suitable disposable syringe.

Phase II
► Rinse/flush the instrument thoroughly (all accessible surfaces) under running water.
► Mobilize non-rigid components, such as set screws and hinges, during rinsing.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	<25/77	3	D–W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH = 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH = 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to the program for cleaning and disinfection device
D–W:	Drinking water				
FD–W:	Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)				
*Recommended:BBraun Helimatic Cleaner alkaline					

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

Inspection, maintenance and checks

- Allow the product to cool down to room temperature.
- Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- Immediately sort out damaged or inoperative products.

Packaging

- Place the product in its holder or on a suitable tray.
- Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Make sure that the packaging will prevent a recontamination of the product whilst in storage.

Sterilization

- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several instruments at the same time in a steam sterilizer, ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

Application

WARNING

Risk of injury caused by incorrect operation of the product!

- Attend appropriate product training before using the product.
- For information about product training, please contact your national B. Braun/ Aesculap agency.

A detailed operation manual for the Targon® nail system is available upon request from Aesculap or from your local Aesculap representative.

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- Information has been obtained from the manufacturer in the event of an ambiguous preoperative situation and if implants are present in the area to be treated

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- Under certain circumstances it can happen that, when a fracture is treated with an implant, humerus anatomy will not be fully restored.
- Following treatment of a fracture, shoulder and elbow joint function may be impaired.
- Pains may occur following the fracture treatment.
- The implants must not be exposed to overload due to excessive strains, hard physical work or sports activities. Overstraining entails the risk of implant loosening or breakage.
- The patient must be informed about the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The risks of transgressing these rules must be explained to the patient.
- Revision surgery may become necessary in cases of implant loosening or pseudarthrosis.
- The patient must undergo regular medical follow-up examinations.
- The attending physician decides when to remove, in part or completely, the osteosynthesis implants through another surgical procedure.

WARNING

Complications due to incorrect use and/or incorrect assessment/treatment of the fracture!

- Only apply nails of appropriate length and diameter.
- Make certain that implants are aligned correctly.

The implantation and implantation site are prepared in the following way:

The application steps for Targon® implant components are described in detail in the operating manuals. The following notes must be observed for all Targon® implantations:

Check for correct mounting of the nail on targeting device: Insert the drill through the drill sleeve and the targeting hole and guide it through the nail hole.

WARNING

Improper use can cause bones to shatter and nails to jam or bend!

- Do not use excessive force, under any circumstances, when handling/applying the instruments and implants.
- Only apply mild hammer taps to Targon® F/T nails and, if necessary, RF nails to insert them into the medullary cavity.
- Check the nail insertion site, and correct if necessary.
- Should a nail be difficult to insert, drill open the medullary cavity to a smaller (for solid nails) or larger (for drilled nailing) extent.

The operating surgeon determines the following, depending upon the indications:

- if implant locking is necessary
- if the implant is to be locked, whether it will be locked statically or dynamically.

Note

For locking an implant without a targeting attachment, Aesculap recommends free-hand locking supported by a C-bow and laser pilot beam, or with a radiolucent angle transmission

WARNING

Avoid drill mistargeting, drill breakage, guide pin breakage and premature implant failure due to crack formation as a result of damaged surfaces in the drill hole regions!

- Use a drill bit of the correct diameter.
- Level the drill insertion site with a flat counterbore or use a center punch.
- Only use sharp drill bits and guide pins and drill carefully and without applying excessive force. Avoid changes of the drilling direction.
- During center punching, flat countersinking and drilling, do not apply force to the targeting device; do not bend or twist the targeting device.

Implant removal requires the following steps:

- Prior to removing the nail, remove all screws.
- Remove the nail with quick, light taps, and with the aid of the impactor, extraction pin and slotted hammer.

Note

Implant removal can give rise to complications due to implants that are stuck, ingrown bone tissue, etc. Such complications can lead to damage to implants and/or instruments. For such cases we recommend using a set of special instruments, which can be ordered from Aesculap. This instrument set should be at hand for every explantation. Other special instruments are available for dealing with broken implants. Note the instructions for use for such case!

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC

3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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