Targon® FN

- GB

Instructions for use
Targon® FN
- D

Gebrauchsanweisung
Targon® FN
- F

Mode d'emploi
Targon® FN
- E

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Intended use

The Targon® FN system is used for splinting, stabilization and fixation of the proximal femur.

System variants:

- Targon® FN femoral plate
- Targon® FN Tiescrew support screw system in various lengths

Materials

The materials used in the implant are listed on the package.

- ISOTAN® titanium forged alloy Ti6Al4V to ISO 5832-3

The titanium implants are coated with a colored oxide. Their color may change slightly, but this has no bearing on the quality of the implant.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for:

- Medial femoral neck fractures

The operating surgeon is responsible for establishing the above indications and indications beyond those, taking into account the specific clinical, biological and biomechanical situation.

Contraindications

Do not use in the presence of:

- Trochanteric femur fractures
- Acute or chronic infections
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Anticipated overloading of the implant
- Inadequate patient compliance
- Foreign body sensitivity to the implant materials

Side-effects and adverse interactions

- Changes in position, loosening, wear and tear or fracture of implant components
- Position change and loosening of fragments
- Delayed healing, or non-union and development of pseudoarthrosis
- Primary and secondary infections
- Venous thrombosis, lung embolism, and cardiac arrest
- Tissue reaction to implant materials
- Trauma to nerves, tendons or vessels
- Hematomas and wound healing disorders
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise
- Pain in the implant entry region and in the region of locking components

Safety notes

- The implant components were tested and approved for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- Always observe the instructions for use of the individual Aesculap implant components.
- Follow the OP manual.
- Use implant components only in the appropriate plate drill holes.
- The operating surgeon must be thoroughly familiar with the bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- Prior to the operation, the operating surgeon must familiarize him/herself with the Targon® implants and associated operating technique.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- The Targon® implants may only be implanted with the Aesculap implantation instruments designed for this purpose.
- Aesculap is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- Under no circumstances may implant components from different suppliers be combined.
- Implants that have already been used must not be reused.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Considerably higher strain due to increased leverage is to be expected if there are locking holes at the height or in the region of the fracture line. This can lead to implant failure. Accordingly, postoperative strain has to be reduced, or full load may be envisaged only after complete callus formation.
- Overloading the implant can cause the implant material to fracture. The dynamic load on the implant must be reduced if healing fails to occur or is slow, if pseudoarthrosis develops, or if the implant is overloaded for an excessively long period. This can be achieved e.g. through dynamizing.
- Always insert the Targon® implants in such a way that forces to be transmitted are slight and can be transmitted by the bone at an early stage.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- In order to promote the earliest possible detection of any problems or complications, the operation results must be followed up at regular intervals with the aid of appropriate examination procedures. A precise diagnosis requires x-rays taken in the directions anterior-posterior and medial-lateral.
- Never use damaged implant components or implant components that have been removed surgically.

- Do not apply strong force when implanting or explanting the implants. Should a problem arise, check the position and orientation of the implants, instruments and fragments and analyze the origin of the problem. If necessary, repeat previous operating steps and check the instruments (e.g. for blocked drill helix).
- Increase the load on the fracture and implant in accordance with the progress of healing.
- The timing of follow-up examinations, the type of postoperative load and the aftercare should be determined for each patient on the basis of their weight, their activities, the type and severity of the fracture, and any coexisting injuries. The dimensions of the implant should also be taken into account.
- Recommended follow-up intervals:
 - prior to the patient's discharge from hospital
 - 10-12 weeks post operationem
 - 6 months post operationem
 - 12 months post operationem

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation (min. dose 25 kGy).
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiry date or whose packaging is damaged.

Resterilization

Risk of contamination of the patient!

- Never resterilize or reuse implants contaminated e.g. with blood or secretions.

- Targon® FN implants to which the above restrictions do not apply can be resterilized through steam sterilization preceded by appropriate pre-cleaning, ultrasound treatment and disinfection.
- Make sure the implant mechanism is not damaged during cleaning and resterilization.

The responsibility for resterilization of implant components lies with the end user. In certain countries (e.g. France) such resterilization is specifically prohibited.

Sterilization method and parameters

- Observe all relevant national regulations and standards with regard to processing.
- For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.
- Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net
- Sterilize with steam, taking note of the following: Carry out steam sterilization through an approved steam sterilization procedure (e.g. sterilizer in conformance with EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated in conformance with EN 554/ISO 13683). For the fractionated vacuum process, sterilization has to be carried out running the 134 °C/2 bar program for a minimum holding time of 5 minutes.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- Operating conditions are highly aseptic.
- Implantation instruments, including the special Aesculap implant system instruments, must be complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant articles by medical specialists from professional literature.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical intervention has been explained to the patient and their agreement concerning the following information has been recorded:

- Treatment of a fracture with implants may not always be successful in completely restoring the original bone anatomy.
- Following the fracture treatment, the functionality of adjacent joints may be restricted.
- Following treatment of a fracture, pain may develop.
- The implants must not be subjected to overload through excessive strains, hard physical work or sports activities. The implant may loosen or suffer material fracture due to excessive strain e.g. suffered in a fall.
- The patient must be informed of the limits to the allowable strain on the implant, and be given appropriate behavior guidelines. The implants will not sustain the load and functionality of a healthy bone. The risks inherent in transgressing the behavior guidelines must be made clear to the patient.
- Corrective surgery may become necessary in cases of implant loosening, fragment migration, pseudarthroses or other complications.
- The patient must undergo a regular medical follow-up examination.
- The attending physician decides the appropriate timing of the surgical removal, in part or completely, of the osteosynthesis implants.
- Check for correct mounting of the plate at the aiming device: Insert the drill through the drill sleeve and the aiming hole and guide it through the plate hole.

Explanation

Note

Implant removal can give rise to complications due to implants that are stuck, ingrown bone tissue, etc. Such complications can lead to damage to implants and/or instruments. For such cases a set of special instruments is recommended, which can be ordered from Aesculap. This instrument set should be at hand for every explantation. Other special instruments are available for dealing with broken implants. Note the instructions for use for such case!

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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