Aesculap Neurosurgery

GB

Instructions for use

USA

YASARGIL aneurysm clip applicator forceps

D

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YASARGIL Aneurysmen-Clip-Anlegezange

F

Mode d'emploi

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
www.aesculap.de

TA-Nr. 013089 07/09 V3

CE marking according to directive 93/42/EEC

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- 1 Jaw piece
- 2 Branches left/right
- 3 Clip applier forceps

Symbols on product

Symbol	Explanation
	Caution: See documentation supplied with the product

Intended use

Aesculap YASARGIL aneurysm clip applier forceps are used for holding and applying Aesculap YASARGIL Phynox and titanium aneurysm clips and Aesculap KOPITNIK AVM microclips in neurosurgical procedures.

Note

In the following, Aesculap YASARGIL aneurysm clips and Aesculap-KOPITNIK AVM microclips are referred to as "clips".

Available sizes

Detailed information on available sizes can be found in the main Neurosurgery Catalog and in the Yasargil Aneurysm Clips brochures.

All clip applier forceps are marked according to size (Mini or Standard, Long or AVM micro) and clip material (Phynox or titanium) so that their correct application with clips of the appropriate size and material is ensured.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

CAUTION

Damage to, imprecise function and incorrect closing force of the clips due to use of incorrect applier forceps!

- Use Aesculap clip applier forceps only.
- Use clips only with the appropriate Aesculap clip applier forceps (see labels).

- Use clips and clip applier forceps only in the following combinations:
 - Titanium Standard clips with the clip applier forceps for Titanium Standard clips
 - Phynox Standard clips with the clip applier forceps for Phynox Standard clips
 - Titanium Mini clips with the clip applier forceps for Titanium Mini clips
 - Phynox Mini clips with the clip applier forceps for Phynox Mini clips
 - Phynox Long clips with the clip applier forceps for Phynox Long clips
 - AVM clips with the applier forceps for AVM clips

Locks

Note

Some clip applier forceps variants are fitted with a lock to engage during clip transfer.

The clip applier forceps may be fitted with one of the following locks:

- Permanent lock
- Button-controlled lock

Operating the permanent lock

Note

With the clip applier forceps with permanent lock, the lock is always active and cannot be switched off.

- Insert the clip into jaw piece 1.
- Compress branches left/right 2 until the lock engages.
The clip is in transfer position.
- Move clip applier forceps 3 with the clip into the operating field.
- To release the lock, slightly compress branches left/right 2.
- Apply the clip.

Operating the button-controlled lock

Note

The clip applier forceps with button-controlled lock can be used with or without the locking function. The locking function of the button-controlled lock needs to be activated deliberately by pressing the lock button.

Transferring the clip and engaging the lock

- Insert the clip into jaw piece 1.
- Push the lock button and simultaneously compress branches left/right 2.
The lock is engaged and the clip is slightly clamped in jaw piece 1.

Releasing the lock

- Compress branches left/right 2.
Jaw piece 1 of clip applier forceps 3 closes and the clip opens.
The lock automatically returns to the unlocked position.

Detaching the clip applier forceps from the clip

- Release the pressure on branches left/right 2.
Jaw piece 1 of clip applier forceps 3 opens and the clip closes.
- Detach clip applier forceps 3 from the clip.

Vario clip applier forceps

The Aesculap Vario clip applier forceps offer an additional rotary function for positioning of the clip in the clip applier forceps.

The maximum range of rotation at the working end of the Aesculap Vario clip applier forceps is 20° in each direction, starting from the middle position, see Fig. 1


Fig. 1

The adjusting plates of the Aesculap Vario clip applier forceps can be used for parallel alignment of the applier jaws, see Fig. 2

Fig. 2

The adjusting plates allow swiveling the clamping jaws by 320°, see Fig. 3


Fig. 3

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of the products.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, the manufacturer's specifications can only serve as an approximate guide for assessing the processing procedures applied by the individual operator/processors.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. To avoid this, the time interval between application and processing should not exceed 6 h, and neither fixating pre-cleaning temperatures >45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process can already cause visible surface changes.
- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use oxidizing process chemicals as these can cause bleaching or layer loss.
- Do not use process chemicals that cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.
- Encrusted debris on HF instruments is dissolved most effectively and conservatively through about 5 minutes of immersion treatment in a 3 % H₂O₂ solution. Subsequently, the debris can be removed by hand, with a medium-hard brush and/or in an ultrasonic cleaning bath. After that treatment, proceed through the normal steps of processing.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, Publications Red Brochure – Proper maintenance of instruments.
- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Open jointed products.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.

Cleaning/Disinfecting

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use a cleaning and disinfecting detergent according to the manufacturer's instructions. The cleaning and disinfecting detergent
 - must be approved for e.g. aluminum, plastics, high-grade steel,
 - must not attack softeners (e.g. silicone).
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55 °C.

➤ Clean and disinfect microsurgical products mechanically, provided they can be securely fixed in machines or storage devices in such a way that they will be thoroughly cleaned.

Manual cleaning/disinfecting

- Clean products with movable links with the links opened or mobilized.
- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Disinfecting	RT (cold)	15	2	D–W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	0.5	–	FD–W	–
V	Drying	RT	–	–	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized)
RT: Room Temperature

Stage I

- Fully immerse the product in the cleaning/disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product under running tap water, using a suitable cleaning brush if necessary, until all visible residues have been removed from the surfaces.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometry, for at least 1 min or until no more residues can be removed. Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
- After cleaning, thoroughly (at least five times) rinse through these components with the cleaning/disinfecting solution, using a disposable syringe (20-ml).
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.

Stage IV

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage V

- Dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting

Note
Categorically, the disinfectant must be of tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note
For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that Ao is >3 000 for the process.

Note
The disinfectant used for processing must be serviced and checked at regular intervals.

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
- Ensure that product links and joints are opened when placing the product on the tray.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfectant without ultrasound

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	BBRAUN HELIMATIC CLEANER alkaline with tensides, application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to disinfectant program

D–W: Drinking water
FD–W: Fully desalinated water (demineralized)

Inspection, maintenance and checks

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® spray JG600 or maintenance oil JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn, or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Check for compatibility with associated products.

Packaging

- Appropriately protect products with fine working tips.
- Store products with ratchet locks fully opened or locked no further than in the first notch.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
- Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.

To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterilcontainer (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F–275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap Sterilcontainer cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Technical Service

Risk of injury and/or malfunction!

- Do not modify the product.

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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