GB

Instructions for use

Metha® expeller instrument for cone 12/14 (ND656R)

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
www.aesculap.de

TA-Nr. 012480 04/09 Änd.-Nr. 32568

CE marking according to directive 93/42/EEC

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Legend

- 1 Tension sleeve
- 2 Sheath
- 3 Threaded hole
- 4 Ledge
- 5 Top of tension sleeve
- 6 Chuck

Symbols on product

	Attention, see instructions for use
--	-------------------------------------

Intended use

The Metha® expeller instrument for cone 12/14 (ND656R) is used for separating Metha® trial collum components from Metha® osteoprofilers and for explanting non-modular Metha® prostheses (in combination with ND655R).

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Clean the new product either manually or mechanically prior to the initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

Separating Metha® trial collum components from Metha® osteoprofilers

- Slide the Metha® expeller instrument for cone 12/14 (ND656R) with chuck 6 over Metha® trial collum component.
- Tighten tension sleeve 1 in the direction of the arrow towards "close".
- Pull off the Metha® trial collum component.

Explanting non-modular Metha® prostheses (in combination with ND655R)

Damage to the Metha® prosthesis/prosthesis cone during explantation!

➤ Do not reuse explanted Metha® prostheses.

WARNING

➤ Slide the combined Metha® expeller instrument for cone 12/14 (ND656R/ND655R) with chuck 6 over the prosthesis cone.

➤ Tighten tension sleeve 1 in the direction of the arrow towards "close".

➤ Explant non-modular Metha® prostheses using a slotted hammer (NF275R).

Disassembling

Separating Metha® trial collum components from Metha® osteoprofilers

- Turn tension sleeve 1 in the direction of the arrow towards "open" until top 5 of tension sleeve 1 sits flush with ledge 4.
- If necessary, remove the Metha® trial collum component.
- Turn tension sleeve 1 further to "open" (as indicated by the arrow) and remove it.

Explanting non-modular Metha® prostheses (in combination with ND655R)

- Unscrew Metha® impactor/expeller instrument ND655R from threaded hole 3 of ND656R.
- Turn tension sleeve 1 in the direction of the arrow towards "open" until top 5 of tension sleeve 1 sits flush with ledge 4.
- If necessary, remove the Metha® prostheses.
- Turn tension sleeve 1 further to "open" (as indicated by the arrow) and remove it.

Assembling

Separating Metha® trial collum components from Metha® osteoprofilers

- Screw tension sleeve 1 onto sheath 2 in the direction of the arrow towards "close" until top 5 of tension sleeve 1 sits flush with ledge 4.

Explanting non-modular Metha® prostheses (in combination with ND655R)

- Screw tension sleeve 1 onto sheath 2 in the direction of the arrow towards "close" until top 5 of tension sleeve 1 sits flush with ledge 4.
- Screw Metha® impactor/expeller instrument ND655R into threaded hole 3 of ND656R and tighten it.

Validated processing procedure

Note
Observe all relevant national regulations and standards with regard to processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the processing of the products.

Note
Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note
Successful processing of this medical product can only be ensured through a validated processing procedure. The user/processor is responsible for the validation.

Due to process tolerances, manufacturer's specifications can only serve as guide values for the processing procedures applied by the individual user.

Note
Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of "stainless" steels. To avoid this, the time interval between application and processing should not exceed 6 h, and neither fixating pre-cleaning temperatures > 45 °C nor any fixating disinfecting agents (active ingredient: aldehyde, alcohol) should be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading of the laser engraving on stainless steel.

Residues containing chlorine or chlorides, as found e.g. in surgical residues, tinctures, medicines, saline solutions, service water used for cleaning, and in cleaning/disinfecting agents, cause corrosion damage (pitting, stress corrosion) and thereby destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals recommended by the chemicals' manufacturers as effective cleaning/disinfecting agents compatible with the materials used in the respective product. Strictly observe all application instructions, e.g. regarding temperatures, concentrations, exposure times, etc. Failure to do so can result in the following problems:

• Optical changes to the material, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH > 8 to cause visible surface changes or

• Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.

➤ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.a-k-i.org, Publications, Red Brochure – Proper maintenance of instruments.

Preparations at the place of use

- Disassemble the product immediately after use, as described in the respective instructions for use.
- Remove visible residues as completely as possible, using a lint-free, single-use cleaning tissue.

Preparation prior to cleaning

- Have the product dry in a disposal container and ready for immediate cleaning and disinfecting within 30 min after use.

Cleaning/Disinfecting

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

➤ Use cleaning and disinfecting agents approved for high-grade steel according to manufacturer's instructions.

➤ Observe specifications regarding concentration, temperature and exposure time.

Manual cleaning/disinfecting

- Check visible surface for residues after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual cleaning with immersion disinfecting and cleaning with brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV free; pH = 9
II	Intermediate rinse	RT (cold)	1	-	D-W	-
III	Final rinse	RT (cold)	0.5	-	FD-W	-
IV	Drying	RT	-	-	-	-

D-W: Drinking water
FD-W: Fully desalinated water (demineralized)
RT: Room Temperature

Stage I

- Fully immerse the product in the disinfecting solution. Make certain that all accessible surfaces are moistened.
- Clean the product under running tap water, using a suitable cleaning brush (toothbrush and tube brush) if necessary, until all visible residues have been removed from the surface.
- Brush through all surfaces that are not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 5 min or until no more residues can be removed.
- After that, thoroughly rinse (at least 5 times) with the cleaning solution these difficult to access areas of the product, using a single-use syringe (20 ml).
- To avoid the risk of corrosion, do not use a metal brush or other abrasives, which would damage the surfaces, for cleaning the product.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Stage III

- Carry out a full rinse of the product (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage IV

- Dry the product completely with compressed air.

Mechanical cleaning/disinfecting

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber cleaning/disinfecting machine without ultrasound

- Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Prerinse	< 25/77	3	D–W	-
II	Cleaning	55/131	10	FD–W	- Concentrate, alkaline: pH = 10.9 < 5 % anionic tensides - 1 % solution: pH = 10.5
III	Neutralization	20/68	2	FD–W	- Concentrate, acid: pH = 2.6 Basis: Citric acid - 0.2 % solution: pH = 3.0
IV	Intermediate rinse	70/158	1	FD–W	-
V	Thermal disinfecting	94/201	10	FD–W	-
VI	Drying	90/194	40	-	-

D–W: Drinking water
FD–W: Fully desalinated water (demineralized)

Control, care and inspection

- Allow the product to cool down to room temperature.
- Lightly lubricate moving parts (e.g. joints and links) with a maintenance oil suitable for the respective sterilization process (e.g. Aesculap STERILIT® spray JG600 or maintenance oil JG598).
- Inspect the product after each cleaning and disinfecting cycle to be sure it is clean, functioning properly, not damaged, has intact insulation and does not have any loose, bent, broken, cracked, worn, or fractured components.
- Check for compatibility with associated products.
- Set aside the product if it is damaged.

Packaging

- Position the product in its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Pack the product in such a way that the packaging will prevent recontamination of the product in the period between reprocessing and reuse.

Sterilization method and parameters

Note
The product can be sterilized either in disassembled or in assembled condition.

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilization acc. to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST 46-1993, validated acc. to EN ISO 17665 or EN 554/ISO 13683
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270–275 °F	4 min	4 min

WARNING for the US market

If this device is/was used in a patient with, or suspected of having Creutzfeldt–Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of cross-contamination!

Storage

- Store processed products in germ-proof packaging in a dust-protected, dry, dark and temperature-controlled room.

Technical Service

For service, maintenance and repairs, please contact your national B. Braun/Aesculap agency.
Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-2700
Fax: +49 7461 16-2887
E-mail: ats@aesculap.de

Or in the US:
Aesculap Implant Systems Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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