Metha®

GB	Instructions for use Metha® short-stem prosthesis
USA	Instructions for use Metha® short-stem prosthesis <i>Note for U.S. users</i> <i>U.S. cleared device instructions for use can be found at</i> <i>www.aesculapimplantsystems.com "Products & Services".</i> <i>If you do not have internet access please call your local Aesculap</i> <i>representative or customer service at 1-800-258-1946.</i>
D	Gebruuchsanweisung Metha® Kurzschaftprothese
F	Mode d'emploi Prothèse Metha® à tige courte
E	Instrucciones de manejo Prótesis de vástago corto Metha®
I	Istruzioni per l'uso Protesi a stelo corto Metha®
P	Instruções de utilização Haste protésica curta Metha®
DK	Brugsanvisning Metha® protese med kort skaft
FIN	Käyttöohjeet Metha® lyhytvartinen proteesi
GR	Οδηγίες χρήσης Metha® πρόθεση κοντού στελέχους
CZ	Návod k použití Metha® protéza s krátkým dílkem
PL	Instrukcja użytkowania Metha® proteza z krótkim trzpieniem
CN	使用说明 组件式短柄假体
KR	사용방법및사용시주의사항 인공영덩이뼈관절 (Metha short hip stem, 형명별 개별기재)

B|BRAUN

SHARING EXPERTISE

Aesculap AG
Am Aesculap-Platz
78532 Tuttlingen
Germany
Phone +49 7461 95-0
Fax +49 7461 95-2600
www.aesculap.de

TA-Nr. 012200 11/09 Änd.-Nr. 32607

CE marking according to directive 93/42/EEC
CE-Kennzeichnung gemäß Richtlinie 93/42/EWG
Marquage CE conforme à la directive 93/42/CEE
Identificación CE en conformidad con la directriz 93/42/CEE
Marchio CE conforme alla direttiva 93/42/CEE
Simbolo CE, em conformidade com a Directiva 93/42/CEE
CE-markering iht. retningslinje 93/42/EEC
93/42/EEC-standardin mukainen CE-hyväksyntä
Ενδειξη CE σύμφωνα με την Οδηγία 93/42/EEC
CE-označení podle směrnice 93/42/EHS
Oznakowanie CE zgodnie z wymaganiami dyrektywy 93/42/WE
根据 93/42/EEC 指令的 CE 标志

Technical alterations reserved
Technische Änderungen vorbehalten
Sous réserve de modifications techniques
Sujeto a modificaciones técnicas
Con riserva di modifiche tecniche
Salvo alterações técnicas
Retten til tekniske ændringer forbeholdes
Oikeusteknisistä syistä johtuviin muutoksiin pidätetään
Επιφύλασσόμεθα για τεχνικές αλλαγές
Technické zmeny vyhradeny
Zmiany techniczne zastrzeżone
Zmiany techniczne zastrzeżone
保留技术变动的权力

Intended use

The implant is used:

- as a component of a human hip endoprosthesis: Hip endoprosthesis stem
- to combine with Aesculap hip endoprosthesis components
- for implantation without bone cement

The short stem is available as either modular or non-modular implant.

The modular short stem consists of:

- Prosthesis stem
- Modular cone adapter

Materials

The materials used in the implant are listed on the package.

- ISOTAN®: titanium forged alloy Ti6Al4V acc. to ISO 5832-3
 - PLASMAPORE®: μ -CaP surface coating made of pure titanium acc. to ISO 5832-2 with additional calcium phosphate layer
 - ISODUR®: cobalt chrome alloy CoCr29Mo acc. to ISO 5832-12
- ISOTAN®, ISODUR® and PLASMAPORE® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Use for acute disorders of the hip that cannot be treated by other therapies:

- Degenerative osteoarthritis
- Rheumatic arthritis
- Femoral head necrosis

Contraindications

Do not use in the presence of:

- Joint diseases that can be treated by reconstructive surgery (e.g. displacement osteotomy)
- Acute or chronic periarticular or systemic infections
- Secondary diseases that could influence joint implant functionality
- Systemic conditions and metabolic disorders influencing the functioning of the joint implant
- Acute osteoporosis or osteomalacia
- Severely damaged bone structures that could prevent stable implantation of implant components
- Bone tumors in the region of implant fixation
- Bone deformities, axis misalignments, or other bone conditions that rule out the implantation of a hip joint prosthesis preserving the collum femoris
- Expected overload on the joint implant, especially due to high patient weight or major physical strains and activities
- Medicinal or other drug abuse or alcohol addiction affecting the functionality of the joint implant, especially through temporary loss of control or impairment of the implant integration in the bone
- Expected lack of patient cooperation in dealing with the joint implant, especially during rehabilitation, and non-compliance with doctors' recommendations
- Foreign body sensitivity to the implant materials
- Skeletal immaturity
- Neuromuscular diseases impairing the affected extremity
- Cups whose implantation, in analogy to ISO 21535:2002, result in a limited range of movement (less than 80° in flexion/extension or less than 60° in abduction/adduction), e.g. Bipolar cup.
- Prosthesis heads with neck length XL in combination with short-stem prosthesis sizes 0 and 1.
- Prosthesis heads with neck length XXL

Side-effects and adverse interactions

- Dislocation, loosening, wear, corrosion and fracture of implant components
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Venous thrombosis, lung embolism, and cardiac arrest
- Tissue reaction to implant materials
- Injury to blood vessels and nerves
- Hematomas and wound healing disorders
- Periarticular calcification
- Reduced joint mobility and flexibility
- Arthralgia and reduced tolerance for exercise

Safety notes

- The implant components were tested and approved for use in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in the present documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be thoroughly familiar with bone anatomy, including the pathways of nerves, blood vessels, muscles and tendons.
- It is the operating surgeon's responsibility to ensure the correct combination of implant components and their implantation.
- Aesculap is not responsible for any complications arising from incorrect establishment of indication, incorrect choice of implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or lack of asepsis.
- The instructions for use of the individual Aesculap implant components must be observed.
- Do not use damaged or surgically excised components under any circumstances.
- Implants that have been used once may not be reused.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.

- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Combine modular implant components only with suitable Aesculap hip endoprostheses.
- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- Avoid damaging the implant, especially in the neck or cone region, by applying instruments (e.g. HF surgical devices) close to the implant.

CAUTION

Risk of implant component breakage due to combination with implant components from other manufacturers!

- Use Aesculap implant components only.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation (min. dose 25 kGy).
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging. Do not use implant components that are past their expiry date or whose packaging is damaged.
- Do not resterilize the prosthesis stem.

Implants made of non-coated metal may be individually resterilized in steam following appropriate preliminary cleaning and disinfecting.

- Make certain that such implants are not damaged during cleaning and sterilization.
- The protective prosthesis cone caps (labeled accordingly) can be resterilized with steam.

WARNING

The responsibility for resterilization of implant components lies with the end-user. Such resterilization is expressly prohibited by law in certain countries (e.g. France).

- Never resterilize and reuse surface-coated implants (PLASMAPORE® μ -CaP).

Sterilization method and parameters

- Sterilize with steam according to the following rules: Sterilization must be carried out through a validated steam sterilization process (e.g. in a sterilizer according to EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated according to EN 554/ISO 13683). For the fractionated vacuum process, sterilization must be carried out using the 134 °C/2 bar program with a minimum holding time of 5 minutes.

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- Implantation instruments, including the special Aesculap implant system instruments, are complete and functional.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is familiar with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical intervention has been explained to the patient and their agreement concerning the following information has been recorded:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
- The prosthetic joint insert can come loose due to excessive load, wear and tear, or infection.
- The lifespan of the prosthetic joint depends on the body weight and the load put on the joint.
- The prosthetic joint insert must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may become necessary if the implant loosens.
- In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
- The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.

The implantation site is prepared in the following way:

- Following a high femoral neck osteotomy, open the femur with the implant-specific awl.
- Introduce the awl until it reaches the lateral corticalis.

WARNING

Bone fractures/perforations in the implant bed will adversely affect the implant fixation!

- Avoid bone fractures by applying careful operative techniques.
- Treat bone fractures with appropriate intra- and postoperative therapies.
- Handle the implant components properly.

- Gradually prepare the implant bed with the implant-specific form rasps (starting with the smallest size).
- Check and, if necessary, correct the implant position, depth and antetorsion.
- Check the stem size against the form rasp that was introduced last in the correct position.
- Attach the rasp trail adapter and perform trial reposition with the appropriate trial heads; change rasp trial adapter or trial heads as appropriate.
- Check the mobility/range of movement and stability of the joint and the length of the leg.
- Remove the trial implants and the rasp.
- When using the modular short stem, couple the prosthesis shaft and the modular cone adapter **prior** to insertion.

WARNING

Implant malfunction due to insufficiently fixated cone adapter!

- Fixate the modular cone adapter in the short stem by applying a hammer blow and using instrument ND401R.
- Securely fixate the modular cone adapter in the short stem.

- Insert the short shaft (modular/non-modular) with specific insertion instruments.

WARNING

Loosening of the cone adapter during insertion caused by impacting on the Metha stem!

- Carry out the insertion procedure by impacting on the cone adapter with the protective cone cape in place.

Optional procedure

Trial reposition with inserted modular short stem

- Remove the trial cone adapter only after the stem component has been impelled.
- Insert the stem component and punch it into the implant bed with well-dosed blows.
- If necessary, change the trial cone adapter, making certain that the arrow points to "MEDIAL".
- Carry out a trial reposition.
- Check the mobility/range of movement and stability of the joint and the length of the leg.
- Remove the trial implants.

WARNING

Fracture and increased wear of the modular implant components, or loosening of the implant connection due to soiled or damaged joining surfaces!

- Carefully rinse the inside cone of the stem and the joining surfaces of the cone adapter; clean with cleaning swabs (ND619 or ND622) and dry.
- Check the joining surfaces for any damage.
- Do not join modular implant components more than once.
- Always rinse, clean and dry the cone prior to joining.

- Select the implant cone adapter according to the fitting trial cone adapter.
- Carefully insert the implant cone adapter in the short stem already implanted, making certain that the arrow points to "MEDIAL".
- Fixate the implant cone adapter with hammer taps on the appropriate insertion instrument.

Concluding steps for applying (modular/non-modular) short stems

- Select a prosthesis ball according to the trial heads.
- Verify that the cone size of the prosthesis stems are the same as the cone size of the prosthesis heads (see cone size on the implant packaging, e.g. 12/14).
- Rinse, clean and dry the outside cone of the stem and, if necessary, the inside cone of the prosthesis heads prior to setting the prosthesis head in place.
- Do not remove the protective cap until immediately before positioning the prosthesis head.
- Attach the prosthesis ball and prosthesis cone to each other at room temperature only. If necessary, allow the implant to cool down to room temperature.
- Put the prosthesis ball in its position.
- In order to forestall abnormal wear and tear on the prosthesis: Remove all obvious bone cement residues and bone chips before closing the wound.

Further information on B. Braun/Aesculap implant systems can be obtained from B. Braun/Aesculap or the relevant B. Braun/Aesculap agency.

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