Rongeur

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Symbols on product and packages

Symbol	Explanation
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The rongeurs are used for the ablation of soft tissue (e.g. intervertebral disk tissue, mucosa, ethmoidal cells, cell tissue, nucleus tissue, sequester material) in neurosurgery and spine surgery.

The rongeurs cannot be used for the ablation of bones or adjacent bone structures, for which purpose Aesculap recommends using bone punches, because of their higher mechanical load capacity.

Available sizes

The rongeurs are available with various:

- jaw lengths
- jaw angulations
- grind orientations
- jaw widths
- jaw shapes
- working lengths of the slider part
- handle designs (branch or ring)

Note

Detailed information about the available sizes/variants can be found in our brochures.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge or experience.
- Read, follow and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Remove the transport packaging and thoroughly clean the new product, either by hand or by a mechanical process, prior to its initial sterilization.
- Store any new or unused products in a dry, clean and safe place.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Replace any damaged components immediately with original spare parts.

Safe operation

WARNING

Risk of injury and/or malfunction

- Do not use instruments with permanently deformed or misaligned jaws or with jaws that do not completely close. Submit such instruments for repair.

CAUTION

Damage or destruction of the jaws due to overloading!

- Use the rongeurs only for to their intended use (i.e. ablation of soft tissue).
- Do not cut or lever bones or adjacent bone structures.
- Avoid overloading the jaws by tearing off, holding or levering tissue; avoid lateral forces on the jaws (do not twist them).
- Do not use rongeurs with blunt edges, but submit them for repair.

Validated processing procedure

Note

Adhere to national statutory regulations, international standards and guidelines, and local, clinical hygiene instructions for sterile processing.

Note

Mechanical processing should be preferred over manual cleaning because of the better and more reliable cleaning results of mechanical processing.

Note

Successful processing of this medical product can only be ensured if processing is performed through a validated processing procedure. The user/processor is responsible for the validation.

Note

Up-to-date information on processing can be found on the Aesculap Extranet at www.aesculap-extra.net

General notes

To avoid unnecessary, excessive contamination of the complete instrument tray during operations, take care that contaminated instruments are collected separately and not put back into the instrument tray.

Encrusted or fixated residues from surgery can make the cleaning process more difficult or ineffective, and can cause corrosion of stainless steels. Therefore the time interval between application and cleaning should not exceed 6 h, and neither potentially fixating pre-cleaning temperatures >45 °C nor any fixating disinfectants (active ingredient: aldehyde, alcohol) be used.

Excessive doses of neutralizers or basic detergents can cause chemical degradation and/or fading and obliteration of laser inscriptions on stainless steel surfaces, regarding visual reading and machine-readability of the inscriptions. Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfecting and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. To remove such residues, the products must be rinsed sufficiently with fully desalinated water and dried thoroughly.

Only use process chemicals that have been tested and approved (e.g. VAH/DGHM or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations. All process parameters specified by the chemical's manufacturer, such as temperatures, concentrations and exposure times, must be strictly observed. Failure to do so can result in the following problems:

- Optical deterioration, e.g. fading or discoloration of titanium or aluminum surfaces. For aluminum, pH >8 in the application/process solution can already cause visible surface changes.
- Material damage, e.g. corrosion, cracks, fracturing, premature aging or swelling.
- Do not use oxidizing chemicals (e.g. H₂O₂) for mechanical cleaning, as these can cause bleaching or layer loss.
- Do not use process chemicals that would cause stress cracking or brittleness of plastics.
- Clean the product immediately after use.

- Encrustations on instruments are softened and dissolved most effectively and gently by immersion treatment in a 3 % H₂O₂ solution for approx. 5 minutes. Subsequently, the debris can be removed by hand, with a medium-hard brush and/or in an ultrasonic cleaning bath. After this pre-treatment, continue through the normal steps of the processing procedure.

For further detailed advice on hygienically safe and material-/value-preserving reprocessing, see www.a-k-i.org

- Use suitable cleaning/disinfecting agents if the product is put away in wet condition. To prevent foam formation and reduced effectiveness of the process chemicals: Prior to mechanical cleaning and disinfecting, rinse the product thoroughly with running water.

Preparations at the place of use

- Rinse surfaces that cannot be visually inspected, e.g. on products with hidden crevices or lumens or products with complex geometries, preferably with distilled water, using e.g. a disposable syringe.
- Remove visible surgical residues as completely as possible, using a lint-free wet wipe.
- Put the dry product into a closed disposal container and have it transferred to cleaning and disinfecting within 6 h.

Cleaning/Disinfecting

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning and disinfecting agents approved for high-grade steel according to manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed the maximum allowable cleaning temperature of 55 °C.

Note

This product cannot be disassembled. It is cleaned and disinfected in fully assembled condition.

- Carry out ultrasound cleaning:
 - as an effective mechanical supplement to manual cleaning/disinfecting.
 - as a pre-cleaning procedure for products with encrusted residues, in preparation for mechanical cleaning/disinfecting.
 - as an integrated mechanical support measure for mechanical cleaning/disinfecting.
 - for additional cleaning of products with residues left after mechanical cleaning/disinfecting.

Manual cleaning/disinfecting

- Clean products with movable links with the links opened or mobilized.
- Inspect visible surfaces for residual contamination after manual cleaning/disinfecting.
- Repeat the cleaning process if necessary.

Manual ultrasonic cleaning and immersion disinfection

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	5	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	1 st Intermediate rinse	RT (cold)	1	-	D-W	-
III	Disinfecting	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
IV	Final rinse	RT (cold)	1	-	FD-W	-
V	Drying	-	-	-	-	-

D-W: Drinking water

FD-W: Fully desalinated water (demineralized)

RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush.
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, use a 20-ml single-use syringe to rinse thoroughly, for at least 5 times, these difficult to access parts of the product.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage III

- Fully immerse the product in the disinfecting solution.

Stage IV

- Rinse the product completely (all accessible surfaces) under running water.
- Allow water to drip off for a sufficient length of time.

Stage V

- Thoroughly dry the product with lint-free tissue or medical-quality filtered compressed air.

Mechanical cleaning/disinfecting with manual pre-cleaning

Note
The disinfectant must be tested and approved effectiveness (e.g. DGHM or FDA approval or CE mark according to DIN EN ISO 15883).

Note
Use appropriate (neutral, enzymatic and mild-alkaline) detergents for this product with aluminum components.

Note
For thermal disinfection, always use fully desalinated (demineralized) water. Ensure that Ao is >3 000 for the process.

Note
The disinfectant used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Stage	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Ultrasonic cleaning	RT (cold)	15	2	D-W	BBraun Stabimed; aldehyde phenol and QAV-free
II	Irrigation	RT (cold)	1	-	D-W	-

D-W: Drinking water
RT: Room Temperature

Stage I

- Clean the product in an ultrasound cleaning bath (frequency 35 kHz). Make sure that all accessible surfaces are moistened and acoustic shadows are avoided.
- Clean the product with a suitable cleaning brush until all visible residues have been removed from the surface.
- Brush through all surfaces not accessible to visual inspection, e.g. in products with hidden crevices, lumens or complex geometries, for at least 1 min or until no more residues can be removed, using a suitable cleaning brush.
- Mobilize non-rigid components, e.g. set screws, links, etc. during cleaning.
- After cleaning, use a 20-ml single-use syringe to rinse thoroughly, for at least 5 times, these difficult to access parts of the product.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.

Stage II

- Rinse the product completely (all accessible surfaces) under running water.

Mechanical alkaline cleaning and thermal disinfecting

Machine type: Single-chamber washer/disinfector without ultrasound
➤ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).
➤ Keep working ends open during cleaning.
➤ Ensure that product links and joints are opened when placing the product on the tray.

Stage	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Prerinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	BBRAUN HELIMATIC CLEANER alkaline with tensides, application solution 0.5 %
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfecting	90/194	5	FD-W	-
V	Drying	-	-	-	According to disinfectant program

D-W: Drinking water
FD-W: Fully desalinated water (demineralized)

Inspection, maintenance and checks

CAUTION

Damage (metal seizure/friction corrosion) to the product caused by insufficient lubrication!

- Prior to function checks, lubricate moving parts (e.g. joints, pusher components and threaded rods) with maintenance oil suitable for the respective sterilization process (e.g. for steam sterilization: Aesculap STERILIT® I oil spray JG600 or STERILIT® I drip lubricator JG598).

- Allow the product to cool down to room temperature.
- After each complete cleaning, disinfecting and drying cycle, check that the product is: dry, clean, operational, and free of damage (e.g. broken insulation or corroded, loose, bent, broken, cracked, worn or fractured components).
- Dry the product if it is wet or moist.
- Repeat cleaning and disinfecting of products that still show impurities or contamination.
- Check the product for proper functioning.
- Immediately sort out damaged or inoperative products and have them sent to Aesculap Technical Service, see Technical Service.
- Send in to have repaired product with blunt cutting edge, see Technical Service.

Packaging

- Appropriately protect products with fine working tips.
- Sort the product into its appropriate storage device or put it on a suitable tray. Make certain that all cutting edges are protected.
- Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- Ensure that the packaging provides sufficient protection against recontamination of the product during storage (DIN EN ISO 11607).

Sterilization method and parameters

- Make certain that all external and internal surfaces will be exposed to the sterilizing agent (e.g. by opening all valves and faucets).
- Validated sterilization process
 - Steam sterilization through fractionated vacuum process
 - Steam sterilizer according to DIN EN 285, validated according to DIN EN ISO 17665
 - Sterilization through fractionated vacuum process at 134 °C/holding time 5 min
- When sterilizing several products at the same time in one steam sterilizer: Make certain that the maximum allowable load capacity of the steam sterilizer, as specified by the manufacturer, is not exceeded.

Sterilization for the US market

- Aesculap does not recommend the device sterilized by flash or chemical sterilization.
 - Sterilization may be accomplished by steam autoclave in a standard prevacuum cycle.
- To achieve a sterility assurance level of 10⁻⁶, Aesculap recommends the following parameters:

Aesculap Orga Tray/sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Pre-vacuum	270 °F–275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these instruments. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.
WARNING for the US market
If this device is/was used in a patient with, or suspected of having Creutzfeldt-Jakob Disease (CJD), the device cannot be reused and must be destroyed due to the inability to reprocess or sterilize to eliminate the risk of crosscontamination.

Storage

- Store sterile products in germ-proof packaging under dust protection in a dry, dark and temperature-controlled room.

Technical Service

- For service and repairs, please contact your national B. Braun/Aesculap agency.
- Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1602
Fax: +49 7461 16-5621
E-Mail: ats@aesculap.de

Or in the US:
Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood, MO 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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