en	Instructions for use/Technical description S4® Spinal System – Augmentation Screw
de	Gebrauchsanweisung/Technische Beschreibung S4® Spinal System – Augmentationsschraube
fr	Mode d'emploi/Description technique S4® Spinal System – vis d'augmentation
es	Instrucciones de manejo/Descripción técnica Tornillo de refuerzo S4 Spinal System
it	Istruzioni per l'uso/Descrizione tecnica S4® Spinal System – vite di augmentatione
pt	Instruções de utilização/Descrição técnica Parafuso de augmentação do S4® Spinal System
nl	Gebruiksaanwijzing/Technische beschrijving S4® Spinal System – augmentatieschroef
sv	Bruksanvisning/Teknisk beskrivning S4® Spinal System – augmentationsskruv
ru	Инструкция по применению/Техническое описание Усилительный винт S4® Spinal System
cs	Návod k použití/Technický popis Augmentační šroub S4® Spinal System
pl	Instrukcja użytkowania/Opis techniczny S4® Spinal System– śruba augmentacyjna
sk	Návod na použitie/Technický opis Augmentačná skrutka systému S4 Spinal System
tr	Kullanım Kılavuzu/Teknik açıklama S4® Spinal System – Ogmentasyon vidası
zh	使用说明 / 技术说明 S4® 胸腰椎后路固定系统 – 骨水泥增强型螺钉

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.aesculap.com

Aesculap® – a B. Braun brand

TA012865 2019-09 V6 Change No. 60104

Aesculap® S⁴® Spinal System – Augmentation Screw

Note
Generally, the S⁴ Spinal System – augmentation screw is covered by the instructions for use of the S⁴ Spinal System – Lumbar/Deformity TA011187. This information on the augmentation screw is intended to complement the respective information contained in the instructions for use of the S⁴ Spinal System – Lumbar/Deformity.

Symbols on product and packages

	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Caution, general warning symbol Caution, see documentation supplied with the product

Intended use

The S⁴ Spinal System implants are used for dorsal monosegmental and multisegmental stabilization of the lumbar and thoracic spine.

The S⁴ Spinal System – augmentation screw can be fixed with bone cement to increase anchoring stability. In this case, the injection cannula is inserted in the S⁴ Spinal System – augmentation screw for application of the bone cement.

The S⁴ Spinal System – augmentation screw comprises:

- S⁴ Monoaxial/polyaxial screws (augmentation screw), supplied in sterile condition
- S⁴ Element monoaxial/polyaxial screws (augmentation screw), supplied in sterile condition
- Cement injection cannula (sterile), see TA013132
- for percutaneous application with S⁴ Element monoaxial/polyaxial screws (augmentation screws): S⁴ Element Augmentation instruments, see TA014315

Note
There are special S⁴ instruments provided for the implantation of these system components and for the augmentation, distraction, compression, and reduction of the lumbar and thoracic spine.

Materials

The materials used in the implant are listed on the packaging:

- ISOTAN®_F titanium forged alloy Ti6Al4V acc. to ISO5832-3
- ISOTAN®_P pure titanium acc. to ISO 5832-2

The titanium implants are coated with a colored oxide. Slight changes in coloration may occur, but do not affect the implant quality.

ISOTAN® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

Indications

Surgically installed implants are designed to support the normal healing process. They are not intended for use either as replacements for natural body parts or to bear loads over the long term if healing does not occur.

Use for acute disorders of the spine that cannot be treated by other therapies::

- Fractures
- Pseudarthrosis or delayed healing
- Dislocation
- Degenerative instability
- Post-discectomy syndrome
- Spondylolisthesis
- Kyphosis
- Scoliosis
- Stenosis
- Post-trauma instability
- Tumors
- Osteoporosis (with the simultaneous use of bone cement and fenestrated screws)

Risk of screw breakage caused by pedicle screws used for spondylolisthesis!
► **Support the stabilization measure with interbody fusion (ALIF, PLIF or TLIF).**

Absolute contraindications

Do not use in the presence of:

- Severe damage to the bone structures of the spine that could prevent the stable implantation of the implant components; for example, osteopenia, severe osteoporosis, Paget's disease, bone tumors etc.
- Metabolic or degenerative metabolic bone diseases that could compromise the stable anchoring of the implant system
- Suspected allergy or sensitivity to the implant materials
- Acute or chronic vertebral infections of a local or systemic nature
- Poor patient compliance or limited ability to follow medical instructions, particularly in the post-op phase, including with regard to the restrictions on range of movement in terms of physical exercise and occupational activity
- Cases not listed under indications

Relative contraindications

In the following circumstances, use of the implant system could represent an increased clinical risk and therefore requires precise, individual assessment by the surgeon:

- Medical or surgical conditions that could negatively impact the success of the implantation, including wound healing disorders
- Conditions that could subject the spine and implants to excessive pressure; for example, pregnancy, obesity, neuromuscular diseases or disorders
- Generally poor condition of the patient; for example, drug or alcohol addiction

Side effects and interactions

The general risks of surgery are assumed known and are therefore not described in these instructions for use. Potential risks and clinical consequences associated with use of the implant system include, but are not limited to:

- Implant failure resulting from excessive load
 - Warping or bending
 - Loosening
 - Breakage
- Poor stabilization, requiring subsequent revision surgery or surgery to remove the implant system
- Failed or delayed fusion
- Cement leakages (cement emerging from the vertebral body)
- Infection
- Fractured vertebral body or bodies
- Injuries to
 - Nerve roots
 - Spinal cord
 - Blood vessels
 - Organs
- Changes of the normal spine lordosis
- Neurological disorders, including impairment of the gastrointestinal, urological and/or reproductive system
- Pain or indisposition
- Bursitis
- Decreased bone density due to load avoidance
- Bone atrophy/fracture above or below the spine section provided for
- Limited physical ability
- Persistence of symptoms that were to be treated by the implantation
- Paresis
- Subsequent surgery to remove the implant system
- Metal/foreign body allergy or hypersensitivity
- Dural tears
- Meningitis

Safety notes

CAUTION
Federal law restricts this device to sale by, or on order of a physician!

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon must be fully conversant with bone anatomy, including the pathways of nerves, blood vessels, muscles, and tendons.
- The operating surgeon is responsible for combining and implanting the implant components, with or without bone cement.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate aseptics.
- The user must follow the instructions that accompany each of the Aesculap implant components.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants already used must not be used again.
- Do not use instruments belonging to another system or made by another manufacturer.
- Delayed healing can cause implant breakage due to metal fatigue.
- The surgeon in charge of the therapy shall make any decision with regard to the removal of implant components that have been used.
- Damage to the weight-bearing bone cement and/or bone structures can cause loosening of the components, fracture of the bone or implant and other grave complications.
- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are gamma-sterilized.
- Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration date or whose packaging is damaged.

Damage to implants caused by processing and resterilization!
► **Do not reprocess or resterilize the implants.**

Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand
- Operating conditions are highly aseptic
- The implantation instruments, including the special Aesculap implant system instruments, are complete and in working condition.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- In the case of delayed or incomplete fusion, the implants can break and loosen due to high loads.
- The life span of the implant depends on the patient's body weight.
- The implant components must not be overloaded by extreme strains, hard physical labor, or sports.
- Corrective surgery may be necessitated by implant loosening, fracture, or loss of correction.
- Smokers present an increased risk of bone fusion failure.
- The patient must undergo medical check-ups of the implant components at regular intervals.
- The patient must be advised that being overweight places such excessive pressure on the implant that the prospect of failure becomes likely.
- The patient must be advised that strenuous physical activity in the form of heavy lifting, excessive muscle strain, exaggerated twisting, repeated bending and bending over or running should be avoided until the bones have healed fully.
- The patient must be instructed to keep pressure off the implant post surgery in order to prevent the onset of clinical problems associated with implant failure.

To implant the S⁴ Spinal System – augmentation screw, proceed as follows:

- Select the appropriate S⁴ Spinal System variant and anchoring implants according to the individual indication, preoperative planning, and bone situation found intraoperatively.
- When selecting S⁴ Spinal System screw types, make certain to select the correct diameter, length, and canal orientation.

WARNING

Incorrect application can result in trauma to the spinal cord and nerve roots!

- **Align and screw in the instruments, the pedicle screw, and the polyaxial screw only under X-ray control or using a navigation system.**

- To avoid internal stress on, and weakening of, the implant: avoid scoring or scratching of the implant components.
- Only allow axial load on the screws during screw insertion. Avoid induction of lateral forces during screw insertion.

WARNING

Danger of screw breakage due to hard bone or contact with lateral cortex!

- **Prepare the vertebral body over the entire screw length with thread cutter.**

WARNING

Risk of connecting rod breakage!

- **Do not bend back or bend excessively the connecting rods and cross connectors.**

CAUTION

Damage to the hexagon socket due to incorrect application of the screwdriver tip!

- **Make certain the screwdriver tip is fully inserted in the hexagon socket of the screw head when turning the screw.**

- If the screw head loosens during revision of the larger polyaxial screws (Ø 7.5 mm/Ø 8.5 mm), remove the remaining screw component with the SW 3.5 screwdriver.

S⁴ Spinal System – fixing the S4 augmentation screw with bone cement

To increase anchoring stability, the S⁴ Spinal System – augmentation screw can be fixed with bone cement. The following must be observed when cementing:

Note
Cementing is carried out only after all screws have been aligned and the connecting rods have been adapted, and before locking the S⁴ Spinal System – augmentation screw!

WARNING

Risk of breaking the injection cannula by inducing excessive shear forces!

- **Do not subject the injection cannula to transverse loads.**
- **Do not bend the injection cannula.**

Note
Once the S⁴ Spinal System – augmentation screw has been cemented, it will be impossible to change the implant position, use the reduction lever, or make corrections!

Note
With the S⁴ Spinal System – augmentation screw finally cemented, explanting the screw can become extremely difficult, especially in cases of osteopenia!

- Prior to applying the bone cement, fully screw in the S⁴ injection cannula at the thread in the bone screw head, without applying perceptible force.
- Check the bone cement for proper texture and viscosity before positioning the applier on the injection cannula. Observe the manufacturer's instructions for the bone cement!
- Cement S⁴ Spinal System – augmentation screw with a cement suitable for screw augmentation.

WARNING

Risk of injury from cement leakages!

- **Always apply imaging techniques for control during the preparation and positioning of the implants.**
- **Make certain the slots of the S⁴ Spinal System – augmentation screw are positioned fully inside the bone.**
- **Make sure that the pedicle cortex is not perforated.**
- **Do not penetrate the opposite cortex and ensure that the opposite cortex is undamaged.**
- **Do not insert or cement the S⁴ Spinal System – augmentation screw in a fractured or partially resected vertebral body. Perform a thin-slice CT scan preoperatively.**
- **Screw the S⁴ injection cannula into the screw head without applying any perceptible force. Make certain the axes of the pedicle screw and the injection cannula are aligned with each other.**
- **Apply cement only under high-quality radiographic control (max. volume of cement per polyaxial screw: 2 ml).**

WARNING

Risk of injury due to material incompatibility!

- **Only use application systems that have been approved for use with bone cement.**

WARNING

Contamination of the locking thread/Insufficient fixation when removing the application system, due to incorrect viscosity of the cement!

- **Make certain the cement has the correct viscosity.**

- Do not allow the cement to cure completely before removing the injection cannula. Remove the cement application system incl. cannula as soon as the cement shows a doughy texture (allow for a certain amount of cement for reference purposes when preparing the cement).

WARNING

Wound contamination due to insufficient curing of cement! Risk of caulking of the locking screw thread, resulting in system instability!

- **Observe the correct curing period after removing the injection cannula.**
- **For monitoring the cement texture, prepare a reference volume of cement together with the cement to be applied.**
- **Make certain the system is not exposed to any load before the cement has reached its final strength.**

Further information on Aesculap implant systems is always available from B. Braun/Aesculap or the appropriate B. Braun/Aesculap office.

Disposal

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Implant Systems LLC
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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