
en Instructions for use/Technical description

ELAN 4 electro remote foot control GA810

USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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



Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
	Manufacturer
	Date of manufacture
	Manufacturer serial number
	Manufacturer's batch name
	Manufacturer's article number
	Delivery quantity
	Non-sterile medical product
	CE mark according to Regulation (EU) 2017/745
	Medical device
	Temperature limits during transport and storage
	Atmospheric pressure limits during transport and storage
	Air humidity limits during transport and storage
	Preselection of motor rotation direction
	Pump selector

Changing the maximum speed

Class AP equipment according to IEC/DIN EN 60601-1

Rx only

US federal law restricts the supply of the product exclusively to or by order of a physician

Labeling of electrical and electronic devices according to directive 2012/19/EU (WEEE), see Disposal

Contents

1.	About this document.	3
1.1	Scope	3
1.2	Safety instructions.	3
2.	General information	3
2.1	Intended use.	3
2.2	Application Environment.	3
3.	Safe handling.	3
4.	Product description	3
4.1	Scope of supply	3
4.2	Components required for operation.	3
4.3	Operating principle	3
5.	Preparation and setup	4
5.1	First use	4
6.	Working with the product.	4
7.	Disassembling the pedal for cleaning	4
8.	Mounting the pedal (after cleaning)	4
9.	Processing procedure.	4
9.1	General safety notes	4
9.2	Reusable products	4
9.3	Preparation before cleaning	4
9.4	Cleaning/disinfection	4
9.5	Wiper disinfection for electrical equipment	5
9.6	Manual cleaning including immersion disinfection (recommended for heavily soiled objects)	5
9.7	Inspection, maintenance and checks	5
10.	Maintenance	6
11.	Troubleshooting list.	6
12.	Technical Service	6
13.	Accessories/Spare parts.	6
14.	Technical data	6
14.1	Classification according to Regulation (EU) 2017/745	6
14.2	Basic data, information about standards.	6
14.3	Ambient conditions	6
15.	Disposal	7

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply to the ELAN 4 electro remote foot control GA810.

- For article specific instructions for use and material compatibility, see B. Braun eIFU at eifu.bbraun.com

1.2 Safety instructions

Safety instructions specify the dangers to patient, user and/or product that can occur during the use of the product. Safety instructions are labeled as follows:

DANGER

Indicates a potential hazard. If not avoided, death or serious injury may result.

WARNING

Indicates a potential hazard. If not avoided, minor or moderate injury may result.

CAUTION

Indicates a potential threat of material damage. If not avoided, the product may be damaged.

2. General information

2.1 Intended use

The ELAN 4 electro remote foot control GA810 is an accessory for the ELAN 4 electro motor system.

With the remote foot control:

- the user activates the motor of the applied part and selects the rotational direction,
- the user activates/deactivates the coolant and irrigation pump,
- the user can change the maximum speed of the motor.

Note

The ELAN 4 electro remote foot control GA810 cannot be used if a manually controlled applied part or a cable-connected foot control is used.

2.2 Application Environment

Application in the unsterile area

Remote foot controls are Class AP devices. The foot control circuit is ignition-safe and approved for operation in medical environments according to IEC /DIN EN 60601-1. The housing is constructed according to Protection Type IPX8.

Set-up location

Floor

3. Safe handling

WARNING

Risk of injury and material damage if this product is not used as intended!

- Use the product only according to its intended use.

WARNING

Risk of injury and damage to property due to improper handling of the product!

This product is an accessory for the ELAN 4 electro control unit GA800.

- Follow the instructions for use of the ELAN 4 electro control unit GA800 (TA014401).
- Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- Prior to use, check that the product is in good working order.
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.

Note

The user agrees to report all serious incidents occurring in connection with the product to the manufacturer and to the competent authority of the country in which the user is registered.

4. Product description

4.1 Scope of supply

Art. no.	Designation
GA810	ELAN 4 Electro remote foot control
TA014411	Instructions for use of remote foot control GA810 (brochure)

4.2 Components required for operation

- ELAN 4 electro control unit GA800
- ELAN 4 electro motor cable for foot control GA806
- ELAN 4 electro applied part
- Tool

4.3 Operating principle

With the ELAN 4 electro remote foot control GA810 the active applied part is connected to the control unit without manual control. The coolant and irrigation pump is activated or deactivated and the maximum speed and direction of the applied part is changed.

Function button 2:

Short actuation	Activation/deactivation of the coolant and irrigation pump. The pump only operates when the motor is running.
Prolonged Pressing	Flush function The pump will operate until the actuation is stopped.

Rotational direction button for motor 4: Selection of right or left rotation.
Pedal 3: The pedal serves as the actuator for the speed setting of the control unit. Depending on the actuating path, the speed of the motor can be selected continuously from 0 to the maximum speed set.

Maximum speed button 5: The maximum speed button is for step-by-step changes to the set maximum speed on the control device GA800 with full operation of the pedal.

5. Preparation and setup

Non-compliance with the following instructions will preclude all responsibility and liability in this respect on the part of Aesculap.

- ▶ When setting up and operating the product, adhere to
 - national regulations for installation and operation,
 - national regulations on fire and explosion protection.

5.1 First use

Assembling/disassembling the transport brace

- ▶ Insert transport brace 1 on one side into the holding receptacle for transport brace 6, see Fig. A2.
- ▶ Widen the transport brace 1 and insert it into the other holding receptacle for transport brace 6. Make certain that the transport brace is correctly mounted.

Note

Disassembly is carried out in the reverse order.

6. Working with the product

How to work with the ELAN 4 electro system and its accessory components is described in the instructions for use of the ELAN 4 electro control unit GA800, see TA014401.

7. Disassembling the pedal for cleaning

- ▶ Lift the pedal 3 upwards, see Fig. B.
- ▶ Remove pedal 3 from the guide bolts for pedal 7.
- ▶ Remove the pedal 3 in an upward direction, see Fig. C.

8. Mounting the pedal (after cleaning)

- ▶ Insert pedal 3 onto guide bolts for pedal 7, see Fig. B.
- ▶ Fold pedal 3 down.

9. Processing procedure

9.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

9.2 Reusable products

The life time of the product is limited due to damage, normal wear and tear, type and duration of use, in addition to the handling, storing and transporting of the product.

Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product.

9.3 Preparation before cleaning

- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Carry out non-fixating/NaCl-free pre-cleaning immediately after use.

Note

The pedal can be removed for cleaning, see Disassembling the pedal for cleaning.

9.4 Cleaning/disinfection

Product-specific safety instructions for the reprocessing procedure

⚠ CAUTION

Damage to, or destruction of the product caused by mechanical cleaning/disinfection!

- ▶ Only clean and disinfect the product manually.
- ▶ Do not sterilize the product under any circumstances.

⚠ CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents!

- ▶ Only use cleaning/disinfecting agents approved for surface cleaning. Follow the manufacturer's instructions for the respective cleaning/disinfecting agent.

Note

If the product is heavily soiled, Aesculap recommends manual cleaning including immersion disinfection.

9.5 Wiper disinfection for electrical equipment

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT	1	-	-	17 % Propan-1-ol, 0.23 % Didecyldimethylammonium chloride*
II	Wipe disinfection	RT	≥ 1	-	-	17 % Propan-1-ol, 0.23 % Didecyldimethylammonium chloride*

RT: Room temperature

*Recommended: Meliseptol® Wipes sensitive (B. Braun)

- Disassemble pedal for cleaning, see Disassembling the pedal for cleaning.

Phase I

- Remove any visible residues with a disposable disinfectant wipe.

Phase II

- Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- Observe contact time requirements (at least 1 min).

9.6 Manual cleaning including immersion disinfection (recommended for heavily soiled objects)

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Disinfecting cleaning	RT (cold)	>15	2	P-W	Concentrate free from aldehydes, phenols and QAV, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	-	P-W	-
III	Disinfection	RT (cold)	15	2	P-W	Concentrate free from aldehydes, phenols and QAV, pH ~ 9*
IV	Final rinse	RT (cold)	1	-	DI-W	-
V	Drying	RT	-	-	-	-

P-W: Potable water

DI-W: de-ionized water (demineralized, microbiological at least in compliance with potable water regulations)

RT: Room temperature

*Recommended: BBraun Stabimed

- Disassemble pedal for cleaning, see Disassembling the pedal for cleaning.

Phase I

- Immerse product at least 15min completely in the actively cleaning disinfecting solution. Ensure that all accessible surfaces are wetted.
- Clean the product in the solution using a suitable cleaning brush until no residue is visible on the surface.

Phase II

- Rinse product completely (all accessible surfaces) under running water.
- Allow sufficient time for water to drip off.

Phase III

- Completely immerse the product in the disinfectant solution.

Phase IV

- Rinse/flush product completely (all accessible surfaces).
- Allow sufficient time for water to drip off.

Phase V

- During the drying phase dry the product using the appropriate equipment (e.g. wipes, compressed air).

9.7 Inspection, maintenance and checks

- Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- Set aside the product if it is damaged.
- Change battery when the energy level is low:
 - Remove the 4 screws on the bottom of the device with a screwdriver (Torx TX20).
 - Remove the cover of the battery compartment. Replace batteries (3x) (battery type D / LR20 / Mono).
 - Note the polarity of the batteries.
 - Replace the battery compartment cover and secure with 4 screws.

10. Maintenance

To ensure reliable operation, the product must be maintained at least once a year.

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

11. Troubleshooting list

Note

For further information, see instructions for use of ELAN 4 electro control unit GA800 (TA014401).

12. Technical Service

⚠ DANGER

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- ▶ Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.
- ▶ Do not modify the product.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

13. Accessories/Spare parts

Art. no.	Name
TA014411	Instructions for use for ELAN 4 electro remote foot control GA810 (brochure)

14. Technical data

14.1 Classification according to Regulation (EU) 2017/745

Art. no.	Name	Class
GA810	ELAN 4 electro remote foot control	I

14.2 Basic data, information about standards

For connection to	ELAN 4 electro control unit GA800
Degrees of protection in accordance with IEC/ DIN EN 60601-1	AP
Degrees of protection provided by enclosure according to IEC/ EN 60529	IPX8
Weight	3810 g ±10 %
Dimensions (L x W x H) without transport bracket	271 mm x 266 mm x 85 mm ±10 %
Dimensions (L x W x H) with transport bracket	271 mm x 310 mm x 200 mm ±10 %
Compliance with standards	IEC/DIN EN 60601-1
EMC	IEC/DIN EN 60601-1-2
Battery type	D / LR20 / Mono (3 x)
Radio standards	ETSI EN 300 440 ETSI EN 301 489-1 ETSI EN 301 489-3 IEC/DIN EN 62311
ISM band	2 403 MHz to 2 480 MHz
CISPR11	Class A

14.3 Ambient conditions

	Operation	Transport and storage
Temperature	10 °C to 40 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

15. Disposal

WARNING

Risk of infection due to contaminated products!

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

Note

The user institution is obliged to process the product before its disposal, see Processing procedure.

The recycling passport can be downloaded from the extranet as a PDF document using the individual part number. (The recycling passport is a disassembly instruction for the equipment containing information on the proper disposal of environmentally harmful components.)

A product marked with this symbol must be disposed separately with electrical and electronic equipment. Within the European Union the manufacturer performs the disposal free of charge.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

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14.1 (AB) 2017/745 yönetmeliği uyarınca sınıflandırma