
en Instructions for use/Technical description

ELAN 4 electro foot control GA808

USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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



Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
	Manufacturer combined with date of manufacture (year)
	Manufacturer's serial number
	Preselection of motor rotation direction
	Pump selector
	Class AP equipment in accordance with IEC/DIN EN 60601-1
	Labeling of electrical and electronic equipment in accordance with Directive 2012/19/EC (WEEE), see Disposal
	Date of manufacture
	Manufacturer's batch designation
	Manufacturer's order number
	Delivery quantity
	Non-sterile medical device
	According to federal law, this product may only be sold by or by the order of a physician
	CE marking in accordance with Regulation (EC) 2017/745
	Medical device

	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage

Contents

1.	About this document.	3
1.1	Scope	3
1.2	Warnings	3
2.	General information	3
2.1	Intended purpose	3
2.2	Application Environment	3
3.	Safe handling	3
4.	Product description	4
4.1	Scope of supply	4
4.2	Components required for operation	4
4.3	Operating principle	4
5.	Preparation and setup	4
5.1	First use	4
6.	Working with the ELAN 4 electro foot control	4
7.	Disassembling the pedal for cleaning	4
8.	Mounting the pedal (after cleaning)	4
9.	Processing procedure	4
9.1	General safety notes	4
9.2	Reusable products	4
9.3	Preparation before cleaning	4
9.4	Cleaning/disinfection	5
9.5	Wiper disinfection for electrical equipment	5
9.6	Inspection, maintenance and checks	5
10.	Maintenance	6
11.	Troubleshooting list	6
12.	Technical Service	6
13.	Accessories/Spare parts	6
14.	Technical data	6
14.1	Classification in accordance with Regulation (EC) 2017/745 ..	6
14.2	Basic data, information about standards	6
14.3	Ambient conditions	6
15.	Disposal	7

1. About this document

Note

General risks of surgical interventions are not described in these instructions for use.

1.1 Scope

These instructions for use apply to the ELAN 4 electro foot control GA808.

- For article specific instructions for use and material compatibility, see B. Braun eFU at eifu.bbraun.com

1.2 Warnings

Warnings indicate hazards to patients, users and/or products that may emerge during the use of the product. Warnings are marked as follows:

DANGER

Indicates a potential danger. If it is not prevented, death or severe injuries may result.

WARNING

Indicates a potential danger. If it is not prevented, minor or moderate injuries may result.

CAUTION

Indicates a potential damage to property. If it is not avoided, the product may be damaged.

2. General information

2.1 Intended purpose

The ELAN 4 electro foot control GA808 is an accessory of the ELAN 4 electro motor system.

With the foot control:

- the user activates the motor of the applied part and selects the rotational direction,
- the user activates or deactivates the cooling or flushing liquid pump, provided that a hand controlled applied part is not being used.

Cable length

5 m

2.2 Application Environment

Application in the unsterile area

Foot controls are Class AP devices. The foot control circuit is ignition-safe and approved for operation in medical environments according to IEC /DIN EN 60601-1. The housing is constructed according to Protection Type IPX8.

Set-up location

Floors

3. Safe handling

WARNING

Risk of injury and material damage due to inappropriate use of the product!

- Use the product only in accordance with the intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

This product is an accessory of the ELAN 4 electro control unit GA800.

- Follow the instructions for use of the ELAN 4 electro control unit GA800 (TA014401).
- Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- Prior to use, check that the product is in good working order.
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.

Note

The user agrees to report all serious incidents occurring in connection with the product to the manufacturer and to the competent authority of the country in which the user is registered.

4. Product description

4.1 Scope of supply

Art. no.	Designation
GA808	ELAN 4 electro foot control
TA014409	Instructions for use of foot control GA808 (leaflet)

4.2 Components required for operation

- ELAN 4 electro control unit GA800
- ELAN 4 electro motor cable for foot control GA806
- ELAN 4 electro applied part
- Tool

4.3 Operating principle

The ELAN 4 electro foot control GA808 is used to activate the active application part connected to the control unit without manual control. The cooling or flushing liquid pump activates or deactivates.

Function button 2:

Short actuation	Turns the cooling or flushing liquid pump on and off. The pump only works when a motor is running
Longer actuation	Flush function The pump operates until activation is ended.

Rotational direction button for motor 4: Selection of right or left rotation.

Pedal 3: The pedal serves as the actuator for the speed setting of the control unit. Depending on the actuating path, the speed of the motor can be selected continuously from 0 to the maximum speed set.

5. Preparation and setup

Non-compliance with the following instructions will preclude all responsibility and liability in this respect on the part of Aesculap.

- ▶ When setting up and operating the product, adhere to
 - national regulations for installation and operation,
 - national regulations on fire and explosion protection.

5.1 First use

Assembling/disassembling the transport brace

- ▶ Insert transport brace 1 on one side into the holding receptacle for transport brace 5, see Fig. A.
- ▶ Widen the transport brace 1 and insert it into the other holding receptacle for transport brace 5. Make certain that the transport brace is correctly mounted.

Note

Disassembly is carried out in the reverse order.

6. Working with the ELAN 4 electro foot control

How to work with the ELAN 4 electro system and its accessory components is described in the instructions for use of the ELAN 4 electro control unit GA800, see TA014401.

7. Disassembling the pedal for cleaning

- ▶ Lift the pedal 3 upwards, see Fig. B.
- ▶ Remove pedal 3 from the guide bolts for pedal 6.
- ▶ Remove the pedal 3 in an upward direction, see Fig. C.

8. Mounting the pedal (after cleaning)

- ▶ Insert pedal 3 onto guide bolts for pedal 6, see Fig. B.
- ▶ Fold pedal 3 down.

9. Processing procedure

9.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

9.2 Reusable products

The service life of the product is limited by damage, normal wear, the nature and duration of the application, and handling, storage and transportation of the product.

Careful visual and functional testing prior to the next use is the best way to identify a malfunctioning product.

9.3 Preparation before cleaning

- ▶ Separate the products from each other immediately after use.
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Carry out non-fixating/NaCl-free pre-cleaning immediately after use.

Note

The pedal can be removed for cleaning, see Disassembling the pedal for cleaning.

9.4 Cleaning/disinfection

Product-specific safety instructions for the reprocessing procedure

⚠ CAUTION

Damage to, or destruction of the product caused by mechanical cleaning/disinfection!

- ▶ Only clean and disinfect the product manually.
- ▶ Do not sterilize the product under any circumstances.

⚠ CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents!

- ▶ Only use cleaning/disinfecting agents approved for surface cleaning. Follow the manufacturer's instructions for the respective cleaning/disinfecting agent.

9.5 Wiper disinfection for electrical equipment

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Cleaning	RT	1	-	-	17 % Propan-1-ol, 0.23 % Didecyldimethylammonium chloride*
II	Wipe disinfection	RT	≥ 1	-	-	17 % Propan-1-ol, 0.23 % Didecyldimethylammonium chloride*

RT: Room temperature

* Validation was performed using Meliseptol® wipes sensitive (B.Braun)

Phase I

- ▶ Remove any visible residues with a disposable disinfectant wipe.

Phase II

- ▶ Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- ▶ Observe contact time (at least 1min).

9.6 Inspection, maintenance and checks

- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Set aside the product if it is damaged.

10. Maintenance

To ensure reliable operation, the product must be maintained at least once a year.

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

11. Troubleshooting list

Note

For further information, see instructions for use of ELAN 4 electro control unit GA800 (TA014401).

12. Technical Service

DANGER

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- ▶ Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.
- ▶ Do not modify the product.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

13. Accessories/Spare parts

Art. no.	Name
TA014408	Instructions for use for ELAN 4 electro foot control GA808 (A4 for ring-binder)
TA014409	Instructions for use of ELAN 4 electro foot control GA808 (leaflet)

14. Technical data

14.1 Classification in accordance with Regulation (EC) 2017/745

Art. no.	Name	Class
GA808	ELAN 4 electro foot control	I

14.2 Basic data, information about standards

For connection to	ELAN 4 electro control unit GA800
Protection class according to IEC/ DIN EN 60601-1	AP
Enclosure protection class according to IEC/ DIN EN 60529	IPX8
Weight	2240 g ±10 %
Dimensions (L x W x H) without transport brace	245 mm x 143.7 mm x 75.5 mm ±10 %
Dimensions (L x W x H) with transport brace	245 mm x 193.5 mm x 196.4 mm ±10 %
Cable length (without connector)	5.0 m ±5 %
Conforming to standard	IEC/EN 60601-1
EMC	IEC/DIN EN 60601-1-2

14.3 Ambient conditions

	Operation	Transport and storage
Temperature	10 °C to 40 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

15. Disposal

WARNING

Risk of infection from contaminated products!

- **Observe national regulations when disposing or recycling the product, its components and their packaging.**

Note

The user institution is obliged to process the product before its disposal, see Processing procedure.

The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.)

Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

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14.1 (AB) 2017/745 yönergesi uyarınca sınıflandırma