
(GB) Instructions for use/Technical description

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









E

Aesculap®

ELAN 4 electro MIS hand piece with shafts and cutters

Legend

- 1 ELAN 4 MIS shaft coupling
- 2 Unlocking sheath for hand piece
- 3 Connector for motor cable
- 4 Slider for tool release
- 5 ELAN 4 MIS hand piece shaft (L10 or L13, straight or curved)
- 6 Length marking (L10 or L13)
- 7 ELAN 4 MIS shaft connection
- 8 ELAN 4 MIS tool (L10 or L13)
- 9 Symbol of the applied part type in display of ELAN 4 electro control unit
- 10 ELAN 4 single spray nozzle for L10 or L13
- 11 ELAN 4 electro oil spray adapter GB600860
- 12 ELAN 4 oil spray adapter for MIS shafts GB600870

Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
 YYYY-MM	Maintenance label Indication of the next maintenance appointment (Date: Year-Month)
	Machine-readable, two-dimensional code The code contains a unique serial number which can be used for electronic tracking of the individual instrument. The serial number is based on the global standard sGTIN (GS1).
	Manufacturer
	Date of manufacture
	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Manufacturer's batch designation
	Manufacturer's serial number
	Manufacturer's article number
	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage
L10/L13	Length coding of shafts and tools

Contents

1.	Applicable to	3
2.	General information	3
2.1	Intended use	3
2.2	Main functions and design characteristics	3
2.3	Indications	4
2.4	Absolute contraindications	4
2.5	Relative contraindications	4
3.	Safe handling	4
3.1	Sterile ELAN 4 MIS tools	5
3.2	Sterile ELAN 4 single-use spray nozzles	5
4.	Product description	5
4.1	Package includes	5
4.2	Components required for operation	5
4.3	Operating principle	5
5.	Preparation	5
6.	Working with the applied part, the shafts and the tools	6
6.1	System set-up	6
6.2	Function check	7
6.3	Operation	7
7.	Validated reprocessing procedure	8
7.1	General safety notes	8
7.2	General information	8
7.3	Preparations at the place of use	8
7.4	Preparation before cleaning	8
7.5	Cleaning/disinfection	9
7.6	Manual cleaning with immersion disinfection – applied part/attachments	10
7.7	Mechanical cleaning/disinfecting with manual pre-cleaning – applied part and shafts	11
7.8	Inspection, maintenance and checks	12
7.9	Packaging	12
7.10	Steam sterilization	12
7.11	Sterilization for the US market	13
7.12	Storage	13
8.	Maintenance	13
9.	Troubleshooting	13
10.	Technical Service	13
11.	Accessories/Spare parts	14
11.1	MIS hand piece shafts	14
11.2	MIS tools	14
12.	Technical Data	14
12.1	Classification acc. to Directive 93/42/EEC	14
12.2	Performance data, information about standards	14
12.3	Tool speed/Setting	15
12.4	Operating mode	15
12.5	Ambient conditions	15

13.	Disposal	15
14.	Distributor in the US/Contact in Canada for product information and complaints	15

1. Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

2. General information

2.1 Intended use

Task/Function

The ELAN 4 electro MIS hand piece GA860 with shafts and tools is an accessory of the ELAN 4 electro motor system.

The MIS hand piece is connected by an ELAN 4 electro motor cable to the control unit.

The MIS hand piece is used in the operation of the ELAN 4 MIS tools in combination with the ELAN 4 MIS hand piece shafts.

Application Environment

Application in the sterile area

2.2 Main functions and design characteristics

Speed

min. 0 rpm to max. 80 000 rpm

Rotational direction

Right and left rotation

Operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

■ 30 s application, 30 s pause

■ ∞ repetitions at 60,000 min⁻¹

■ Max. Temperature 48 °C

Electrical systems generally heat up during continual operation. It is advised to give the system a break after use to cool down, as indicated here.

Aesculap®

ELAN 4 electro MIS hand piece with shafts and cutters

Heating depends on the tool used and the load. After a certain number of repetitions, the system should cool down. This procedure prevents the system overheating as well as possible injury to the patient or user.

The user is responsible for the use and adherence to the pause sequence described.

2.3 Indications

Modes of application	Separation, removal and modeling of hard tissue, cartilage, etc. and bone replacement material
Surgical discipline/areas of application	Neurosurgery, ENT surgery, orthopaedics

Note

The type and area of application depend on the hand piece shafts and tools selected.

2.4 Absolute contraindications

The product is not licensed for use on the central nervous system or central circulatory system.

2.5 Relative contraindications

The safe and effective use of the product greatly depends on influences which can only be controlled by the user. Therefore the specifications provided represent framework conditions only.

Clinically successful use of the product is dependent on the knowledge and experience of the surgeon. The surgeon must decide which structures it is sensible to treat and take into account the safety and warning information contained in these instructions for use.

3. Safe handling

CAUTION

Federal law restricts this device to sale by or on order of a physician!

WARNING

Risk of injury and material damage if this product is not used as intended!

- Use the product only for its intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

This product is an accessory of the ELAN 4 electro control unit GA800.

- Follow the instructions for use of the ELAN 4 electro control unit GA800 (TA014401).
- Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
 - Remove the transport packaging and clean the new applied parts and shafts, either manually or mechanically, prior to its initial sterilization.
 - Prior to use, check that the product is in good working order.
 - Observe "Notes on Electromagnetic Compatibility (EMC)", see TA022130.
 - To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.
- Ensure that the electrical installations at the premises where the device is to be used meet all requirements acc. to IEC/DIN EN.
- Do not use the control unit and applied part in potentially explosive areas.
- Carry out sterile processing for any reusable applied parts, shafts and motor hoses prior to use.
- When handling, observe instructions for use that are relevant to Aesculap holder systems TA009721, see Aesculap Extranet at <https://extranet.bbraun.com>

3.1 Sterile ELAN 4 MIS tools

DANGER

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!
► **Do not reprocess the product.**

The product is gamma-sterilized and supplied in sterile packaging.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see General information.
- Do not use products from open or damaged sterile packaging.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Do not use the product after its use-by date.

3.2 Sterile ELAN 4 single-use spray nozzles

DANGER

Risk of infection for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product!
► **Do not reprocess the product.**

The product is gamma-sterilized and supplied in sterile packaging.

The product must not be reused.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see General information.
- Do not use products from open or damaged sterile packaging.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Do not use the product after its use-by date.

4. Product description

4.1 Package includes

Art. no.	Designation
GA860	ELAN 4 electro MIS hand piece
GB600870	ELAN 4 oil spray adapter for MIS shafts
TA014441	Operating instructions for ELAN 4 electro MIS hand piece with shafts and cutters (brochure)

4.2 Components required for operation

- Operation-ready ELAN 4 electro control unit GA800, see TA014401
- ELAN 4 MIS hand piece shaft
- ELAN 4 MIS-tool

4.3 Operating principle

The ELAN 4 electro MIS hand piece is equipped with an ELAN 4 MIS shaft coupling and an ELAN 4 MIS tool coupling.

The MIS hand piece is used with ELAN 4 electro foot controls or hand controls.

With the MIS hand piece, the coupled tool rotates at the set motor speed. The motor speed of the MIS hand piece can be adjusted continuously by the hand/foot control.

The MIS hand piece can be operated in right or left (clockwise or counter-clockwise) rotation mode.

Note

The length marking L10/L13 shows the shaft length of the tool to be utilized. Only ELAN 4 MIS tools may be attached to the applied parts with the markings (L10/L13) of shaft and tool matching.

5. Preparation

Non-compliance with these rules will result in complete exclusion of liability on the part of Aesculap.

- Do not use products from open or damaged sterile packaging.
- Prior to use, inspect the product and its accessories for any visible damage.
- Use the products and their accessories only if they are in perfect technical condition.

Aesculap®

ELAN 4 electro MIS hand piece with shafts and cutters

6. Working with the applied part, the shafts and the tools

WARNING

Risk of infection and contamination!
Applied parts, shafts, and motor cable are delivered unsterile!

- Carry out sterile processing of the applied parts, shafts and motor cable prior to use, in accordance with the instructions for use.

WARNING

Risk of injury and material damage due to inadvertent activation of the applied part!

- Secure applied parts that are not actually used against inadvertent activation ("Off" position), see instructions for use for ELAN 4 electro control unit GA800 (TA014401).

WARNING

Damage to the product if dropped!

- Use the products only if they are in perfect technical condition, see Function check

WARNING

Risk of burns to skin and tissue due to blunt tools/ lack of maintenance of the applied part!

- Use tools only if they are in perfect condition.
- Replace blunt tools.
- Ensure correct maintenance of the applied part, see Maintenance.

6.1 System set-up

Note

Operating elements on the system components of the ELAN 4 electro motor system are labeled with a gold marking.

Connecting the accessories

Do not use accessories in combinations that are not mentioned in these instructions for use.

Coupling the tool to the applied part

- Lock applied part ("Off" position), see instructions for use of ELAN 4 electro control unit GA800 (TA014401).
- Push the ELAN 4 MIS hand piece shaft 5 up to the stop into ELAN 4 MIS shaft coupling 1.
The ELAN 4 MIS hand piece shaft clicks into place.
- Pull on the ELAN 4 MIS hand piece shaft 5 to make sure it is securely attached.

Detaching the shaft from the applied part

- Pull back the unlocking sheath 2 on the applied part and pull the ELAN 4 MIS hand piece shaft 5 out of the ELAN 4 MIS shaft coupling 1.

Coupling the tool to the applied part

WARNING

Injury and damage to equipment by wrong combinations of shafts and tools!

- Do not use L13 tool in L10 shaft.

WARNING

Risk of injury when coupling/uncoupling tools in the "On" position due to accidental activation of the applied part!

- Only couple/uncouple tools in the "Off" position.

Note

The tool is difficult to insert into/remove from the bent hand piece shafts due to the bend.

- Lock applied part ("Off" position), see instructions for use of ELAN 4 electro control unit GA800 (TA014401).
- Push tool 8 up to the stop in ELAN 4 MIS hand piece shaft 5, see Fig. A. The tool engages.
- Pull at tool 8 to check that it is securely coupled.

Uncoupling the tool from the applied part

WARNING

Risk of injury when coupling/uncoupling tools in the "On" position due to accidental activation of the applied part!

- Only couple/uncouple tools in the "Off" position.

- Lock the applied part (Off position), see instructions for use of ELAN 4 electro control unit GA800 (TA014401).
- Pull the slide for tool release 4 on the applied part and pull tool 8 out of ELAN 4 MIS hand piece shaft 5, see Fig. A.

Mounting the spray nozzle and tube set

- Connect the hose set GA395SU to the pump of the ELAN 4 electro control unit GA800, see the instructions for use TA014401.
- Push application-side tube end onto tube connection of the single-use spray nozzle 10.
- Push the single-use spray nozzle 10 onto ELAN 4 MIS hand piece shaft 5 and position it by pushing/turning, see Fig. A.

Note

The spray tube of the single-use spray nozzle for is bendable by hand and can be adapted to the requirements of application (e.g. aligned to tool head).

6.2 Function check

- ▶ Prior to each use, check that all products to be used are in good working order.
- ▶ Check the secure connection of all products to be used.
- ▶ Pull on the hand piece shaft to check that it is securely coupled.
- ▶ Pull on the tool to check that it is securely coupled.
- ▶ Briefly run the applied part at maximum speed.
- ▶ Check for damage to, and irregular running noises, excessive vibrations and excessive heat-up of applied part and the shaft.
- ▶ Check that the applied part or shaft is only warm to touch after a short period of use, no hotter.
- ▶ Do not use the product if it is damaged or defective. Set aside the product if it is damaged.

6.3 Operation

WARNING

Coagulation of patient tissue or risk of burns for patients and users from hot applied part/tool!

- ▶ Cool the tool during operation.
- ▶ Put down the applied part/tool beyond reach of the patient.
- ▶ Allow the applied part/tool to cool down.
- ▶ Use a cloth to protect against burns when changing the tool.

WARNING

Risk of infection from aerosol formation!

Risk of injuries caused by particles coming loose from the tool!

- ▶ Take appropriate protective measures, e. g. watertight protective clothing, face mask and protective goggles, suction extraction.

WARNING

Risk of injury and/or malfunction!

- ▶ Always carry out a function test prior to each use of the product.

WARNING

Risk of injury due to application of the product outside the field of view!

- ▶ Only use the product under vision.

WARNING

Risk of injury and damage to the tool/system!

The rotating tool may catch cover drapes (textile or other).

- ▶ Never allow the tool to touch cover drapes (textile or other) during operation.

WARNING

Risk of injury from bent or snapped tools!

- ▶ Only use perfectly straight tools.
- ▶ Only apply mild pressure when using a tool.
- ▶ Avoid overstraining (e. g. bending) of tools during operation.

WARNING

Risk of injury due to small parts that could enter the operation site!

- ▶ Do not perform any change of reamer over the operating field.

WARNING

Danger of injury from foreign body reactions (e. g. inflammation, encapsulation) if particles remain in the body!

When using diamond tools, it is always possible for diamond grains/nickel particles to break away/off.

- ▶ Carefully rinse and vacuum while and after working with diamond tools.

WARNING

Risk of injury if the tool comes loose accidentally!

- ▶ Check for secure fixation of the tool after any tool change.

WARNING

Risk of injury and infection if surgical gloves are damaged by sharp cutting edges!

- ▶ Avoid contact with cutting tools.

WARNING

Risk of injury and material damage due to inappropriate use of tools!

- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ When coupling/uncoupling, handle tools with cutting edges with care.

WARNING

Risk of damage to the tool due to removal or dissection of other materials than bone or when removing hard or soft tissue (e. g. sawing into instruments or implants). Risk of injuries and infections caused by the abrasion debris produced in this way!

- ▶ When operating, be careful to avoid contact with tools, implants and similar objects.

CAUTION

Damage to micro instruments caused by incorrect handling!

- ▶ Protect products against overstraining.
- ▶ Do not bend the hand piece shaft.
- ▶ Do not use the shaft as a lever.

ELAN 4 electro MIS hand piece with shafts and cutters

Note

The hand piece shaft has a limited service life. When used as intended, a minimum of 30 uses are guaranteed. Correct oiling of the shafts is prerequisite for this; see *Inspection, maintenance and checks*. The end of the service life will manifest itself through a relatively quick increase in the temperature of the shaft.

It is only possible to operate the applied part and change the setting parameters on the control unit if:

- the applied part is connected to the control unit,
- there is no second applied part released ("On" position) at the same time and
- the applied part type 9 is shown on the display of the control unit.

Note

For further information, see *instructions for use of ELAN 4 electro control unit GA800 (TA014401)*.

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

The recommended chemical was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bb Braun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

Residues containing chlorine or chlorides, e. g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result damage to metallic products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e. g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

7.3 Preparations at the place of use

- ▶ Separate the products from each other immediately after use.
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

7.4 Preparation before cleaning

Applied parts and shafts

- ▶ Before the first machine cleaning/disinfection: mount the ECCOS holder GB084R /GB723R in the suitable sieve (such as JF222R).
- ▶ Insert the application part and shafts in the correct position in the ECCOS holder GB084R/GB723R, see Fig. B or Fig. C.

7.5 Cleaning/disinfection

Product-specific safety guidelines on the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents which
 - be approved for plastic material and high-grade steel,
 - do not attack softeners (e. g., in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C during chemical cleaning and/or disinfection.
- ▶ Do not exceed the maximum temperature of 96 °C during thermal disinfection with FD water.
- ▶ Dry the product for at least 10 min at a maximum of 120 °C.

Note

The indicated drying time is a guide time only. It must be checked taking into account the specific conditions (e. g. load) and if applicable adjusted.

Validated cleaning and disinfection procedure

Product	Validated procedure	Reference
Applied parts and shafts	Manual cleaning with immersion disinfection	see Manual cleaning with immersion disinfection – applied part/attachments
	Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection	see Mechanical cleaning/disinfecting with manual pre-cleaning – applied part and shafts

7.6 Manual cleaning with immersion disinfection – applied part/attachments

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Pre-cleaning	RT (cold)	-	-	D-W	-
II	Cleaning	RT (cold)	>5	1	D-W	pH-neutral, pH ~ 9*
III	Intermediate rinse	RT (cold)	-	-	D-W	-
IV	Disinfection	RT (cold)	>15	2	D-W	Concentrate-free, aldehyde-free, phenol-free, QUAT-free, pH ~ 9**
V	Final rinse	RT (cold)	-	-	D-W	-
VI	Drying	RT	-	-	-	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Helizyme

**Recommended: BBraun Stabimed

- Follow the instructions for use of the ELAN 4 rinsing devices and rinsing adapters TA014447 and TA014448.
- Do not clean the applied part in an ultrasound bath.

Phase I

- Clean the applied part/attachments under running water, using a suitable cleaning brush until all visible residues have been removed from the surfaces. Do not brush the inside of the shafts.
- Mobilize non-rigid components, such as release button, unlocking sleeve, etc. during cleaning.
- Flush the applied part and the flush adapter GB698R/GB699R at least 3 times for 5 s with the pressurized water gun.

Phase II

- Fill the interior of the applied part and the shafts through the flush adapter GB698R/GB699R with single-use syringe filled with cleaning solution.
- Place the applied part and shaft, with rinsing adapter connected, immersed in enzymatic cleaning solution for 5 min.

Note

The shafts can be put in the cleaning solution without the flush adapter as well.

Phase III

- Rinse the entire applied part and shafts (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as release button, unlocking sleeve, etc. during rinsing.
- Flush the applied part and the flush adapter GB698R/GB699R at least 3 times for 5 s with the pressurized water gun.

Phase IV

- Prior to manual disinfection, allow rinsing water to drain off the applied part and shafts for a sufficient length of time and use compressed air to blow all liquid out of the product through the rinsing adapter GB698R/GB699R in order to prevent dilution of the disinfecting solution.
- Fill the interior of the applied part and the shafts through the flush adapter GB698R/GB699R with single-use syringe filled with disinfection solution.
- Place the applied part and shaft, with rinsing adapter connected, in enzymatic cleaning solution for 15 min.

Phase V

- Rinse the entire applied part and shafts (all accessible surfaces) under running water.
- Mobilize non-rigid components, such as release button, unlocking sleeve, etc. during cleaning.
- Flush the applied part and the flush adapter GB698R/GB699R at least 3 times for 5 s with the pressurized water gun.

Phase VI

- Dry the applied part and shafts in the drying phase with suitable equipment (e. g. cloth, compressed air).
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfection process if necessary.

7.7 Mechanical cleaning/disinfecting with manual pre-cleaning – applied part and shafts

Note

The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Brushes	RT (cold)	–	–	D–W	–
II	Rinsing	RT (cold)	5	–	D–W	–

D–W: Drinking water

RT: Room temperature

- ▶ Follow the instructions for use of the ELAN 4 rinsing devices and rinsing adapters TA014447 and TA014448.
- ▶ Do not clean the applied part in an ultrasound bath.

Phase I

- ▶ Clean the applied part and shafts with a suitable cleaning brush until all discernible residues have been removed. Do not brush the inside of the shafts.
- ▶ Mobilize non-rigid components, such as release button, unlocking sleeve, etc. during cleaning.

Phase II

- ▶ Connect the connector for motor cable 3 with ELAN 4 electro rinsing device GB692R.
- ▶ Connect the shaft connection 7 to the ELAN 4 flush system for four MIS hand piece shafts GB679R.
- ▶ Flush the applied part and shafts:
 - 3 times for 5 s with a water gun or
 - for 5 min with a tap/hose

Aesculap®

ELAN 4 electro MIS hand piece with shafts and cutters

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH ~ 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Connect the connector for motor cable **3** with ELAN 4 electro rinsing device GB692R.
- ▶ Connect the shaft connection **7** to the ELAN 4 flush system for four MIS hand piece shafts GB679R.
- ▶ Place the rinsing device on a wire basket suitable for cleaning.
- ▶ Connect the rinsing device to the rinsing connection of the rinsing cart.
- ▶ After the mechanical cleaning/disinfection process:
 - Remove the residual rinse water from the applied parts and shafts with a compressed air gun, see the instructions for use for ELAN 4 flush system and flush adapter TA014447 or TA014448.
 - Check visible surfaces for residues.

7.8 Inspection, maintenance and checks

- ▶ Let the applied part and shafts cool to room temperature.
- ▶ Inspect the applied part and shafts after each cleaning and disinfecting cycle to be sure it is clean, functional and undamaged.
- ▶ Connect ELAN 4 electro oil spray adapter GB600860 (gray) to connector for motor cable **3** and spray through the applied part for approx. 2 s with Aesculap-STERILIT oil spray GB600, see Fig. D.
- ▶ Place the ELAN 4 oil spray adapter for MIS shafts GB600870 (white) on the shaft connection **7** and spray through for about 2 s with Aesculap-STERILIT oil spray GB600, see Fig. E.
- ▶ Check the applied part and shafts for any damage, abnormal running noise, overheating or excessive vibration.
- ▶ Set aside the product if it is damaged.

7.9 Packaging

- ▶ Always follow the instructions for use of the respective packaging and storage devices (e. g. instructions for use TA009721 for Aesculap-ECCOS holder system).
- ▶ Insert the application part and shafts in the correct position in the ECCOS holder GB084R/GB723R, see Fig. B or Fig. C.
- ▶ Pack trays appropriately for the intended sterilization process (e. g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

7.10 Steam sterilization

Note

For further information, see *Technical Data*.

- ▶ Make certain that all external and internal surfaces of the product will be exposed to the sterilizing agent.
- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer: ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

7.11 Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min

*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

7.12 Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

8. Maintenance

Note

The following instructions only apply to the ELAN 4 electro MIS hand piece GA860. No maintenance necessary for shafts and tools.

- Do not change or repair shafts.
- Do not modify, re-sharpen or repair the tools.

To ensure reliable operation, the product must be maintained in accordance with the maintenance labeling or at least once a year.

e. g. 2016-07

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

9. Troubleshooting

Note

For further information, see instructions for use of ELAN 4 electro control unit GA800 (TA014401).

10. Technical Service

DANGER

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.
- Do not modify the product.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen /Germany

Phone: +49 7461 95-1601

Fax: +49 7461 14-939

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042

Aesculap Repair Hotline

Phone: +1 800 214-3392

Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

ELAN 4 electro MIS hand piece with shafts and cutters

11. Accessories/Spare parts

Art. no.	Designation
GB796SU	Single-use spray nozzle for MIS hand piece shafts L10
GB797SU	Single-use spray nozzle for MIS hand piece shafts L13
GA395SU	ELAN 4 electro single-use hose set
GB084R	ELAN 4 electro ECCOS holder for MIS hand piece GA860
GB600	STERILIT Power Systems oil spray
GB600860	ELAN 4 electro oil spray adapter
GB600870	ELAN 4 oil spray adapter for MIS shafts
GB692R	ELAN 4 electro flush system for MIS hand piece
GB698R	ELAN 4 electro flush adapter for MIS hand piece
GB679R	ELAN 4 flush system for four MIS hand piece shafts
GB699R	ELAN 4 flush system for MIS hand piece shafts
GB723R	ELAN 4 ECCOS holder for two MIS hand piece shafts
TA014440	Operating instructions for ELAN 4 electro MIS hand piece with shafts and cutters (brochure)
TA014441	Operating instructions for ELAN 4 electro MIS hand piece with shafts and cutters (brochure)

11.1 MIS hand piece shafts

Art. no.	Designation
GB920R	ELAN 4 MIS hand piece shaft L10 straight
GB921R	ELAN 4 MIS hand piece shaft L10 curved
GB922R	ELAN 4 MIS hand piece shaft L10 heavy curved
GB925R	ELAN 4 MIS hand piece shaft L13 straight
GB926R	ELAN 4 MIS hand piece shaft L13 curved
GB927R	ELAN 4 MIS hand piece shaft L13 heavy curved

11.2 MIS tools

Art. no.	Designation
GP420SU-GP437SU	ELAN 4 MIS tools L10, single-use
GP440SU-GP457SU	ELAN 4 MIS tools L13, single-use

12. Technical Data

12.1 Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
GA860	ELAN 4 electro MIS hand piece	Ila
GB920R	ELAN 4 MIS hand piece shaft L10 straight	Ila
GB921R	ELAN 4 MIS hand piece shaft L10 curved	Ila
GB922R	ELAN 4 MIS hand piece shaft L10 heavy curved	Ila
GB925R	ELAN 4 MIS hand piece shaft L13 straight	Ila
GB926R	ELAN 4 MIS hand piece shaft L13 curved	Ila
GB927R	ELAN 4 MIS hand piece shaft L13 heavy curved	Ila
GP420SU-GP437SU	ELAN 4 MIS tools L10, single-use	Ila
GP440SU-GP457SU	ELAN 4 MIS tools L13, single-use	Ila

12.2 Performance data, information about standards

ELAN 4 electro MIS hand piece GA860

Max. power	approx. 140 W
Max. torque	approx. 2.2 Ncm
Max. motor speed	80 000 rpm
Weight	100 g ±10 %
Dimensions	17 mm x 118 mm ±5 %
Tool connection	ELAN 4 MIS shaft
Applied part	Type BF
EMC	IEC/DIN EN 60601-1-2
Conforming to standard	IEC/DIN EN 60601-1

The product can be reprocessed 350 times. This is certified by a high voltage test after 350 cycles (1.6 kV).

MIS hand piece shafts

Art. no.	Designation	Dimen- sions (L)	Weight
GB920R	ELAN 4 MIS hand piece shaft L10 straight	111 mm ±5 %	14 g ± 10 %
GB921R	ELAN 4 MIS hand piece shaft L10 curved	111 mm ± 5 %	14 g ± 10 %
GB922R	ELAN 4 MIS hand piece shaft L10 heavy curved	111 mm ±5 %	14 g ± 10 %
GB925R	ELAN 4 MIS hand piece shaft L13 straight	141 mm ±5 %	19 g ± 10 %
GB926R	ELAN 4 MIS hand piece shaft L13 curved	141 mm ±5 %	19 g ± 10 %
GB927R	ELAN 4 MIS hand piece shaft L13 heavy curved	141 mm ± 5 %	19 g ± 10 %

12.3 Tool speed/Setting

Adjustable motor direction	Clockwise/counterclockwise rotation
Adjustable maximum speed setting	10 000 rpm to 80 000 rpm
Step width for maximum speed setting	5 000 rpm
Factory pre-setting for maximum speed setting	75 000 rpm, right rotation

12.4 Operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

- 30 s application, 30 s pause
- ∞ repetitions at 60,000 min⁻¹
- Max. Temperature 48 °C

12.5 Ambient conditions

	Operation	Storage and transport
Temperature	10 °C to 27 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

13. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

Adhere to national regulations when disposing of or recycling the product, its components and its packaging! The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.) Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

14. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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