




- (GB)** Instructions for use/Technical description
- (USA)** ELAN 4 electro control unit GA800
- (D)** Gebrauchsanweisung/Technische Beschreibung
- ELAN 4 electro Steuereinheit GA800
- (F)** Mode d'emploi/Description technique
- Unité de commande ELAN 4 electro GA800
- (E)** Instrucciones de manejo/Descripción técnica
- Unidad de control ELAN 4 electro GA800
- (I)** Istruzioni per l'uso/Descrizione tecnica
- Centralina ELAN 4 electro GA800
- (P)** Instruções de utilização/Descrição técnica
- Unidade de controlo ELAN 4 electro GA800
- (NL)** Gebruiksaanwijzing/Technische beschrijving
- ELAN 4 electro bedieningssysteem GA800
- (DK)** Brugsanvisning/Teknisk beskrivelse
- ELAN 4 electro styreenhed GA800

- (S)** Bruksanvisning/Teknisk beskrivning
- ELAN 4 electro styrenhet GA800
- (FIN)** Käyttöohje/Tekninen kuvaus
- ELAN 4 electro-ohjausyksikkö GA800
- (EST)** Kasutusjuhend/Tehniline kirjeldus
- ELAN 4 electro, juhtseadis GA800
- (RUS)** Инструкция по применению/Техническое описание
- Блок управления ELAN 4 electro GA800
- (CZ)** Návod k použití/Technický popis
- Řídící jednotka ELAN 4 electro GA800
- (PL)** Instrukcja użytkowania/Opis techniczny
- Jednostka sterująca ELAN 4 electro GA800
- (SK)** Návod na použitie/Technický opis
- Ovládacia jednotka ELAN 4 electro GA800
- (TR)** Kullanım Kılavuzu/Teknik açıklama
- ELAN 4 electro kumanda ünitesi GA800

Aesculap®

ELAN 4 electro control unit GA800

Legend

- 1 ELAN 4 electro control unit GA800
- 2 Display with touch control panel
- 3 Coolant pump
- 4 Flap
- 5 Indicator light
- 6 Connection socket for foot control
- 7 Connection sockets for motor cable
- 8 Power OFF switch
- 9 Power ON indicator
- 10 Power ON switch
- 11 Ventilation grate
- 12 Fuse holder
- 13 Power socket
- 14 Terminal for equipotential bonding conductor
- 15 Bottle holder mount
- 16 Bottle holder
- 17 USB interface: Only intended for use by the manufacturer or by Service technicians authorized by Aesculap.
- 18 Interface RS232: Only intended for use by the manufacturer.

Motor cable/applied part

- 19 Connection for applied part
- 20 Visual field "Off"
- 21 Visual field "On"
- 22 Release button
- 23 Lug
- 24 Connection for motor cable on applied part
- 25 Slider for tool release
- 26 Pull back unlocking sleeve
- 27 Connection for control unit

Foot control

- 28 Function button
- 29 Pedal
- 30 Rotational direction button for motor

Symbols on product and packaging

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
	Follow the instructions for use

	"OFF" (power supply)
	"ON" (power supply)
	Type BF applied part
	Foot control
	Terminal for equipotential connector according to IEC/DIN EN 60601-1
	Fuse
	Alternating current
	Manufacturer combined with date of manufacture (year)
	Labeling of electrical and electronic devices according to directive 2012/19/EU (WEEE), see Disposal
	Date of manufacture
	Manufacturer's batch designation
	Manufacturer's serial number
	Manufacturer's article number
	Delivery quantity
	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage

Applied part types

Symbol	Text	Art. no.	Designation
	Perforator	GA822	ELAN 4 electro perforator driver
	Intra	GA824	ELAN 4 electro Low-speed motor with intra coupling
	Saw	GA831	ELAN 4 electro sagittal saw
		GA832	ELAN 4 electro reciprocating saw
		GA836	ELAN 4 electro micro sagittal saw
		GA837	ELAN 4 electro micro reciprocating saw
	Highspeed	GA849	ELAN 4 electro craniotome and multifunctional handpiece (2-ring)
		GA860	ELAN 4 electro MIS handpiece
		GA861	ELAN 4 electro Standard handpiece (1-ring) L4
		GA862	ELAN 4 electro Standard handpiece (1-ring) L7
		GA863	ELAN 4 electro Standard handpiece (1-ring) L10
		GA864	ELAN 4 electro Standard handpiece (1-ring) L13
		GA865	ELAN 4 electro Standard handpiece (2-ring) L7
		GA866	ELAN 4 electro Standard handpiece (2-ring) L10
		GA867	ELAN 4 electro Standard handpiece (2-ring) L13
		GA868	ELAN 4 electro Standard handpiece (2-ring) L17
		GA869	ELAN 4 electro Standard handpiece (2-ring) L22

Display/control elements in the applied part control panel

Symbol	Designation
	Activation mode of motor via foot control GA808
	Activation mode of motor via hand control GA804/GA805
	Motor direction indicator for clockwise rotation preselected The display depends on whether the motor is activated or disabled.
	Motor direction indicator for counterclockwise rotation preselected The display depends on whether the motor is activated or disabled.
	Reduce the maximum speed setting
	Increase the maximum speed setting

Display/control elements in the pump control panel

Symbol	Designation
	Identification of the control panel for the coolant pump Display in active applied part control panel
	Switch on the pump
	Switch off the pump
	Activate the "Flush" function (continuous flushing)
	Reduce the delivery rate
	Increase the delivery rate

Aesculap®

ELAN 4 electro control unit GA800

Display/control elements in the system settings menu

Symbol	Designation
	Call up the system settings menu
	Leave the system settings menu
	Applied part settings
	Control unit settings
	Control unit information
	Maximum speed/stroke rate
	Speed up
	Slow down
	Flow rate
	Increase the value
	Reduce the value
	System volume
	Display brightness
	System language
	Reset device to factory settings
	Go back to menu structure

Symbol	Designation
	Call up the submenu
	Move down list
	Move back up list

Contents

1.	About this document	5
2.	Intended use	5
2.1	Task/function within system	5
2.2	Application Environment	5
2.3	Indications.	5
2.4	Contraindications	5
3.	Safe handling	6
4.	Product description	6
4.1	Scope of supply.	6
4.2	Components required for operation	6
4.3	Operating principle.	7
5.	Preparation and setup	8
5.1	Set-up environment/set-up location	8
5.2	Stacking of units.	8
6.	Working with the ELAN 4 electro system.	8
6.1	System set-up.	8
6.2	Function checks	11
6.3	Safe operation	12
6.4	System settings menu	14
7.	Validated reprocessing procedure	16
7.1	General safety notes.	16
7.2	Preparation before cleaning	16
7.3	Cleaning/disinfection	16
7.4	Wipe disinfection for electrical devices without sterilization .	17
7.5	Inspection, maintenance and checks	17
8.	Maintenance.	17
9.	Troubleshooting list	18
9.1	Display error messages.	18
9.2	Fuse replacement	22
10.	Technical Service	22
11.	Accessories/Spare parts	23
11.1	ELAN 4 electro motor cables, applied parts and foot controls .	23
11.2	Coolant pump	23

11.3	Power cord.....	23
11.4	Equipotential bonding leads.....	23
11.5	Spare parts.....	23
12.	Technical data.....	24
12.1	Classification acc. to Directive 93/42/EEC.....	24
12.2	Performance data, information about standards.....	24
12.3	Ambient conditions.....	24
13.	Disposal.....	25
14.	Distributor in the US/Contact in Canada for product information and complaints.....	25

1. About this document

This document describes all necessary instructions and steps for the provision, setting and safe operation of the ELAN 4 electro system and its accessory components.

Additional instructions and steps regarding the accessory components, in particular on their connection and reprocessing, can be found in the respective instructions for use or supplement for the components.

2. Intended use

2.1 Task/function within system

The ELAN 4 electro control unit GA800 together with the accessory forms an electrically operated motor system.

The ELAN 4 electro control unit GA800 supplies the power and monitors the motors in the ELAN 4 electro applied parts. The speed commands are sent to the control unit via the hand or foot control. The rotational direction and activation of the pump are controlled via buttons on the foot control.

The purpose of the integrated coolant pump is to deliver coolant or irrigation fluid to the operating field for cooling the tool and tissue or irrigating the operating field.

Maximum delivery rate of the pump	65 ml/min
-----------------------------------	-----------

2.2 Application Environment

The motor system fulfills the requirements of Type BF in accordance with IEC/DIN EN 60601-1.

For use in operating rooms, outside of the explosive hazard zone (e.g. areas with highly pure oxygen or anesthetic gases).

Control unit	
Application Environment	In the unsterile area
Set-up location	Table, ceiling supply unit, equipment cart, or similar

2.3 Indications

Modes of application	Separating, removing and shaping hard tissue, cartilage and similar, as well as for drilling holes in bone and bone replacement materials.
Surgical discipline/areas of application	Neurosurgery, ENT surgery, OMF surgery, orthopedics and trauma surgery

Note

The type and area of application depend on the applied parts and tools selected.

2.4 Contraindications

The ELAN 4 electro motor system is not authorized for use in the central nervous system or central circulatory system.

Note

The safe and effective use of electrically driven applied parts depends heavily on influences that only the user can control him or herself. Therefore the specifications provided represent framework conditions only.

Note

The clinical success of the use of the ELAN 4 electro motor system is dependent on the knowledge and experience of the surgeon. The surgeon must decide which structures it is sensible to treat and take into account the safety and warning information contained in these instructions for use.

3. Safe handling

CAUTION

Federal law restricts this device to sale by, or on order of a physician!

DANGER

Risk of death by electric shock!

- ▶ Do not open the product.
- ▶ Connect the product only to a grounded power supply

WARNING

Risk of injury and material damage due to inappropriate use of the product!

- ▶ Use the product only in accordance with the intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

The ELAN 4 electro control unit GA800 together with the accessory forms an electronically operated motor system.

- ▶ Follow the instructions for use of ELAN 4 electro accessories.
- ▶ Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- The ELAN 4 electro control unit GA800 fulfills the requirements of CISPR11, Class A.
- ▶ Inspect the new product after removing its transport packaging and prior to first use to ensure it is in good working order.
- ▶ Observe "Notes on Electromagnetic Compatibility (EMC)", see TA022130.
- ▶ To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- ▶ Keep the instructions for use accessible for the user.
- ▶ Always adhere to applicable standards.
- ▶ Do not pull any of the cables by the actual cable but only by the connector.

4. Product description

4.1 Scope of supply

Art. no.	Designation
GA800	ELAN 4 electro control unit
GD412804	Bottle holder
TA014401	Instructions for use of ELAN 4 electro control unit
TA014482	Supplement for ELAN 4 electro motor system
TA022130	Notes on Electromagnetic Compatibility

4.2 Components required for operation

- ELAN 4 electro control unit GA800
- Power cord, see Accessories/Spare parts
- ELAN 4 electro motor cable for foot control GA806 and ELAN 4 electro foot control GA808
 - or –
- ELAN 4 electro motor cable with hand lever GA804 (from software version 2.00)
 - or –
- ELAN 4 electro motor cable with hand control GA805 (from software version 2.00)
- ELAN 4 electro applied part

Note

The installed software version is displayed in the system settings menu in the device information sub-item.

For upgrading the control unit to the current software version, contact the national B. Braun/Aesculap agency, see Technical Service.

Note

The term "applied part" covers all handpieces and motors of the ELAN 4 electro motor system, see Accessories/Spare parts.

When using the coolant pump:

- Bottle holder GD412804
- Single-use tube set GA395SU
- Coolant or irrigation fluid: Physiological saline solutions up to 1 000 ml

Note

The coolant/irrigation fluid is not an Aesculap accessory.

4.3 Operating principle

Control unit

The ELAN 4 electro control unit **1** is designed for main voltage ranges of 100 V~ to 120 V~ and 220 V~ to 240 V~ from 50 Hz to 60 Hz. The mains voltage is converted into safety extra low voltage for supplying the micro motors in the applied parts.

The control unit features two connection sockets for applied parts for connecting two different applied parts, and has one connection socket for a foot control. Only one applied part can be operated at a time.

Note

The essential performance features are defined by the "speed" and "rotational direction" specifications. The exception is a conditionally defined motor stop to a safe state caused by a detected error condition

Display/control concept

The display **2** indicates the current device status (operating and error status) at all times. The display is divided into the applied part control panel and the pump control panel.

The display shows the group of the applied part currently connected.

The display shows information about the speed, rotational direction, activation and delivery rate of the irrigation pump. By activating the relevant panel, the control elements appear. Settings can then be changed. If the control elements are not activated, they will disappear again after a short time.

Explanation with an example

Note

If two applied parts are connected to the control unit, the display is divided into 2/3 for the active applied part and 1/3 for the inactive applied part.


Fig. 1 Example of control concept

Legend

- a Pump control panel
- Selected delivery rate % (here 30 %)
- b Applied part control panel
- c Applied part type (here ELAN 4 electro high speed motor)
- d Activation mode: (for foot control)
- e Qualitative bar display:
Display of the maximum set speed (white bar).
The actual current speed in the range from 0 to the maximum set speed (green bar)
Difference between the maximum set speed and the upper limit of the maximum speed (gray bar)
- f Maximum set speed: (here 75 000 rpm)
- g Maximum speed (uppermost limit) with this group: max. 80 000 rpm
- h Call up the system settings menu
- i Rotational direction (here clockwise rotation)

Applied part types

The control unit recognizes the different types of applied parts (motors and handpieces). These are shown in the display as a combination of symbol and text. The size and display depend on which socket the applied part is connected to and whether the applied part is active or locked.

Identification of control elements on system components

Control elements on the system components of the ELAN 4 electro motor system are labeled with a gold marking.

Recognition of connected motor cable and applied part

The settings (upper speed limit, pump status and delivery volume) last selected with this application part type on this bushing are retrieved upon inserting the same application part type.

Overload protection

To protect the micro motors in the applied parts against damage from overheating, the motor temperature is monitored. If the temperature is too high a warning signal sounds and a thermometer symbol appears in the display **2**.

If the temperature continues to be too high, the applied part is switched off. The following message appears in the display **2**: "The motor currently being used is overheating. Please let the motor cool down or use another suitable motor".

After a cooling phase the applied part can be reactivated.

It is recommended to have a second applied part ready at hand.

Aesculap®

ELAN 4 electro control unit GA800

Coolant pump

The control unit is equipped with a coolant pump 3.

The coolant pump can be activated either in the corresponding pump control panel or via the corresponding button on the foot control. It can be started by activating the motor or via the "Flush" function (continuous flushing). The delivery rate can only be adjusted via the pump control panel.

5. Preparation and setup

Non-compliance with the following instructions will preclude all responsibility and liability in this respect on the part of Aesculap.

- ▶ When setting up and operating the product, adhere to
 - national regulations for installation and operation,
 - national regulations on fire and explosion protection.

Note

For the safety of patients and users it is essential that the mains power cord and, especially, the protective earth connection are intact. In many cases defective or missing protective earth connections are not registered immediately.

- ▶ Connect the device via the potential equalization terminal at the rear panel of the device to the potential equalization system of the room used for medical purposes.

Note

The potential equalization lead can be ordered from the manufacturer as art. no. GK535 (4 m length) or TA008205 (0.8 m length).

5.1 Set-up environment/set-up location

Fire and explosion hazard!

- ▶ Use the product outside of explosion hazard zones (e.g. areas with highly pure oxygen or anesthetic gases).

The ELAN 4 electro control unit is authorized for operation in the operating room.

Note

After being installed and put into operation, the control unit must not be transported or moved to a different set-up location.

Note

The control unit must not be placed on an Aesculap mobile stand (GA415, GA416 and GD416M).

- ▶ Make certain that the ventilation slots in the bottom of the housing and in the rear panel of the control unit are not covered, e.g. by an OR cloth.
- ▶ Make certain that the control elements, mains power switch and power socket 13 are freely accessible for the user.
- ▶ Make certain the system is set up on a sufficiently stable support (e.g. a table, ceiling support, equipment cart, etc.).
- ▶ Follow the instructions for use of the support.

5.2 Stacking of units

- ▶ Do not exceed a stack height of 415 mm.
- ▶ Place the units in a stable position.
- ▶ Aesculap devices must be stacked on top of each other facing identically.
- ▶ Never move stack devices.

6. Working with the ELAN 4 electro system

6.1 System set-up

WARNING

Risk of injury and material damage due to inadvertent activation of the applied part!

- ▶ Secure applied parts that are not actually used against inadvertent activation ("Off" position), see section on locking applied part ("Off" position).

Connecting the accessories

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

Also note that any equipment connected at the interfaces must demonstrably meet the respective IEC standards (e.g. IEC 60950 for data processing equipment, IEC/DIN EN 60601-1 for electromedical devices).

All configurations must fulfill the basic standard IEC/DIN EN 60601-1. Any individual connecting devices with one another is responsible for such configuration and must ensure compliance with basic standard IEC/DIN EN 60601-1 or applicable national standards.

- ▶ Follow the instructions for use of individual accessories.
- ▶ Please address your B. Braun/Aesculap partner or Aesculap Technical Service with any inquiries in this respect; for a contact address, see Technical Service.

Connecting the power supply

Risk of death by electric shock!

- ▶ **Connect the product only to a grounded power supply.**

Note

The mains voltage must correspond to the voltage indicated on the type plate at the back of the unit.

- ▶ Plug in the power cord at power cord socket 13.
- ▶ Connect power supply by inserting into a wall socket.

Switching on the control unit

- ▶ Activate the power ON switch 10.
The power ON indicator 9 and indicator light 5 light up.
The control unit 1 always performs a power-on self-test when the unit is switched on.
If a malfunction is detected, an error message will appear in the display 2, see system errors.

Switching off the control unit

- ▶ Activate the power OFF switch 8.
The power ON indicator 9, the indicator light 5 and the display with touch control panel 2 will fade out.

Putting out of operation

Note

The safe and all-pole disconnection of the product from the main power supply is only guaranteed when the power cord is unplugged.

- ▶ Switching off the product: Activate the power OFF switch 8.
- ▶ Unplug the power cord from the power socket 13.
The operation of the device is safely terminated.

Connecting the ELAN 4 electro foot control GA808 to the control unit

Note

The plug connection of the foot control has a yellow coded ring and a filled-in dot.

- ▶ Position connector of foot control c so that marking b on the connector is aligned with marking a on the connection socket for foot control 6, see Fig. 2.
- ▶ Plug the connector for foot control c firmly into the connection socket for foot control 6.


Fig. 2 Connecting the foot control

Legend

- a Marking for connection socket
- b Marking for connector
- c Connector for foot control

Connecting the ELAN 4 electro motor cable GA804/GA805/GA806 to the control unit

Note

The motor cable is sterile. Sterile separation occurs at the motor cable.

Note

The plug connection of the motor cable has a blue coded ring and a hollow dot.

Aesculap®

ELAN 4 electro control unit GA800

- Position connector for control unit **c** on the motor cable so that marking **b** on the connector is aligned with marking **a** on the connection socket for applied parts **7**, see Fig. 3.
- Plug connector for control unit **c** on the motor cable firmly into one of the two connection sockets for applied parts **7**.


Fig. 3 Connecting the motor cable

Legend

- a** Marking for connection socket
- b** Marking for connector
- c** Connector for control unit

Note

The motor cable must be connected to the control unit without an applied part or with the applied part locked ("Off" position).

Otherwise the applied part will not be recognized by the control unit and a message will appear in the display.

- If the applied part is not recognized;
 - Lock the applied part, see Locking the applied part ("Off" position).
 - Release the applied part again, see Releasing the applied part for operation ("On" position).

Connecting the ELAN 4 electro single-use tube set GA395SU

Note

The tube set is sterile. Sterile separation occurs at the tube set.

- Open flap **a** of the coolant pump **3**, see Fig. 4.
- Insert single-use tube set **c**:
 - Place pump tube with loop over the roller wheel **b**.
 - Push nib **e** of the single-use tube set under the pull slide **d** until the nib engages.

- Close flap of the coolant pump **a**. When doing this, make certain the pump tube does not get clamped in.
- Insert the bottle holder **16** for the sterile fluid bottle in bottle holder mount **15**.
- Insert the punching spike into the sterile fluid bottle.
- When using sterile fluid bottles made from glass: Open the ventilation flap on the punching spike.
- Hang the sterile fluid bottle from the bottle holder **16**.
- Fasten the tube to the motor cable using fixation clamps.
- Shorten the length of the tubes to correspond with the applied part being used and connect to the spray nozzle.

Fig. 4 Connecting the single-use tube set

Legend

- a** Flap of coolant pump
- b** Roller wheel
- c** Single-use tube set
- d** Pull slide
- e** Nib of the single-use tube set

Connecting the applied part to the motor cable

Note

All ELAN 4 electro motor cables (GA804/GA805/GA806) can be used in combination with all applied parts that do not have a fixed connection with a cable.

- Connect connector for motor cable **24** to connector for application part **19** of the motor cable. Ensure that the nose **23** on the motor cable is aligned with the groove on the coupling of the application part.
The applied part engages. When the visual field is "Off" **20** on the motor cable, a gold-colored marking is visible.
The control unit **1** recognizes the applied part type and displays this type in the corresponding applied part control panel in the display **2**.
The most recent settings for this applied part type at the same motor connection socket are shown on the display **2**.

Note

The applied part connected to this motor cable is only ready for operation if a gold-colored marking is visible in the visual field "On" **21** on the motor cable.

Releasing the applied part for operation ("On" position)

- Activate the release button **22** on the motor cable and push the applied part further onto the motor cable.
The applied part engages. When the visual field is "On" **21** on the motor cable, a gold-colored marking is visible.

Note

For applied parts with slider for tool release **25**, in the "On" position **21** the nib **23** on the motor cable blocks the slider. Uncoupling of the tools is therefore only possible in "Off" position **20**.

For applied parts with unlocking sleeve **26**, in the "On" position **21** the connection for applied part **19** on the motor cable blocks the unlocking sleeve. Uncoupling of the attachments is therefore only possible in "Off" position **20**.

For applied parts without slider for tool release, it is possible to couple/uncouple tools in the "On" position, but this must not be done due to the risk of injury from unintentional activation of the applied part.

Locking the applied part ("Off" position)

- Activate the release button **22** on the motor cable and detach the motor cable from the applied part.
The applied part engages. When the visual field is "Off" **20** on the motor cable, a gold-colored marking is visible.

Disconnecting the applied part from the motor cable

- Activate the release button **22** on the motor cable and pull at the connection for applied part **19** to disconnect the motor cable from the applied part.

Disconnecting the ELAN 4 electro foot control GA808 from the control unit

- Unplug connector for foot control **c** from the connection socket for foot control **6**, see Fig. 2.

Disconnecting the ELAN 4 electro motor cable GA804/GA805/GA806 from the control unit

- Unplug the motor cable on the connector for control unit **c** from the connection socket for applied parts **7**, see Fig. 3.

6.2 Function checks

- Prior to each use and after every change of applied part, check that all products to be used are in good working order.
- Check the secure connection of all products to be used.
- Ensure that the parameters are set and the unit is operated according to the instructions for use and the safety information for the applied parts or tools.
- Check that the cutting edges of the tools are not showing any mechanical damage.
- Check that the correct applied part type for the applied parts connected is showing in the corresponding applied part control panel in the display.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Release the applied part for operation.
- Activate the foot or hand control down to the limit stop.
The applied part starts running in the preselected rotational direction and reaches the maximum speed as indicated in the applied part control panel in the display.
The applied part runs quietly at a constant speed.
The qualitative bar display for the actual current applied part speed in the applied part control panel is fully illuminated.
- If necessary, release the coolant pump in the corresponding applied part control panel or by pressing the function button on the foot control.
The coolant pump starts as soon as the applied part is operated.

Aesculap®

ELAN 4 electro control unit GA800

6.3 Safe operation

Risk of injury and/or malfunction!

- Always carry out a function check prior to using the product.

Risk of injury and material damage due to inadvertent activation when repositioning/moving the foot control!

- To reposition the foot control: Use a transport brace.
- Prior to repositioning: Secure the applied part against inadvertent activation ("Off" position).

It is only possible to operate the applied part and change the setting parameters on the control unit if:

- the applied part is connected to the control unit by a motor cable,
- a foot control is connected to the connection socket for foot control 6 (when using the ELAN 4 electro motor cable GA806),
- there is no second applied part released at the same time and
- the applied part type is shown on the display of the control unit.

The following motor settings of the applied part cannot be changed when the applied part is running:

- Rotational direction
- Maximum speed

Activating the motor/pump control panel

Note

The motor settings of the applied part cannot be changed when the applied part is running.

- Changing the setting parameters of the applied part: Activate the applied part control panel **b** in the display 2, see Fig. 5.
- Changing the setting parameters of the irrigation pump: Activate the pump control panel **a** in the display 2.

The activated control panel switches to the Settings mode. Now the setting parameters listed below can be changed.


Fig. 5 Activating the control panel

Legend

- a Pump control panel
- b Applied part control panel

Changing the maximum speed setting/rotational direction of the applied part

- Activate the applied part control panel, see Activating the motor/pump control panel.
- Change the rotational direction: Activate the inactive, gray rotational direction button for clockwise/counterclockwise rotation **a/b**, see Fig. 6.
The rotational direction is changed from clockwise rotation to counterclockwise rotation and vice versa.
- Changing the maximum speed setting: Activate the buttons for reducing/increasing the maximum speed setting **c/d**.
The maximum speed setting is reduced/increased by increments.

Note

The increments for changing the maximum speed setting depend on the applied part connected.


Fig. 6 Changing the maximum speed setting/rotational direction

Legend

- a Rotational direction button for counterclockwise rotation
- b Rotational direction button for clockwise rotation
- c Button for reducing the maximum speed setting
- d Button for increasing the maximum speed setting

Activating the irrigation pump/"Flush" function (continuous flushing)

- Activate the pump control panel, see Activating the motor/pump control panel.
- Activate the irrigation pump: Activate the button for switching on the irrigation pump **a**, see Fig. 7.
The irrigation pump is active and delivers coolant at the currently set delivery rate.
- Activate the "Flush" function: Activate and hold down the "Flush" button **b**.
The "Flush" function is now active. The irrigation pump delivers coolant at the maximum delivery rate until the "Flush" button **b** is released.

Note

The activation of the "Flush" function does not depend on whether an applied part is connected or activated.


Fig. 7 Activating the irrigation pump/"Flush" function

Legend

- a Button for switching the irrigation pump on/off
- b "Flush" button
- c Button for reducing the delivery rate
- d Button for increasing the delivery rate

Deactivating the irrigation pump and changing the delivery rate

Note

The delivery rate of the irrigation pump can only be changed when the irrigation pump is activated.

- Activate the pump control panel, see Activating the motor/pump control panel.
- Deactivate the irrigation pump: Activate the button for switching off the irrigation pump **a**, see Fig. 7.
The irrigation pump is inactive and no longer delivers coolant.
- Change the delivery rate: Activate the buttons for reducing/increasing the delivery rate **c/d**.
The delivery rate of the irrigation pump is reduced/increased by increments.

The delivery rate can be adjusted in the following increments:

- 1 % to 5 %: 1 % increments
- 5 % to 100 %: 5 % increments

Activating the applied part via foot control GA808

Note

When an applied part is plugged in, clockwise rotation is always pre-set.

To activate clockwise rotation:

- Set the rotational direction button for motor **30** to clockwise rotation.
The motor direction indicator for clockwise rotation is illuminated in the applied part control panel.
- Activate the pedal **29**.
The applied part runs clockwise.

To activate counterclockwise rotation:

- Set the rotational direction button for motor **30** to counterclockwise rotation.
The motor direction indicator for counterclockwise rotation is illuminated in the applied part control panel.
- Activate the pedal **29**.
The applied part rotates counterclockwise.
The control unit emits an acoustic signal.

Activating the applied part via hand control GA804/GA805

Note

When an applied part is plugged in, clockwise rotation is always pre-set.

To activate clockwise rotation:

- Set the motor direction to clockwise rotation on the control unit.
The motor direction indicator for clockwise rotation is illuminated in the applied part control panel.
- Activate the control element of the hand control.
The applied part rotates clockwise.

To activate counterclockwise:

- Set the motor direction to counterclockwise rotation on the control unit.
The motor direction indicator for counterclockwise rotation is illuminated in the applied part control panel.
- Activate the control element of the hand control.
The applied part runs counterclockwise.
The control unit emits an acoustic signal.

Activating the coolant pump via foot control

- Switch the coolant pump on/off: Briefly activate function button **28**.
- Activate the "Flush" function (continuous flushing): Activate the function button **28** for a longer period and keep pressed down.

6.4 System settings menu

Note

The system settings menu can only be called up when no applied part is being operated.

While the system settings menu is active, the operation of the applied part is blocked.

- Call up the system settings menu: Activate the button for the system settings menu **h**, see Example of control concept.
The system settings menu opens, see Fig. 8.


Fig. 8 System settings menu

Legend

- a "Device information" button
- b "Device settings" button
- c "Motor settings" button
- d "Leave system settings" button

Menu	Description
Motor settings	Display and change the settings for the individual applied part types
Device settings	Display and change the basic settings for the control unit
Device Information	Display information about the control unit

- Call up menu: Activate the menu button.
- Leave the system settings menu: Activate the "Leave system settings" **d** button.

Motor settings

The motor settings menu displays the applied part types, see Fig. 9.


Fig. 9 Motor settings – Overview of applied part types

Legend

- a Buttons for applied part types
- b Leave the menu

- Leave the menu: Activate the "Leave menu" **b** button.
 - Display/change the motor settings of an applied part type: Activate the button for the applied part type **a**.
- The submenu for the applied part type selected opens, see Fig. 10.


Fig. 10 Motor settings – Selected applied part type

Legend

- a** Pump Flow
- b** Slow down
- c** Speed up
- d** Maximum speed
- e** Leave the menu
- f** Reduce the value
- g** Increase the value

Setting	Description
Maximum speed	Maximum speed/stroke rate
Speed up	Acceleration rate of the applied part
Slow down	Braking rate the applied part
Pump Flow	Delivery rate of the coolant pump

- Leave the menu: Activate the "Leave menu" **e** button.
 - Change the motor setting: Activate the "Increase value" button **g** or "Decrease value" button **f**.
- The change is implemented immediately. The set value is displayed.

Device settings


Fig. 11 Device settings

Legend

- a** Factory settings
- b** Language
- c** Brightness
- d** Volume
- e** Leave the menu
- f** Reduce the value
- g** Increase the value

Setting	Description
Volume	Adjust system volume
Brightness	Adjust brightness of display
Language	Set system language
Factory settings	Reset device to factory settings (includes motor settings)

- Leave the menu: Activate the "Leave menu" **e** button.
 - Change the system volume/display brightness: Activate the "Increase value" button **g** or "Decrease value" button **f**.
- The change is implemented immediately. The set value is displayed.
- Change the system language:
 - Select the "Language" button **b**.
 - Select the required language.
 - Restore the default settings:
 - Select "Reset device to factory settings" **a**.
 - Confirm the message.

Aesculap®

ELAN 4 electro control unit GA800

Device information

The "device information" menu displays general information about the device and the device software.

- Leave the menu: Activate the "Leave menu" e button.

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

The recommended chemical was used for validation.

7.2 Preparation before cleaning

- Separate the products from each other immediately after use.
- Shut down control unit 1 immediately after use, see Putting out of operation.
- Remove any visible surgical residues as much as possible with a damp, lint-free cloth.

7.3 Cleaning/disinfection

Product-specific safety instructions for the reprocessing procedure

DANGER

Risk of electric shock and fire hazard!

- Unplug the device before cleaning.
- Do not use flammable or explosive cleaning or disinfecting solutions.
- Ensure that no fluids will penetrate the product.

CAUTION

Damage to, or destruction of the product caused by mechanical cleaning/disinfection!

- Clean the product with wipe disinfection only.
- Do not sterilize the product under any circumstances.

CAUTION

Damage to the product due to inappropriate cleaning/disinfecting agents!

- Only use cleaning/disinfecting agents approved for surface cleaning. Follow the manufacturer's instructions for the respective cleaning/disinfecting agent.

- Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids.

7.4 Wipe disinfection for electrical devices without sterilization

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Wipe disinfection	RT	≥1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol

RT: Room temperature

Phase I

- ▶ Remove any visible residues with a disposable disinfectant wipe.
- ▶ Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- ▶ Observe the specified application time (1 min minimum).

7.5 Inspection, maintenance and checks

- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Set aside the product if it is damaged.

8. Maintenance

To ensure reliable operation, the product must be maintained at least once a year.

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

9. Troubleshooting list

- Have defective products repaired by Aesculap Technical Service, see Technical Service.

9.1 Display error messages

Faults that are recognized by the control unit are shown as error messages in the display.

There are three types of error messages:

- System error (text in red field): It is not possible to work with the control unit or system.
- Accessory error (text in yellow field): It is possible to work with a different component.
- Operating error (text in blue field): It is possible to work with the system after rectifying the cause of the error.

Note

Many faults cannot be clearly categorized. They can be caused by an operating error or an accessory error. In such cases, an operating error is assumed first, to avoid any unnecessary changing or sending in of new products.

System error

Message displayed	Cause	Remedy
System error Switch the control unit off and on. If the error re-occurs, replace the control unit.	The internal monitoring of the control unit recognizes an error or fault.	Switch the control unit off and on again. If the message appears again: Replace the control unit.

Accessory error

Message displayed	Cause	Remedy
Motor or motor cable could be damaged	Motor cable or applied part defective	Replace motor cable or applied part.
Malfunction of foot control Critical error in foot control module	Foot control defective	Replace foot control.
Manual control could be damaged Critical error in manual control module	Hand control defective	Replace hand control.
Prohibited combination of accessories A prohibited combination of motor and cable coding was found on socket 1 / socket 2.	A prohibited combination of application part and motor cable is connected to the control unit.	Check inserted products and replace with a permissible combination of application part and motor cable as applicable.
Unknown motor type Please verify that the software version of the device is suitable for this motor type.	Software version of the control unit is incompatible with the application part used.	Verify software version of the control unit. Have control unit software updated by Technical Service.

Operating error

Message displayed	Cause	Remedy
Motor overheated The motor currently being used is overheating. Please let the motor cool down or use another suitable motor.	Motor of the applied part has overheated	Allow the applied part to cool down. If the applied part overheats excessively: Replace the applied part.
Motor blocked Stop motor actuation and release blockage If the error re-occurs, replace the product.	Applied part blocked	Stop applied part actuation and release blockage. If the error occurs when the applied part is running idle: Replace the applied part.
Motor not recognized Switch the motor to "Off" position and then back to "On" position	Applied part connected to control unit in "On" position	Lock the applied part ("Off" position). The control unit recognizes the applied part type. In order to work: Release the applied part ("On" position).
Two motors in ON position simultaneously Please move one of the motors to the Off position.	Foot control is activated while two applied parts are released on the motor cable ("On" position) Note: Only one applied part can be operated at a time.	Only release the applied part intended to be operated ("On" position). Lock the applied part which isn't to be operated ("Off" position).
Activation in OFF position Activation of the motor in Off position Switch the motor into the On position before activation.	Foot control is activated while the applied part is locked on the motor cable ("Off" position)	Release the applied part ("On" position).
Activation of foot control without motor Connect a motor to the control unit.	Foot control is activated while the applied part is not connected to the control unit	Connect the motor cable to the control unit. Connect the applied part to the motor cable.
Activation of foot control during initialization The pedal or a foot control button was activated during initialization. Please release pedal and buttons.	Pedal or button on the foot control is activated during automatic self-test.	Release the activation. If the pedal or button is not activated, the foot control is defective. If necessary, replace foot control.
Multiple control elements of the foot control activated Multiple control elements were activated simultaneously. Please release all control elements. If this message does not disappear, please have the foot control inspected.	Multiple control elements were activated simultaneously.	Release all foot control operating elements. If this message continues to appear: foot control must undergo a technical inspection.
Extended button activation on foot control A foot control button was activated for longer than 20 seconds without interruption. Please release all buttons. If this message does not disappear, please have the foot control inspected.	A foot control operating element was activated for longer than 20 seconds without interruption.	Release all foot control operating elements. If this message continues to appear: foot control must undergo a technical inspection.
Activation of manual control without motor Connect a motor to the control unit.	Manual control was activated without an application part being connected to the control unit.	Connect application part to control unit.
Manual control activated upon connecting Please move the manual control lever/slider to an unactivated position one time.	Manual control lever/slider was activated during connection.	Move manual control level/slider to an unactivated position one time.

Message displayed	Cause	Remedy
Touchscreen timeout The touchscreen was activated for longer than 20 seconds without interruption. Please release the touchscreen. If this message continues to display, the device may be damaged.	The touchscreen was activated for longer than 20 seconds without interruption.	Release touchscreen. If this message continues to display, the control unit may be damaged. Replace the control unit.

Other faults in control unit

Malfunction	Cause	Detection	Remedy
Control unit cannot be switched on.	No power to control unit	Control unit not connected to mains power or not switched on (power ON indicator dark, black display)	Connect control unit to the main power. Switch on the control unit.
	Fuses blown	Power ON indicator does not light up, black display	Replace the fuses.
No coolant flow.	Coolant container empty	Coolant container empty	Replace coolant container.
	Tubing inserted incorrectly	Tubing inserted incorrectly	Insert tubing correctly.
	Tubing leaking	Coolant is leaking out	Replace tubing.
	Spray nozzle blocked	Coolant pump running. Coolant is not being delivered.	Replace spray nozzle.
	Coolant pump motor defective	Coolant pump not running.	Replace the control unit.

Faults when working with the applied part

Malfunction	Cause	Detection	Remedy
Disconnecting/connecting the applied part from/with the motor cable not possible (change from "On" to "Off" position and vice versa)	Release button on the motor cable is not released between connecting/disconnecting and releasing/locking the applied part	Coupling processes cannot be performed	Release the release button on the motor cable between connecting/disconnecting and releasing/locking the applied part and activate again.
	Motor coupling defective		Change the applied part or motor cable.
Tool cannot be coupled	ELAN 4 electro High speed: Automatic tool lock blocked	High speed tool doesn't engage	Activate slider for tool release and hold, then couple tool.
	ELAN 4 electro High speed: Slider for tool release locked	Applied part released for operation ("On" position)	Lock the applied part ("Off" position).
	Incompatible tool	Wrong tool	Select an appropriate tool according to the instructions for use of the applied part.
	Tool connection or coupling deformed, defective	Tool is difficult to couple or cannot be coupled	Use a new tool. Change applied part.

Malfunction	Cause	Detection	Remedy
Attachment for GA849 (craniotome) cannot be coupled or uncoupled	Automatic attachment lock blocked	Attachment doesn't engage	Pull back unlocking sleeve and hold, then couple the attachment.
	Unlocking sleeve blocked	Applied part released for operation ("On" position)	Lock the applied part ("Off" position).
	Attachment coupling defective	Attachment is difficult to couple or cannot be coupled	Use a new attachment.
	Attachment coupling is dirty		Change applied part. Clean the attachment or use a new attachment. Clean the applied part.
Pivotable dura protection GB947R is difficult to turn	Pivotable dura protection is difficult to move	Bearing dirty or worn	Follow the instructions for use (TA014438/TA014439) (reprocessing, maintenance).
			Change turnable dura guard.
Loud noise coming from applied part	Gearings/ball bearings defective	Loud, irregular noise	Change applied part. Preventive measure: Ensure that the applied part is lubricated at regular intervals.
High speed applied part vibrates excessively	Excessive running noises, vibrations	Shaft of applied part bent	Change applied part.
		Applied part defective	Follow instructions for use of the applied part (reprocessing, maintenance)
		Reprocessing carried out incorrectly	
Applied part is getting too hot	Blunt tool	Tool overheats	Change tool.
	Applied part defective	Cutting edges of the tool are sharp, but the applied part is heating up	Change applied part. Preventive measure: Ensure that the applied part is lubricated at regular intervals.
	Excessive use	Applied part heating up Loud running noise Vibration	Follow instructions for use of the applied part (interval operation).
	Reprocessing/maintenance carried out incorrectly		Follow instructions for use of the applied part (reprocessing, maintenance)
	Shaft of applied part bent		Change applied part.
	Applied part defective		
Applied part not running	Applied part defective	Tool not moving	Change applied part.
	Foot control defective	Pedal not moving	Change foot control.

Malfunction	Cause	Detection	Remedy
Insufficient power	Blunt tool	Cutting edge worn	Change tool.
		Cutting edges coated with substances, e.g. from too little rinsing	
	Applied part operated in counterclockwise rotation	Serrated tool run in counterclockwise rotation mode	Operate serrated tool in clockwise rotation.
	No central alignment of the craniotomy reamer to the dura protector bracket	Dura protector bracket bent Poor progress during craniotomy	Follow instructions for use (TA014438/TA014439). Change dura guard
	Applied part defective	Insufficient power of applied part Large amount of heat generated after a short time	Follow instructions for use of the applied part (reprocessing, maintenance) Change applied part.

9.2 Fuse replacement

DANGER

Risk of death by electric shock!

- Unplug the device before changing the fuses!

Specified fuse set: 2 units IEC 127 – T 6.3 A switching capacity H (1 500 A at 250 V/50 Hz)

- Use a small screwdriver to release the clip on the fuse holder **12**.
- Pull out the fuse holder **12**.
- Replace both fuses.
- Reinsert the fuse holder **12** in such a way that it audibly snaps into place.

Note

If the fuses burn out frequently, the device is faulty and should be repaired, see Technical Service.

10. Technical Service

DANGER

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.
- Do not modify the product.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: ats@aesculap.de

Or in the US:

Aesculap Inc.

Attn. Aesculap Technical Services

615 Lambert Pointe Drive

Hazelwood

MO, 63042 USA

Aesculap Repelectro Hotline

Phone: +1 (800) 214 -3392

Fax: +1 (314) 895 -4420

Other service addresses can be obtained from the address indicated above.

11. Accessories/Spare parts

11.1 ELAN 4 electro motor cables, applied parts and foot controls

Art. no.	Designation
GA804	ELAN 4 electro motor cable with hand lever (from software version 2.00)
GA805	ELAN 4 electro motor cable with hand control (from software version 2.00)
GA806	ELAN 4 electro motor cable for foot control
GA808	ELAN 4 electro foot control
GA822	ELAN 4 electro perforator driver
GA824	ELAN 4 electro Low-speed motor with intra coupling
GA831	ELAN 4 electro sagittal saw
GA832	ELAN 4 electro reciprocating saw
GA836	ELAN 4 electro micro sagittal saw
GA837	ELAN 4 electro micro reciprocating saw
GA849	ELAN 4 electro craniotome and multifunctional handpiece (2-ring)
GA860	ELAN 4 electro MIS hand piece
GA861	ELAN 4 electro Standard handpiece (1-ring) L4
GA862	ELAN 4 electro Standard handpiece (1-ring) L7
GA863	ELAN 4 electro Standard handpiece (1-ring) L10
GA864	ELAN 4 electro Standard handpiece (1-ring) L13
GA865	ELAN 4 electro Standard handpiece (2-ring) L7
GA866	ELAN 4 electro Standard handpiece (2-ring) L10
GA867	ELAN 4 electro Standard handpiece (2-ring) L13
GA868	ELAN 4 electro Standard handpiece (2-ring) L17
GA869	ELAN 4 electro Standard handpiece (2-ring) L22
TA014401	Instructions for use of ELAN 4 electro control unit GA800 (A4 for ring binder)

11.2 Coolant pump

Art. no.	Designation
GA395SU	ELAN 4 electro single-use tube set
GD412804	Bottle holder
–	Physiological saline solutions up to 1 000 ml Note: No Aesculap accessories

11.3 Power cord

Art. no.	Approvals	Color	Length
TE780	Europe	black	1.5 m
TE730	Europe	black	5 m
TE734	Great Britain	black	5 m
TE735	USA, Canada, Japan	gray	3.5 m

11.4 Equipotential bonding leads

Art. no.	Designation
GK535	Equipotential bonding lead (4 m)
TA008205	Equipotential bonding lead (0.8 m)

11.5 Spare parts

Art. no.	Designation
TA021473	Fuse: Fusible cartridge T 6.3 AH

12. Technical data

12.1 Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
GA800	ELAN 4 electro control unit	Ila
GA804	ELAN 4 electro motor cable with hand lever	I
GA805	ELAN 4 electro motor cable with hand control	I
GA806	ELAN 4 electro motor cable for foot control	I
GA808	ELAN 4 electro foot control	I
GA822	ELAN 4 electro perforator driver	Ila
GA824	ELAN 4 electro Low-speed motor with intra coupling	Ila
GA831	ELAN 4 electro sagittal saw	Ila
GA832	ELAN 4 electro reciprocating saw	Ila
GA836	ELAN 4 electro micro sagittal saw	Ila
GA837	ELAN 4 electro micro reciprocating saw	Ila
GA849	ELAN 4 electro craniotome and multifunctional handpiece (2-ring)	Ila
GA860	ELAN 4 electro MIS hand piece	Ila
GA861	ELAN 4 electro Standard handpiece (1-ring) L4	Ila
GA862	ELAN 4 electro Standard handpiece (1-ring) L7	Ila
GA863	ELAN 4 electro Standard handpiece (1-ring) L10	Ila
GA864	ELAN 4 electro Standard handpiece (1-ring) L13	Ila
GA865	ELAN 4 electro Standard handpiece (2-ring) L7	Ila
GA866	ELAN 4 electro Standard handpiece (2-ring) L10	Ila
GA867	ELAN 4 electro Standard handpiece (2-ring) L13	Ila
GA868	ELAN 4 electro Standard handpiece (2-ring) L17	Ila
GA869	ELAN 4 electro Standard handpiece (2-ring) L22	Ila
GA395SU	ELAN 4 electro single-use tube set	Ila

12.2 Performance data, information about standards

Protection class (acc. to IEC/DIN EN 60601-1)	I
Housing protection rating according to IEC/DIN EN 60529	IP20
Applied part	Type BF
Mains voltage ranges	100 V~ to 120 V~ ±10 % 220 V~ to 240 V~ ±10 %
Current consumption (ready for operation)	0.2 A (at 100 V~ to 120 V~) 0.3 A (at 220 V~ to 240 V~)
Current consumption (maximum load)	5.4 A to 4.4 A (at 100 V~ to 120 V~) 2.3 A to 2.2 A (at 220 V~ to 240 V~)
Frequency	50 Hz to 60 Hz
Operating mode	Continuous operation
Device fuse (acc. to IEC 60127-1)	T 6.3 AH, 250 V Form factor: 5 x 20 mm
Maximum delivery rate of the coolant pump	65 ml/min ± 15 %
Weight	9.5 kg ± 10 %
Dimensions (L x W x H)	380 mm x 330 mm x 201 mm ± 5 %
Dimensions (L x W x H) with bottle holder	380 mm x 379 mm x 427 mm ± 5 %
Conforming to standard	IEC/DIN EN 60601-1
EMC	IEC/DIN EN 60601-1-2
CISPR11	Class A

12.3 Ambient conditions

	Operation	Storage and transport
Temperature	10 °C to 40 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

13. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

Adhere to national regulations when disposing of or recycling the product, its components and its packaging!

The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number.

(The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.)

Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

14. Distributor in the US/Contact in Canada for product information and complaints

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