
(GB) Instructions for use/Technical description

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



E

Aesculap®

ELAN 4 air Standard hand piece (1-ring) GA761 to GA764 / ELAN 4 air handpiece HD (2-ring) GA765 to GA769

Legend

- 1 ELAN 4 Tool coupling (1-ring or 2-ring)
- 2 Ring marking
- 3 Connection for motor hose
- 4 Slider for tool release
- 5 ELAN 4 Tool (1-ring or 2-ring)

ELAN 4 ECCOS holder for 12 burrs GB718R

- 6 Side handles
- 7 Holding lugs
- 8 Tool coupling
- 9 Tab
- 10 Connection for applied part
- 11 Field of vision "On"
- 12 Field of vision "Off"
- 13 Release button
- 14 Nib on motor hose
- 15 Groove on connection for motor hose

Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
 YYYY-MM	Maintenance label Recommended date of next service (format: YYYY-MM) by international B. Braun-/Aesculap representation, see Technical Service
	Machine-readable, two-dimensional code The code contains a unique serial number which can be used for electronic tracking of the individual instrument. The serial number is based on the global standard sGTIN (GS1).
	Date of manufacture
	Sterilization using irradiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Manufacturer's batch designation
	Manufacturer's serial number
	Manufacturer's article number
	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage

Contents

1.	Applicable to	3
2.	Intended use	3
2.1	Task/function within system	3
2.2	Application Environment	3
2.3	Indications	3
2.4	Contraindications	3
3.	Safe handling	4
3.1	Sterile ELAN 4 tools	4
4.	Product description	4
4.1	Available sizes	4
4.2	Package includes	5
4.3	Components required for operation	5
4.4	Operating principle	5
5.	Preparation	5
6.	Working with the applied part	5
6.1	System set-up	6
6.2	Function check	6
6.3	Operation	7
7.	Validated reprocessing procedure	7
7.1	General safety notes	7
7.2	General information	7
7.3	Preparations at the place of use	8
7.4	Preparation before cleaning	8
7.5	Cleaning/disinfection	8
7.6	Manual cleaning with immersion disinfection – applied part ..	9
7.7	Mechanical cleaning/disinfecting with manual pre-cleaning –	10
	applied part	
7.8	Mechanical cleaning/disinfecting with manual pre-cleaning –	11
	tools	
7.9	Inspection, maintenance and checks	13
7.10	Packaging	13
7.11	Steam sterilization	13
7.12	Sterilization for the US market	13
7.13	Storage	13
8.	Maintenance	13
9.	Troubleshooting	14
10.	Technical Service	14
11.	Accessories/Spare parts	14
12.	Technical Data	14
12.1	Classification acc. to Directive 93/42/EEC	14
12.2	Performance Data	15
12.3	Short interval operation	15
12.4	Ambient conditions	15
13.	Disposal	15
14.	Distributor in the US/Contact in Canada for product	
	information and complaints	15

1. Applicable to

- For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

2. Intended use

2.1 Task/function within system

The hand-held highspeed ELAN 4 air Standard handpiece (1-ring) GA761 to GA764 or ELAN 4 air handpiece HD (2-ring) GA765 to GA769 is an accessory of the ELAN 4 air motor system.

The highspeed handpiece is connected to an ELAN 4 air motor hose with hand control or by an ELAN 4 air motor hose to a foot control.

The highspeed handpiece is used to drive Aesculap ELAN 4 tools (1-ring or 2-ring).

2.2 Application Environment

Application in the sterile area

2.3 Indications

Modes of application	Separating, removing and shaping hard tissue, cartilage and similar, as well as for drilling holes in bone and bone replacement materials.
Surgical discipline/areas of application	Neurosurgery, ENT surgery, OMF surgery, orthopedics and trauma surgery

Note

The type and area of application depend on the tools selected.

2.4 Contraindications

The ELAN 4 air motor system is not authorized for use in the central nervous system or central circulatory system.

Note

The clinical success of the use of the ELAN 4 air motor system is dependent on the knowledge and experience of the surgeon. The surgeon must decide which structures it is sensible to treat and take into account the safety and warning information contained in these instructions for use.

Aesculap®

ELAN 4 air Standard hand piece (1-ring) GA761 to GA764 / ELAN 4 air handpiece HD (2-ring) GA765 to GA769

3. Safe handling

CAUTION

Federal law restricts this device to sale by or on order of a physician!

WARNING

Risk of injury and material damage due to inappropriate use of the product!

- Use the product only in accordance with the intended purpose.

WARNING

Risk of injury and damage to property due to improper handling of the product!

This product is part of the ELAN 4 air motor system.

- Follow the general safety advice given in the instructions for use of ELAN 4 air foot control GA708, ELAN 4 air wall hoses, ELAN 4 air wall adapters (TA014450) and/or instructions for use of ELAN 4 air motor hose with hand control (TA014481).
- Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Prior to use, check that the product is in good working order.
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
 - Use the product only according to these instructions for use.
 - Follow the safety and maintenance instructions.
 - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Always adhere to applicable standards.
- Carry out sterile processing in any reusable applied parts, motor hoses and tools prior to use.
- When handling, observe instructions for use that are relevant to Aesculap holder systems TA009721, see Aesculap Extranet at <https://extranet.bbraun.com>
- Use the applied parts and tools within the specified pressure range only.

3.1 Sterile ELAN 4 tools

The product is gamma-sterilized and supplied in sterile packaging.

- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- read, follow, and keep the instructions for use.
- Use the product only in accordance with its intended use, see Intended use.
- Do not use products from open or damaged sterile packaging.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.
- Do not use the product after its use-by date.

4. Product description

4.1 Available sizes

Note

The shaft length can be seen on the shaft of the drill handpiece.

ELAN 4 air Standard handpiece (1-ring)

Art. no.	Length of shaft
GA761	L4 (approx. 4 cm)
GA762	L7 (approx. 7 cm)
GA763	L10 (approx. 10 cm)
GA764	L13 (approx. 13 cm)

ELAN 4 air HD handpiece (2-ring)

Art. no.	Length of shaft
GA765	L7 (approx. 7 cm)
GA766	L10 (approx. 10 cm)
GA767	L13 (approx. 13 cm)
GA768	L17 (approx. 17 cm)
GA769	L22 (approx. 22 cm)

4.2 Package includes

ELAN 4 air Standard handpiece (1-ring)

Art. no.	Designation
GA761–GA764	ELAN 4 air Standard handpiece (1 ring)
TA014478	Instructions for use of ELAN 4 air Standard handpiece (1-ring) GA761 to GA764/ELAN 4 air handpiece HD (2-ring) GA765 to GA769

ELAN 4 air HD handpiece (2-ring)

Art. no.	Designation
GA765–GA769	ELAN 4 air HD handpiece (2-ring)
TA014478	Instructions for use of ELAN 4 air Standard handpiece (1-ring) GA761 to GA764/ELAN 4 air handpiece HD (2-ring) GA765 to GA769

4.3 Components required for operation

- ELAN 4 air wall adapter GA710R, GA712R, GA713R or GA715R
- ELAN 4 air wall hose GA702R or GA703R, ELAN 4 air foot control GA708 and ELAN 4 air motor hose GA705R
 - or –
- ELAN 4 air motor hose with hand control GA706R or GA707R
- ELAN 4 Tool (1-ring or 2-ring)

4.4 Operating principle

Note

The ring marking 2 shows the shaft diameter of the tool to be used. ELAN 4 tools can be only coupled to the applied part if the ring markings (1 ring/2 ring) of the applied part and tool match.

The reamer handpiece is fitted with an ELAN 4 tool coupling (1-ring or 2-ring). It is used with ELAN 4 air foot controls/hand controls.

The drill handpiece is operated in a clockwise direction.

5. Preparation

Non-compliance with these rules will result in complete exclusion of liability on the part of Aesculap.

- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Prior to use, inspect the product and its accessories for any visible damage.
- ▶ Use the products and their accessories only if they are in perfect technical condition.

6. Working with the applied part

WARNING

Risk of infection and contamination!

Applied parts and motor hose are delivered unsterile!

- ▶ Carry out sterile processing of the applied parts and motor hose prior to use, in accordance with the instructions for use.

WARNING

Risk of injury and material damage due to inadvertent activation of the applied part!

- ▶ Secure applied parts that are not actually used against inadvertent activation (off position).

WARNING

Risk of injury and material damage due to inappropriate use of tools!

- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ When coupling/uncoupling, handle tools with cutting edges with care.

WARNING

Damage to the product if dropped!

- ▶ Use the products only if they are in perfect technical condition, see Function check

WARNING

Risk of burns to skin and tissue due to blunt tools/lack of maintenance of the applied part!

- ▶ Use tools only if they are in perfect condition.
- ▶ Replace blunt tools.
- ▶ Ensure correct maintenance of the applied part, see Maintenance.

Aesculap®

ELAN 4 air Standard hand piece (1-ring) GA761 to GA764 / ELAN 4 air handpiece HD (2-ring) GA765 to GA769

6.1 System set-up

Note

Operating elements on the system components of the ELAN 4 air motor system are labeled with a gold marking.

Connecting the accessories

The reamer handpiece is fitted with an ELAN 4 tool coupling (1-ring or 2-ring).

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

- Follow the instructions for use of individual accessories.
- Please contact your B. Braun/Aesculap Partner or Aesculap Technical Service (address: see Technical Service with any inquiries in this respect).

Connecting the applied part to the motor hose

- Attach the connection for motor hose of the applied part **3** to the connection for applied part **10** of the motor hose, see Fig. E. When doing so, ensure that the nib **14** on the motor hose is aligned with the groove on the connection for motor hose **15**.

The applied part engages. When the field of vision is "Off" **12** on the motor hose, a gold-colored marking is visible.

Note

The applied part attached to this motor hose is only ready for operation when the release button is in "On" **11** position.

Removing ELAN 4 ECCOS holder for 12 burrs GB718R from the tray

The top part of the ELAN 4 ECCOS holder for 12 burrs GB718R can be removed from the tray together with the tools.

- Pull the tab **9** and remove the top part of the holder, see Fig. C.

Coupling the tool to the applied part

WARNING

Risk of injury when coupling/uncoupling tools in the "On" position due to accidental activation of the applied part!

- Only couple/uncouple tools in the "Off" position.

- Lock the applied part ("Off" position), see Locking the applied part ("Off" position).
- Push tool shaft **5** as far as it will go into ELAN 4 tool coupling (1-ring or 2-ring) **1**, see Fig. A.
The tool engages. The tool is correctly coupled when the ring marking **2** on the tool is no longer visible.
- Pull at tool **5** to check that it is securely coupled.

Releasing the applied part for operation ("On" position)

- Activate the release button **13** on the motor hose and push the connection for applied part **10** onto the applied part, see Fig. E.

The applied part engages. When the field of vision is "On" **11** on the motor hose, a gold-colored marking is visible.

Locking the applied part ("Off" position)

- Activate the release button **13** on the motor hose and detach the motor hose from the applied part, see Fig. E.

The applied part engages. When the field of vision is "Off" **12** on the motor hose, a gold-colored marking is visible.

Uncoupling the tool from the applied part

WARNING

Risk of injury when coupling/uncoupling tools in the "On" position due to accidental activation of the applied part!

- Only couple/uncouple tools in the "Off" position.

- Lock the applied part ("Off" position), see Locking the applied part ("Off" position).
- Pull back slider for tool release **4** on the applied part and remove tool **5** from ELAN 4 tool coupling (1-ring or 2-ring) **1**, see Fig. A.

Disconnecting the applied part from the motor hose

- Lock the applied part ("Off" position), see Locking the applied part ("Off" position).
- Activate the release button **13** on the motor hose and pull at the connection for applied part **10** to disconnect the motor hose from the applied part, see Fig. E.

6.2 Function check

- Prior to each use, check that all products to be used are in good working order.
- Check the secure connection of all products to be used.
- Check that the tool is coupled properly. To do this, pull on the tool.
- Check that the cutting edges of the tools are not showing any mechanical damage.
- Briefly run the applied parts at maximum speed.
- Do not use the product if it is damaged or defective. Set aside the product if it is damaged.

6.3 Operation

WARNING

Coagulation of patient tissue or risk of burns for patients and users from hot applied part/tool!

- Cool the tool during operation.
- Put down the applied part/tool beyond reach of the patient.
- Allow the applied part/tool to cool down.
- Use a cloth to protect against burns when changing the tool.

WARNING

Risk of infection from aerosol formation!

Risk of injuries caused by particles coming loose from the tool!

- Take appropriate protective measures, e.g. watertight protective clothing, face mask and protective goggles, suction extraction.

WARNING

Risk of injury and/or malfunction!

- Always carry out a function test prior to each use of the product.

WARNING

Risk of injury due to application of the product outside the field of view!

- Only use the product under vision.

WARNING

Risk of injury and damage to the tool/system!

The rotating tool may catch cover drapes (textile or other).

- Never allow the tool to touch cover drapes (textile or other) during operation.

Operation of the applied part is only possible if:

- the applied part is released ("On" position)
- a foot control or hand control is connected.

7. Validated reprocessing procedure

7.1 General safety notes

Note

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.

Note

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note

It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.

The recommended chemistry was used for validation.

Note

If there is no final sterilization, then a virucidal disinfectant must be used.

Note

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

Residues containing chlorine or chlorides, e.g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization, will cause corrosion damage (pitting, stress corrosion) and result damage to metallic products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Aesculap®

ELAN 4 air Standard hand piece (1-ring) GA761 to GA764 / ELAN 4 air handpiece HD (2-ring) GA765 to GA769

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Material damage such as corrosion, cracks, fracturing, premature aging, or swelling.
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- ▶ For further detailed information on hygienically safe and material-preserving/value-preserving reprocessing, see www.a-k-i.org, link to Publications, Red Brochure – Proper maintenance of instruments.

7.3 Preparations at the place of use

- ▶ Separate the products from each other immediately after use.
- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

7.4 Preparation before cleaning

Applied part

- ▶ Prior to first mechanical cleaning/disinfection process: Mount ECCOS holder GB062R/GB063R in a suitable tray (e.g. JF222R).
- ▶ Place the product in the correct position in the ECCOS holder GB062R/GB063R, see Fig. B.

Tools

- ▶ Thoroughly rinse the product with running, cold water.
- ▶ Prior to first mechanical cleaning/disinfection process: Mount the bottom part of the ECCOS holder GB718R in a suitable tray (e.g. JF222R).
- ▶ Place the tools in the correct position in the tool holder 8 of the ECCOS holder GB718R, see Fig. C.
- ▶ Engage the top part of the holder in the bottom part. When doing so, ensure that the holding lugs 7 on the top part are inserted in the recesses in the bottom part.

7.5 Cleaning/disinfection

Product-specific safety guidelines on the reprocessing procedure

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Following the manufacturer's instructions, use cleaning and disinfecting agents which
 - are approved for plastic material and high-grade steel,
 - do not attack softeners (e.g., in silicone).
- ▶ Do not use cleaning agents that contain acetone.
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum temperature of 60 °C during chemical cleaning and/or disinfection.
- ▶ Do not exceed the maximum temperature of 90 °C during thermal disinfection with DI water.
- ▶ Dry the product for at least 10 minutes at a maximum of 120 °C.

Note

The indicated drying time is a guide time only. It must be checked taking into account the specific conditions (e.g. load) and if applicable adjusted.

Validated cleaning and disinfection procedure

Product	Validated procedure	Reference
Applied part	Manual cleaning with immersion disinfection	see Manual cleaning with immersion disinfection - applied part
	Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection	see Mechanical cleaning/disinfecting with manual pre-cleaning - applied part
Tools	Manual pre-cleaning with ultrasound and brush, and subsequent mechanical alkaline cleaning and thermal disinfection	see Mechanical cleaning/disinfecting with manual pre-cleaning - tools

7.6 Manual cleaning with immersion disinfection – applied part

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Pre-cleaning	RT (cold)	-	-	D-W	-
II	Cleaning	RT (cold)	>5	1	D-W	pH-neutral, pH ~ 9*
III	Intermediate rinse	RT (cold)	-	-	D-W	-
IV	Disinfection	RT (cold)	>15	2	D-W	Concentrate-free, aldehyde-free, phenol-free, QUAT-free, pH ~ 9**
V	Final rinse	RT (cold)	-	-	D-W	-
VI	Drying	RT	-	-	-	-

D-W: Drinking water

RT: Room temperature

*Recommended: BBraun Helizyme

**Recommended: BBraun Stabimed

- ▶ Follow the instructions for use of the ELAN 4 rinsing devices and rinsing adapters TA014447 and TA014448.
- ▶ Do not clean the applied part in an ultrasound bath.
- ▶ Never place the applied part in liquid without an ELAN 4 air rinsing adapter GB697R connected. To avoid the risk of corrosion/malfunctioning, allow any fluid that has entered the product to drain out immediately.
- ▶ Leave the rinsing adapter GB697R on the applied part for the entire cleaning process.

Phase I

- ▶ Clean the product under running water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.
- ▶ Mobilize non-rigid components, such as release button, unlocking sleeve, etc. during cleaning.
- ▶ Rinse the applied part through the rinsing adapter GB697R at least 3 times for 5 seconds with a water gun.

Phase II

- ▶ Using a disposable syringe, fill the inside of the applied part with cleaning solution through the rinsing adapter GB697R.
- ▶ Place the applied part, with rinsing adapter connected, fully in enzymatic cleaning solution for 5 minutes.

Phase III

- ▶ Rinse the entire product (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as release button, unlocking sleeve, etc. during rinsing.
- ▶ Rinse the applied part through the rinsing adapter GB697R at least 3 times for 5 seconds with a water gun.

Phase IV

- ▶ Prior to manual disinfection, allow rinsing water to drain off the product for a sufficient length of time and use compressed air to blow all liquid out of the product through the rinsing adapter GB697R in order to prevent dilution of the disinfecting solution.
- ▶ Using a disposable syringe, fill the inside of the applied part with disinfectant solution through the rinsing adapter GB697R.
- ▶ Place the applied part, with rinsing adapter connected, fully in disinfectant solution for 15 minutes.

Phase V

- ▶ Rinse the entire product (all accessible surfaces) under running water.
- ▶ Mobilize non-rigid components, such as release button, unlocking sleeve, etc. during cleaning.
- ▶ Rinse the applied part through the rinsing adapter GB697R at least 3 times for 5 seconds with a water gun.

Phase VI

- ▶ Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air).
- ▶ After manual cleaning/disinfection, check visible surfaces visually for residues.
- ▶ Repeat the cleaning/disinfection process if necessary.

Aesculap®

ELAN 4 air Standard hand piece (1-ring) GA761 to GA764 / ELAN 4 air handpiece HD (2-ring) GA765 to GA769

7.7 Mechanical cleaning/disinfecting with manual pre-cleaning – applied part

Note

The cleaning and disinfecting machine must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Brushes	RT (cold)	–	–	D–W	–
II	Rinsing	RT (cold)	5	–	D–W	–

D–W: Drinking water

RT: Room temperature

- ▶ Follow the instructions for use of the ELAN 4 rinsing devices and rinsing adapters TA014447 and TA014448.
- ▶ Do not clean the applied part in an ultrasound bath.
- ▶ Never place the applied part in liquid without an ELAN 4 air rinsing adapter GB697R connected. To avoid the risk of corrosion/malfunctioning, allow any fluid that has entered the product to drain out immediately.

Phase I

- ▶ Clean the product with a suitable cleaning brush until all discernible residues have been removed.

- ▶ Mobilize non-rigid components, such as release button, unlocking sleeve, etc. during cleaning.

Phase II

- ▶ Connect the connection for motor hose 3 with ELAN 4 air rinsing device GB691R.
- ▶ Rinse the applied part:
 - For 5 min with a tap/hose or
 - 3 times for 5 seconds with a water gun

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	■ Concentrate, alkaline: – pH ~ 13 – <5 % anionic surfactant ■ 0.5 % working solution – pH ~ 11*
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfection	90/194	5	FD–W	–
V	Drying	–	–	–	In accordance with the program for the cleaning and disinfecting machine

D–W: Drinking water

FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Connect the connection for motor hose **3** with ELAN 4 air rinsing device GB691R.
- ▶ Place the rinsing device on a wire basket suitable for cleaning.
- ▶ Connect the rinsing device to the rinsing connection of the rinsing cart.
- ▶ After the mechanical cleaning/disinfection process:
 - Remove any residual rinsing water from the applied part using a compressed-air gun, see instructions for use ELAN 4 rinsing device and rinsing adapter TA014447 and TA014448.
 - Check visible surfaces for residues.

7.8 Mechanical cleaning/disinfecting with manual pre-cleaning – tools

Note

The cleaning and disinfecting machine must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

Manual pre-cleaning with ultrasound and brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemicals
I	Immersion	RT (cold)	>30	50	D–W	Enzymatic detergent*
II	Rinsing	RT (cold)	–	–	D–W	–
III	Ultrasonic cleaning	55/131	>15	2	D–W	Concentrate-free, aldehyde-free, phenol-free, QUAT-free, pH ~ 9**
IV	Cleaning with brush	RT (cold)	–	–	D–W	–

D–W: Drinking water

RT: Room temperature

*Recommended: BBraun Helizyme

**Recommended: BBraun Stabimed

- ▶ Only perform ultrasound cleaning in the ECCOS holder GB718R to prevent any damage to the cutting edges.

Phase I

- ▶ Fully immerse the product in enzymatic cleaning solution for at least 30 minutes. Ensure that all accessible surfaces are moistened.

Phase II

- ▶ Rinse the entire product (all accessible surfaces) under running water.

Phase III

- ▶ Tilt up the tool holder **8** of the ECCOS holder GB718R by the side handles **6**, see Fig. C.
- ▶ Clean the product in an ultrasonic cleaning bath (frequency 35 kHz) for at least 15 minutes. Ensure that all accessible surfaces are immersed and acoustic shadows are avoided.

Phase IV

- ▶ Clean the product with a suitable cleaning brush until all discernible residues have been removed.

Aesculap®

ELAN 4 air Standard hand piece (1-ring) GA761 to GA764 / ELAN 4 air handpiece HD (2-ring) GA765 to GA769

Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	■ Concentrate, alkaline: - pH ~ 13 - <5 % anionic surfactant ■ 0.5 % working solution - pH ~ 11*
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	-	-	-	In accordance with the program for the cleaning and disinfecting machine

D-W: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

*Recommended: BBraun Helimatic Cleaner alkaline

- Place the tools in the correct position in the tool holder **8** of the ECCOS holder GB718R, see Fig. C.
- Check visible surfaces for residues after mechanical cleaning/disinfecting.

7.9 Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.
- ▶ Inspect the product after each cleaning and disinfecting cycle to be sure it is: clean, functional, and undamaged.
- ▶ Connect ELAN 4 air oil spray adapter GB600850 (purple) to connection for motor hose 3 and spray through the applied part for approx. 2 seconds with Aesculap-STERILIT oil spray GB600, see Fig. E.
- ▶ Check the product for any damage, abnormal running noise, overheating or excessive vibration.
- ▶ Inspect tools for broken, damaged or blunt edges.
- ▶ Set aside the product if it is damaged.

7.10 Packaging

- ▶ Always follow the instructions for use of the respective packaging and storage devices (e.g. instructions for use TA009721 for Aesculap-ECCOS holder system).
- ▶ Place the applied part in the correct position in the ECCOS holder GB062R/GB063R, see Fig. B.
- ▶ Place the tools in the correct position in the ECCOS holder GB718R, see Fig. C.
- ▶ Pack trays appropriately for the intended sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

7.11 Steam sterilization

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
- ▶ Validated sterilization process
 - Steam sterilization using fractional vacuum process
 - Steam sterilizer according to DIN EN 285 and validated according to DIN EN ISO 17665
 - Sterilization using fractionated vacuum process at 134 °C/holding time 5 min
- ▶ When sterilizing several instruments at the same time in a steam sterilizer: ensure that the maximum load capacity of the steam sterilizer specified by the manufacturer is not exceeded.

7.12 Sterilization for the US market

- Aesculap advises against sterilizing the device by flash sterilization or chemical sterilization.
- Sterilization may be accomplished by a standard prevacuum cycle in a steam autoclave.

To achieve a sterility assurance level of 10^{-6} , Aesculap recommends the following parameters:

Aesculap Orga Tray/Sterile container (perforated bottom) Minimum cycle parameters*			
Sterilization method	Temp.	Time	Minimum drying time
Prevacuum	270 °F/275 °F	4 min	20 min


*Aesculap has validated the above sterilization cycle and has the data on file. The validation was accomplished in an Aesculap sterile container cleared by FDA for the sterilization and storage of these products. Other sterilization cycles may also be suitable, however individuals or hospitals not using the recommended method are advised to validate any alternative method using appropriate laboratory techniques. Use an FDA cleared accessory to maintain sterility after processing, such as a wrap, pouch, etc.

7.13 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

8. Maintenance

To ensure reliable operation, the product must be maintained in accordance with the maintenance labeling or at least once a year.

e.g. 2016-07

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

ELAN 4 air Standard hand piece (1-ring) GA761 to GA764 / ELAN 4 air handpiece HD (2-ring) GA765 to GA769

9. Troubleshooting

Note

For further information, see the instructions for use of ELAN 4 air foot control GA708, ELAN 4 air wall hoses, ELAN 4 air wall adapters (TA014450) and/or instructions for use of ELAN 4 air motor hose with hand control (TA014481).

10. Technical Service

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.
- Do not modify the product.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses

Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen /Germany
Phone: +49 7461 95-1601
Fax: +49 7461 14-939
E-Mail: ats@aesculap.de

Or in the US:

Aesculap Inc.
Attn. Aesculap Technical Services
615 Lambert Pointe Drive
Hazelwood
MO, 63042
Aesculap Repair Hotline
Phone: +1 800 214-3392
Fax: +1 314 895-4420

Other service addresses can be obtained from the address indicated above.

11. Accessories/Spare parts

Art. no.	Designation
GB062R	ELAN 4 air ECCOS holder for reamer handpiece GA761 to GA769, 1 handpiece
GB063R	ELAN 4 air ECCOS holder for reamer handpiece GA761 to GA769, 3 handpieces
GB600	STERILIT Power Systems oil spray
GB600850	ELAN 4 air oil spray adapter
GB691R	ELAN 4 air rinsing device
GB697R	ELAN 4 air rinsing adapter
GB718R	ELAN 4 ECCOS holder for 12 burs
GP111R-GP243R	ELAN 4 tools (1 ring), reusable
GP301R-GP397R	ELAN 4 tools (2 ring), reusable
GP398SU GP399TC-SU	ELAN 4 tools (2 ring), single-use
TA014477	Instructions for use of ELAN 4 air Standard hand-piece (1 ring) GA761 to GA764 / ELAN 4 air hand-piece HD (2 ring) GA765 to GA769 (A4 for ring binder)
TA014478	Instructions for use of ELAN 4 air Standard hand-piece (1 ring) GA761 to GA764 / ELAN 4 air hand-piece HD (2 ring) GA765 to GA769 (leaflet)

12. Technical Data

12.1 Classification acc. to Directive 93/42/EEC

Art. no.	Designation	Class
GA761–GA764	ELAN 4 air Standard handpiece (1 ring)	Ila
GA765–GA769	ELAN 4 air handpiece HD (1 ring)	Ila
GP111R–GP243R	ELAN 4 tools (1 ring), reusable	Ila
GP301R–GP397R	ELAN 4 tools (2 ring), reusable	Ila
GP398-SU GP399TC-SU	ELAN 4 tools (2 ring), single-use	Ila

12.2 Performance Data

Max. power	approx. 80 W at 8 bar
Max. torque	approx. 3,5 Ncm at 8 bar
Max. motor speed	approx. 80 000 rpm at 8 bar
Tool connection	ELAN 4 1-ring (GA761 to GA764) ELAN 4 2-ring (GA765 to GA769)
Pressure range	8 bar $\pm$ 2 bar
Rotational direction	Clockwise rotation

Art. no.	Weight	Dimensions ($\varnothing$ x L)
GA761	68 g $\pm$ 10 %	18 mm x 125 mm $\pm$ 5 %
GA762	71 g $\pm$ 10 %	18 mm x 150 mm $\pm$ 5 %
GA763	75 g $\pm$ 10 %	18 mm x 180 mm $\pm$ 5 %
GA764	79 g $\pm$ 10 %	18 mm x 210 mm $\pm$ 5 %
GA765	76 g $\pm$ 10 %	18 mm x 150 mm $\pm$ 5 %
GA766	82 g $\pm$ 10 %	18 mm x 180 mm $\pm$ 5 %
GA767	88 g $\pm$ 10 %	18 mm x 210 mm $\pm$ 5 %
GA768	95 g $\pm$ 10 %	18 mm x 250 mm $\pm$ 5 %
GA769	105 g $\pm$ 10 %	18 mm x 300 mm $\pm$ 5 %

12.3 Short interval operation

Ambient temperature	10 °C to 27 °C
Application cycle	30 second application 30 second pause
Max. application cycles	∞

12.4 Ambient conditions

	Operation	Storage and transport
Temperature	10 °C to 27 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

13. Disposal

Note

The user institution is obliged to process the product before its disposal, see Validated reprocessing procedure.

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.
- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical Service.

14. Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA, 18034,
USA

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