- GB** **Instructions for use/Technical description**  
Modular ISOCER® hip endoprosthesis head
- USA** **Instructions for use/Technical description**  
Modular ISOCER® hip endoprosthesis head  
Note for U.S. users  
This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapimplantsystems.com](http://www.aesculapimplantsystems.com) and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
- D** **Gebrauchsanweisung/Technische Beschreibung**  
Modularer ISOCER® Hüftendoprothesenkopf
- F** **Mode d'emploi/Description technique**  
Tête d'endoprothèse de hanche ISOCER® modulaire
- E** **Instrucciones de manejo/Descripción técnica**  
Cabeza de endoprótesis de cadera modular ISOCER®
- I** **Istruzioni per l'uso/Descrizione tecnica**  
Testa di endoprotesi coxo-femorale modulare ISOCER®
- P** **Instruções de utilização/Descrição técnica**  
Esfera modular ISOCER® para endoprótese de quadril
- NL** **Gebruiksaanwijzing/Technische beschrijving**  
Modulaire ISOCER® heupprothesekop
- S** **Bruksanvisning/Teknisk beskrivning**  
Modulärt ISOCER® höftendoproteshuvud
- RUS** **Инструкция по применению/Техническое описание**  
Модульная головка тазобедренного эндопротеза ISOCER®
- CZ** **Návod k použití/Technický popis**  
Modulární hlava kyčelní endoprotézy ISOCER®
- PL** **Instrukcja użytkowania/Opis techniczny**  
Modułowa główka endoprotezy stawu biodrowego ISOCER®
- SK** **Návod na použitie/Technický opis**  
Modulárna bedrová protetická hlavica ISOCER®
- TR** **Kullanım Kılavuzu/Teknik açıklama**  
Modüler ISOCER® kalça protezi başı
- KR** **사용 설명서 / 기술 설명**  
모듈식 ISOCER® 둔부 보철물 헤드

**B | BRAUN**  
SHARING EXPERTISE

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

TA-Nr. 014606 2018-06 V6 Änd.-Nr. 58293

**CE** 0123 - DIR 93/42/EEC

## Aesculap® Modular ISOCER® hip endoprosthesis head

### Intended use

The implant is used

- as a component of a human hip endoprosthesis: Prosthesis head
- to combine with Aesculap hip endoprosthesis components
- as a bearing for combination with Aesculap cup components made from UHMWPE and Vitelene®
- not for combination with TrendHip®, HistorHip® and cemented SLA hip endoprosthesis stems

### Materials

The materials used in the implant are listed on the packaging:

- ISOCER® Ceramic Al2O3/ZrO2 according to ISO 6474-2
- ISOCER® is a registered trademark of Aesculap AG, 78532 Tuttlingen / Germany.

### Indications

Use for acute disorders of the hip that cannot be treated by other therapies:

- Degenerative osteoarthritis
- Rheumatic arthritis
- Fractures of the joint
- Femoral head necrosis

### Contraindications

#### Absolute contraindications

Do not use in the presence of:

- Joint diseases that can be treated by reconstructive surgery (e.g. displacement osteotomy)
- Acute or chronic infections near the joint, or systemic infections
- Secondary diseases that could influence joint implant functionality
- Severe osteoporosis or osteomalacia
- Insufficient quality of the bone or osseous deformities, disorders in the implant fixation area which could impair the stability of the fixation of the joint replacement primarily or over time
- Known hypersensitivity to the implant materials
- All areas of application not mentioned under 'Indications'

#### Relative contraindications

The following conditions, individual or combined, can lead to delayed healing or compromise the success of the operation:

- Anticipated excessive load on the joint implant
- Pronounced osseous deformities
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Expected lack of patient cooperation

Use of the product in these cases requires the individual, critical assessment of the surgeon.

The revision of an endoprosthetic implant is also a difficult procedure with individual requirements that must be appropriately assessed by the surgeon.

### Side effects and interactions

Currently known possible risks relating to the use of the product, including resulting clinical consequences, include:

- Changes in position, loosening, wear and tear of, or fracture of implant components
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Thromboses, embolisms
- Tissue reaction to implant materials
- Injury to the surrounding tissue, including nerve and vascular damage
- Hematomas and wound healing disorders
- Periarticular calcification
- Reduced joint mobility and flexibility
- Restricted load-bearing capacity of the joint and joint pain

### Safety notes

- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed properly.
- General risk factors associated with surgical procedures are not described in this documentation.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- The operating surgeon is responsible for combining and implanting the implant components, with or without bone cement.
- Aesculap is not responsible for any complications arising from erroneous indication, wrong choice of implant, incorrect combination of implant components and operating technique, the limitations of the treatment method, or inadequate asepsis.
- The user instructions for individual Aesculap implant components must be followed.
- The implant components were tested and approved in combination with Aesculap components. If other combinations are used, the responsibility for such action lies with the operating surgeon.
- Do not, under any circumstances, combine implant components from different manufacturers.
- Do not, under any circumstances, use damaged or surgically excised components.
- Implants that have been used before must not be reused.
- Damage to load-bearing bone structures can result in loosening of the components, bone or implant fractures or other acute complications.
- To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.
- Combine modular implant components only with suitable Aesculap hip endoprostheses.
- Observe material, friction coupling diameter and cone specifications.
- Observe further restrictions for combined implants.
- Avoid damage to the implant, especially in the neck or cone region, caused by instruments (e.g. HF surgical devices) applied close to the implant.

**The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.**

- The implant components applied, along with their article numbers, the name of the implant, as well as the batch number and serial number (if available) must be documented in all patient records.
- Postoperatively, individual patient information, as well as mobility and muscle training, is of particular importance.

CAUTION

**Risk of head breakage and increased wear and tear due to a combination with implant components from other manufacturers or in combination with ceramic inserts!**

- ▶ Use Aesculap implant components only.
- ▶ Use only with Aesculap polyethylene cups and inserts.

### Safety advisory for ceramic implant components

Any load on the joint continued after a breakage of a ceramic implant in vivo causes damage to all implants, the socket articulation surface and the prosthesis cone.

As soon as a breakage event is noticed or suspected, the patient must be examined and treated at a specialist endoprosthetic clinic without delay.

- ▶ Prior to the operation, inform the patient that breakage of a ceramic implant component can occur in rare cases and that the sudden breakage of a ceramic implant would usually be noticed by the patient.
- ▶ Prior to the operation, inform the patient that no load must be put on the affected hip joint anymore as soon as a breakage of a ceramic implant component is suspected, in order to prevent further destruction of endoprosthetic components.

### Sterility

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components have been sterilized by irradiation.
- ▶ Store implant components in their original packaging. Remove them from their original protective packaging only just prior to application.
- ▶ Prior to use, check the product expiry date and verify the integrity of the sterile packaging.
- ▶ Do not use implant components that are past their expiration date or whose packaging is damaged.

WARNING

**Damage to implants caused by processing and reesterilization!**

- ▶ Do not reprocess or reesterilize the implants.

### Application

The operating surgeon shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

The following conditions must be fulfilled prior to application:

- All requisite implant components are ready to hand.
- All requisite implantation instruments must be available and in working order, including specialized Aesculap implantation systems.
- The operating surgeon and operating room team are thoroughly conversant with the operating technique and with the available range of implants and instruments; information materials on these subjects must be complete and ready to hand.
- The operating surgeon is fully conversant with the rules governing medical practice, the current state of scientific knowledge, and the contents of relevant scientific articles by medical authors.
- The manufacturer has been consulted if the preoperative situation was unclear and if implants were found in the area operated on.

The surgical procedure and following information has been explained to the patient, and the patient's consent has been documented:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic joint can bring about only limited improvement in the patient's condition vis-à-vis their condition prior to the operation.
- The prosthetic joint can come loose due to excessive load, wear and tear, or infection.
- The lifespan of the prosthetic joint depends on the patient's body weight and the strain to which the joint is subjected.
- The prosthetic joint must not be overstrained by extreme loads, hard physical labor or sporting activities.
- Corrective surgery may become necessary if the implant loosens.
- In the event that corrective surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
- The revision of a hip endoprosthesis stem is a complex joint replacement operation.
- The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.

Steps of application and intraoperative handling:

- Modular prosthesis heads are positioned on the conical connecting piece of Aesculap-hip endoprosthesis stems.
- The nominal diameter, cone size and femoral neck length are noted down explicitly on the packages and the implants.
- ▶ Make absolutely sure the inside cone size of the hip endoprosthesis head matches the outside cone size of the Aesculap hip endoprosthesis shaft.
  - The Aesculap prosthesis shaft and prosthesis head cones, respectively, are marked with the designations 12/14 or 8/10.
- ▶ Verify that the nominal diameter of the prosthesis head matches the nominal diameter of the acetabular implant.
- ▶ When implanting hip endoprosthesis heads, always take care that:
  - the protective cap on the prosthesis cone is not removed until immediately before putting the prosthesis head into position
  - both the outside cone of the prosthesis stem and the inside cone of the prosthesis head are clean and dry prior to installing the prosthesis head
  - all foreign bodies, such as tissue fragments and bone or cement particles, are removed
  - the prosthesis head is fixed with a slight turn under axial pressure and with a plastic (not metal) head punch
  - the modular hip end prosthesis heads are combined with Aesculap implant components only Aesculap-implant components only
  - all exposed bone cement and bone residues are removed prior to closing the wound, since bone cement particles as well as bone fragment caught between the sliding surfaces can cause abnormal wear of the artificial joint.
- ▶ When using ceramic ISOCER® prosthesis heads, always take care that:
  - under no circumstances a damaged prosthesis cone is used
  - prior to installing the prosthesis head, the prosthesis cone is carefully inspected for damage
  - visibly and, possibly, invisibly damaged ISOCER® prosthesis heads are not reused under any circumstances. For instance, an ISOCER® prosthesis head that has dropped to the floor, or one that has already been punched onto the prosthesis shaft cone and has been removed again, must not be reused. Unintentional reuse must also be excluded.
  - a surgically removed ISOCER® prosthesis head (e.g. replacement of incorrect neck length, or in revision cases) is only replaced by a prosthesis head made of metal or by a prosthesis head from the BILOX® option system. In such cases the prosthesis cone must be undamaged, or a revision of the prosthesis stem must be performed. A damaged prosthesis shaft cone presents the risk that the new prosthesis head is not seated firmly. This can lead to metal attrition between the prosthesis cone and the new prosthesis head. If this occurs, replace the prosthesis stem.

\*BILOX® is a registered trademark of CeramTec GmbH, 73207 Plochingen / Germany.

## Procedure in case of revision

In case of a revision, perform the following steps:

- ▶ Remove any ceramic splinters carefully and completely.
- ▶ If possible, clean the site with a pulsed irrigation/suction system (jet lavage).
- ▶ The polyethylene or ceramic friction partner must be replaced even if the cup is jammed.

### Note

Follow the relevant instructions for use when replacing the components!

WARNING

**Risk of wear and tear due to ceramic fragments inside the joint following revision surgery that involved breaking ceramic components!**

- ▶ Do not implant a metal head (risk of three-body wear).

- ▶ Exchange prosthesis shaft or treat joint with a modular BIOLOX®-option system-hip endoprosthesis head.
  - ▶ Replace the inlay and the jammed PE cup.
  - ▶ Ceramic-ceramic articulation is preferable to ceramic-PE.
- Further information on Aesculap implant systems is available at any time from B. Braun/Aesculap or from the responsible B. Braun/Aesculap-branch.

## Disposal

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

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