


TABLES OF DIMENSIONS STERILE CONTAINERS

TABLE OF DIMENSIONS FOR STERILE CONTAINERS SIZE 1/1
WITH AND WITHOUT EDS²

Item no.	Lid - outer dimensions (LxWxH mm) 	Sterile container bottom - outer dimensions (LxWxH mm) 	Sterile container bottom - inner dimensions [LxWxH mm] (loading height) 	Sterile container with lid - outer dimensions (LxWxH mm) 
Lid				
JJ410	598 x 297 x 28			
Sterile container bottom				
JJ110/ JJ111*		598 x 297 x 86	548 x 258 x 64	598 x 297 x 100
JJ120/ JJ121*		598 x 297 x 136	548 x 258 x 114	598 x 297 x 150
JJ130/ JJ131*		598 x 297 x 186	548 x 258 x 164	598 x 297 x 200
JJ140/ JJ141*		598 x 297 x 236	548 x 258 x 214	598 x 297 x 250

TABLE OF DIMENSIONS FOR STERILE CONTAINERS SIZE 1/2
WITH AND WITHOUT EDS²


Item no.	Lid - outer dimensions (LxWxH mm) 	Sterile container bottom - outer dimensions (LxWxH mm) 	Sterile container bottom - inner dimensions [LxWxH mm] (loading height) 	Sterile container with lid - outer dimensions (LxWxH mm) 
Lid				
JJ420	300 x 297 x 28			
Sterile container bottom				
JJ210/ JJ211*		300 x 297 x 86	250 x 258 x 64	300 x 297 x 100
JJ220/ JJ221*		300 x 297 x 136	250 x 258 x 114	300 x 297 x 150
JJ230/ JJ231*		300 x 297 x 186	250 x 258 x 164	300 x 297 x 200
JJ240/ JJ241*		300 x 297 x 236	250 x 258 x 214	300 x 297 x 250

² EDS = Enhanced Drying System

* Container bottoms with EDS (Enhanced Drying System)

TABLE OF DIMENSIONS FOR STERILE CONTAINERS SIZE 3/4
WITH AND WITHOUT EDS²

Item no.	Lid - outer dimensions (L x W x H mm) 	Sterile container bottom - outer dimensions (L x W x H mm) 	Sterile container bottom - inner dimensions [L x W x H mm] (loading height) 	Sterile container with lid - outer dimensions (L x W x H mm) 
Lid				
JJ430	464 x 297 x 28			
Sterile container bottom				
JJ310 / JJ311*		464 x 297 x 86	414 x 258 x 64	464 x 297 x 100
JJ320 / JJ321*		464 x 297 x 136	414 x 258 x 114	464 x 297 x 150
JJ330 / JJ331*		464 x 297 x 186	414 x 258 x 164	464 x 297 x 200
JJ340 / JJ341*		464 x 297 x 236	414 x 258 x 214	464 x 297 x 250



See catalog C40402 for additional sterile containers and baskets

Scan the QR code or type in the link below.

LINK: <http://bbraun.info/ij>

NOTES:

² EDS = Enhanced Drying System

* Container bottoms with EDS (Enhanced Drying System)

AESCU LAP Aicon®

TABLES OF DIMENSIONS BASKETS

TABLE OF DIMENSIONS FOR
BASKETS SIZE 1/1

Item no.	Basket lid outer dimensions (LxWxH mm)	Basket outer dimensions (LxWxH mm)	Basket inner dimensions (LxWxH mm)
			
Basket lid			
JJ167R	528 x 253 x 18		
Basket			
JJ161R		537 x 250 x 35	535 x 248 x 32
JJ162R		537 x 250 x 55	535 x 248 x 52
JJ163R		537 x 250 x 80	535 x 248 x 77
JJ164R		537 x 250 x 105	535 x 248 x 102

TABLE OF DIMENSIONS FOR
BASKETS SIZE 1/2

Item no.	Basket lid outer dimensions (LxWxH mm)	Basket outer dimensions (LxWxH mm)	Basket inner dimensions (LxWxH mm)
			
Basket lid			
JJ267R	231 x 253 x 18		
Basket			
JJ261R		240 x 250 x 35	238 x 248 x 32
JJ262R		240 x 250 x 55	238 x 248 x 52
JJ263R		240 x 250 x 80	238 x 248 x 77
JJ264R		240 x 250 x 105	238 x 248 x 102

TABLE OF DIMENSIONS FOR
BASKETS SIZE 3/4

Item no.	Basket lid outer dimensions (LxWxH mm)	Basket outer dimensions (LxWxH mm)	Basket inner dimensions (LxWxH mm)
			
Basket lid			
JJ367R	394 x 253 x 18		
Basket			
JJ361R		403 x 250 x 35	401 x 248 x 32
JJ362R		403 x 250 x 55	401 x 248 x 52
JJ363R		403 x 250 x 80	401 x 248 x 77
JJ364R		403 x 250 x 105	401 x 248 x 102

TABLE OF DIMENSIONS FOR
BASKETS SIZE DIN

Item no.	Basket lid outer dimensions (LxWxH mm)	Basket outer dimensions (LxWxH mm)	Basket inner dimensions (LxWxH mm)
			
Basket lid			
JJ467R	476 x 253 x 18		
Basket			
JJ461R		485 x 250 x 35	483 x 248 x 32
JJ462R		485 x 250 x 55	483 x 248 x 52
JJ463R		485 x 250 x 80	483 x 248 x 77
JJ464R		485 x 250 x 105	483 x 248 x 102