























Aesculap Orthopaedics

e.motion® PS/Revision

Surgical Procedure Operationsablauf

	Procedure / Operationsschritt		
1	Option 1: 	Option 2: 	
Tibia Resection / Tibiaresektion			
2			
Tibia Preparation / Tibiapräparation			
3			
Positioning of the tibia trial implants / Positionierung der Tibia-Probe- komponenten			
4			
Measurement of the extension gap / flexion gab / Messung des Streckspaltes / des Beugespaltes			
5			
Femoral Resection / Femurresektion			
6			
Patella Preparation / Präparation der Patella			

	Procedure / Operationsschritt	
7		
Trial Reposition / Probereposition		
8		
Femur Box Preparation / Femurbox- präparation		
9		
Femur Stem Preparation / Femurschaft präparation		
10		
Final Trial Reposition / Endgültige Probereposition		
11		
Assembling the final implants / Montage der endgültigen Implantate		
12		
Final Implantation / Endgültige Implantation		