

Monosyn® Quick

Surgeon's quick choice where a reliable suture is needed



Monosyn® Quick

Recommended codes

Codes	Description
C0025047	Monosyn® Quick 2/0 70 cm FR26
C0025052	Monosyn® Quick 4/0 70 cm HR17
C0025058	Monosyn® Quick 3/0 70 cm HR26
C0025059	Monosyn® Quick 2/0 70 cm HR26
C0025101	Monosyn® Quick 4/0 45 cm HSMP15
C0025024	Monosyn® Quick 5/0 45 cm DSMP13
C0025043	Monosyn® Quick 3/0 45 cm DSMP24

Try it and send your comments to:
monofilaments@bbraun.com

C Box of 36 units



Urology & Paediatric Surgery

B. Braun Surgical SA | Carretera de Terrassa, 121 | 08191 Rubí | Spain
Phone +34 93 5 86 62 00 | Fax +34 93 6 99 73 03 | www.bbraun.es

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 74 61 95-0 | Fax +49 74 61 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

B | BRAUN
SHARING EXPERTISE

Monosyn® Quick

The short-term absorption profile
monofilament suture of B. Braun

...where a reliable suture is
needed



Monosyn® Quick provides:

- Short-term absorption profile for fast healing tissues
- Superior handling and flexibility to work in deep areas
- Monofilament structure avoids the stone formation and the spread of infection over the suture
- Comfortable and a quick recovery

Monosyn® Quick, Surgeon's quick choice ...

Main indications:

- Perineal closure
- Phimosis
- Hypospadias

Product features

Structure	Monofilament
Chemical composition	Glyconate
Size	USP 1 - USP 6/0
Knot tensile strength retention	50 % 6 - 7 days post implantation 0 % 14 - 21 days post implantation
Mass Absorption	Degradation by hydrolysis after 56 days
Sterilization	Ethylene Oxide (EO)