


# B. Braun Meshes

## Optilene® Mesh LP

Mesh for hernia repair



Super-light and thin-walled mesh  
with full transparency

# Optilene® Mesh LP

Maximum patient comfort – Lightweight Polypropylene Mesh

super-light



With its weight of 36 g/m<sup>2</sup> Optilene Mesh LP is the alternative lightweight mesh for all open hernia repair procedures.

## **Super-light Polypropylene**

for surgeons who prefer to implant less foreign material

## **Soft and pliable**

for maximum patient comfort

## **Full transparency**

to see the tissue below

## **Thinner, more conformable**

than traditional polypropylene meshes

## **Flexible, strong and secure**


mesh offers several times the strength of the human abdominal wall

## **Ideal for inguinal hernia repair**

performed in open mesh techniques

## **Easy placement**

without the need of absorbable material to improve handling



# maximum comfort

The softer, microfilament construction minimizes inflammation and pain.

It helps to improve the patient's comfort.

Due to its excellent flexibility the mesh adapts to the movements of the patient resulting from body motion.




Product range

Article No.	Dimensions			Contents
1064735	5	x	10 cm	5 pieces
1064715	7.5	x	15 cm	5 pieces
1064725	10	x	15 cm	5 pieces
1064705	15	x	15 cm	5 pieces

# Optilene® Mesh LP

Lighter, soft and pliable, thinner, flexible and strong

## Weight comparison



With the lightweight polypropylene **Optilene® Mesh LP** the material implanted is considerably reduced in comparison to traditional meshes.

## Product features

Material	Monofilament Polypropylene
Construction	Knitted
Weight	36 g/m²
Pore size	1 mm
Thickness	0.39 mm
Indications	Inguinal Hernia Incisional Hernia Reconstruction of chest wall
Contraindications	Optilene® Mesh LP should not be implanted in contaminated and infected areas as well as in children during the growth phase
Sterilisation	Ethylene oxide



AESCULAP®

**B | BRAUN**  
SHARING EXPERTISE

**Aesculap AG & Co. KG**

Am Aesculap-Platz  
78532 Tuttlingen  
Telefon +49 7461 95-0  
Fax +49 7461 95-2600  
[www.aesculap.de](http://www.aesculap.de)

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.