

Aesculap Sterile Technology

PrimeLine Clear Lid



Think Innovation. Think Aesculap.

To keep things in view.

PrimeLine Clear Lid made
of transparent plastic

- impact-resistant,
scratch-proof and
thermostable
- quick identification
of instrument trays
- colour change of the
indicator visible from
the outside



 Standard 1/1	Color	 Standard 3/4	Color	 Standard 1/2	Color	 Flat 1/2 (JK339)	Color
JP001	red	JP011	red	JP021	red	JP031	red
JP002	blue	JP012	blue	JP022	blue	JP032	blue
JP003	green	JP013	green	JP023	green	JP033	green
JP004	ocher	JP014	ocher	JP024	ocher	JP034	ocher
JP005	gray	JP015	gray	JP025	gray	JP035	gray
JP006	anthracite	JP016	anthracite	JP026	anthracite	JP036	anthracite
JP007	transparent	JP017	transparent	JP027	transparent	JP037	transparent



JP007
592 x 285 x 46 mm
fitting for the adjoining container bodies



Container bodies without base perforation = JK

JK	External dimensions Length x Width x Height (mm)
JK440	592 x 274 x 90
JK441	592 x 274 x 120
JK442	592 x 274 x 135
JK444	592 x 274 x 187
JK446	592 x 274 x 247



JP017
463 x 285 x 46 mm
fitting for the adjoining container bodies



Container bodies without base perforation = JK

JK	External dimensions Length x Width x Height (mm)
JK740	470 x 274 x 90
JK741	470 x 274 x 120
JK742	470 x 274 x 135



JP027
298 x 285 x 46 mm
fitting for JK340 – JK346



Container bodies without base perforation = JK

JK	External dimensions Length x Width x Height (mm)
JK339	300 x 274 x 57
without carrying handle	
JK340	300 x 274 x 90
JK341	300 x 274 x 120
JK342	300 x 274 x 135
JK344	300 x 274 x 187
JK346	300 x 274 x 247



JP037
298 x 285 x 46 mm
fitting for JK339



AESCULAP®

B | BRAUN
SHARING EXPERTISE

Aesculap AG & Co. KG

Am Aesculap-Platz
78532 Tuttlingen
Germany

Phone +49 7461 95-0
Fax +49 7461 95-2600

www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.