



AESCULAP® PRIME*Line*® PRO

THE AESCULAP® STERILE CONTAINER SYSTEM

THE AESCULAP® STERILE CONTAINER SYSTEM



The right INSTRUMENTS

At the right TIME

In sterile CONDITION

At the right PLACE OF USE

At a high level of PROFITABILITY

In a reliable PACKAGING SYSTEM





Fast, reliable and
efficiently organized

THE AESCULAP® STERILE CONTAINER SYSTEM

- Integrated permanent germ barrier system made of PTFE designed for 5,000 sterilisation cycles. No filter changes required - saves time and processing costs
- Stable, robust aluminium lid - reliable and long-lasting
- Outstanding cleaning results
- Good drying properties
- Available in five colours (red, blue, green, gold, silver) to help you organize your instrument sets





Article number	Size / colour	Dimensions mm (length x width x height)
JP101	1/1 Lid, red	593 x 294 x 37
JP102	1/1 Lid, blue	593 x 294 x 37
JP103	1/1 Lid, green	593 x 294 x 37
JP104	1/1 Lid, gold	593 x 294 x 37
JP105	1/1 Lid, silver	593 x 294 x 37
JP111	3/4 Lid, red	468 x 294 x 37
JP112	3/4 Lid, blue	468 x 294 x 37
JP113	3/4 Lid, green	468 x 294 x 37
JP114	3/4 Lid, gold	468 x 294 x 37
JP115	3/4 Lid, silver	468 x 294 x 37
JP121	1/2 Lid, red	303 x 294 x 37
JP122	1/2 Lid, blue	303 x 294 x 37
JP123	1/2 Lid, green	303 x 294 x 37
JP124	1/2 Lid, gold	303 x 294 x 37
JP125	1/2 Lid, silver	303 x 294 x 37
JP131	1/2 Lid for JK339, red	303 x 294 x 37
JP132	1/2 Lid for JK339, blue	303 x 294 x 37
JP133	1/2 Lid for JK339, green	303 x 294 x 37
JP134	1/2 Lid for JK339, gold	303 x 294 x 37
JP135	1/2 Lid for JK339, silver	303 x 294 x 37

Compatible with all bottom parts (1/1, 1/2, 3/4) manufactured since 1988
Consumables and accessories can be found in catalogue no. C40402



Fast, reliable and
efficiently organized
**THE AESCULAP®
STERILE CONTAINER
SYSTEM**

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "Aesculap" and the product trademarks "PrimeLine" are registered trademarks of Aesculap AG. Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.