

THE WHOLE PORTFOLIO AT YOUR FINGERTIPS


B. Braun Barcode App

Design Features	Clinical Benefit
Access the whole B. Braun Portfolio	<ul style="list-style-type: none">App is not limited to Ennovate® products.Find the correct and up-to-date product information immediately.Access to details on material and reprocessing, as well as related products.
Independent setup from client internal system	<ul style="list-style-type: none">Quick and easy to implement.No interference with client internal systems.Client benefits from B. Braun IT infrastructure.



BETTER TOGETHER

Ennovate®
SOLUTIONS BEYOND FUSION



Ennovate®
SOLUTIONS BEYOND FUSION

B. BRAUN BARCODE APP

Ennovate® PROCESS MANAGEMENT SUPPORT

AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "Aesculap" and the product trademarks "Ennovate", "PentaCore" and "PolyLock" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

THERE IS AN APP FOR EVERYTHING

B. Braun Barcode App

“Solutions beyond fusion” is what Ennovate® stands for and this means offering client support beyond the scope of surgery itself.

Based on experiences made at our client organizations we have developed an intelligent, informative and intuitive barcode app supporting our clients' processes according to their needs.

Currently many of our client organizations apply individual process solutions, sometimes based on handwritten notes, to tackle the reordering of inventory or commissioning repair orders. These processes often lack standardization and open up unnecessary error sources as well as additional costs once mistakes are made.

Aiming to resolve these error sources the barcode app has been developed offering quick access based on the B. Braun IT infrastructure, helping to standardize basic processes.


Additional services will be added to the application in near future, allowing direct reordering of implants and commissioning of repairs by a tap of a finger. This will increase the safety of the reordering process and allow further inventory optimization.

How does the app add value?	
OR	CSSD
■ Process support	■ Quick access to reprocessing guidelines
■ Quick access to product and system information	■ Material information

HOW TO USE IT

B. Braun Barcode App

1. Download the “B. Braun Barcode App” from your BBraun App Catalog
2. Open the App
3. Tap “Barcode Scan”



HOW TO USE IT

B. Braun Barcode App

4. Place camera over barcode or data matrix code to scan and access further product information. The app will scan automatically.
5. Tap “Product Detail Information” to access BBraun Product Centre
6. Enter “Product Centre” to find information on material, reprocessing and related products

