

# B. Braun Vascular Systems

## Angiodyn<sup>®</sup>

Diagnostic Catheters for Angiography



# Angiodyn® – Torque

# Torque



## High performance design

- Maximal control
- Precise placement
- Optimal stability



## Ergonomic hub

- Safe connection
- Facilitated handling
- Easy identification of French sizes

4 French **Red**

5 French **Grey**

6 French **Green**

## Symmetrically braided shaft

- Superior torque
- 1:1 force transmission

# Angiodyn<sup>®</sup> – Push

# Push

## Non-braided tip

- Gentle navigation due to ultimate tip flexibility
- Atraumatic rounded distal tip design
- Soft tip on all selective tip configurations

## Unique material composition

- Excellent push
- Optimal kink resistance

## Technical Data

### 4 French

|                           | Selective       | Pigtail         |
|---------------------------|-----------------|-----------------|
| Outer Diameter Body & Tip | 0.055"/ 1.40 mm | 0.055"/ 1.40 mm |
| Inner Diameter Body       | 0.044"/ 1.12 mm | 0.044"/ 1.12 mm |
| Inner Diameter Tip        | 0.042"/ 1.07 mm | 0.042"/ 1.07 mm |

### 5 French

|                           | Selective       | Pigtail         |
|---------------------------|-----------------|-----------------|
| Outer Diameter Body & Tip | 0.065"/ 1.65 mm | 0.065"/ 1.65 mm |
| Inner Diameter Body       | 0.050"/ 1.27 mm | 0.050"/ 1.27 mm |
| Inner Diameter Tip        | 0.047"/ 1.18 mm | 0.042"/ 1.07 mm |

### 6 French

|                           | Selective       | Pigtail         |
|---------------------------|-----------------|-----------------|
| Outer Diameter Body & Tip | 0.079"/ 2.00 mm | 0.079"/ 2.00 mm |
| Inner Diameter Body       | 0.057"/ 1.44 mm | 0.057"/ 1.44 mm |
| Inner Diameter Tip        | 0.057"/ 1.44 mm | 0.042"/ 1.07 mm |

# Angiodyn® – Order information


## 4F – 5F – 6F

5 units per box. Minimum order quantity 20 units.


|                                                                       |             | 4 French     | 5 French     | 6 French     |
|-----------------------------------------------------------------------|-------------|--------------|--------------|--------------|
| Shape description                                                     | Length (cm) | Product code | Product code | Product code |
| <b>Angiodyn Pigtail – catheter for ventriculography – 8 sideholes</b> |             |              |              |              |
| PIG                                                                   | 110         | 5011573      | 5011536      | 5011521      |
| PIG                                                                   | 130         | 5011592      | 5011527      | 5011526      |
| PIG 145                                                               | 110         | 5011595      | 5011572      | 5011525      |
| PIG 145                                                               | 130         | 5011577      | 5011503      | 5011569      |
| PIG 155                                                               | 110         | 5011534      | 5011541      | 5011528      |
| PIG 155                                                               | 130         | 5011532      | 5011557      | 5011502      |
| <b>Angiodyn Judkins – catheter for the left coronary artery</b>       |             |              |              |              |
| JL3,5                                                                 | 100         | 5011594      | 5011512      | 5011511      |
| JL4                                                                   | 100         | 5011563      | 5011545      | 5011517      |
| JL5                                                                   | 100         | 5011514      | 5011586      | 5011518      |
| <b>Angiodyn Judkins – catheter for the right coronary artery</b>      |             |              |              |              |
| JR3,5                                                                 | 100         | 5011588      | 5011535      | 5011589      |
| JR4                                                                   | 100         | 5011500      | 5011546      | 5011508      |
| JR5                                                                   | 100         | 5011576      | 5011558      | 5011509      |
| <b>Angiodyn Amplatz – catheter for the left coronary artery</b>       |             |              |              |              |
| AL1                                                                   | 100         | 5011552      | 5011582      | 5011548      |
| AL2                                                                   | 100         | 5011599      | 5011565      | 5011550      |
| AL3                                                                   | 100         | 5011506      | 5011520      | 5011553      |
| <b>Angiodyn Amplatz – catheter for the right coronary artery</b>      |             |              |              |              |
| AR 1 MOD                                                              | 100         | 5011513      | 5011551      | 5011530      |
| AR 2 MOD                                                              | 100         | 5011537      | 5011516      | 5011542      |
| AR MOD                                                                | 100         | 5011522      | 5011529      | 5011547      |
| <b>Angiodyn Multi Purpose catheter</b>                                |             |              |              |              |
| MPA                                                                   | 100         | 5011501      | 5011561      | 5011564      |
| MPA 2 – 2 sideholes                                                   | 100         | 5011580      | 5011540      | 5011556      |
| <b>Angiodyn Mammaria Interna catheter</b>                             |             |              |              |              |
| IM                                                                    | 100         | 5011533      | 5011505      | 5011562      |
| <b>Angiodyn Sones catheter</b>                                        |             |              |              |              |
| SONES 1 – 2 sideholes                                                 | 100         | 5011539      | 5011570      | 5011554      |
| SONES 2 – 2 sideholes                                                 | 100         | 5011111      | 5011593      | 5011559      |
| <b>Angiodyn bypass catheter</b>                                       |             |              |              |              |
| 3DS Williams right                                                    | 100         | 5011596      | 5011597      | 5011598      |
| LCB                                                                   | 100         | 5012222      | 5011591      | 5011510      |
| RCB                                                                   | 100         | 5013333      | 5011560      | 5011519      |
| RADIAL                                                                | 100         | 5014444      | 5011543      | 5011544      |
| <b>Angiodyn angiography catheter kits</b>                             |             |              |              |              |
| JL4/JR4/PIG                                                           | 100/100/110 | 5011538      | 5011549      | 5011567      |
| JL4/JR4/PIG 145                                                       | 100/100/110 | 5014892      | 5011524      | 5011568      |

# Angiodyn® – Flow


# Flow



Left Coronary Artery





Right Coronary Artery



Ventricle


## Large inner lumen

- High quality angiographic results
- Superior flow rates
- Catheter stability during injection
- Easy injection



## Radio-opaque tip

- Optimal flow dynamics
- Smooth inner and outer surfaces



## High concentration of x-ray filler

- Excellent radiopacity



For product information  
please contact:

**B. Braun Melsungen AG**  
Vascular Systems

Sieversufer 8  
12359 Berlin

Phone +49 30 689897-0

Fax +49 30 689897-30

[www.bbraun.com](http://www.bbraun.com)

**B | BRAUN**  
SHARING EXPERTISE

**Aesculap AG & Co. KG**

Am Aesculap-Platz  
78532 Tuttlingen  
Germany

Phone +49 7461 95-0

Fax +49 7461 95-2600

[www.aesculap.de](http://www.aesculap.de)

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.