

Aesculap Orthopaedics

BiCONTACT[®] H

Hip endoprosthesis stem
with high offset



Bone Preservation. For Years to Come.



BiCONTACT® H

Hip endoprosthesis stem with high offset




BiCONTACT® S
Standard Offset

BiCONTACT® H
High offset

The BiCONTACT® S standard implant program has been expanded to include the BiCONTACT® H series.

BiCONTACT® H stems have a 6 mm higher offset and are designed with a CCD angle of 128°.

The S and H prosthesis stems are implanted using the same A and B osteoprofilers. For trial reduction of the H implants, special trial heads are available that simulate the reduced CCD angle and the increased offset situation.



Offset characteristic of the BiCONTACT® S and H hip endoprosthesis stems

Size	BiCONTACT® H stems cementless	BiCONTACT® H trial prosthesis heads	
10	NK110T	ø28 S	NG186
11	NK111T	ø28 M	NG187
12	NK112T	ø28 L	NG188
13	NK113T	ø28 XL	NG189
14	NK114T	ø28 XXL	NG190
15	NK115T		
16	NK116T	ø32 S	NG196
17	NK117T	ø32 M	NG197
18	NK118T	ø32 L	NG198
19	NK119T	ø32 XL	NG199
21	NK121T	ø32 XXL	NG200

ND746

x-ray templates

BiCONTACT® H cementless

B | BRAUN
SHARING EXPERTISE

Aesculap AG & Co. KG

Am Aesculap-Platz
78532 Tuttlingen

Phone +49 7461 95-0

Fax +49 7461 95-2600

www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.