























Aesculap e.motion® PS/Revision

Surgical Procedure – Operationsablauf

Procedure / Operationsschritt			Procedure / Operationsschritt			
1	Option 1: 	Option 2: 		7		
Tibia Resection / Tibiaresektion				Trial Reposition / Probereposition		
2				8		
Tibia Preparation / Tibiapräparation				Femur Box Preparation / Femurbox- präparation		
3				9		
Positioning of the tibia trial implants / Positionierung der Tibia-Probe- komponenten				Femur Stem Preparation / Femurschaft präparation		
4				10		
Measurement of the extension gap / flexion gab / Messung des Streckspaltes / des Beugespaltes				Final Trial Reposition / Endgültige Probereposition		
5				11		
Femoral Resection / Femurresektion				Assembling the final implants / Montage der endgültigen Implantate		
6				12		
Patella Preparation / Präparation der Patella				Final Implantation / Endgültige Implantation		