

# Monosyn® Quick

Surgeon's quick choice for excellent cosmetic results



## Monosyn® Quick

Recommended codes

Codes	Description
C0025014	Monosyn® Quick 3/0 70 cm DS30
C0025015	Monosyn® Quick 2/0 70 cm DS30
C0025020	Monosyn® Quick 6/0 45 cm DSMP11
C0025021	Monosyn® Quick 5/0 45 cm DSMP11
C0025022	Monosyn® Quick 4/0 45 cm DSMP11
C0025024	Monosyn® Quick 5/0 45 cm DSMP13
C0025029	Monosyn® Quick 4/0 45 cm DSMP16
C0025036	Monosyn® Quick 4/0 70 cm DSMP19
C0025043	Monosyn® Quick 3/0 70 cm DSMP24

Try it and send your comments to:  
[monofilaments@bbraun.com](mailto:monofilaments@bbraun.com)

C Box of 36 units



Plastic Surgery

B. Braun Surgical SA | Carretera de Terrassa, 121 | 08191 Rubí | Spain  
Phone +34 93 5 86 62 00 | Fax +34 93 6 99 73 03 | [www.bbraun.es](http://www.bbraun.es)

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany  
Phone +49 74 61 95-0 | Fax +49 74 61 95-26 00 | [www.aesculap.com](http://www.aesculap.com)

Aesculap – a B. Braun company

**B | BRAUN**  
SHARING EXPERTISE

## Monosyn® Quick

The short-term absorption profile  
monofilament suture of B. Braun

...for excellent cosmetic results



### Monosyn® Quick provides:

- Short-term absorption profile
- Quick mass absorption
- Atraumatic penetration in the skin
- Superior handling and good knot security
- Outsourcing knot run down
- Perfect range of cutting needles (DSMP needles):
  - Good penetration performance
  - Sharpness of the point
  - Advanced technology, providing a new level of smoothness during tissue penetration

## Monosyn® Quick, the surgeon's quick choice ...

### Main indications:

- Skin closure
- Blepharoplastias
- Lip cleft
- Periocular and oral surgery

### Product features

Structure	Monofilament
Chemical composition	Glyconate
Size	USP 1 - USP 6/0
Knot tensile strength retention	50 % 6 - 7 days post implantation 0 % 14 - 21 days post implantation
Mass Absorption	Degradation by hydrolysis after 56 days
Sterilization	Ethylene Oxide (EO)