

# Aesculap Single Use Surgical Instruments

## Pack per case

Procedure sets





## PACK PER CASE

### Pack per case – Procedure sets


„All inclusive“: – high quality single use instrumentation  
– all required consumables

At your finger tips! Ready for use when required:

- instantly available
- always complete
- always sterile (sterilised in acc. to validated sterilisation process)
- 100% reliable
- procedure related cost transparency
- no additional service cost (e.g. repair cost, etc..)
- always up-to-date



The extensive B Braun range offers you new opportunities and establish a new quality standard in the area of procedure sets.



### Pack per case – The special set

B Braun Aesculap procedure sets make a difference, independently if used for small procedures or for treatment of large wounds. Our "Pack per case" range of procedure sets include all the single use products which hospitals, clinics and practitioners rely upon every day. However the surgical instruments in the sets make the difference. They are high performance precision instruments, although designed for single use. No wonder, Aesculap instruments represent quality, security and reliability since generations.

### SUSI (Single Use Surgical Instruments) – The difference

SUSI Instrumentation is a completely new and unique range of single use surgical instruments. Utilising proven established design patterns and the use of light high performance polymers for the production provides you with high precision instruments with excellent mechanical properties. The material used allows simple, secure and environmentally friendly incineration utilising existing disposal systems in hospitals.

SUSI instruments are also available individually sterile packed in a dispenser box.

### The choice is yours – Establish a new standard

"Pack per case" procedure sets have proven their suitability on hospital wards, in intensive care units, A&E Departments, Outpatient Departments, Nursing homes, GP surgeries & practices and can also be used by ambulance services. Should you have specific needs and requirements (i.e custom made sets) please do not hesitate to contact us.



### 5061000 Suture Set 1



Contents	pieces
Defix kidney tray, white	1
Bowl, transparent, 120 ml	1
Gauze swabs, plumsized	5
SUSI forceps Adson, sharp, 120 mm	1
SUSI needle holder Mayo-Hegar, 150 mm	1
SUSI scissors, sharp/blunt, 145 mm	1
SUSI haemostat forceps Pean, straight, 145 mm	1
Order unit: PAK = Pack of 20 sets	

### 5061010 Surgical Dressing Set 1



Contents	pieces
Defix kidney tray, white	1
Bowl, transparent, 120 ml	1
Gauze compress, 7.5 x 7.5 cm, 8-fold (ply)	5
Gauze swabs, plumsized	5
SUSI forceps Adson, sharp, 120 mm	1
SUSI needle holder Mayo-Hegar, 150 mm	1
SUSI scissors, sharp/blunt, 145 mm	1
Order unit: PAK = Pack of 20 sets	



### 5061020 Surgical Dressing Set 2



Contents	pieces
Defix kidney tray, white	1
Bowl, transparent, 120 ml	1
Drape, 75 x 90 cm	1
Hole towel, 75 x 75 cm, adhesive area	1
Gauze compress, 7.5 x 7.5 cm, 12-fold (ply)	10
Gauze swabs, plumsize	5
SUSI forceps Adson, sharp, 120 mm	1
SUSI needle holder Mayo-Hegar, 150 mm	1
SUSI scissors, sharp/blunt, 145 mm	1
SUSI haemostat forceps Pean, straight, 145 mm	1
Order unit: PAK = Pack of 20 sets	

### 5061030 Surgical Dressing Set 3



Contents	pieces
Defix kidney tray, white	1
Bowl, transparent, 120 ml	1
Gauze compress, 7.5 x 7.5 cm, 8-fold (ply)	5
Gauze swabs, plumsize	5
SUSI tissue forceps, sharp, 145 mm	1
SUSI scissors, sharp/blunt, 145 mm	1
Order unit: PAK = Pack of 20 sets	



### 5061040 Suture Removal Set 1



Contents	pieces
Cardboard tray, 180 x 90 x 17 mm	1
Mini compress, 25 x 25 mm	17
SUSI forceps, sharp, 120 mm	1
SUSI scissors, sharp/sharp, 110 mm	1
Order unit: PAK = Pack of 20 sets	

### 5062000 Thorax Drainage Set, small



Contents	pieces
Plastic tray	1
Vasco surgical gloves, size 7.5	1
Drape, 45 x 75 cm	1
Hole towel, 75 x 110 cm, blue	1
Gauze compress with slit, 10 x 10 cm, 12-fold (ply)	2
Gauze compress, 10 x 10 cm, 12-fold (ply)	5
Disposable scalpel size 11	1
SUSI Metzenbaum dissecting scissors, curved, 180 mm	1
SUSI needle holder Mayo-Hegar, 150 mm	1
SUSI tubing clamp, 135 mm	1
Suture Miralene 2/0, 75 cm, 1 x DS39	2
Order unit: PAK = Pack of 20 sets	



## 5062010 Thorax Drainage Set, large



Contents	pieces
Plastic tray, 64 x 38 x 3 cm	1
Drape, 75 x 90 cm	1
Surgical gloves, size 7.5, powder free	1
Gauze compress, 10 x 10 cm, 12-fold (ply)	8
Slit compress, 10 x 10 cm, 12-fold (ply)	2
Bowl, round, 250 ml	1
SUSI dressing forceps, curved, 265 mm	1
Drape, 75 x 100 cm, self adhesive	2
Injection syringe, Luer, 10 ml	1
Sterican cannula size 1, G 20 x 1 1/2, Ø 0.90 mm, 40 mm, yellow	1
Sterican cannula, LL, G 18 x 1 1/2, 1.20 x 40 mm, pink	1
Sterican cannula, LL, G 20 x 2 3/4, Ø 0.90 mm, 70 mm, yellow	1

Contents	pieces
Sterican Kanüle size 12, G 22 x 1 1/4, Ø 0.70 mm, 30 mm, black	1
Disposable scalpel size 11	1
SUSI tissue forceps, sharp, 145 mm	1
SUSI Metzenbaum dissecting scissors, 180 mm	1
SUSI needle holder Mayo-Hegar, 150 mm	1
SUSI tubing clamp, 135 mm	1
Suture Synthofil green 0 (3.5), 75 cm, 1 x DS 30	1
Suture Dafilon blue 3/0 (2), 75 cm, 1 x DS 19	1
Askina fixation plaster, 15 x 20 cm	2
Order unit: PAK = Pack of 20 sets	
Developed in cooperation with: Dr. med. Thomas Kiefer, Klinikum Offenburg, Ebertplatz 12, 77654 Offenburg, Germany	



AESCULAP®

**B | BRAUN**  
SHARING EXPERTISE

**Aesculap AG & Co. KG**

Am Aesculap-Platz  
78532 Tuttlingen  
Germany

Phone +49 7461 95-0  
Fax +49 7461 95-2600

[www.aesculap.de](http://www.aesculap.de)

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.