



DEGENERATIVE
SPINAL DISORDERS

AESCULAP® Ennovate®

STERILE PACKED IMPLANT CONCEPT

SOLUTIONS BEYOND FUSION

AESCULAP® Ennovate®

STERILE PACKED IMPLANTS





CHAPTERS

A	GS1 CONFORMITY
B	IMPLANT OVERVIEW
C	SAFETY AND HANDLING ASPECTS
D	MODULES / SUPPLY CART
E	ORDERING INFORMATION ■ E1 Ennovate® Screws Poly Solid ■ E2 Ennovate® Screws Mono Solid ■ E3 Ennovate® Rods Open Ti/Connectors ■ E4 Ennovate® Screws Poly Canulated ■ E5 Ennovate® Screws Mono Canulated ■ E6 Ennovate® Rods MIS Ti/CoCr ■ E7 Ennovate® Screws Poly Fenestrated / Cannula ■ E8 Ennovate® Screws Mono Fenestrated / Cannula
F	Ennovate® BROCHURES

SOLUTIONS BEYOND FUSION

Designed around your needs –
We provide added value as your
spine partner fulfilling global
standards with GS1 conformity.

IMPROVED PATIENT SAFETY

Various regulatory bodies responsible for medical devices such as the US FDA with the UDI (Unique Device Identification¹) regulation and the upcoming harmonization of the regulations in the European Union² as well as recommendations from medical societies³ address the need of an improved transparency and traceability of medical devices.

Maintaining the highest possible level of safety is paramount to Aesculap. The importance of machine-readable identification of healthcare products by using barcodes is increasing continuously. B. Braun is marking the instruments and implants, both sterile packed and non-sterile packed, with a Data Matrix or GS1-128 barcode. Only very few products cannot be marked due to space restriction. GS1 standards are accepted by the UDI as standard coding. It is to be expected that the UDI will be applicable worldwide in the near future. GS1 standards are globally accepted and incorporate a GTIN (Global Trade Item Number).

With this implementation we are setting standards by following our core values of innovation and efficiency. Unambiguous identification is the key to improving patient safety and efficient healthcare business processes (logistics, documentation etc.).

Aesculap's standard for sterile packed items contains:

- Clear information about product including item number and other descriptive information.
- GS1-Datamatrix or GS1-128 barcode containing the GTIN, lot number and expiry date.
- Full traceability through patient ID stickers with GS1 data matrix and instructions for use in every implant component.

¹ Code of Federal Regulation Part 21 801 et seqq.

² European Commission, Revisions of Medical Device Directives, Regulation (EU) 2017/745, Regulation (EU) 2017/746.

³ Amman B, Bröcheler P, Carter A, Diedrich D, Fiedler S, Forster A, et al. Aufbereitung von Implantaten, die unsteril geliefert werden, für die Fachbereiche Orthopädie und Traumatologie: Teil 2. Zentralsterilisation. 2015;5:371-3. Retrieved from http://www.mhp-verlag.de/fileadmin/MHP/Zeitschriften/Zentralsterilisation/DGSV/FA_Qualitaet/deutsch/AKQ_d_ZT_5_2015.pdf.

TRACK AND TRACE

Main targets of the GS1 system are:

- numbering systems, e.g. Global Trade Item Numbers (GTIN) for product identification
- barcodes to supporting automatic data capturing
- standardized electronic messages to allowing efficient data sharing.

The GS1 Data Matrix coding according to the global standard enables you to track forward by monitoring product movements through your hospital pathway.

The GS1 Data Matrix according to the global standard also allows backtracing of the product from the patient to the manufacturer.

BENEFITS:

- Facilitating patient safety
- Improvement of implant management
- Full traceability of implants back to the manufacturer
- Fully compliant with legal regulations and recommendations
- Ensure that the right product reaches the right patient
- Straightforward documentation in electronic health records


DATA MATRIX AND GTIN



Application Identifier

GTIN	#040094098548
LOT	52350528
	2022-08-31

Product Detail Information



AESCULAP® Ennovate®

B | STERILE PACKED IMPLANTS – IMPLANT OVERVIEW

SOLUTIONS BEYOND FUSION

Designed around your needs –
We provide added value as
your spine partner offering
the complete range of implants
sterile packed.



RANGE OF IMPLANTS

Ennovate® provides the complete range of implants in sterile packed boxes with clear color coding elements.

IMPLANT TYPE	TYPE I	MATERIAL	DIAMETER AND LENGTH MM	UNIT	COLOR CODING
Pedicle Screws	Monoaxial	Titanium Alloy	Ø 4.5-10.5	1 pc	Diameter – Screw Head (various)
	Polyaxial				
	Solid		25-80 / 110		
	Canulated		25-80 / 110		Type – Screw Shaft (silver / yellow)
	Fenestrated		35-80 / 110		
Set Screws		Titanium Alloy		10 pcs (in box) 1 pc	Silver
Rods Open / MIS	Curved	Titanium Alloy	30-150	1 pc	Yellow
	Straight	CoCr*	30-500		Silver
R2R Connectors	Axial	Titanium Alloy	19/34	1 pc	Green
	Domino		7 / 11		Dark pink
	Offset		20/35/50		Yellow
Cross Connectors	Transconnector	Titanium Alloy	30-100	1 pc	Silver
	Multiaxial		26-100		



COLOR CODING /EASY HANDLING

All implant components have a clearly legible product description and Ennovate® marking on the front side for easy visual identification.

Additionally, the sterile packaging of the pedicle screws are marked with the same color coding as the screw head of the pedicle screw itself, indicating the diameter.

Further color coding concerning the type of screw, type of rod material or type of connectors ease the differentiation between implants.

All implants are designed around your needs for a fast identification.

BENEFIT:

- Time saving – easy identification of implant
- Improvement of your workflow
- Complete range of implants

AESCULAP® Ennovate®

C | STERILE PACKED IMPLANTS – SAFETY AND HANDLING ASPECTS

SOLUTIONS BEYOND FUSION

Designed around your needs –
We provide added value as your
spine partner providing safety in
sterile packaging and handling

SAFETY

Risks of postoperative infections through contamination or retaining chemical residues on implants, corrosion or surface oxidation and other impacts on the implants after hundreds of reprocessing cycles cause the need for sterile packed implants. Furthermore, the lack of full implant traceability also causes a lot of concerns. Sterile packed implants reduce the risk of contamination before, during and after the implantation¹⁻³.

All Ennovate® implant components are sterile packed in safe transport boxes with primary packaging and single sterile packed blister.

All implant components have an expiration time of 8 years.


BENEFIT:

- Ease of use
- Safe handling
- Reduced risk of contamination
- Long-term sterility

¹ Jesse E, Bible MD, MHSa, Kevin R, O'Neill MDa, Colin G, Crosby MDa, Jonathan G, Schoenecker MD, PhDa, Matthew J, McGirt MDb, Clinton J, Devin MDa. Implant contamination during spine surgery. The Spine Journal. 2013;13:637-40.

² Menekse G, Kuscu F, Suntur BM, Gezercan Y, Ates T, Ozsoy KM, Okten AI. Evaluation of the Time-dependent Contamination of Spinal Implants: Prospective Randomized Trial. Spine (Phila Pa 1976). 2015 Aug 15;40(16):1247-51.

³ Poelstra KA, Barekzi NA, Grainger DW, Gristina AG, Schuler TC. A novel spinal implant infection model in rabbits. Spine (Phila Pa 1976). 2000 Feb 15;25(4):406-10.



PEDICLE SCREW

Safe primary packaging with pedicle screw.



SET SCREWS

Set screws are available single sterile packed or in a box of 10 pcs each sterile packed.



CONNECTORS, RODS

Cross connectors, R2R connectors and rods are packed safely in slide packaging.

SOLUTIONS BEYOND FUSION

Designed around your needs –
Customization of implant
modules according to
specific need or standard
implant modules.

CUSTOMIZATION

You may choose between two options for sterile packed implants:

Option 1:

Choose your desired
implant module with our
recommended quantities.

Option 2:

Configure your own implant
module from our range of products
according to your specific need.

For the implant list and quantities please see chapter
"order information" with its implant modules.

AVAILABLE MODULES

● Ennovate® Open

- Ennovate® Screws POLY Solid
- Ennovate® Screws MONO Solid
- Ennovate® Rods OPEN Ti/Connectors

● Ennovate® MIS

- Ennovate® Screws POLY Canulated
- Ennovate® Screws MONO Canulated
- Ennovate® Rods MIS Ti

● ● Ennovate® Trauma

- Ennovate® Screws POLY Fenestrated/Cannula
- Ennovate® Screws MONO Fenestrated/Cannula

BENEFIT:

- Reduced inventory
- Less storage capacity required
- Cost-savings
- Use of modules for different instrument sets possible
- Ease of use
- Time-saving
- Adaption to specific requirements



Supply cart with sample configuration of 2 recommended modules

- Module Screws POLY Solid
- Rods OPEN Ti/Connectors



SUPPLY CART

To support you in the supply of sterile packed implants we may also offer you a mobile supply cart including especially fitted baskets.

The supply cart has a high standard of hygiene, by surfaces designed for wash disinfection, with a high level of flexibility and positioning. It may easily be moved from the sterile storage to the OR and can be used for various Ennovate® implant modules.

The supply cart offers space for the storage of two of the recommended implant modules or may of course be configured to your special needs.

For further information your local sales partner will support you with the configuration of your desired implant modules.

BENEFIT:

- Simple handling and moving
- Process optimization
- Configuration according to your needs
- Easy adaption for a variety of instrument sets
- Time-saving, ready to use

AESCULAP® Ennovate®


E | STERILE PACKED IMPLANTS – ORDERING INFORMATION CARTS/BASKETS

Product	Article No.	Description	Barcode
SUPPLY CARTS AND BASKETS			




**SCAN-
M-49500L**

ISO-Trolley 500 L,
456x655x1890 mm



**SCAN-
M-40413100**

Scancell Trolley Open
462x706x1205 mm




**SCAN-
M-1464881**

Transport Trolley
625x700x1520 mm






Product	Article No.	Description	Barcode
SUPPLY CARTS AND BASKETS			



QW232

Ennovate® Transport Box Implants Screws
(With Telescopic Slide-Out)
(Without Printing, Baskets and Partitions)



QW233

Ennovate® Transport Box Implants Rods
(With Telescopic Slide-Out)
(Without Printing, Baskets and Partitions)



BASKETS / PARTITIONS / SHELF



**SCAN-
M-2364251**

Basket 600x400x100 mm



**SCAN-
M-2364252**

Basket 600x400x200 mm



**SCAN-
M-2364253**

Partition 400x100 mm



**SCAN-
M-2364254**

Partition 600x100 mm



**SCAN-
M-2364255**

Partition 400x200 mm



**SCAN-
M-2364256**

Partition 600x200 mm



**SCAN-
M-2064266**

T-Piece for Partitions



AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION CARTS/BASKETS

Product	Article No.	Description	Barcode
BASKETS / PARTITIONS / SHELF			
	SCAN-M-2364259	Partition Lock – Optional	
	SCAN-M-2264400	Laminate Shelf 600x400x6 mm	
HANDLES			
	SCAN-M-1961880	Handle Grey with Stop 400 mm	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	
	SCAN-M-1961883	Handle Red with Stop 400 mm	
	SCAN-M-1961885	Handle Light Blue with Stop 400 mm	
ACCESSORIES			
	SCAN-M-1961870	Safety Stop for Baskets – Optional (Per Basket 2 Pieces Necessary)	
	SCAN-M-1869970	Roller Support Wheel – Optional (Per Basket 2 Pieces Necessary)	

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E1 | Ennovate® SCREWS POLY SOLID

Product	Article No.	Description	QTY	Barcode
● Ennovate® SCREWS POLY SOLID				
	BASKET 1	1 Ennovate® Screws POLY Solid Ø 4.5–5.5 mm I Set Screws Basket 600x400x100 mm	1	
	BASKET 2	2 Ennovate® Screws POLY Solid Ø 6.5–7.5 mm Basket 600x400x100 mm	1	
	BASKET 3	3 Ennovate® Screws POLY Solid Ø 7.5–10.5 mm Basket 600x400x100 mm	1	
	BASKET 4	4 Ennovate® Screws POLY Solid Ø 6.5–10.5 mm I 70–110 mm I Set Screws Basket 600x400x200 mm	1	
● BASKET 1 – 1 Ennovate® SCREWS POLY SOLID Ø 4.5–5.5 MM I Set Screws				

	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	40	40	Set Screw	Set Screw
	40	40	50	50
	35	35	50	50
	35	35	50	50
	30	30	45	45
	30	30	45	45
	25	25	45	45
5.5	25	25	40	40
	35	35	50	50
	35	35	50	50
	30	30	45	45
	30	30	45	45
	25	25	40	40
4.5	25	25	40	40

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E1 | Ennovate® SCREWS POLY SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS POLY SOLID Ø 4.5-5.5 MM I Set Screws				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961883	Handle Red with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	2	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 4.5 mm	SY110TS	Polyaxial Screw 4.5x30 mm Solid	4	
	SY111TS	Polyaxial Screw 4.5x30 mm Solid	4	
	SY112TS	Polyaxial Screw 4.5x35 mm Solid	4	
	SY113TS	Polyaxial Screw 4.5x40 mm Solid	4	
	SY114TS	Polyaxial Screw 4.5x45 mm Solid	4	
	SY115TS	Polyaxial Screw 4.5x50 mm Solid	4	
 Ø 5.5 mm	SY120TS	Polyaxial Screw 5.5x30 mm Solid	4	
	SY121TS	Polyaxial Screw 5.5x30 mm Solid	4	
	SY122TS	Polyaxial Screw 5.5x35 mm Solid	4	



Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS POLY SOLID Ø 4.5-5.5 MM I Set Screws				
 Ø 5.5 mm	SY123TS	Polyaxial Screw 5.5x40 mm Solid	6	
	SY124TS	Polyaxial Screw 5.5x45 mm Solid	6	
	SY125TS	Polyaxial Screw 5.5x50 mm Solid	6	
	SY001TS	Set Screw*	10	
● BASKET 2 – 2 Ennovate® SCREWS POLY SOLID Ø 6.5-7.5 MM				

7.5	35	35	45	45
	35	35	45	45
	30	30	45	45
	30	30	40	40
	25	25	40	40
	25	25	40	40
	45	45	60	60
	40	40	60	60
	40	40	55	55
	40	40	55	55
6.5	35	35	55	55
	35	35	50	50
	30	30	50	50
	30	30	50	50
	25	25	45	45
	25	25	45	45

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E1 | Ennovate® SCREWS POLY SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS POLY SOLID Ø 6.5–7.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961883	Handle Red with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 6.5 mm	SY130TS	Polyaxial Screw 6.5x25 mm Solid	4	
	SY131TS	Polyaxial Screw 6.5x30 mm Solid	4	
	SY132TS	Polyaxial Screw 6.5x35 mm Solid	4	
	SY133TS	Polyaxial Screw 6.5x40 mm Solid	6	
	SY134TS	Polyaxial Screw 6.5x45 mm Solid	6	
	SY135TS	Polyaxial Screw 6.5x50 mm Solid	6	
	SY136TS	Polyaxial Screw 6.5x55 mm Solid	6	
	SY137TS	Polyaxial Screw 6.5x60 mm Solid	4	



Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS POLY SOLID Ø 6.5-7.5 MM				



SY140TS	Polyaxial Screw 7.5x25 mm Solid	4	
SY141TS	Polyaxial Screw 7.5x30 mm Solid	4	
SY142TS	Polyaxial Screw 7.5x35 mm Solid	4	
SY143TS	Polyaxial Screw 7.5x40 mm Solid	6	
SY144TS	Polyaxial Screw 7.5x45 mm Solid	6	

● BASKET 3 – 3 Ennovate® SCREWS POLY SOLID Ø 7.5-10.5 MM				
--	--	--	--	--

10.5	45	45	60	60
	40	40	55	55
	35	35	50	50
	45	45	60	60
9.5	40	40	55	55
	35	35	50	50
	45	45	60	60
	40	40	55	55
8.5	40	40	50	50
	35	35	50	50
	30	30	45	45
	55	55	60	60
7.5	50	50	60	60
	50	50	55	55
	50	50	55	55

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E1 | Ennovate® SCREWS POLY SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS POLY SOLID Ø 7.5-10.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961883	Handle Red with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	2	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 7.5 mm	SY145TS	Polyaxial Screw 7.5x50 mm Solid	6	
	SY146TS	Polyaxial Screw 7.5x55 mm Solid	6	
	SY147TS	Polyaxial Screw 7.5x60 mm Solid	4	
 Ø 8.5 mm	SY151TS	Polyaxial Screw 8.5x30 mm Solid	2	
	SY152TS	Polyaxial Screw 8.5x35 mm Solid	2	
	SY153TS	Polyaxial Screw 8.5x40 mm Solid	4	
	SY154TS	Polyaxial Screw 8.5x45 mm Solid	4	
	SY155TS	Polyaxial Screw 8.5x50 mm Solid	4	
	SY156TS	Polyaxial Screw 8.5x55 mm Solid	2	
	SY157TS	Polyaxial Screw 8.5x60 mm Solid	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS POLY SOLID Ø 7.5-10.5 MM				
 Ø 9.5 mm	SY162TS	Polyaxial Screw 9.5x35 mm Solid	2	
	SY163TS	Polyaxial Screw 9.5x40 mm Solid	2	
	SY164TS	Polyaxial Screw 9.5x45 mm Solid	2	
	SY165TS	Polyaxial Screw 9.5x50 mm Solid	2	
	SY166TS	Polyaxial Screw 9.5x55 mm Solid	2	
	SY167TS	Polyaxial Screw 9.5x60 mm Solid	2	
 Ø 10.5 mm	SY172TS	Polyaxial Screw 10.5x35 mm Solid	2	
	SY173TS	Polyaxial Screw 10.5x40 mm Solid	2	
	SY174TS	Polyaxial Screw 10.5x45 mm Solid	2	
	SY175TS	Polyaxial Screw 10.5x50 mm Solid	2	
	SY176TS	Polyaxial Screw 10.5x55 mm Solid	2	
	SY177TS	Polyaxial Screw 10.5x60 mm Solid	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E1 | Ennovate® SCREWS POLY SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® SCREWS POLY SOLID Ø 6.5-10.5 MM 70-110 MM SET SCREWS				

Set Screws		Set Screw	Set Screw	Set Screw	Set Screw
		Set Screw	Set Screw	Set Screw	Set Screw
		Set Screw	Set Screw	Set Screw	Set Screw
		Set Screw	Set Screw	Set Screw	Set Screw
		Set Screw	Set Screw	Set Screw	Set Screw
		X 90	X 90	X 110	X 110
		X 80	X 80	X 100	X 100
	10.5	X 70	X 70	X 110	X 110
		X 80	X 80	X 100	X 100
	9.5	X 70	X 70	X 90	X 90
		X 90	X 90	X 110	X 110
		X 80	X 80	X 100	X 100
	8.5	X 70	X 70	X 110	X 110
		X 80	X 80	X 100	X 100
	7.5	X 70	X 70	X 90	X 90
	6.5	X 70	X 70	X 80	X 80

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



SCAN-
M-2364252

Basket 600x400x200 mm

1



SCAN-
M-1961883

Handle Red with Stop 400 mm

1



SCAN-
M-2364255

Partition 400x200 mm

2



SCAN-
M-2364256

Partition 600x200 mm

1





Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® SCREWS POLY SOLID Ø 6.5-10.5 MM 70-110 MM SET SCREWS				
 Ø 6.5 mm	SY138TS	Polyaxial Screw 6.5x70 mm Solid	2	
	SY139TS	Polyaxial Screw 6.5x80 mm Solid	2	
 Ø 7.5 mm	SY148TS	Polyaxial Screw 7.5x70 mm Solid	2	
	SY149TS	Polyaxial Screw 7.5x80 mm Solid	2	
	SY188TS	Polyaxial Screw 7.5x90 mm Solid	2	
	SY189TS	Polyaxial Screw 7.5x100 mm Solid	2	
	SY190TS	Polyaxial Screw 7.5x110 mm Solid	2	
 Ø 8.5 mm	SY158TS	Polyaxial Screw 8.5x70 mm Solid	2	
	SY159TS	Polyaxial Screw 8.5x80 mm Solid	2	
	SY191TS	Polyaxial Screw 8.5x90 mm Solid	2	
	SY192TS	Polyaxial Screw 8.5x100 mm Solid	2	
	SY193TS	Polyaxial Screw 8.5x110 mm Solid	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E1 | Ennovate® SCREWS POLY SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® SCREWS POLY SOLID Ø 6.5-10.5 MM 70-110 MM SET SCREWS				
 Ø 9.5 mm	SY168TS	Polyaxial Screw 9.5x70 mm Solid	2	
	SY169TS	Polyaxial Screw 9.5x80 mm Solid	2	
	SY194TS	Polyaxial Screw 9.5x90 mm Solid	2	
	SY195TS	Polyaxial Screw 9.5x100 mm Solid	2	
	SY196TS	Polyaxial Screw 9.5x110 mm Solid	2	
 Ø 10.5 mm	SY178TS	Polyaxial Screw 10.5x70 mm Solid	2	
	SY179TS	Polyaxial Screw 10.5x80 mm Solid	2	
	SY197TS	Polyaxial Screw 10.5x90 mm Solid	2	
	SY198TS	Polyaxial Screw 10.5x100 mm Solid	2	
	SY199TS	Polyaxial Screw 10.5x110 mm Solid	2	
	SY001TS	Set Screw*	20	

E2 | Ennovate® SCREWS MONO SOLID




Product	Article No.	Description	QTY	Barcode
● Ennovate® SCREWS MONO SOLID				
	BASKET 1	1 Ennovate® Screws MONO Solid Ø 4.5-5.5 mm I Set Screws Basket 600x400x100 mm	1	
	BASKET 2	2 Ennovate® Screws MONO Solid Ø 6.5-7.5 mm Basket 600x400x100 mm	1	
	BASKET 3	3 Ennovate® Screws MONO Solid Ø 7.5-10.5 mm Basket 600x400x100 mm	1	
	BASKET 4	4 Ennovate® Screws MONO Solid Ø 6.5-10.5 mm I 70-80 mm I Set Screws Basket 600x400x200 mm	1	
● BASKET 1 – 1 Ennovate® SCREWS MONO SOLID Ø 4.5-5.5 MM I Set Screws				

5.5	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	40	40	Set Screw	Set Screw
	40	40	50	50
	35	35	50	50
	35	35	50	50
	30	30	45	45
	30	30	45	45
	25	25	45	45
	25	25	40	40
4.5	35	35	50	50
	35	35	50	50
	30	30	45	45
	30	30	45	45
	25	25	40	40
	25	25	40	40
	25	25	40	40

SCREWS			
Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	
SET SCREW			
B 7.5	T 2.75	H 11.8	
BOX CM			
B 7.8	T 3.0		
H 14.5	X H 19.5		

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E2 | Ennovate® SCREWS MONO SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS MONO SOLID Ø 4.5-5.5 MM I Set Screws				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961883	Handle Red with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	2	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 4.5 mm	SY010TS	Monoaxial Screw 4.5x25 mm Solid	4	
	SY011TS	Monoaxial Screw 4.5x30 mm Solid	4	
	SY012TS	Monoaxial Screw 4.5x35 mm Solid	4	
	SY013TS	Monoaxial Screw 4.5x40 mm Solid	4	
	SY014TS	Monoaxial Screw 4.5x45 mm Solid	4	
	SY015TS	Monoaxial Screw 4.5x50 mm Solid	4	
 Ø 5.5 mm	SY020TS	Monoaxial Screw 5.5x25 mm Solid	4	
	SY021TS	Monoaxial Screw 5.5x30 mm Solid	4	
	SY022TS	Monoaxial Screw 5.5x35 mm Solid	4	



Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS MONO SOLID Ø 4.5-5.5 MM I Set Screws				



Ø 5.5 mm

SY023TS

Monoaxial Screw 5.5x40 mm Solid

6



SY024TS

Monoaxial Screw 5.5x45 mm Solid

6



SY025TS

Monoaxial Screw 5.5x50 mm Solid

6



SY001TS

Set Screw*

10



● BASKET 2 – 2 Ennovate® SCREWS MONO SOLID Ø 6.5-7.5 MM				
---	--	--	--	--

7.5	35	35	45	45
	35	35	45	45
	30	30	45	45
	30	30	40	40
	25	25	40	40
	25	25	40	40
	45	45	60	60
	40	40	60	60
	40	40	55	55
	40	40	55	55
6.5	35	35	55	55
	35	35	50	50
	30	30	50	50
	30	30	50	50
	25	25	45	45
	25	25	45	45
	25	25	45	45
	25	25	45	45

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E2 | Ennovate® SCREWS MONO SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS MONO SOLID Ø 6.5–7.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961883	Handle Red with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 6.5 mm	SY030TS	Monoaxial Screw 6.5x25 mm Solid	4	
	SY031TS	Monoaxial Screw 6.5x30 mm Solid	4	
	SY032TS	Monoaxial Screw 6.5x35 mm Solid	4	
	SY033TS	Monoaxial Screw 6.5x40 mm Solid	6	
	SY034TS	Monoaxial Screw 6.5x45 mm Solid	6	
	SY035TS	Monoaxial Screw 6.5x50 mm Solid	6	
	SY036TS	Monoaxial Screw 6.5x55 mm Solid	6	
	SY037TS	Monoaxial Screw 6.5x60 mm Solid	4	



Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS MONO SOLID Ø 6.5-7.5 MM				



SY040TS	Monoaxial Screw 7.5x25 mm Solid	4	
SY041TS	Monoaxial Screw 7.5x30 mm Solid	4	
SY042TS	Monoaxial Screw 7.5x35 mm Solid	4	
SY043TS	Monoaxial Screw 7.5x40 mm Solid	6	
SY044TS	Monoaxial Screw 7.5x45 mm Solid	6	

● BASKET 3 – 3 Ennovate® SCREWS MONO SOLID Ø 7.5-10.5 MM				
--	--	--	--	--

10.5	45	45	60	60
	40	40	55	55
	35	35	50	50
	45	45	60	60
9.5	40	40	55	55
	35	35	50	50
	45	45	60	60
	40	40	55	55
8.5	40	40	50	50
	35	35	50	50
	30	30	45	45
	55	55	60	60
7.5	50	50	60	60
	50	50	55	55
	50	50	55	55

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E2 | Ennovate® SCREWS MONO SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS MONO SOLID Ø 7.5-10.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961883	Handle Red with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	2	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 7.5 mm	SY045TS	Monoaxial Screw 7.5x50 mm Solid	6	
	SY046TS	Monoaxial Screw 7.5x55 mm Solid	6	
	SY047TS	Monoaxial Screw 7.5x60 mm Solid	4	
 Ø 8.5 mm	SY051TS	Monoaxial Screw 8.5x30 mm Solid	2	
	SY052TS	Monoaxial Screw 8.5x35 mm Solid	2	
	SY053TS	Monoaxial Screw 8.5x40 mm Solid	4	
	SY054TS	Monoaxial Screw 8.5x45 mm Solid	4	
	SY055TS	Monoaxial Screw 8.5x50 mm Solid	4	
	SY056TS	Monoaxial Screw 8.5x55 mm Solid	2	
	SY057TS	Monoaxial Screw 8.5x60 mm Solid	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS MONO SOLID Ø 7.5-10.5 MM				
 Ø 9.5 mm	SY062TS	Monoaxial Screw 9.5x35 mm Solid	2	
	SY063TS	Monoaxial Screw 9.5x40 mm Solid	2	
	SY064TS	Monoaxial Screw 9.5x45 mm Solid	2	
	SY065TS	Monoaxial Screw 9.5x50 mm Solid	2	
	SY066TS	Monoaxial Screw 9.5x55 mm Solid	2	
	SY067TS	Monoaxial Screw 9.5x60 mm Solid	2	
 Ø 10.5 mm	SY072TS	Monoaxial Screw 10.5x35 mm Solid	2	
	SY073TS	Monoaxial Screw 10.5x40 mm Solid	2	
	SY074TS	Monoaxial Screw 10.5x45 mm Solid	2	
	SY075TS	Monoaxial Screw 10.5x50 mm Solid	2	
	SY076TS	Monoaxial Screw 10.5x55 mm Solid	2	
	SY077TS	Monoaxial Screw 10.5x60 mm Solid	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E2 | Ennovate® SCREWS MONO SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® SCREWS MONO SOLID Ø 6.5-10.5 MM 70-80 MM SET SCREWS				

Set Screws				
	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	10.5	X 70	X 70	X 80
	9.5	X 70	X 70	X 80
	8.5	X 70	X 70	X 80
	7.5	X 70	X 70	X 80
	6.5	X 70	X 70	X 80

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



SCAN-
M-2364252

Basket 600x400x200 mm

1



SCAN-
M-1961883

Handle Red with Stop 400 mm

1



SCAN-
M-2364255

Partition 400x200 mm

2



SCAN-
M-2364256

Partition 600x200 mm

1





Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® SCREWS MONO SOLID Ø 6.5-10.5 MM 70-80 MM SET SCREWS				
 Ø 6.5 mm	SY038TS	Monoaxial Screw 6.5x70 mm Solid	2	
	SY039TS	Monoaxial Screw 6.5x80 mm Solid	2	
 Ø 7.5 mm	SY048TS	Monoaxial Screw 7.5x70 mm Solid	2	
	SY049TS	Monoaxial Screw 7.5x80 mm Solid	2	
 Ø 8.5 mm	SY058TS	Monoaxial Screw 8.5x70 mm Solid	2	
	SY059TS	Monoaxial Screw 8.5x80 mm Solid	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS


E2 | Ennovate® SCREWS MONO SOLID

Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® SCREWS MONO SOLID Ø 6.5-10.5 MM 70-80 MM SET SCREWS				
 Ø 9.5 mm	SY068TS	Monoaxial Screw 9.5x70 mm Solid	2	
	SY069TS	Monoaxial Screw 9.5x80 mm Solid	2	
 Ø 10.5 mm	SY078TS	Monoaxial Screw 10.5x70 mm Solid	2	
	SY079TS	Monoaxial Screw 10.5x80 mm Solid	2	
	SY001TS	Set Screw*	20	

E3 | Ennovate® RODS OPEN TI/CONNECTORS



Product	Article No.	Description	QTY	Barcode
● Ennovate® RODS OPEN TI / CONNECTORS				
	SHELF 1	1 Ennovate® Rods OPEN Ti Straight 300-500 mm Laminate Shelf 600x400x6 mm	1	
	BASKET 2	2 Ennovate® Rods OPEN Ti Straight 30-200 mm Basket 600x400x200 mm	1	
	BASKET 3	3 Ennovate® Rods OPEN Ti Curved 30-150 mm Connectors Offset Basket 600x400x200 mm	1	
	BASKET 4	4 Ennovate® Connectors Trans Multi Straight Domino Basket 600x400x100 mm	1	
● SHELF 1 – 1 Ennovate® RODS OPEN TI STRAIGHT 300-500 MM				



RODS TI CONNECT			
STR	CVD	Offset	Trans
Multi	Axial	Domino	
BOX CM			
B 7.8	T 3.0		
H 14.5	X H 19.5	XX H 25.5	XXX H 57

● Usable for OPEN module / ● Usable for MIS module / *alternatively SY002TS 2 x Set Screw 10-Pack

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E3 | Ennovate® RODS OPEN TI/CONNECTORS

Product	Article No.	Description	QTY	Barcode
● SHELF 1 – 1 Ennovate® RODS OPEN TI STRAIGHT 300-500 MM				
 	SCAN-M-2264400	Laminate Shelf 600x400x6 mm	1	
	SY421TS	Straight Rod 5.5x300 mm	2	
	SY422TS	Straight Rod 5.5x400 mm	2	
	SY423TS	Straight Rod 5.5x500 mm	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® RODS OPEN TI STRAIGHT 30-200 MM				

140-200	XX 2 x 140 XX 2 x 150	XX 2 x 160 XX 2 x 170	XX 2 x 180 XX 2 x 190 XX 2 x 200	
	60	60	X 130	X 130
	55	55	X 120	X 120
	50	50	X 110	X 110
	45	45	X 100	X 100
30-130	40	40	X 90	X 90
	35	35	80	80
	30	30	70	70

RODS TI CONNECT

STR	CVD	Offset	Trans
Multi	Axial	Domino	

BOX CM

B 7.8	T 3.0		
H 14.5	X H 19.5	XX H 25.5	XXX H 57



**SCAN-
M-2364252**

Basket 600x400x200 mm

1



**SCAN-
M-1961880**

Handle Grey with Stop 400 mm

1



**SCAN-
M-2364255**

Partition 400x200 mm

1



**SCAN-
M-2364256**

Partition 600x200 mm

1



ORDERING
INFORMATION

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E3 | Ennovate® RODS OPEN TI/CONNECTORS

Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® RODS OPEN TI STRAIGHT 30-200 MM				
	SY400TS	Straight Rod 5.5x30 mm	2	
	SY401TS	Straight Rod 5.5x35 mm	2	
	SY402TS	Straight Rod 5.5x40 mm	2	
	SY403TS	Straight Rod 5.5x45 mm	2	
	SY404TS	Straight Rod 5.5x50 mm	2	
	SY405TS	Straight Rod 5.5x55 mm	2	
	SY406TS	Straight Rod 5.5x60 mm	2	
	SY407TS	Straight Rod 5.5x70 mm	2	
	SY408TS	Straight Rod 5.5x80 mm	2	
	SY409TS	Straight Rod 5.5x90 mm	2	
	SY410TS	Straight Rod 5.5x100 mm	2	
	SY411TS	Straight Rod 5.5x110 mm	2	
	SY412TS	Straight Rod 5.5x120 mm	2	
	SY413TS	Straight Rod 5.5x130 mm	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® RODS OPEN TI STRAIGHT 30-200 MM				
	SY414TS	Straight Rod 5.5x 140 mm	2	 SY414TS
	SY415TS	Straight Rod 5.5x 150 mm	2	 SY415TS
	SY416TS	Straight Rod 5.5x 160 mm	2	 SY416TS
	SY417TS	Straight Rod 5.5x 170 mm	2	 SY417TS
	SY418TS	Straight Rod 5.5x 180 mm	2	 SY418TS
	SY419TS	Straight Rod 5.5x 190 mm	2	 SY419TS
	SY420TS	Straight Rod 5.5x 200 mm	2	 SY420TS

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E3 | Ennovate® RODS OPEN TI/CONNECTORS

Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® RODS OPEN TI CURVED 30-150 MM CONNECTORS OFFSET				

Offset	Closed 50	Closed 50	Open 50	Open 50
	Closed 50	Closed 50	Open 50	Open 50
	Closed 35	Closed 35	Open 35	Open 35
	Closed 35	Closed 35	Open 35	Open 35
	Closed 20	Closed 20	Open 20	Open 20
	Closed 20	Closed 20	Open 20	Open 20
			X 130	
	80	80	X 130	
	70	70	X 120	
	60	60	X 120	
30-150	55	55	X 110	
	50	50	X 110	
	45	45	X 100	
	40	40	X 100	
	35	35	X 90	
	30	30	X 90	
			XX 12 x 140	XX 12 x 150

RODS TI CONNECT

STR	CVD	Offset	Trans
Multi	Axial	Domino	

BOX CM

B 7.8	T 3.0		
H 14.5	X H 19.5	XX H 25.5	XXX H 57



SCAN-
M-2364252

Basket 600x400x200 mm

1



SCAN-
M-1961880

Handle Grey with Stop 400 mm

1



SCAN-
M-2364255

Partition 400x200 mm

1



SCAN-
M-2364256

Partition 600x200 mm

1





Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® RODS OPEN TI CURVED 30-150 MM CONNECTORS OFFSET				
	SY430TS	Curved Rod 5.5x30 mm	2	
	SY431TS	Curved Rod 5.5x35 mm	2	
	SY432TS	Curved Rod 5.5x40 mm	2	
	SY433TS	Curved Rod 5.5x45 mm	2	
	SY434TS	Curved Rod 5.5x50 mm	2	
	SY435TS	Curved Rod 5.5x55 mm	2	
	SY436TS	Curved Rod 5.5x60 mm	2	
	SY437TS	Curved Rod 5.5x70 mm	2	
	SY438TS	Curved Rod 5.5x80 mm	2	
	SY439TS	Curved Rod 5.5x90 mm	2	
	SY440TS	Curved Rod 5.5x100 mm	2	
	SY441TS	Curved Rod 5.5x110 mm	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E3 | Ennovate® RODS OPEN TI/CONNECTORS

Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® RODS OPEN TI CURVED 30-150 MM I CONNECTORS OFFSET				
	SY442TS	Curved Rod 5.5x120 mm	2	
	SY443TS	Curved Rod 5.5x130 mm	2	
	SY444TS	Curved Rod 5.5x140 mm	2	
	SY445TS	Curved Rod 5.5x150 mm	2	
	SY711TS	Closed Lateral Offset Connector 20 mm	4	
	SY712TS	Closed Lateral Offset Connector 35 mm	4	
	SY713TS	Closed Lateral Offset Connector 50 mm	4	
	SY716TS	Open Lateral Offset Connector 20 mm	4	
	SY717TS	Open Lateral Offset Connector 35 mm	4	
	SY718TS	Open Lateral Offset Connector 50 mm	4	



Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® CONNECTORS TRANS I MULTI I STRAIGHT I DOMINO				

Axial	Domino	Closed 11	Closed 11	Cl/Op 11	Cl/Op 11
		Closed 11	Closed 11	Cl/Op 11	Cl/Op 11
		Closed 7	Closed 7	Cl/Op 7	Cl/Op 7
		Closed 7	Closed 7	Cl/Op 7	Cl/Op 7
		Str 19	Str 19	Str 34	Str 34
		Str 19	Str 19	Str 34	Str 34
Multi		31-37	2x80-100		
		28-34	34-40	45-58	58-80
		28-34	34-40	39-45	58-80
		26-30	31-37	39-45	45-58
		26-30	2x80-100		
		38	38	60-80	60-80
Trans		36	36	50-60	50-60
		34	34	45-50	45-50
		20	32	42-45	42-45
		20	32	40	40

RODS TI CONNECT

STR	CVD	Offset	Trans
Multi	Axial	Domino	

BOX CM

B 7.8	T 3.0		
H 14.5	X H 19.5	XX H 25.5	XXX H 57

	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961880	Handle Grey with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	
	SCAN-M-2064266	T-Piece for Partitions	1	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E3 | Ennovate® RODS OPEN TI/CONNECTORS

Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® CONNECTORS TRANS I MULTI I STRAIGHT I DOMINO				
	SX942TS	Transconnector 30 mm	2	
	SX943TS	Transconnector 32 mm	2	
	SX944TS	Transconnector 34 mm	2	
	SX945TS	Transconnector 36 mm	2	
	SX946TS	Transconnector 38 mm	2	
	SX947TS	Transconnector 40 mm	2	
	SX948TS	Transconnector 42-45 mm	2	
	SX949TS	Transconnector 45-50 mm	2	
	SX950TS	Transconnector 50-60 mm	2	
	SX951TS	Transconnector 60-80 mm	2	
	SX952TS	Transconnector 80-100 mm	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® CONNECTORS TRANS I MULTI I STRAIGHT I DOMINO				
	SX932TS	Multiaxial Cross-Connector, 26-30 mm	2	
	SX933TS	Multiaxial Cross-Connector, 28-34 mm	2	
	SX934TS	Multiaxial Cross-Connector, 31-37 mm	2	
	SX935TS	Multiaxial Cross-Connector, 34-40 mm	2	
	SX936TS	Multiaxial Cross-Connector, 39-45 mm	2	
	SX937TS	Multiaxial Cross-Connector, 45-58 mm	2	
	SX938TS	Multiaxial Cross-Connector, 58-80 mm	2	
	SX939TS	Multiaxial Cross-Connector, 80-100 mm	2	
	SY701TS	Axial Connector 19 mm	4	
	SY702TS	Axial Connector 34 mm	4	
	SY703TS	Closed Domino Connector 7 mm	4	
	SY704TS	Closed Domino Connector 11 mm	4	
	SY705TS	Closed/Open Domino Connector 7 mm	4	
	SY706TS	Closed/Open Domino Connector 11 mm	4	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E4 | Ennovate® SCREWS POLY CANULATED

Product	Article No.	Description	QTY	Barcode
● ● Ennovate® SCREWS POLY CANULATED				
	BASKET 1	1 Ennovate® Screws POLY Canulated Ø 4.5-5.5 mm I Set Screws Basket 600x400x100 mm	1	
	BASKET 2	2 Ennovate® Screws POLY Canulated Ø 6.5-7.5 mm Basket 600x400x100 mm	1	
	BASKET 3	3 Ennovate® Screws POLY Canulated Ø 7.5-10.5 mm Basket 600x400x100 mm	1	
	BASKET 4	4 Ennovate® Screws POLY Canulated Ø 6.5-10.5 mm I 70-110 mm I Set Screws Basket 600x400x200 mm	1	
● ● BASKET 1 – 1 Ennovate® SCREWS POLY CANULATED Ø 4.5-5.5 MM I Set Screws				

	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	40	40	Set Screw	Set Screw
	40	40	50	50
	35	35	50	50
	35	35	50	50
	30	30	45	45
	30	30	45	45
	25	25	45	45
5.5	25	25	40	40
	35	35	50	50
	35	35	50	50
	30	30	45	45
	30	30	45	45
	25	25	40	40
4.5	25	25	40	40

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS POLY CANULATED Ø 4.5-5.5 MM				
 Ø 4.5 mm	 SCAN-M-2364251	Basket 600x400x100 mm	1	
	 SCAN-M-1961885	Handle Light Blue with Stop 400 mm	1	
	 SCAN-M-2364253	Partition 400x100 mm	2	
	 SCAN-M-2364254	Partition 600x100 mm	1	
	SY610TS	Polyaxial Screw 4.5x25 mm Canulated	4	
	SY611TS	Polyaxial Screw 4.5x30 mm Canulated	4	
	SY612TS	Polyaxial Screw 4.5x35 mm Canulated	4	
	SY613TS	Polyaxial Screw 4.5x40 mm Canulated	4	
	SY614TS	Polyaxial Screw 4.5x45 mm Canulated	4	
	SY615TS	Polyaxial Screw 4.5x50 mm Canulated	4	
 Ø 5.5 mm	SY620TS	Polyaxial Screw 5.5x25 mm Canulated	4	
	SY621TS	Polyaxial Screw 5.5x30 mm Canulated	4	
	SY622TS	Polyaxial Screw 5.5x35 mm Canulated	4	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E4 | Ennovate® SCREWS POLY CANULATED

Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS POLY CANULATED Ø 4.5–5.5 MM I Set Screws				
 Ø 5.5 mm	SY623TS	Polyaxial Screw 5.5x40 mm Canulated	6	
	SY624TS	Polyaxial Screw 5.5x45 mm Canulated	6	
	SY625TS	Polyaxial Screw 5.5x50 mm Canulated	6	
	SY001TS	Set Screw*	10	
● BASKET 2 – 2 Ennovate® SCREWS POLY CANULATED Ø 6.5–7.5 MM				

7.5	35	35	45	45
	35	35	45	45
	30	30	45	45
	30	30	40	40
	25	25	40	40
	25	25	40	40
	45	45	60	60
	40	40	60	60
	40	40	55	55
	40	40	55	55
6.5	35	35	55	55
	35	35	50	50
	30	30	50	50
	30	30	50	50
	25	25	45	45
	25	25	45	45

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS POLY CANULATED Ø 6.5-7.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961885	Handle Light Blue with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 6.5 mm	SY630TS	Polyaxial Screw 6.5x25 mm Canulated	4	
	SY631TS	Polyaxial Screw 6.5x30 mm Canulated	4	
	SY632TS	Polyaxial Screw 6.5x35 mm Canulated	4	
	SY633TS	Polyaxial Screw 6.5x40 mm Canulated	6	
	SY634TS	Polyaxial Screw 6.5x45 mm Canulated	6	
	SY635TS	Polyaxial Screw 6.5x50 mm Canulated	6	
	SY636TS	Polyaxial Screw 6.5x55 mm Canulated	6	
	SY637TS	Polyaxial Screw 6.5x60 mm Canulated	4	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E4 | Ennovate® SCREWS POLY CANULATED

Product	Article No.	Description	QTY	Barcode
● ● BASKET 2 – 2 Ennovate® SCREWS POLY CANULATED Ø 6.5-7.5 MM				
 Ø 7.5 mm	SY640TS	Polyaxial Screw 7.5x25 mm Canulated	4	
	SY641TS	Polyaxial Screw 7.5x30 mm Canulated	4	
	SY642TS	Polyaxial Screw 7.5x35 mm Canulated	4	
	SY643TS	Polyaxial Screw 7.5x40 mm Canulated	6	
	SY644TS	Polyaxial Screw 7.5x45 mm Canulated	6	
● ● BASKET 3 – 3 Ennovate® SCREWS POLY CANULATED Ø 7.5-10.5 MM				

10.5	45	45	60	60
	40	40	55	55
	35	35	50	50
	45	45	60	60
9.5	40	40	55	55
	35	35	50	50
	45	45	60	60
	40	40	55	55
8.5	40	40	50	50
	35	35	50	50
	30	30	45	45
	55	55	60	60
7.5	50	50	60	60
	50	50	55	55
	50	50	55	55

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS POLY CANULATED Ø 7.5-10.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961885	Handle Light Blue with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	2	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 7.5 mm	SY645TS	Polyaxial Screw 7.5x50 mm Canulated	6	
	SY646TS	Polyaxial Screw 7.5x55 mm Canulated	6	
	SY647TS	Polyaxial Screw 7.5x60 mm Canulated	4	
 Ø 8.5 mm	SY651TS	Polyaxial Screw 8.5x30 mm Canulated	2	
	SY652TS	Polyaxial Screw 8.5x35 mm Canulated	2	
	SY653TS	Polyaxial Screw 8.5x40 mm Canulated	4	
	SY654TS	Polyaxial Screw 8.5x45 mm Canulated	4	
	SY655TS	Polyaxial Screw 8.5x50 mm Canulated	4	
	SY656TS	Polyaxial Screw 8.5x55 mm Canulated	2	
	SY657TS	Polyaxial Screw 8.5x60 mm Canulated	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E4 | Ennovate® SCREWS POLY CANULATED

Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS POLY CANULATED Ø 7.5-10.5 MM				
 Ø 9.5 mm	SY662TS	Polyaxial Screw 9.5x35 mm Canulated	2	
	SY663TS	Polyaxial Screw 9.5x40 mm Canulated	2	
	SY664TS	Polyaxial Screw 9.5x45 mm Canulated	2	
	SY665TS	Polyaxial Screw 9.5x50 mm Canulated	2	
	SY666TS	Polyaxial Screw 9.5x55 mm Canulated	2	
	SY667TS	Polyaxial Screw 9.5x60 mm Canulated	2	
 Ø 10.5 mm	SY672TS	Polyaxial Screw 10.5x35 mm Canulated	2	
	SY673TS	Polyaxial Screw 10.5x40 mm Canulated	2	
	SY674TS	Polyaxial Screw 10.5x45 mm Canulated	2	
	SY675TS	Polyaxial Screw 10.5x50 mm Canulated	2	
	SY676TS	Polyaxial Screw 10.5x55 mm Canulated	2	
	SY677TS	Polyaxial Screw 10.5x60 mm Canulated	2	



Product	Article No.	Description	QTY	Barcode
● ● BASKET 4 – 4 Ennovate® SCREWS POLY CANULATED Ø 6.5-10.5 MM 70-110 MM SET SCREWS				

Set Screws	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	X 90	X 90	X 110	X 110
	X 80	X 80	X 100	X 100
	X 70	X 70	X 110	X 110
	X 80	X 80	X 100	X 100
	X 70	X 70	X 90	X 90
10.5	X 90	X 90	X 110	X 110
	X 80	X 80	X 100	X 100
	X 70	X 70	X 110	X 110
9.5	X 80	X 80	X 100	X 100
	X 70	X 70	X 90	X 90
	X 90	X 90	X 110	X 110
8.5	X 80	X 80	X 100	X 100
	X 70	X 70	X 110	X 110
	X 80	X 80	X 100	X 100
7.5	X 70	X 70	X 90	X 90
	X 70	X 70	X 80	X 80
6.5	X 70	X 70	X 80	X 80

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



SCAN-
M-2364252

Basket 600x400x200 mm

1



SCAN-M-2364252



SCAN-
M-1961885

Handle Light Blue with Stop 400 mm

1



SCAN-M-1961885



SCAN-
M-2364255

Partition 400x200 mm

2



SCAN-M-2364255

SCAN-
M-2364256

Partition 600x200 mm

1



SCAN-M-2364256

ORDERING
INFORMATION

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E4 | Ennovate® SCREWS POLY CANULATED

Product	Article No.	Description	QTY	Barcode
● ● BASKET 4 – 4 Ennovate® SCREWS POLY CANULATED Ø 6.5-10.5 MM 70-110 MM SET SCREWS				
 Ø 6.5 mm	SY638TS	Polyaxial Screw 6.5x70 mm Canulated	2	
	SY639TS	Polyaxial Screw 6.5x80 mm Canulated	2	
 Ø 7.5 mm	SY648TS	Polyaxial Screw 7.5x70 mm Canulated	2	
	SY649TS	Polyaxial Screw 7.5x80 mm Canulated	2	
	SY688TS	Polyaxial Screw 7.5x90 mm Canulated	2	
	SY689TS	Polyaxial Screw 7.5x100 mm Canulated	2	
	SY690TS	Polyaxial Screw 7.5x110 mm Canulated	2	
 Ø 8.5 mm	SY658TS	Polyaxial Screw 8.5x70 mm Canulated	2	
	SY659TS	Polyaxial Screw 8.5x80 mm Canulated	2	
	SY691TS	Polyaxial Screw 8.5x90 mm Canulated	2	
	SY692TS	Polyaxial Screw 8.5x100 mm Canulated	2	
	SY693TS	Polyaxial Screw 8.5x110 mm Canulated	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® SCREWS POLY CANULATED Ø 6.5-10.5 MM 70-110 MM SET SCREWS				
 Ø 9.5 mm	SY668TS	Polyaxial Screw 9.5x70 mm Canulated	2	
	SY669TS	Polyaxial Screw 9.5x80 mm Canulated	2	
	SY694TS	Polyaxial Screw 9.5x90 mm Canulated	2	
	SY695TS	Polyaxial Screw 9.5x100 mm Canulated	2	
	SY696TS	Polyaxial Screw 9.5x110 mm Canulated	2	
 Ø 10.5 mm	SY678TS	Polyaxial Screw 10.5x70 mm Canulated	2	
	SY679TS	Polyaxial Screw 10.5x80 mm Canulated	2	
	SY697TS	Polyaxial Screw 10.5x90 mm Canulated	2	
	SY698TS	Polyaxial Screw 10.5x100 mm Canulated	2	
	SY699TS	Polyaxial Screw 10.5x110 mm Canulated	2	
	SY001TS	Set Screw*	20	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E5 | Ennovate® SCREWS MONO CANULATED

Product	Article No.	Description	QTY	Barcode
● ● Ennovate® SCREWS MONO CANULATED				
	BASKET 1	1 Ennovate® Screws MONO Canulated Ø 4.5–5.5 mm I Set Screws Basket 600x400x100 mm	1	
	BASKET 2	2 Ennovate® Screws MONO Canulated Ø 6.5–7.5 mm Basket 600x400x100 mm	1	
	BASKET 3	3 Ennovate® Screws MONO Canulated Ø 7.5–10.5 mm Basket 600x400x100 mm	1	
	BASKET 4	4 Ennovate® Screws MONO Canulated Ø 6.5–10.5 mm I 70–80 mm I Set Screws Basket 600x400x200 mm	1	
● ● BASKET 1 – 1 Ennovate® SCREWS MONO CANULATED Ø 4.5–5.5 MM I Set Screws				

	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
	40	40	Set Screw	Set Screw
	40	40	50	50
	35	35	50	50
	35	35	50	50
	30	30	45	45
	30	30	45	45
	25	25	45	45
5.5	25	25	40	40
	35	35	50	50
	35	35	50	50
	30	30	45	45
	30	30	45	45
	25	25	40	40
4.5	25	25	40	40

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS MONO CANULATED Ø 4.5-5.5 MM I Set Screws				
 Ø 4.5 mm	 SCAN-M-2364251	Basket 600x400x100 mm	1	
	 SCAN-M-1961885	Handle Light Blue with Stop 400 mm	1	
	 SCAN-M-2364253	Partition 400x100 mm	2	
		Partition 600x100 mm	1	
	SY510TS	Monoaxial Screw 4.5x25 mm Canulated	4	
	SY511TS	Monoaxial Screw 4.5x30 mm Canulated	4	
	SY512TS	Monoaxial Screw 4.5x35 mm Canulated	4	
	SY513TS	Monoaxial Screw 4.5x40 mm Canulated	4	
	SY514TS	Monoaxial Screw 4.5x45 mm Canulated	4	
	SY515TS	Monoaxial Screw 4.5x50 mm Canulated	4	
 Ø 5.5 mm	SY520TS	Monoaxial Screw 5.5x25 mm Canulated	4	
	SY521TS	Monoaxial Screw 5.5x30 mm Canulated	4	
	SY522TS	Monoaxial Screw 5.5x35 mm Canulated	4	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E5 | Ennovate® SCREWS MONO CANULATED

Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS MONO CANULATED Ø 4.5-5.5 MM I Set Screws				
 Ø 5.5 mm	SY523TS	Monoaxial Screw 5.5x40 mm Canulated	6	
	SY524TS	Monoaxial Screw 5.5x45 mm Canulated	6	
	SY525TS	Monoaxial Screw 5.5x50 mm Canulated	6	
	SY001TS	Set Screw*	10	
● BASKET 2 – 2 Ennovate® SCREWS MONO CANULATED Ø 6.5-7.5 MM				

7.5	35	35	45	45
	35	35	45	45
	30	30	45	45
	30	30	40	40
	25	25	40	40
	25	25	40	40
	45	45	60	60
	40	40	60	60
	40	40	55	55
	40	40	55	55
6.5	35	35	55	55
	35	35	50	50
	30	30	50	50
	30	30	50	50
	25	25	45	45
	25	25	45	45

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



Product	Article No.	Description	QTY	Barcode
● ● BASKET 2 – 2 Ennovate® SCREWS MONO CANULATED Ø 6.5-7.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961885	Handle Light Blue with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 6.5 mm	SY530TS	Monoaxial Screw 6.5x25 mm Canulated	4	
	SY531TS	Monoaxial Screw 6.5x30 mm Canulated	4	
	SY532TS	Monoaxial Screw 6.5x35 mm Canulated	4	
	SY533TS	Monoaxial Screw 6.5x40 mm Canulated	6	
	SY534TS	Monoaxial Screw 6.5x45 mm Canulated	6	
	SY535TS	Monoaxial Screw 6.5x50 mm Canulated	6	
	SY536TS	Monoaxial Screw 6.5x55 mm Canulated	6	
	SY537TS	Monoaxial Screw 6.5x60 mm Canulated	4	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E5 | Ennovate® SCREWS MONO CANULATED

Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS MONO CANULATED Ø 6.5-7.5 MM				
 Ø 7.5 mm	SY540TS	Monoaxial Screw 7.5x25 mm Canulated	4	
	SY541TS	Monoaxial Screw 7.5x30 mm Canulated	4	
	SY542TS	Monoaxial Screw 7.5x35 mm Canulated	4	
	SY543TS	Monoaxial Screw 7.5x40 mm Canulated	6	
	SY544TS	Monoaxial Screw 7.5x45 mm Canulated	6	
● BASKET 3 – 3 Ennovate® SCREWS MONO CANULATED Ø 7.5-10.5 MM				

10.5	45	45	60	60
	40	40	55	55
	35	35	50	50
	45	45	60	60
9.5	40	40	55	55
	35	35	50	50
	45	45	60	60
	40	40	55	55
8.5	40	40	50	50
	35	35	50	50
	30	30	45	45
	55	55	60	60
7.5	50	50	60	60
	50	50	55	55
	50	50	55	55

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS MONO CANULATED Ø 7.5-10.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961885	Handle Light Blue with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	2	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 7.5 mm	SY545TS	Monoaxial Screw 7.5x50 mm Canulated	6	
	SY546TS	Monoaxial Screw 7.5x55 mm Canulated	6	
	SY547TS	Monoaxial Screw 7.5x60 mm Canulated	4	
 Ø 8.5 mm	SY551TS	Monoaxial Screw 8.5x30 mm Canulated	2	
	SY552TS	Monoaxial Screw 8.5x35 mm Canulated	2	
	SY553TS	Monoaxial Screw 8.5x40 mm Canulated	4	
	SY554TS	Monoaxial Screw 8.5x45 mm Canulated	4	
	SY555TS	Monoaxial Screw 8.5x50 mm Canulated	4	
	SY556TS	Monoaxial Screw 8.5x55 mm Canulated	2	
	SY557TS	Monoaxial Screw 8.5x60 mm Canulated	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E5 | Ennovate® SCREWS MONO CANULATED

Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS MONO CANULATED Ø 7.5-10.5 MM				
 Ø 9.5 mm	SY562TS	Monoaxial Screw 9.5x35 mm Canulated	2	
	SY563TS	Monoaxial Screw 9.5x40 mm Canulated	2	
	SY564TS	Monoaxial Screw 9.5x45 mm Canulated	2	
	SY565TS	Monoaxial Screw 9.5x50 mm Canulated	2	
	SY566TS	Monoaxial Screw 9.5x55 mm Canulated	2	
	SY567TS	Monoaxial Screw 9.5x60 mm Canulated	2	
 Ø 10.5 mm	SY572TS	Monoaxial Screw 10.5x35 mm Canulated	2	
	SY573TS	Monoaxial Screw 10.5x40 mm Canulated	2	
	SY574TS	Monoaxial Screw 10.5x45 mm Canulated	2	
	SY575TS	Monoaxial Screw 10.5x50 mm Canulated	2	
	SY576TS	Monoaxial Screw 10.5x55 mm Canulated	2	
	SY577TS	Monoaxial Screw 10.5x60 mm Canulated	2	



Product	Article No.	Description	QTY	Barcode
● ● BASKET 4 – 4 Ennovate® SCREWS MONO CANULATED Ø 6.5–10.5 MM 70–80 MM SET SCREWS				

Set Screws	Set Screw			
	Set Screw			
	Set Screw			
	Set Screw			
	Set Screw			
	Set Screw			
	Set Screw			
	Set Screw			
	Set Screw			
	Set Screw			
10.5	X 70	X 70	X 80	X 80
9.5	X 70	X 70	X 80	X 80
8.5	X 70	X 70	X 80	X 80
7.5	X 70	X 70	X 80	X 80
6.5	X 70	X 70	X 80	X 80

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



**SCAN-
M-2364252**

Basket 600x400x200 mm

1



**SCAN-
M-1961885**

Handle Light Blue with Stop 400 mm

1



**SCAN-
M-2364255**

Partition 400x200 mm

2



**SCAN-
M-2364256**

Partition 600x200 mm

1



ORDERING
INFORMATION

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E5 | Ennovate® SCREWS MONO CANULATED

Product	Article No.	Description	QTY	Barcode
● ● BASKET 4 – 4 Ennovate® SCREWS MONO CANULATED Ø 6.5–10.5 MM 70–80 MM SET SCREWS				
 Ø 6.5 mm	SY538TS	Monoaxial Screw 6.5x70 mm Canulated	2	
	SY539TS	Monoaxial Screw 6.5x80 mm Canulated	2	
 Ø 7.5 mm	SY548TS	Monoaxial Screw 7.5x70 mm Canulated	2	
	SY549TS	Monoaxial Screw 7.5x80 mm Canulated	2	
 Ø 8.5 mm	SY558TS	Monoaxial Screw 8.5x70 mm Canulated	2	
	SY559TS	Monoaxial Screw 8.5x80 mm Canulated	2	




Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® SCREWS MONO CANULATED Ø 6.5–10.5 MM 70–80 MM SET SCREWS				
 Ø 9.5 mm	SY568TS	Monoaxial Screw 9.5x70 mm Canulated	2	
	SY569TS	Monoaxial Screw 9.5x80 mm Canulated	2	
 Ø 10.5 mm	SY578TS	Monoaxial Screw 10.5x70 mm Canulated	2	
	SY579TS	Monoaxial Screw 10.5x80 mm Canulated	2	
	SY001TS	Set Screw*	20	

AESCULAP® Ennovate®


E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E6 | Ennovate® RODS MIS TI/COCR

Product	Article No.	Description	QTY	Barcode
● ● Ennovate® RODS MIS RODS TI/COCR				
	SHELF 1	1 Ennovate® Rods MIS Ti I CoCr Straight 300-500 mm* Laminate Shelf 600x400x6 mm	1	
	BASKET 2	2 Ennovate® Rods MIS Ti Straight 30-120 mm I Curved 30-150 mm Basket 600x400x200 mm	1	
	BASKET 3	3 Ennovate® Rods MIS CoCr Straight 30-120 mm I Curved 30-150 mm* Basket 600x400x200 mm	1	
	BASKET 4	4 Ennovate® Rods MIS Ti I CoCr Straight 130-200 mm* Basket 600x400x100 mm	1	
● ● SHELF 1 – 1 Ennovate® RODS MIS TI I COCR STRAIGHT 300-500 MM*				



RODS TI CONNECT			
STR	CVD	Offset	Trans
Multi	Axial	Domino	
BOX CM			
B 7.8	T 3.0		
H 14.5	X H 19.5	XX H 25.5	XXX H 57



Product	Article No.	Description	QTY	Barcode
● ● SHELF 1 – 1 Ennovate® RODS MIS TI I COCR STRAIGHT 300-500 MM*				
	SCAN-M-2264400	Laminate Shelf 600x400x6 mm	1	
	SY421TS	Straight Rod 5.5x300 mm	2	
	SY422TS	Straight Rod 5.5x400 mm	2	
	SY423TS	Straight Rod 5.5x500 mm	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E6 | Ennovate® RODS MIS TI/COCR

Product	Article No.	Description	QTY	Barcode
● ● BASKET 2 – 2 Ennovate® RODS MIS TI STRAIGHT 30-120 MM I CURVED 30-150 MM				

XX I 2 x 150			
XX I 2 x 140			
XX I 2 x 130			
X 120	X 120	X 120	X 120
X 110	X 110	X 110	X 110
X 100	X 100	X 100	X 100
X 90	X 90	X 90	X 90
X 80	X 80	X 80	X 80
70	70	70	70
60	60	60	60
55	55	55	55
50	50	50	50
45	45	45	45
40	40	40	40
35	35	35	35
30	30	30	30

Ti Straight 30-120

Ti Curved 30-150

RODS TI CONNECT

STR	CVD	Offset	Trans
Multi	Axial	Domino	

BOX CM

B 7.8	T 3.0		
H 14.5	X H 19.5	XX H 25.5	XXX H 57



SCAN-
M-2364252

Basket 600x400x200 mm

1



SCAN-
M-1961880

Handle Grey with Stop 400 mm

1



SCAN-
M-2364255

Partition 400x200 mm

1



SCAN-
M-2364256

Partition 600x200 mm

1



SCAN-
M-2064266

T-Piece for Partitions

1





Product	Article No.	Description	QTY	Barcode
● ● BASKET 2 – 2 Ennovate® RODS MIS TI STRAIGHT 30-120 MM CURVED 30-150 MM				
	SY900TS	MIS Straight Rod 5.5x30 mm	2	
	SY901TS	MIS Straight Rod 5.5x35 mm	2	
	SY902TS	MIS Straight Rod 5.5x40 mm	2	
	SY903TS	MIS Straight Rod 5.5x45 mm	2	
	SY904TS	MIS Straight Rod 5.5x50 mm	2	
	SY905TS	MIS Straight Rod 5.5x55 mm	2	
	SY906TS	MIS Straight Rod 5.5x60 mm	2	
	SY907TS	MIS Straight Rod 5.5x70 mm	2	
	SY908TS	MIS Straight Rod 5.5x80 mm	2	
	SY909TS	MIS Straight Rod 5.5x90 mm	2	
	SY910TS	MIS Straight Rod 5.5x100 mm	2	
	SY911TS	MIS Straight Rod 5.5x110 mm	2	
	SY912TS	MIS Straight Rod 5.5x120 mm	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E6 | Ennovate® RODS MIS TI/COCR

Product	Article No.	Description	QTY	Barcode
● ● BASKET 2 – 2 Ennovate® RODS MIS TI STRAIGHT 30-120 MM I CURVED 30-150 MM				
	SY930TS	MIS Curved Rod 5.5x30 mm	2	
	SY931TS	MIS Curved Rod 5.5x35 mm	2	
	SY932TS	MIS Curved Rod 5.5x40 mm	2	
	SY933TS	MIS Curved Rod 5.5x45 mm	2	
	SY934TS	MIS Curved Rod 5.5x50 mm	2	
	SY935TS	MIS Curved Rod 5.5x55 mm	2	
	SY936TS	MIS Curved Rod 5.5x60 mm	2	
	SY937TS	MIS Curved Rod 5.5x70 mm	2	
	SY938TS	MIS Curved Rod 5.5x80 mm	2	
	SY939TS	MIS Curved Rod 5.5x90 mm	2	
	SY940TS	MIS Curved Rod 5.5x100 mm	2	
	SY941TS	MIS Curved Rod 5.5x110 mm	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® RODS MIS TI STRAIGHT 30-120 MM I CURVED 30-150 MM				



SY942TS

MIS Curved Rod 5.5x120 mm

2



SY943TS

MIS Curved Rod 5.5x130 mm

2



SY944TS

MIS Curved Rod 5.5x140 mm

2



SY945TS

MIS Curved Rod 5.5x150 mm

2



● BASKET 3 – 3 Ennovate® RODS MIS COCR STRAIGHT 30-120 MM I CURVED 30-150 MM*				
---	--	--	--	--



SCAN-M-2364252

Basket 600x400x200 mm

1



SCAN-M-1961880

Handle Grey with Stop 400 mm

1



SCAN-M-2364255

Partition 400x200 mm

1



SCAN-M-2364256

Partition 600x200 mm

1



SCAN-M-2064266

T-Piece for Partitions

1




NOTE: Rods CoCr only available in 2018

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E6 | Ennovate® RODS MIS TI/COCR

Product	Article No.	Description	QTY	Barcode
● ● BASKET 4 – 4 Ennovate® RODS MIS TI I COCR STRAIGHT 130-200 MM*				



RODS TI CONNECT

STR	CVD	Offset	Trans
Multi	Axial	Domino	

BOX CM

B 7.8	T 3.0		
H 14.5	X H 19.5	XX H 25.5	XXX H 57

	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961880	Handle Grey with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	



Product	Article No.	Description	QTY	Barcode
● ● BASKET 4 – 4 Ennovate® RODS MIS TI I COCR STRAIGHT 130-200 MM*				
	SY913TS	MIS Straight Rod 5.5x 130 mm	2	
	SY914TS	MIS Straight Rod 5.5x 140 mm	2	
	SY915TS	MIS Straight Rod 5.5x 150 mm	2	
	SY916TS	MIS Straight Rod 5.5x 160 mm	2	
	SY917TS	MIS Straight Rod 5.5x 170 mm	2	
	SY918TS	MIS Straight Rod 5.5x 180 mm	2	
	SY919TS	MIS Straight Rod 5.5x 190 mm	2	
	SY920TS	MIS Straight Rod 5.5x 200 mm	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E7 | Ennovate® SCREWS POLY FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● ● Ennovate® SCREWS POLY FENESTRATED / CANNULAS				
	BASKET 1	1 Ennovate® Screws POLY Fenestrated Ø 5.5-6.5 mm Basket 600x400x100 mm	1	
	BASKET 2	2 Ennovate® Screws POLY Fenestrated Ø 7.5-8.5 mm Basket 600x400x100 mm	1	
	BASKET 3	3 Ennovate® Screws POLY Fenestrated Ø 9.5-10.5 mm I Set Screws Basket 600x400x100 mm	1	
	BASKET 4	4 Ennovate® Cannula Short I Long Basket 600x400x200 mm	1	
● ● BASKET 1 – 1 Ennovate® SCREWS POLY FENESTRATED Ø 5.5-6.5 MM				

6.5	50	50	60	X 2 x 70 X 2 x 80
	50	50	60	
	45	45	60	
	45	45	60	
	45	45	55	
	40	40	55	
	40	40	55	
	40	40	55	
	35	35	50	
	35	35	50	
5.5	45	45		
	40	40	50	
	40	40	50	
	40	40	50	
	35	35	45	
	35	35	45	

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS POLY FENESTRATED Ø 5.5-6.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 5.5 mm	SY822TS	Polyaxial Screw 5.5x35 mm Fenestrated	4	
	SY823TS	Polyaxial Screw 5.5x40 mm Fenestrated	6	
	SY824TS	Polyaxial Screw 5.5x45 mm Fenestrated	6	
	SY825TS	Polyaxial Screw 5.5x50 mm Fenestrated	6	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E7 | Ennovate® SCREWS POLY FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS POLY FENESTRATED Ø 5.5-6.5 MM				
 Ø 6.5 mm	SY832TS	Polyaxial Screw 6.5x35 mm Fenestrated	4	
	SY833TS	Polyaxial Screw 6.5x40 mm Fenestrated	6	
	SY834TS	Polyaxial Screw 6.5x45 mm Fenestrated	6	
	SY835TS	Polyaxial Screw 6.5x50 mm Fenestrated	6	
	SY836TS	Polyaxial Screw 6.5x55 mm Fenestrated	4	
	SY837TS	Polyaxial Screw 6.5x60 mm Fenestrated	4	
	SY838TS	Polyaxial Screw 6.5x70 mm Fenestrated	2	
	SY839TS	Polyaxial Screw 6.5x80 mm Fenestrated	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS POLY FENESTRATED Ø 7.5-8.5 MM				

8.5	50	50	60	X 2x70 X 2x80
	45	45	60	
	45	45	55	
	40	40	55	
	40	40	50	
	35	35	50	
7.5	50	50	60	X 2x70 X 2x80
	45	45	60	
	45	45	60	
	45	45	60	
	40	40	55	
	40	40	55	
	40	40	55	
	35	35	50	
	35	35	50	50

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E7 | Ennovate® SCREWS POLY FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● ● BASKET 2 – 2 Ennovate® SCREWS POLY FENESTRATED Ø 7.5-8.5 MM				
 Ø 7.5 mm	SY842TS	Polyaxial Screw 7.5x35 mm Fenestrated	4	
	SY843TS	Polyaxial Screw 7.5x40 mm Fenestrated	6	
	SY844TS	Polyaxial Screw 7.5x45 mm Fenestrated	6	
	SY845TS	Polyaxial Screw 7.5x50 mm Fenestrated	6	
	SY846TS	Polyaxial Screw 7.5x55 mm Fenestrated	4	
	SY847TS	Polyaxial Screw 7.5x60 mm Fenestrated	4	
	SY848TS	Polyaxial Screw 7.5x70 mm Fenestrated	2	
	SY849TS	Polyaxial Screw 7.5x80 mm Fenestrated	2	



Ø 8.5 mm

Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS POLY FENESTRATED Ø 7.5-8.5 MM				
	SY852TS	Polyaxial Screw 8.5x35 mm Fenestrated	2	
	SY853TS	Polyaxial Screw 8.5x40 mm Fenestrated	4	
	SY854TS	Polyaxial Screw 8.5x45 mm Fenestrated	4	
	SY855TS	Polyaxial Screw 8.5x50 mm Fenestrated	4	
	SY856TS	Polyaxial Screw 8.5x55 mm Fenestrated	2	
	SY857TS	Polyaxial Screw 8.5x60 mm Fenestrated	2	
	SY858TS	Polyaxial Screw 8.5x70 mm Fenestrated	2	
	SY859TS	Polyaxial Screw 8.5x80 mm Fenestrated	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E7 | Ennovate® SCREWS POLY FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● ● BASKET 3 – 3 Ennovate® SCREWS POLY FENESTRATED Ø 9.5-10.5 MM I SET SCREWS				

Set Screws	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
				Set Screw
	60	60	X 2 x 70 X 2 x 80	Set Screw
	55	55		Set Screw
	50	50		Set Screw
	45	45		Set Screw
	40	40		Set Screw
	35	35		Set Screw
				Set Screw
	60	60	X 2 x 70 X 2 x 80	Set Screw
	55	55		Set Screw
	50	50		Set Screw
	45	45		Set Screw
	40	40		Set Screw
	35	35		Set Screw

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	2	
	SCAN-M-2364254	Partition 600x100 mm	1	



Ø 9.5 mm

Product	Article No.	Description	QTY	Barcode
● ● BASKET 3 – 3 Ennovate® SCREWS POLY FENESTRATED Ø 9.5-10.5 MM SET SCREWS				
	SY862TS	Polyaxial Screw 9.5x35 mm Fenestrated	2	
	SY863TS	Polyaxial Screw 9.5x40 mm Fenestrated	2	
	SY864TS	Polyaxial Screw 9.5x45 mm Fenestrated	2	
	SY865TS	Polyaxial Screw 9.5x50 mm Fenestrated	2	
	SY866TS	Polyaxial Screw 9.5x55 mm Fenestrated	2	
	SY867TS	Polyaxial Screw 9.5x60 mm Fenestrated	2	
	SY868TS	Polyaxial Screw 9.5x70 mm Fenestrated	2	
	SY869TS	Polyaxial Screw 9.5x80 mm Fenestrated	2	

AESCULAP® Ennovate®


E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E7 | Ennovate® SCREWS POLY FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● BASKET 3 – 3 Ennovate® SCREWS POLY FENESTRATED Ø 9.5-10.5 MM I SET SCREWS				
 Ø 10.5 mm	SY872TS	Polyaxial Screw 10.5x35 mm Fenestrated	2	
	SY873TS	Polyaxial Screw 10.5x40 mm Fenestrated	2	
	SY874TS	Polyaxial Screw 10.5x45 mm Fenestrated	2	
	SY875TS	Polyaxial Screw 10.5x50 mm Fenestrated	2	
	SY876TS	Polyaxial Screw 10.5x55 mm Fenestrated	2	
	SY877TS	Polyaxial Screw 10.5x60 mm Fenestrated	2	
	SY878TS	Polyaxial Screw 10.5x70 mm Fenestrated	2	
	SY879TS	Polyaxial Screw 10.5x80 mm Fenestrated	2	
	SY001TS	Set Screw*	22	



Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® CANNULA SHORT LONG				



BOX CM		
B 7.7	T 1.6	H 28.5
X B 11.5	X T 2.6	X H 38

	SCAN-M-2364252	Basket 600x400x200 mm	1	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	1	
	SCAN-M-2364255	Partition 400x200 mm	1	
	SCAN-M-2364256	Partition 600x200 mm	1	
	SCAN-M-2064266	T-Piece for Partitions	1	
	SR136SU	Ennovate® Open Injection Cannula	12	
	SR138SU	Ennovate® MIS Injection Cannula	12	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E8 | Ennovate® SCREWS MONO FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● ● Ennovate® SCREWS MONO FENESTRATED / CANNULAS				
	BASKET 1	1 Ennovate® Screws MONO Fenestrated Ø 5.5-6.5 mm Basket 600x400x100 mm	1	
	BASKET 2	2 Ennovate® Screws MONO Fenestrated Ø 7.5-8.5 mm Basket 600x400x100 mm	1	
	BASKET 3	3 Ennovate® Screws MONO Fenestrated Ø 9.5-10.5 mm I Set Screws Basket 600x400x100 mm	1	
	BASKET 4	4 Ennovate® Cannula Short I Long Basket 600x400x200 mm	1	
● ● BASKET 1 – 1 Ennovate® SCREWS MONO FENESTRATED Ø 5.5-6.5 MM				

6.5	50	50	60	X 2 x 70 X 2 x 80
	50	50	60	
	45	45	60	
	45	45	60	
	45	45	55	
	40	40	55	
	40	40	55	
	40	40	55	
	35	35	50	
	35	35	50	
5.5	45	45		
	40	40	50	
	40	40	50	
	40	40	50	
	35	35	45	
	35	35	45	

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5



Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS MONO FENESTRATED Ø 5.5-6.5 MM				
	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	
 Ø 5.5 mm	SY722TS	Monoaxial Screw 5.5x35 mm Fenestrated	4	
	SY723TS	Monoaxial Screw 5.5x40 mm Fenestrated	6	
	SY724TS	Monoaxial Screw 5.5x45 mm Fenestrated	6	
	SY725TS	Monoaxial Screw 5.5x50 mm Fenestrated	6	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E8 | Ennovate® SCREWS MONO FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● BASKET 1 – 1 Ennovate® SCREWS MONO FENESTRATED Ø 5.5-6.5 MM				
 Ø 6.5 mm	SY732TS	Monoaxial Screw 6.5x35 mm Fenestrated	4	
	SY733TS	Monoaxial Screw 6.5x40 mm Fenestrated	6	
	SY734TS	Monoaxial Screw 6.5x45 mm Fenestrated	6	
	SY735TS	Monoaxial Screw 6.5x50 mm Fenestrated	6	
	SY736TS	Monoaxial Screw 6.5x55 mm Fenestrated	4	
	SY737TS	Monoaxial Screw 6.5x60 mm Fenestrated	4	
	SY738TS	Monoaxial Screw 6.5x70 mm Fenestrated	2	
	SY739TS	Monoaxial Screw 6.5x80 mm Fenestrated	2	



Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS MONO FENESTRATED Ø 7.5-8.5 MM				

8.5	50	50	60	X 2x70 X 2x80
	45	45	60	
	45	45	55	
	40	40	55	
	40	40	50	
	35	35	50	
7.5	50	50	60	X 2x70 X 2x80
	45	45	60	
	45	45	60	
	45	45	60	
	40	40	55	
	40	40	55	
	40	40	55	
	35	35	50	
	35	35	50	50
	35	35	50	50

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	1	
	SCAN-M-2364254	Partition 600x100 mm	1	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E8 | Ennovate® SCREWS MONO FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS MONO FENESTRATED Ø 7.5-8.5 MM				
 Ø 7.5 mm	SY742TS	Monoaxial Screw 7.5x35 mm Fenestrated	4	
	SY743TS	Monoaxial Screw 7.5x40 mm Fenestrated	6	
	SY744TS	Monoaxial Screw 7.5x45 mm Fenestrated	6	
	SY745TS	Monoaxial Screw 7.5x50 mm Fenestrated	6	
	SY746TS	Monoaxial Screw 7.5x55 mm Fenestrated	4	
	SY747TS	Monoaxial Screw 7.5x60 mm Fenestrated	4	
	SY748TS	Monoaxial Screw 7.5x70 mm Fenestrated	2	
	SY749TS	Monoaxial Screw 7.5x80 mm Fenestrated	2	



Ø 8.5 mm

Product	Article No.	Description	QTY	Barcode
● BASKET 2 – 2 Ennovate® SCREWS MONO FENESTRATED Ø 7.5-8.5 MM				
	SY752TS	Monoaxial Screw 8.5x35 mm Fenestrated	2	
	SY753TS	Monoaxial Screw 8.5x40 mm Fenestrated	4	
	SY754TS	Monoaxial Screw 8.5x45 mm Fenestrated	4	
	SY755TS	Monoaxial Screw 8.5x50 mm Fenestrated	4	
	SY756TS	Monoaxial Screw 8.5x55 mm Fenestrated	2	
	SY757TS	Monoaxial Screw 8.5x60 mm Fenestrated	2	
	SY758TS	Monoaxial Screw 8.5x70 mm Fenestrated	2	
	SY759TS	Monoaxial Screw 8.5x80 mm Fenestrated	2	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS

E8 | Ennovate® SCREWS MONO FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● ● BASKET 3 – 3 Ennovate® SCREWS MONO FENESTRATED Ø 9.5–10.5 MM I SET SCREWS				

Set Screws	Set Screw	Set Screw	Set Screw	Set Screw
	Set Screw	Set Screw	Set Screw	Set Screw
				Set Screw
	60	60	X 2 x 70 X 2 x 80	Set Screw
	55	55		Set Screw
	50	50		Set Screw
	45	45		Set Screw
	40	40		Set Screw
	35	35		Set Screw
				Set Screw
	60	60	X 2 x 70 X 2 x 80	Set Screw
	55	55		Set Screw
	50	50		Set Screw
	45	45		Set Screw
	40	40		Set Screw
	35	35		Set Screw

SCREWS

Ø 4.5	Ø 5.5	Ø 6.5	Ø 7.5
Ø 8.5	Ø 9.5	Ø 10.5	

SET SCREW

B 7.5	T 2.75	H 11.8
-------	--------	--------

BOX CM

B 7.8	T 3.0
H 14.5	X H 19.5

	SCAN-M-2364251	Basket 600x400x100 mm	1	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	1	
	SCAN-M-2364253	Partition 400x100 mm	2	
	SCAN-M-2364254	Partition 600x100 mm	1	



Ø 9.5 mm

Product	Article No.	Description	QTY	Barcode
● ● BASKET 3 – 3 Ennovate® SCREWS MONO FENESTRATED Ø 9.5-10.5 MM SET SCREWS				
	SY762TS	Monoaxial Screw 9.5x35 mm Fenestrated	2	
	SY763TS	Monoaxial Screw 9.5x40 mm Fenestrated	2	
	SY764TS	Monoaxial Screw 9.5x45 mm Fenestrated	2	
	SY765TS	Monoaxial Screw 9.5x50 mm Fenestrated	2	
	SY766TS	Monoaxial Screw 9.5x55 mm Fenestrated	2	
	SY767TS	Monoaxial Screw 9.5x60 mm Fenestrated	2	
	SY768TS	Monoaxial Screw 9.5x70 mm Fenestrated	2	
	SY769TS	Monoaxial Screw 9.5x80 mm Fenestrated	2	


AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E8 | Ennovate® SCREWS MONO FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● ● BASKET 3 – 3 Ennovate® SCREWS MONO FENESTRATED Ø 9.5–10.5 MM SET SCREWS				
 Ø 10.5 mm	SY772TS	Monoaxial Screw 10.5x35 mm Fenestrated	2	
	SY773TS	Monoaxial Screw 10.5x40 mm Fenestrated	2	
	SY774TS	Monoaxial Screw 10.5x45 mm Fenestrated	2	
	SY775TS	Monoaxial Screw 10.5x50 mm Fenestrated	2	
	SY776TS	Monoaxial Screw 10.5x55 mm Fenestrated	2	
	SY777TS	Monoaxial Screw 10.5x60 mm Fenestrated	2	
	SY778TS	Monoaxial Screw 10.5x70 mm Fenestrated	2	
	SY779TS	Monoaxial Screw 10.5x80 mm Fenestrated	2	
	SY001TS	Set Screw*	22	



Product	Article No.	Description	QTY	Barcode
● BASKET 4 – 4 Ennovate® CANNULA SHORT LONG				



BOX CM		
B 7.7	T 1.6	H 28.5
X B 11.5	X T 2.6	X H 38

	SCAN-M-2364252	Basket 600x400x200 mm	1	
	SCAN-M-1961882	Handle Yellow with Stop 400 mm	1	
	SCAN-M-2364255	Partition 400x200 mm	1	
	SCAN-M-2364256	Partition 600x200 mm	1	
	SCAN-M-2064266	T-Piece for Partitions	1	

AESCULAP® Ennovate®

E | STERILE PACKED IMPLANTS – ORDERING INFORMATION IMPLANTS
E8 | Ennovate® SCREWS MONO FENESTRATED / CANNULA

Product	Article No.	Description	QTY	Barcode
● ● BASKET 4 – 4 Ennovate® CANNULA SHORT I LONG				
	SR136SU	Ennovate® Open Injection Cannula	12	
	SR138SU	Ennovate® MIS Injection Cannula	12	



DESIGNED AROUND YOUR NEEDS

We provide added value as your spine partner providing solutions for the improvement in your hospital pathway

Ennovate® BROCHURES

**Ennovate® OPEN –
Surgical Technique**
048102



**Ennovate® MIS –
Surgical Technique**
000702



**Ennovate® ELSA –
European Loan Service
of Aesculap 048502**



**Ennovate® Instrument
and Tray Concept**
048402



**Ennovate® Listing
Sterile Packed Implants**
001002



**Ennovate® Sterile
Packed Implant Concept**
048602



**Ennovate® Non-Sterile
Packed Implant Concept**
000602



Ennovate® Health Economics
001602



AESCULAP® – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "Aesculap" and the product trademark "Ennovate" are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.